
(IT) MANIPOLI, CAVI BIPOLARI E MONOPOLARI PER ELETTROCHIRURGIA
(EN) HANDLES, MONOPOLAR AND BIPOLAR CABLE FOR ELECTROSURGERY
(ES) MANGOS, CABLES BIPOLAR Y MONOPOLAR PARA ELECTROCIRUGIA
(FR) MANCHES, CABLE BIPOLAR ET MONOPOLAR POUR ELECTROCHIRURGIE

	- PENCIL -	- PENCIL S -	- CONNECTION -
REF	00206.00	00205.00	00402.00/00402.01/00402.02 00411.00/00412.00/00413.00

ISTRUZIONI PER L'USO.....2
INSTRUCTION FOR USE.....7
INSTRUCCIONES DE USO.....12
MODE D'EMPLOI.....18

LED SpA Via Selciatella 40 04011 Aprilia (LT) ITALY

CE
0051

ITALIANO

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ENGLISH

USE: Through the connection to a high-frequency generator (ESU) compatible, cutting and coagulation during an electrosurgical procedure.

WARNINGS:

- Check the integrity of the package and the product, if damaged do not use and return it to the dealer or directly to the LED SpA.
- The handpieces accept electrodes with a nominal diameter **2.4mm (3/32")**. The spindle system allows you to securely attach the electrodes. Check that the electrode or attached accessory insulation is intact and when inserted there is no part of electrical discoveries between the product and the electrode or accessory.
- Make sure that the product connections with the ESU are stable and the metal part of the connectors are fully inserted.
- Mount electrodes or accessories connected to the cables, such as cleaning, should only be made with the handpiece or cable disconnected from the ESU.
- Make some performance tests before use. Discard the defective products.
- When removing the ESU products do not pull on the cord itself but to act on the intended slot.
- When not in use store the products in containers and / or suitable supports.
- The product must be used with ESU and / or accessories comply with current regulations.
- To use the ESU and other necessary accessories refer to the instructions of the individual elements.
- Do not let the wires come into direct contact with the patient's skin, a capacitive effect could cause burns and / or burns.
- The products must be connected and used by qualified personnel.
- The products are not supplied sterile, **STERILISE** before use (see **RECONDITIONING**).

PREPARATION / USE:

Handpieces: to fix the electrode enlarge (loosen) slightly present in the spindle tip, fully insert the electrode, tighten acting (screwed) on the spindle.

Handpieces with buttons: press the yellow button to activate the function of cutting (CUT) and press the blue button to activate the coagulation (COAG) function. First use a test correct activation by selecting the minimum on ESU

output power, taking away the handpiece by people and objects, provide alternately CUT and COAG checking against the corresponding reporting unit.

Handpieces without buttons: must be controlled by foot pedal connected with ESU.

Cables: Connect to a compatible accessory cable, connect the cable to the unit. To brew you have to act on the foot controller ESU.

RECONDITIONING (ISO17664)

Instruments must not be used without previously being cleaned, disinfected and sterilized. Effective cleaning and disinfection is a precondition for an effective sterilization of instruments. Remove the electrodes before cleaning. Observe that only sufficient instruments and product-specified valid methods for cleaning, disinfection and sterilization are used and that the valid parameters maintained with each cycle.

Observe the valid legal regulations of the country and the regulations for hygiene of the hospital/clinic.

Do not sterilize electrodes and handle bunched together! The amount for reapplication of the individual components is, with a sterilization time of 20 minutes and with a sterilization temperature of 134°C, as follows:

Handle / cables: up **100 times**, # Electrodes: up to **30 times**.

Every exceeding further application falls within the responsibility of the user. Surgical handles and electrodes can be cleaned in the instrument washing machine.

It is imperative that the data of the disinfection agent manufacturer be adhered to. Mechanical pre-cleaning of the knife electrodes and the needle electrodes can be made the same as with conventional surgical electrodes.

CAUTION Do not sterilize in hot air!

1. Soaking

Soaking must be made immediately after application (within a maximum of 2 hours). For this only use aldehyde-free disinfection agents (otherwise fixing of blood contamination) with tested effectiveness (e.g. DGHM or FDA Approval, respectively CE certification mark) which are suitable for instrument disinfection and which are compatible with the instruments.

For preliminary removal of dirt residue it is possible, if necessary, to use plastic fleece material or a plastic cleaning brush. Do not prewash the HF cable and the handle in an ultrasonic bath.

IMPORTANT: The applied disinfection agent used for soaking is only for personnel protection and cannot, after implemented cleaning, replace the disinfection sequence carried out later.

2. Disassembly / Dismantling

Dismountable instruments are to be disassembled for further cleaning according to the operating instructions.

3. Cleaning / Disinfecting

A cleaning machine method should be used for cleaning and disinfection (cleaning and disinfection equipment). Cleaning by hand is not recommended due to the distinctly lower effectiveness. LED recommends the use of a neutral to light alkaline cleaning agent, respectively cleaning and disinfection agent, which does not contain any critical ingredients (depending of concentration). Preferably alcoholic and/or aldehyde contents should be used because the materials used are less likely to be affected. When alkaline cleaning agents are used (pH 9.5 - 11.5) it is possible that discolouring will take place on the metallic surface. This does not cause the product to be inoperable. Strong alkaline cleaning agents (pH > 11.5) are to be avoided. Only cleaning agents which have a proven effectiveness are permitted to be used and which are comparable with the applied system (disinfector and cleaning agent) in the available examination.

When selecting the disinfector it is to be observed

that the disinfector has a certified effectiveness (e.g. DGHM or FDA Approval, resp. CE certification mark according to EN ISO 15883),

that as far as possible a tested program for thermal disinfection (at least 5 min at 90°C or A0-value > 3000) is used (with chemical disinfection there is a danger of disinfectant residue on the instruments),

that the program used is suitable for the instruments and contains sufficient rinsing cycles,

that only sterile or germ-free (max. 10 germs/ml) and endotoxinfree (max. 0.25 endotoxin units/ml) rinsing water is used,

#that the air used for drying is filtered,

#that the disinfector is regularly maintained and inspected.

When selecting the applied cleaning agent system it is to be observed

#that this is basically suitable for cleaning instruments,

#that, as long as no thermal disinfection is used, additionally a suitable disinfection agent with certified effectiveness is used (e.g. DGHM or FDA Approval, resp. CE certification mark) and that this is compatible with the applied cleaning agent,

#that the chemicals used are compatible with the instruments.

It is imperative that the data given by the manufacturer of the cleaning agent and disinfection agent concerning the specified concentrations and effectiveness times be adhered to.

Sequence:

1. Place the instruments into the disinfectant. Ensure that the cable is not bent or squeezed.
 2. Start the program.
 3. Remove the instruments from the disinfectant after the program has finished.
 4. Check and pack as far as possible immediately after removal (after additional drying at a clean location).
- If necessary blow out with filtered compressed air. Only dry the products with compressed air (< 3 bar) in order to avoid possible damage.

4. Inspection

Check the products for operation as well as visual damage according to the corresponding operating instructions. LED recommend that with moveable or swivelling parts (e.g. forceps) instrument oil should not be used (exception: surgical instruments made of stainless steel) due to the fact that certain plastics tend to swell and owing to the oil even hinder movement.

5. Packing

Before sterilization, the disassembled instruments must be packed into a suitable one-way sterilization packing (single layer or double layer packing) and/or a suitable sterilization container:

according to EN 868/ANSI AAMI ISO 11607

suitable for steam sterilization (temperature resistive up to 137°C, sufficient steam porosity)

regular maintenance (sterilization container)

IMPORTANT: It is not permitted to sterilize in the transport packing.

6. Autoclaving

The instruments must only be sterilized in the disassembled condition. Only the following sterilization methods are to be used for sterilization, other sterilization methods are not permitted.

Steam sterilization

Fractionated vacuum method* (with sufficient product drying)

Steam sterilizer according to EN 13060 resp. EN 285

According to EN 554/ANSI AAMI ISO 11134 validated (valid commissioning and specific product service assessment)

Maximum sterilization temperature 134°C (plus tolerance according to EN 554/ANSI AAMI ISO 11134)

Sterilization time min. 20 min (at 121°C 1.0 bar) resp. 5 to 20 min at 132°/134 °C 2.0 bar.

The suitability of the product for effective sterilization was proven with the fractionate vacuum method with the above sterilization times/temperature. The hot air method for sterilization is basically not to be used (destruction of the instruments).

The use of other methods of sterilization (e.g. ethylene oxide, formaldehyde, beam and low temperature plasma sterilization) is outside the responsibility of the manufacturer. In this case the corresponding valid standards are to be observed (EN ISO 14937/ ANSI AAMI ISO 14937, resp. standards for specific methods) and the suitability and principle effectiveness of the method proved (if necessary, including residue inspection of the sterilization material) taking the specific product geometry into consideration within the scope of the validation.

7. Storing for use

Once sterilised, the products must be stored dry and dust-free in their sterilisation packing. The products must also be protected against sun and heat.

8. Monitoring/inspection/function testing of parts in the operating theatre

The disassembled instruments are to be assembled according to the corresponding operating instructions for the application. A visual inspection and function test is to be carried out before every operation.

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