**Gebrauchsanweisung  
Infrarot-Multifunktions-  
Thermometer**  
Instructions  
**Infrared Multifunction  
Thermometer**  
Mode d'emploi  
**Thermomètre à infrarouge  
multifonctions**  
Instrucciones para el uso  
**Termómetro multifuncional  
infrarrojo**  
Инструкция по эксплуатации  
**Инфракрасный  
многофункциональный термометр**  
Istruzioni per l'uso  
**Termometro multifunzionale  
a infrarossi**

## Sommaire

## Page

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Guidance and manufacturer's declaration – electromagnetic emission –  
for all EQUIPMENT AND SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic emission  |            |  |
|---|------------|--|
| The IRIDE1 Infrared Ear Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the IRIDE1 Infrared Ear Thermometer should assure that it is used in such an environment. |            |  |
| Emissions test  | Compliance | Electromagnetic environment - guidance   |
| RF emissions<br>CISPR 11  | Group 1    | The IRIDE1 Infrared Ear Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions<br>CISPR 11  | Class B    |  |
| Harmonic emissions<br>IEC 61000-3-2   | N/A        |  |
| Voltage fluctuations / flicker emissions<br>IEC 61000-3-3   | N/A        |  |

Guidance and manufacturer's declaration – electromagnetic immunity –  
for all EQUIPMENT AND SYSTEMS

| Guidance and manufacturer's declaration – electromagnetic immunity  |  |                              |  |
|---|--|------------------------------|--|
| The IRIDE1 Infrared Ear Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the IRIDE1 Infrared Ear Thermometer should assure that it is used in such an environment. |  |                              |  |
| Immunity test   | IEC 60601 test level   | Compliance level             | Electromagnetic environment - guidance   |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2  | ± 6 kV contact<br>± 8 kV air   | ± 6 kV contact<br>± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.   |
| Electrostatic transient / burst<br>IEC 61000-4-4  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines   | N/A                          | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5  | ± 1 kV differential mode<br>± 2 kV common mode   | N/A                          | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11  | < 5 % $U_T$<br>(>95 % dip in $U_T$ )<br>for 0.5 cycle<br>40 % $U_T$<br>(60 % dip in $U_T$ )<br>for 5 cycles<br>70 % $U_T$<br>(30 % dip in $U_T$ )<br>for 25 cycles | N/A                          | Mains power quality should be that of a typical commercial or hospital environment. If the user of the IRIDE1 Infrared Ear Thermometer requires continued operation during power mains interruptions, it is recommended that the IRIDE1 Infrared Ear Thermometer be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8  | 3 A/m  | 3 A/m                        | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |
| NOTE  | $U_T$ is the a. c. mains voltage prior to application of the test level.   |                              |  |

**Guidance and manufacturer's declaration – electromagnetic immunity –  
for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING**

| <b>Guidance and manufacturer's declaration – electromagnetic immunity</b>   |                             |                         |  |
|---|-----------------------------|-------------------------|--|
| The IRIDE1 Infrared Ear Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the IRIDE1 Infrared Ear Thermometer should ensure that it is used in such an environment.   |                             |                         |  |
| <b>Immunity test</b>  | <b>IEC 60601 test level</b> | <b>Compliance level</b> | <b>Electromagnetic environment - guidance</b>  |
| Conducted RF<br>IEC 61010-4-6   | 3 Vrms                      | N/A                     | Portable and mobile RF communications equipment should be used no closer to any part of the IRIDE1 Infrared Ear Thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><br><b>Recommended separation distance</b><br><br>$d = \left[ \frac{3.5}{F_1} \right] \sqrt{P}$  |
|   | 150 kHz to 80 MHz           |                         |  |
| Radiated RF<br>IEC 61010-4-7  | 3 V/m                       | 3 V/m                   | $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz<br><br>$d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz<br><br>where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m). <sup>1</sup><br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>2</sup> should be less than the compliance level in each frequency range. <sup>3</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |
|   | 80 MHz to 2.5 GHz           |                         |  |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.   |                             |                         |  |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.  |                             |                         |  |
| <sup>1</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IRIDE1 Infrared Ear Thermometer is used exceeds the applicable RF compliance level above, the IRIDE1 Infrared Ear Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the IRIDE1 Infrared Ear Thermometer. |                             |                         |  |
| <sup>2</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 1V/m.  |                             |                         |  |

**Recommended separation distances between portable and mobile  
RF communications equipment and the EQUIPMENT or SYSTEM -  
for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING**

| <b>Recommended separation distances between<br/>portable and mobile RF communications equipment and the IRIDE1 Infrared Ear Thermometer</b>   |  |   |   |
|---|--|---|---|
| The IRIDE1 Infrared Ear Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the IRIDE1 Infrared Ear Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IRIDE1 Infrared Ear Thermometer as recommended below, according to the maximum output power of the communications equipment. |  |   |   |
| Rated maximum<br>output of transmitter<br><br>W   | <b>Separation distance according to frequency of transmitter</b><br><b>m</b> |   |   |
|   | 150 kHz to 80 MHz  | 80 MHz to 800 MHz                             | 800 MHz to 2.5 GHz                          |
|   | $d = \left[ \frac{3.5}{F_1} \right] \sqrt{P}$                                | $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ | $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ |
| 0.01  | /  | 0.12  | 0.23  |
| 0.1   | /  | 0.38  | 0.73  |
| 1   | /  | 1.2   | 2.3   |
| 10  | /  | 3.8   | 7.3   |
| 100   | /  | 12  | 23  |
| For transmitters rated at a maximum output power not listed above the recommended separation distance $d$ in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.   |  |   |   |
| NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.   |  |   |   |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.  |  |   |   |

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→ **Riester offers a large selection of products in the areas of**

Blood pressure measuring devices | Instruments for ENT, Ophthalmological instruments | Dermatological instruments | Thermometers | Stethoscopes | Head mirrors, Head lights, Examination lights | Laryngoscopes | Gynaecological instruments | Percussion hammers | Tuning forks | Products for blood stasis | Pulmonary pressure measuring devices | Dynamometers | Pressure infusion instruments | Veterinary instruments | Doctor's cases and bags

**Detailed descriptions of the products can be found in the respective sections of the omnibus edition catalogue (Order No. 51232-50). Or online under [www.riester.de](http://www.riester.de).**

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**Rudolf Flester GmbH**
→ **Rudolf Riester GmbH**