|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use</b><br>Spring-handle instruments                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| USA | Note for U.S. users<br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusa.com">www.aesculapusa.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung</b><br>Instrumente mit Federgriff                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| fr  | <b>Mode d'emploi</b><br>Instruments médicaux à poignée avec ressort                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| es  | <b>Instrucciones de manejo</b><br>Instrumentos con mango de resorte                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| it  | <b>Istruzioni per l'uso</b><br>Strumenti con impugnatura a molla                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| pt  | <b>Instruções de utilização</b><br>Instrumentos com pega de mola                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| nl  | <b>Gebruiksaanwijzing</b><br>Instrumenten met veergreep                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| da  | <b>Brugsanvisning</b><br>Instrumenter med fjederhåndtag                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| sv  | <b>Bruksanvisning</b><br>Instrument med fjäderhandtag                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| fi  | <b>Käyttöohje</b><br>Jousikahvalla varustetut instrumentit                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| lv  | <b>Lietošanas instrukcijas</b><br>Instrumenti ar atsperes rokturi                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| lt  | <b>Naudojimo instrukcija</b><br>Instrumentai su spyruokline rankenėle                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| ru  | <b>Инструкция по применению</b><br>Инструменты с пружинной рукояткой                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| cs  | <b>Návod k použití</b><br>Nástroje s pružinovou rukojetí                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| pl  | <b>Instrukcja użytkowania</b><br>Przyrządy z uchwytem sprężynowym                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| sk  | <b>Návod na použitie</b><br>Nástroje s pružinovou rúkoväťou                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| hu  | <b>Használati útmutató</b><br>Rugós markolatú műszerek                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| sl  | <b>Navodila za uporabo</b><br>Instrumenti z vzmetnim ročajem                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| hr  | <b>Upute za uporabu</b><br>Instrumenti s opružnom drškom                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| ro  | <b>Manual de utilizare</b><br>Instrumente cu mâner cu arc                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| bg  | <b>Упътване за употреба</b><br>Инструменти с пружинен захват                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| tr  | <b>Kullanım Kılavuzu</b><br>Yaylı sapı olan aletler                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| el  | <b>Οδηγίες χρήσης</b><br>Όργανα με ελατηριωτή λαβή                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

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