 ***EURODEFIPADS[®]*** 

 ***FIAB[®]*** 

CE
2797

- IT** Elettrodi multifunzione monouso Eurodefipads®
- EN** Disposable multifunction electrodes Eurodefipads®
- FR** Électrodes multifonctions à usage unique Eurodefipads®
- DE** Multifunktions-Einwegelektroden Eurodefipads®
- NL** Multifunctionele wegwerpelektroden Eurodefipads®
- ES** Electrodos multifuncionales desechable Eurodefipads®
- PT** Eléktrodos multifunção descartáveis Eurodefipads®
- RU** Многофункциональные одноразовые электроды Eurodefipads®
- EL** Ηλεκτροδια πολλαπλών λειτουργιών μιας mono χρήσης Eurodefipads®
- SV** Multifunktionella elektroder för engångsbruk Eurodefipads®
- HR** Višenamjenske elektrode za jednokratnu upotrebu Eurodefipads®
- PL** Elektrody wielofunkcyjne jednorazowego użytku Eurodefipads®
- CS** Jednorázové multifunkční elektrody Eurodefipads®
- HU** Egyszerhasználatos többfunkciós elektródák Eurodefipads®

ISTRUZIONI D'USO
OPERATING INSTRUCTION
INSTRUCTIONS D'UTILISATION
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZINGEN
INSTRUCCIONES DE USO
INSTRUÇÕES
РУКОВОДСТВО ПО ЭКСПЛУАТАЦИИ
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
BRUKSANVISNING
UPUTE ZA UPOTREBU
INSTRUKCJE UŻYTKOWANIA
POKYNÝ K POUŽITÍ
HASZNÁLATI UTASÍTÁS

 **EURODEFIPADS®** 

- [1] Ian Jacobs, Kjetil Sunde, et al. "Part 6: Defibrillation. 2010 International Consensus on Cardiopulmonary Resuscitation and Emergency Cardiovascular Care Science With Treatment Recommendations" *Circulation* 2010;122:S325-S337
- [2] Jerry Nolan et al. "European Resuscitation Council Guidelines for Resuscitation 2010" *Resuscitation* Vol. 81, 2010, Editor Jerry Nolan, Bath, UK
- [3] A. Luciano et al., "BODY-MASS INDEX REFERENCE CURVES FOR CHILDREN AGED 3-19 YEARS FROM VERONA, ITALY", *European journal of clinical nutrition*, 51(1), 1997, pp. 6-10
- [4] Edward J. Sondik, et al. "2000 CDC Growth Charts for the United States: Methods and Development" *Vital and Health Statistics. Series 11, Number 246*: May 2002
- [5] John M. Field et al. "2010 American Heart Association Guidelines for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care" *Circulation* 2010, 122:S640-S946 *Journal of the American Heart Association* Editor Mary Fran Hazinski, RN, MSN
- [6] Luis A. Pagan-Carlo, Kirk T. Spencer, et al. "Transthoracic defibrillation: importance of avoiding electrode placement directly on the female breast" *Journal of American College of Cardiology* 1996;27:449-452
- [7] Edward A. Panacek, Mark A. Munger et al. "Report of nitropatch explosions complicating defibrillation" *American Journal of Emergency Medicine* 1992;10:128-9
- [8] Keith Wrenn "The hazards of defibrillation through nitroglycerin patches" *Annals of Emergency Medicine* 1990;19:1327-8
- [9] Johannes C. Manegold, Carsten W. Israel et al. "External cardioversion of atrial fibrillation in patients with implanted pacemaker or cardioverter-defibrillator systems: a randomized comparison of monophasic and biphasic shock energy application" *European Heart Journal* 2007;28:1731-8

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