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ED2007 REV02_IE

multiCare[®] GLC

Reagent strips for the determination of glucose concentration in whole blood or venous blood with lithium heparin or EDTA

To be used only with the multiCare meter
Instructions for use

Introduction. The self-checking of the glycaemia allows the diabetic patient to have immediate information about his/her glucose level, but it must be performed only after accurate training with a physician or with a qualified technician. It is important to point out that the self-checking cannot replace a physician, and that it cannot be performed without constant follow up by the physician (to make up a timetable, to evaluate the results and to define a therapy).

Test principle. Colorimetric reaction with glucose oxidase and peroxidase.

Reagents. The reactive area of each strip contains the following reagents: Chromogen min. 0,012 mg; Glucose-oxidase min. 1,8 U; Peroxidase 0,57 U.

Performing the test

- Open the box of multiCare Glucose strips of blue colour and take the data-chip of the same colour contained inside the box.

- Insert the data-chip into the appropriate space behind the multiCare Reflectometer with the arrow turned upward (Picture 1).


- Open the vial, take out the strip and close it immediately.
- Insert the strip into the reflectometer, with the arrow turned towards the reflectometer and the parameter symbol on the up side of the strip (Picture 2) until you hear a click sound. The instrument switches on automatically - Check on the display if the GLC symbol appears and check that the number displayed is the same of the one printed on the strip vial label (Picture 3).

- Execute the puncture of the finger using the Lancing Device contained in the multiCare Reflectometer set and a sterile lancet.

- Remove the first drop of blood and gently press the finger tip to obtain a second drop (minimum 3 microliters).
- When blood is sufficient, place the drop of blood into the white rounded part of the strip as indicated in the picture 4 and after some seconds the result will be displayed.

Only one measurement has to be carried out for each finger pricking!
Carefully follow the instructions contained in the multiCare user manual.

All components of the pack has to be disposed properly, avoiding contaminations.

Correct test execution. Before each test, verify that the reactive area of the strip is uniformly white coloured. Should the colour be altered, do not use the strip for measurement. After each testing, turn the strip and verify that the blue colour is homogeneous. If this is not true, it is necessary to repeat the test with a new strip. You can perform only one measurement with each finger prick.

Checking the system. To test the user's skill, the multiCare system performance and its reactive strips, a periodic check should be carried out using the multiCare control solutions as described in the instructions contained in the package. The glucose concentration intervals are printed on the labels of the reagent strips' vial.

Storage and expiry date. Keep the reactive strips at temperatures between 10°C (50°F) and 35°C (95°F) and always keep them in the original vial. Do not expose to direct sun light. Close the vial with the original cap immediately after pull out the strip. In this way, the strips will last up to 3 months after the first opening. The expiry date on the package of strips refers to a new unopened vial. The use of expired strips can cause wrong results. Do not expose to excessive humidity. Do not freeze!

Normal Values at fasting (in a person that is not a diabetic) 70-110 mg/dl (3,9-6,1 mmol/l), if the glucose concentration is considered in the plasma. The multiCare system allows to obtain glycaemia values calibrated on reference method that measures the glucose concentration on serum/plasma. The glucose concentration measured on serum/plasma is 12% higher than the one measured on whole blood or hemolysate. The diabetic patient has to agree with the physician the optimal values that are appropriate for a well-controlled metabolism.

Table of results :

Glucose below 25 mg/dL (1,4 mmol/l) = LO

Glucose between 25 mg/dl (1,4 mmol/l) and 500 mg/dL (28 mmol/l) = the numerical result

Glucose above 500 mg/dL (28 mmol/l) = HI*

*Such results mean that the glucose concentration in the sample is very elevated and it indicates a serious risk. Contact a physician immediately.

The analytical results have to be communicated to the physician.

The self-checking can be useful for the physician in order to formulate a diagnosis.

Warnings. Keep the multiCare strips away from children. Do not perform tests at temperatures below 5°C (41°F) or above 45°C (113°F). Use the strip immediately after removing it from the vial.

Note for the physician and laboratory technicians. To perform a glycaemia test with the multiCare system, it is possible to use either fresh capillary blood or venous blood treated with heparin or K-EDTA. Anticoagulants as fluoride or oxalate must not be used.

Hematocrit range: 35-55%. Abnormal hematocrit values (greater than 55% or less than 35%) can alter the result obtained with the multiCare strips.

Limitations

* The system is optimised for capillary or venous blood analysis. Do not use serum or plasma.

* Repeat the test if the glucose value obtained does not correspond to the sensations felt by the patient or to the basic knowledge of health conditions and symptoms. Consult the user manual to verify the correct procedure of determination, otherwise consult your physician.

* Capillary blood glucose values may be significantly lower than "true glucose levels" in the hyperglycemic-hyperosmolar state, with or without ketosis. Very ill patients should not be tested by this system or at least tested with extreme caution.

* If glucose values below 50 mg/dl (2,7 mmol/l) or above 250 mg/dl (13,8 mmol/l) are obtained, consult a physician as soon as possible.

* In the case of severe dehydration and excessive water loss the tested values can be significantly inferior in respect to the real values. Reasons for dehydration include vomiting, diarrhoea and the use of certain drugs such as diuretics..

* Reducing agents like for example, ascorbic acid at normal concentrations do not interfere with the results, but if high quantities are taken, abnormally low results could be obtained.

* Marked lipemic samples may interfere with the results.

* Inaccurate results may be obtained at high altitudes due to the lower oxygen concentration.

* The multiCare system is not suitable for blood samples from babies.

Accuracy. The multiCare system average error compared to the hexokinase reference method (bibliography) is <10%. In a series of measurements performed in antidiabetic centre, the following linear regression has been obtained: y[mg/dl]= 1.0016x-0.5635.

Precision. Repeatability: average imprecision is <5%. In a series of measurements performed in laboratory an average coefficient of variation 2.3% was obtained. Reproducibility: average imprecision is <5%. In a series of measurements performed in laboratory an average coefficient of variation 2.3% was obtained.

Minimum blood sample volume: 3 µL.

Packaging:	25 strips pack	50 (2x25) strips pack
Manufacturer:	Biochemical Systems International srl	
Customer Service:	Phone +39 0575 984164	e-mail: multicare@biosys.it

CE 0123 CE mark

IVD In vitro Diagnosticum

Store at...

Use until...

LOT Batch number

Follow instructions

Control solutions

Low

Mid

High