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Measurement range: 34.0-42.5°C (93.2-108.5°F)

Room temperature working range: 10/40°C  
55/104.0°F

Accuracy	from 34.0 to 35.9°C: ±0.3°C
°C	from 36.0 to 39.0°C: ±0.2°C (1)
	from 39.1 to 42.5°C: ±0.3°C

Accuracy	from 93.2 to 96.7°F: ±0.5°F
°F	from 96.8 to 102.2°F: ±0.4°F (1)
	from 102.3 to 108.5°F: ±0.5°F

### Other readings (2)

Measurement range: 1.0/55.0°C (33.8-131°F)

Room temperature working range: 10/40°C  
55/104.0°F (3)

Accuracy	from 1.0 to 19.9°C: ±1.0°C
°C	from 20.0 to 35.9°C: ±0.3°C
	from 36.0 to 39.9°C: ±0.2°C
	from 39.1 to 42.5°C: ±0.3°C
	from 42.5 to 55°C: ±1.0°C

Accuracy	from 33.8 to 67.9°F: ±1.8°F
°F	from 68.0 to 96.7°F: ±0.5°F
	from 96.8 to 102.2°F: ±0.4°F
	from 102.3 to 108.5°F: ±0.5°F
	from 108.7 to 131°F: ±1.8°F

(1) The ASTM E1965-98-2009 standard requires ±0.2°C (±0.4°F) accuracy for infrared thermometers in the 37-39°C (98.6-102.2°F) temperature range while the ASTM E667-86 and E1112-86 standards require ±0.1°C (±0.2°F) precision for mercury and electronic thermometers in that same temperature range.

(2) Only with 0700A2 model: subtract 2°C (3°F) from the temperature shown on the display to get a rough value. If the thermometer is set in rectal or internal temperature, subtract 3°C (5°F).

(3) Thermofocus can also take surface temperatures in rooms where the temperature is lower than 10°C (55°F) although accuracy and the operating range are not guaranteed and the display shows "Lo.5" and the value alternately.

Power supply:	4 AAA (LR03) alkaline batteries - 1.5V
Life of high quality batteries:	up to 3 years or 10,000 readings (depending on use)
Dimensions:	m 165 x 40 x 22 (6.5 x 1.45 x 0.86 inches) - including cap
Weight:	gr. 92 (3.17 oz.) - batteries included

Distance from the subject: calculated using an optical aiming system (approximately 3 cm / 1.2 inches).

Large, clearly visible display.

Atmospheric pressure range of operating conditions: from 700 hPa to 1,060 hPa.

Relative humidity range of operating conditions: from 15 % to 93 %, non-condensing.

Keep in a clean, dry place, preferably at a temperature ranging between +16 and +40°C (60.8 and 104°F). Store at a temperature included between -10 and +60°C (14 and 140°F) and in any case not lower than -18°C (0.4°F) or higher than +70°C (158°F).

Expected life: 10 years.

Thermofocus is a class II medical device (IIa as per 93/42/EEC Directive and subsequent amendments and integrations) tested in hospitals, private clinics and medical offices.

Contact the manufacturer for clinical tests.

Tecnimed srl, P.le Coccoi, 12 - 21040 Vedano Olona (VA) - Italy undertakes full responsibility for this product's compliance with the reference standards.

Thermofocus is compliant with the pertinent ASTM standard (E1965-98-2009) and CEI EN 60601-1, CEI EN 60601-11 standards.

Internally powered equipment for continuous operation.

The aiming LEDs of Thermofocus emit low luminous radiation largely under the level foreseen from IEC 60825-1:1993+A1:1997+A2:2001 Standard for Class 1 lights (class I, wavelength: 610 nm; maximum power: Pmax<1mW).

	Follow instructions for use
	ATTENTION: read the warnings
	Device designed for continuous use
	Applied part: type BF
	Recyclable material

Dispose of the device in compliance with applicable legislation regarding electrical equipment and batteries.

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EN

## USER MANUAL

Read carefully the instruction before using the thermometer

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Manufacturing Ideas

## INDICATION FOR USE

Thermofocus is an infrared thermometer intended for intermittent measurement of human body temperature in people of all ages.

## 1. INTRODUCTION

Dear Customer, thank you for buying Thermofocus®, the world's first non-contact thermometer.

can be affected by profuse sweating, superficial wounds or head injury.

**8. Do not use the thermometer on a sweaty forehead**, since the temperature reading would be unreliable. Read the par.#4.2.

9. The gill waveguide (see figure 1) is the most delicate part of the device. It is composed of a gold-plated concave mirror that must be kept clean, crystal clear and intact. Any damage, dust or dirt will alter the temperature reading.

10. Do not handle the thermometer, and in particular the tip (fig. 1), for longer than strictly necessary before taking the reading.

11. Do not use the thermometer in direct contact with the ear or other parts of the body.

12. Do not use the thermometer in direct contact with objects or liquids.

Do not submerge the thermometer in water or other liquids and keep it away from sources of heat and out of direct sunlight. If water seeps into the thermometer, contact your Dealer immediately for Technical Service.

13. To prevent dust or dirt infiltrations in the tip, always close the cap.

14. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

15. Avoid knocking and dropping it,

and do not use it if damaged or if not functioning properly.

To correctly judge a fever, you need to know the normal temperatures of your family members when they are in good health condition and at various times of the day.

To help you with this, you can fill out the table at [www.thermofocus.com](http://www.thermofocus.com), FAQ no. 7. and keep it as a reference.

16. Prevent dust or dirt infiltrations in the tip, always close the cap.

17. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

18. Avoid knocking and dropping it,

and do not use it if damaged or if not functioning properly.

19. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

20. Avoid knocking and dropping it,

and do not use it if damaged or if not functioning properly.

21. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

22. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

23. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

24. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

25. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

26. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

27. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

28. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

29. Do not use the device on a subject

making a call with a mobile or

cordless telephone or in the presence of strong electromagnetic fields.

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