GIMA

PROFESSIONAL MEDICAL PRODUCTS

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GB - USER'S GUIDE FOR REUSABLE NON-SURGICAL INSTRUMENTS

Read safe use instructions carefully before use. These instruments are intended to be used by medical practitioners who are specially trained on how to use and care of them. The incorrect use, poor or inappropriate maintenance can rapidly lead to deterioration of the instruments. The first time, and in general before and after every use, it is recommended to clean, dry and sterilise, if required, the instruments. The device must be cleaned with warm water, disinfected and rinsed if possible with distilled water. For automated cleaning no particular requirement is requested. For manual cleaning remove visible dust/spots with disposable cloth/paper, use approved and certified detergents. Always follow instructions on how to use the detergent. Also make sure that the product

does not come into contact with acids or other aggressive disinfectants that could corrode it. Disinfectants must conform to the applicable national regulation for disinfectants (such as FDA Approval, and CE Mark).

Instruments that are not dried could suffer damage by corrosion. Always dry the instruments with lint-free disposable cloths (do not use recycled cloths). If sterilization is required, European Standard ISO 17665 i.e. Steam Sterilization for Medical Devices is to be followed. Sterilization is achieved by exposing products to saturated steam at high temperatures (121°C to 134°C). Product(s) are placed in device called the autoclave and heated through pressurized steam to kill all microorganisms including spores. The device's exposure time to steam would be anywhere between 3 to 15 minutes.

Examples of minimum temperatures and times established for adequate level of microbial lethality in sterilization processes (refer table – 1 of ISO/TS 17665-2).

Temperature °C	Time Minutes
121	15
126	10
134	03

In the context of validating the sterilisation process, check the suitability of the specific measures for drying. The humidity in the container can cause the instruments to rust. Often bad, and insufficient drying, is due to the incorrect positioning of the load and the use of unsuitable types of cloths for drying.

There is no advice about maximum number of sterilisation cycles, this depends largely on the state of the product. Repeated processing has minimal effect on the devices; instruments that show signs of corrosion, damage or crack must be immediately replaced and not to be used. Always perform a visual inspection for damage or signs of wear; check for smooth functioning of medical device. It is recommended to store the device in clean and dry condition. Whilst the instruments are under warranty, repairs and/or replacement of all defective parts, for reasons that are accepted by the manufacturer, are free, with the exception of costs for labour, travel, transport, packing, etc. Damage caused by improper use of the product is not covered by the warranty.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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