# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**BILANCIA ASTRA**  
**ASTRA SCALE**  
**BALANCE ASTRA**  
**WAAGE ASTRA**  
**BÁSCULA ASTRA**  
**BALANÇA ASTRA**  
**ΖΥΓΑΡΙΑ ASTRA**  
**ASTRA ميزان**


**REF** 27310

Gima S.p.A.  
Via Marconi, 1 - 20060 Gessate (MI) Italy  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)  
Made in China

CE 0476 CE M24 0474

## CARATTERISTICHE

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## FEATURES

The GIMA Astra scale is a device designed to measure and weigh people. The broad measuring range, high precision and sensitivity make it particularly suitable for hospitals, clinics, schools, family doctors' offices, gyms, etc.

### Technical Specifications

#### Weight measurements

Maximum capacity: 200 kg

Precision: 100 g

#### Height measurements

Measuring range: 75÷200 cm

Precision: 0,5 cm

Platform dimensions (L x W): 37,5x27,5 cm

Overall dimensions (L x W x H): 56,5x43,5x148,5 cm

Weight of scale: 15,3 kg

#### Box contents

- Column and balance beam
- Platform
- Height measuring rod
- Wheels
- Screws and washers

## PRESCRIPTIONS

*Do not use the equipment in case it is damaged. Apply to your retailer.*

*Avoid precarious repairs. Repairs shall be carried out with original spare parts only, which shall be installed according to the intended use.*

Since the product is made of corrosion-proof materials suitable for the environmental conditions foreseen for its normal use, does not require special care, however it is necessary to store it in a closed place making sure that is protected from dust and dirt to assure its hygienic conditions. Moreover, it is recommended to store the product in a place which can be reached easily by the personnel in case of necessity.

## UNPACKING

*Always remember that packing elements (paper, cellophane, stitches, adhesive tape, etc.) can cut and/or hurt if they are not carefully handled. They shall be removed with adequate means and shall not be left at the mercy of irresponsible persons; the same is valid for tools used to remove packages (scissors, knives, etc.).*

After opening the packages, first of all it is necessary to check all pieces and parts composing the product. Check that they are all present and in perfect conditions.

## INSTALLATION

The device is mostly factory pre-assembled, therefore assembly is very easy (see illustration):

- Fix the column to the base using the supplied screws and washers
- Lay the scale on the floor and hook the transmission rod onto its seat
- Fix the wheels to the base with the supplied screws and washers
- Fix the measuring ruler to the two hexagonal-head screws pre-assembled on the column.
- Fix the ruler by tightening the hexagonal-head screws using the relative wrench. Don't tighten too much!

To calibrate the scale, after setting both weights to zero, work on the adjustment screws (see illustration) until the indicator matches the measurement tab.

The base must be positioned on a perfectly flat floor to calibrate the unit correctly.

### OPERATION

Before the person steps onto the platform, bring the height cursor up above the person's height. After the person is on the platform, slowly lower the cursor until it rests on the person's head.

The graduated measuring beam can be extended on both sides: the lower part can be used to measure heights up to 101,5 cm. The height value is indicated right by the cursor. The upper ruler is used to measure heights up to 200 cm. The height value is indicated right where the two graduated measuring beams intersect, indicated by the READ arrow (see illustration).

Note: to correctly use the upper measuring beam block the cursor on the upper part of the sliding measuring ruler.

*Attention: Do not keep the cursor horizontally because it may be dangerous.*

### MAINTENANCE

The device was designed to last a lifetime. It is resistant to chemical attack, except for M.E.K. ketons, diluting solutions and polish removers.

*Periodically verify the calibration and the leveling of the scale. Perform this operation every time it is moved (see installation).*

### GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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