PLICOMETRO FAT-1
SKINFOLD CALIPER FAT-1
PLICOMÈTRE FAT-1
HAUTFALTENDICKENMESSER FAT-1
PLICÓMETRO FAT-1
ADIPÔMETRO FAT-1
ПАХΥΜΕΤΡΟ FAT-1
FAT-1 جهاز قياس ثنية البشرة
PLICOMETRU FAT-1
SUWMIRAKA FAT-1
KALIPER FAT-1
HUIDPLOOIDIKTEMETER FAT-1
HUDVECKSMÄTAREN FAT-1
TOLÓMÉRŐT FAT-1
HUDFOLDSMÅLEREN FAT-1
Инструментът FAT-1
KALIBRS FAT-1

REF 27344

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

CE 0476


Fig. 1


Fig. 2


Fig. 3

Fig. 4

ITALIANO

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TECHNICAL CHARACTERISTICS

The skinfold caliper is an instrument created by Tanner and Whitehouse in 1955 which permits to assess the nutritional state of a subject thanks to the measurement of the thickness of the skinfolds. It must therefore meet two precise and now universally adopted characteristics, which are:

- 1) the two arms of the instrument must exercise on the skinfold a pressure of approximately 10 g./mm²
- 2) the grasp on the skinfold must correspond to an area of approximately 6 x 11 millimetres.

This instrument, in compliance with the above technical characteristics, thanks to the reading scale in mm on the upper edge and indicated by the black arrow, permits to measure with an approximation of +/- 2 mm the skinfold.

MODE OF USE

In order to make the measurement it is necessary to lift a skinfold between the fingers and measure its thickness with the skinfold caliper.

For the tricep measurement (Fig. 1) lift the skinfold in the back area of the arm, exactly halfway between the acromion and the proximity end of the olecranon. For the subscapula measurement (Fig. 2) lift the skinfold in the region immediately under one of the two scapulas. For the bicep measurement (Fig. 3) lift the skinfold at mid-length of the bicep. For the suprailiac measurement (Fig. 4) lift the skinfold found on the side.

Once completed this operation, proceed with the measurement using the instrument and immediately you can read the dimension in mm of the skinfold thickness.

The following page shows some values of skinfolds divided by sex, adults, athletes and children.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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