
ALZAWATER CON VITE DI BLOCCAGGIO FRONTALE
RAISED TOILET SEAT WITH FRONTAL LOCKING SCREW
REHAUSSEUR WC AVEC VIS DE BLOCAGE FRONTAL
ALZAWATER CON TORNILLO DE BLOQUEO FRONTAL
TOILETTENSITZERHÖHUNG MIT FRONTALER
SICHERUNGSSCHRAUBE
ΑΝΥΨΩΤΙΚΟ ΛΕΚΑΝΗΣ ΜΕ ΜΠΡΟΣΤΙΝΗ ΒΙΔΑ ΚΛΕΙΔΩΜΑΤΟΣ

**ITALIANO**

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cal authority beforehand for the correct and environmentally friendly way to dispose it.

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REF	Product code	LOT	Lot number
CE	Medical Device compliant with Regulation (EU) 2017/745		Manufacturer
SN	Serial number		Date of manufacture
MD	Medical Device		Consult instructions for use
UDI	Unique device identifier		

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GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

ENGLISH

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

INTENDED USAGE

The products are designed to increase the height of the toilet bowl, in order to facilitate sitting down and getting up and to promote independence and comfort. The frontal locking system and the armrests make the use of the product safer.

INDICATIONS AND CONTRAINDICATIONS

The Raised Toilet Seat is aimed to elderly people, wheelchair users, people with reduced mobility of lumbar, knee or hip and to all those who might have difficulty in sitting down and getting up from the toilet bowl. The version provided with folding armrests makes this Raised Toilet Seat particularly suitable for wheelchair users who resort to lateral transfer to the toilet seat. The ability to sit independently is required during the use. It is not recommended to use the Raised Toilet Seat in case of perception disorders or balance problems.

SAFETY WARNINGS AND INDICATIONS

- The maximum user weight indicated in the fitting instructions must not be exceeded.
- Before every use, be sure the Raised Toilet Seat is in stable position on the WC and correctly fixed with the special fixing accessories supplied.
- The armrests must be properly assembled and, during use, the moving part must be correctly hooked as shown in the assembly pictures of the User Manual.
- During the assembly, strictly adhere to what is shown in the User Manual.
- Do not modify the Raised Toilet Seat and its components in any way.
- Do not use the Raised Toilet Seat if damaged or altered in any way. Cracks and deformations may strongly affect the performance and safety of the device.
- Be sure the height of the Raised Toilet Seat is such as to make the use of it easy and safe for the user.
- Do not lean on the Raised Toilet Seat while sitting or standing.
- During use, avoid positions that keep the body weight too much forward. Wrong positions may lead the user to slip and fall.
- The Raised Toilet Seat must not be stored outdoors.

MAINTENANCE INFO

The products maintenance consists of the cleaning and sanitation procedures. In order to guarantee a safe use and adequate standard of hygiene, the user should perform these procedures before every use.

The user has to make certain of the structural integrity of the device and its components. For cleaning and sanitation procedures, follow the steps below:

1. Gloves should be worn.
2. Remove the evident dirt before to proceed with the most accurate disinfection.
3. Use water or non aggressive detergent for the cleaning procedure.
4. Wipe with a dry and clean cloth.
Alternatively, the Raised Toilet Seat can be cleaned in domestic dish-washers at the maximum temperature.

REUSE REQUIREMENTS

The product is suitable for reuse by other users. Before reuse by others, the product must be correctly cleaned and sanitized. Before the transfer of the product to others, the product must be checked by qualified dealers.

CORRECT DISPOSAL PROCEDURE

Do not dispose of the product in the normal household waste. Ask your lo-

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Gima S.p.A.