MAXIASPEED

ISTRUZIONI D'USO

INSTRUCTION MANUAL

MANUEL D'INSTRUCTIONS

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MANUAL DE INSTRUCCIONES


Fig - Abb. 1


Fig - Abb. 1a


Fig - Abb. 1b

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EN

1. Aspirator
 - 1.1 Telescopic pulling handle
 - 1.2 Vacuum gauge
 - 1.3 Vacuum regulator
 - 1.4 Vessel selector
 - 1.5 Main switch
 - 1.6 Operating switch with pedal
(for versions with the relevant socket only)
 - 1.7 Anti-static wheels with brake
 - 1.8 Power supply socket with fuse (inspection type)
 - 1.9 Cable and pedal switch compartment
 - 1.10 Socket for pedal switch (where foreseen).
2. Antibacterial filter
3. Short tube Ø 8x14mm
4. Vessels (2000cc or 4000cc)
 - 4.1 VACUUM connection
 - 4.2 PATIENT connection
5. Long tube Ø 8x14mm
6. Aspiration fitting
7. Tool for handle

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The **MAXIASPEED** aspirator is a professional device suitable for outpatient use, specifically for the aspiration of secretions. It is also suitable for use in hospital wards, on tracheotomy patients, or for surgical applications. It is equipped with 4 anti-static wheels (1.7) with a brake device, a telescopic pulling handle (1.1), vacuum regulator (1.3), vacuum gauge (1.2) and two vessels of 2 L. or 4 L. (depending on version). The vessels feature a protection device against liquid entering the vacuum pump, which interrupts the aspiration flow. It is lubrication free, easy to handle, simple to use, reliable, strong and quiet. The **MAXIASPEED** aspirator also has a switch (1.4) for switching the aspiration flow from one vessel to the other, very easily and without having to touch any possibly contaminated parts. The **MAXIASPEED** aspirator has the following accessories: 2 Vessels of 2 L. or 4 L. with protection device, power supply cable, 2 sterilisable silicone connection tubes Ø 8x14 - short, 2 sterilisable silicone connection tubes Ø 8x14 - long, 2 aspiration fittings, 2 antibacteria filters, pedal switch (for versions with the relevant socket only). **N.B.:** Only use genuine accessories supplied by 3A Health Care; the accessories supplied with the device by the manufacturer have been tested and proven compliant with the latest current safety standards. **WARNING!** Correct operation of the device is not guaranteed in the event of use of accessories other than those supplied.

INTENDED USE

Mains-operated medical /surgical aspirator intended for outpatient use.

Medical purposes: This product is intended for use for the aspiration of body fluids.

Intended users of the product:

- Legally certified medical personnel (doctors, nurses and therapists).
- Users must also be able to understand the operation of the medical device, and the contents of the instruction manual, in general terms.

Intended patients for the product: People who need to remove body fluids (saliva, blood, etc.).

Environment: This product is intended for outpatient use.

Expected duration: Duration may vary based on the operating environment. The lifetime of the device is 5 years and that of the collection vessel and the silicone tubes 1 year or 30 sterilisation cycles. The antibacterial filter is disposable device and must therefore be replaced after each application. Frequent use of the product may shorten the duration.

Precautions for use: The warnings and cautions described in the instruction manual must be observed.

IMPORTANT WARNINGS

This is a medical device and it must be used by authorised personnel. It must be operated as instructed in this use manual. It is important for the operator to read and understand the use and maintenance information of the device. Contact your retailer for any questions you may have. MICROBIAL CONTAMINATION: where pathologies are present with microbial infection and contamination risks, perform the cleaning and disinfecting operations scrupulously after every use on the patient.

The manufacturer ensures that everything necessary is done to ensure that every device is of the highest quality and safety. Nevertheless, basic safety rules must always be followed, as for every electrical device.

- The device must only be used by specialised personnel.
- This device must be exclusively used for the purpose for which it was designed; in particular as a surgical aspirator. All other uses are improper and dangerous, and the manufacturer cannot be held responsible for any consequences deriving from improper use.
- Never use adaptors for power supply voltages different to those shown on the plate data label on the back of the device.
- Keep the cable far away from hot surfaces.
- The device is not suitable for use in the presence of anaesthetic mixes inflammable with air, or oxygen, or nitrogen oxide.
- Never handle the power supply plug with wet hands.

 Never leave the device close to water; do not submerge it in liquids, do not wet it. If this should occur, pull the plug out of the power supply socket immediately, before touching the device. Do not use the device if the plug or power supply cable are damaged, or wet (send it to your retailer immediately).

- The device casing is not protected against the entry of liquids.
- Always take out the power supply plug immediately after use.
- Maintenance and/or repair must only be performed by authorised personnel. Unauthorised repairs invalidate the guarantee.
- Ensure that the connections and vessel closing is performed carefully to prevent aspiration leakages.
- Do not upset the vessels while they are connected to the device when running, because the liquid comes into direct contact with the hydrophobic antibacteria filter and immediately blocks aspiration; if this should happen, empty the vessel and replace the antibacteria filter.
- In the event of aspiration without vessel and/or antibacterial filter, or if you suspect that solid or liquid substances have entered the aspiration circuit, send the device to customer service.


- In view of their length, the power cord and connection hose could constitute a strangulation hazard.
- The device is equipped with a safety fuse, easily inspected in the event of a fault, located on the back of the device in the panel socket. **Pull out the power supply plug before inspection.**
- The antibacteria filter is a **disposable** element which must be replaced after every application.
- The device comes into contact with the patient via the single-use cannula (**NOT supplied with the device**): any aspiration cannulas, purchased separately from the device, must comply with the requirements of the ISO 10993-1 standard on the biocompatibility of materials.
- **Do not lift the device using the transportation handle in any circumstance.**
- Before each use, follow the cleaning and disinfection operations as shown in the **"CLEANING AND DISINFECTION OPERATIONS"** paragraph of this instruction manual.

INSTRUCTIONS FOR USE

Before each use, ensure that all accessories are perfectly clean according to the instructions in the **"CLEANING AND DISINFECTION OPERATIONS"** section.

1. Connect the device as shown in figure 1 (page 2):

1.1 Connect the antibacterial filter (2.) to the "INLET" connector of the aspirator, ensuring that the side marked "IN" is facing towards the connector. **Do not use the aspirator without the antibacterial filter, because from a bacteriological point of view, it becomes dangerous for the patient.**

1.2 Connect one end of the short silicone tube to the hose connection of the antibacterial filter and the other end to the "VACUUM" connection (4.1) of the vessel. The collection vessel is complete with overflow valve, lid and vessel in a clear material (polycarbonate). The antibacterial filter also protects the suction circuit from possible contaminants sucked in during use.

1.3 Connect the long tube to the "PATIENT" connector (4.2) of the selected plug.

- 2. Operation of the telescopic handle:** for space reasons the aspirator is supplied with the telescopic handle lowered; to raise it, pull it up until the travel end hooks engage. To put the telescopic handle down again, press the travel end hooks on the external sides of the telescopic tubes with the tool provided (Fig. 2) and push the handle downwards until it closes on the aspirator casing.
- 3.** Using the vacuum regulator (1.4) it is possible to set the desired vacuum value (bar). By turning the knob towards "+" a higher vacuum is obtained, while turning it towards "-" a lower vacuum is obtained: these values can be read on the vacuum gauge (1.2).
- 4.** Start the device by switching the main switch to the "I" (ON) position (1.5).

Version with pedal switch: If the device is equipped with the pedal switch, to use it connect the plug on the cable to the socket provided on the panel (1.10) on the back of the device and turn the switch (1.6) next to the main switch (1.5) to "I" (ON), allowing the device to be switched on/off using the pedal switch.

- 5.** After the application, switch of the device, remove the power supply cable from the current socket and perform the cleaning operations as illustrated in the **"CLEANING AND DISINFECTION OPERATIONS"** paragraph.
- 6.** There is a very capacious compartment (1.9) on the back of the device for putting the power supply cable back in after the application, and if the device is equipped with the pedal switch, also for this latter.

CLEANING AND DISINFECTION OPERATIONS

- Detach the long sterilisable silicone tubes Ø 8 x 14 - long (5.) from the vessel (4.). Disconnect the sterilisable silicone tubes Ø 8 x 14 - short (3.) from the respective vessels and from the respective antibacterial filters and empty the vessel/s. **The liquids collected in the vessels must be disposed of in the biological waste container.**
- In order to clean the accessories, rinsing them under hot running water is sufficient. Before performing the vessel cleaning operations empty their content into the biological material vessels and clean them by sterilisation in an autoclave with a cycle at 121°C. The silicone aspiration tubes and aspiration fitting can be sterilised in an autoclave with a cycle at 121°C. The complete vessel should be replaced every 30 sterilisation cycles. During the sterilisation operations the operator must use latex gloves and suitable clothing in order not to come into contact with any possibly contaminating substances.
- The disposable antibacteria filter (2.) must be replaced after every application.
- Never wash the device under water or by immersion; clean the outside of the device using a cloth moistened with detergent (non abrasive).

PERIODICAL CHECK FOR THE SAFETY OF THE DEVICE

MAXIASPEED does not require maintenance and/or lubrication. Nevertheless, several simple checks must be performed before every use:

- Check the integrity of the casing, power supply cable and the pedal switch, if supplied.
- Block the aspiration connector with a finger and check that the vacuum level reaches -0.9 bar.
- Check that there are no disturbing noises suspicious of malfunctioning.
- Check that the cage is positioned correctly inside its housing. It must be in axis with the plug aspiration hole so that the floating valve can block it in the circumstance in which the aspirated liquid exceeds the maximum volume the vessel can receive. **Do not connect the cage as shown in figure 3b; it must be connected as shown in figure 3a.**
- Check that the float is fitted in the correct position and that it is free to move inside the cage (dirt or encrustation may hinder its movement). Figure 4 shows the correct insertion.
- Check that the plug seal is inserted with the edging as illustrated in figure 5.

PROBLEMS, CAUSES and SOLUTIONS

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Excessive noise	Damaged pump or obstructions in the internal suction duct	Send to the customer care service
The unit turns on but it does not suck	- Damaged pump - Vacuum regulator fully open. Connection tubes disconnected and/or badly connected, broken connection tubes. Vessel not in a vertical position, full, or defective overflow valve. Possible blockage of the hydraulic circuit inside the unit	- Send to the customer care service. - Check the position of the vacuum regulator. Check the connections and the integrity of the tubes. Position the vessel in a vertical position, check the overflow valve (blocked) and/or replace the silicon tubes
The vacuum value cannot be adjusted	Damage to the internal hydraulic circuit or obstruction of the connection tubes to the aspiration unit	Send to the customer care service
The protection fuse is activated any time the device is turned on	Pump damaged or in short circuit	Send to the customer care service
The vacuum gauge does not work	Liquid penetrating the pneumatic circuit	Send to the customer care service

Note: if you experience faults or malfunctioning problems different to those listed above, always and exclusively contact authorised assistance centres.

TECHNICAL SPECIFICATIONS

Piston electric-compressor,
with thermal protection

	Single pump version	Double pump version
		High flow, high vacuum device
Supply voltage	230V ~ 50/60Hz 300VA	230V ~ 50/60Hz 330VA
Fuse	T2, 5A - 250V	
Adjustable vacuum level	0 ÷ -0,90 bar (-90 kPa)	
Air flow	approx. 60 l/min.	approx. 90 l/min.
Protection device – EN 60529	IP20	
Noise level	approx. 60 dBA (1 m)	approx. 50 dBA (1 m)
Operation cycle	Continuous use	
Dimensions	470 x 580 x 560(H) mm	
Weight with 2 vessels of 2000cc	approx. 12 kg	approx. 16 kg
Weight with 2 vessels of 4000cc	approx. 13 kg	approx. 17 kg
Class of risk according to the 93/42/EEC directive	IIa	
Operating conditions	Temperature min. 10° C max 40° C - Humidity: min. 10% max 95%	
Storage conditions	Temperature min. -25° C max 70° C  - Humidity: min. 10% max 95% 	
Operating-storage air pressure	min. 690 hPa - max 1060 hPa 	

SYMBOLS

 Typ B device

 Class II device

 Fuse

 Main switch off

 Main switch on

 Alternate current

 Never use the device when taking a bath or a shower

 It is compulsory to carefully read the instructions before using this device

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 0051 Compliant with Medical Devices Directive 93/42/EEC

Electromagnetic Compatibility Compliance levels according to EN 60601-1-2:2015 standard

- ESD immunity: 15 kV air, 8 kV contact (EN 61000-4-2)
- Burst immunity: 2 kV/100 kHz (EN 61000-4-4)
- Surge immunity (EN 61000-4-5): 1 kV common mode /2 kV differential mode
- Magnetic field (EN 61000-4-8): 30 A/m
- Immunity to rf currents in the 150 kHz-80 MHz range (EN 61000-4-6) 3 V modulation 80% 1 kHz
- RF emissions, CISPR 11: Class B
- Harmonics emissions, EN 61000-3-2: Class A

Rf field immunity (EN 61000-4-3):

Field (V/m)	Frequency	Modulation
3	80MHz 2700MHz	1kHz AM 80%
27	380MHz 390MHz	18Hz PM 50%
28	430MHz 470MHz	18Hz PM 50%
9	704MHz 787MHz	217Hz PM 50%
28	800MHz 960MHz	18Hz PM 50%
28	1700MHz 1990MHz	217Hz PM 50%
28	2400MHz 2570MHz	217Hz PM 50%
9	5100MHz 5800MHz	217Hz PM 50%

Warnings:

Although compliant with the EN 60601-1-2 standard, the **MAXIASPEED** medical device may interfere with other devices in the vicinity. The device must not be used in proximity to or stacked on top of other equipment. Install the device well away from other equipment that emits high frequencies (short waves, microwaves, electric scalpels, cell phones).

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are under control. The customer or user can help prevent electromagnetic interference by maintaining a minimum distances between mobile and portable RF communication equipment (transmitters) and the medical device as recommended below, according to the maximum output power of the radio communication equipment.

Rated maximum output power of transmitter (W)	Separation distance (m) in relation to transmitter frequency		
	from 150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	from 80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	from 800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with rated maximum output power not listed above, the recommended separation distance d in metres (m) may be determined using the equation adopted for the transmitter frequency, where P is the maximum rated output power of the transmitter in Watts (W) stated by the transmitter manufacturer.

Notes:

- (1) At 80 MHz and 800 MHz the highest frequency range applies.
 - (2) These guidelines might not apply in all situations.
- Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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DISPOSAL PROCEDURE (Dir. 2012/19/Ue-WEEE) The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment. At the end of life of the device, do not dispose it as mixed solid municipal waste, but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions. This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same. **Caution!** The wrong disposal of electric and electronic equipment may involve sanctions.

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3A HEALTH CARE S.r.l.