**DINAMOMETRO SMEDLEY PER LA MANO -
alluminio
SMEDLEY HAND DYNAMOMETER -
aluminium
DYNAMOMÈTRE MANUEL SMEDLEY -
aluminium
DINAMÓMETRO DE MANO SMEDLEY-
aluminio**

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ENGLISH

About Smedley hand dynamometer

The Smedley spring-type dynamometer has been used successfully for many years to measure grip strength. Adjustable handle comfortably accommodates any hand size. Unit is calibrated in both pounds and kilograms.

Force indication pointer remains at maximal force until reset. To reset, simply rotate pointer counter-clockwise until it rests against the zero pin.

Force Range: Measures from zero up to 220 pounds (100 kilograms).

Use the child size (110 lb. maximal force reading) for smaller forces and higher resolution. Resolution and Precision Measurement gradations are every five pounds (within black outer ring) and are every 2 kilograms (within white inner ring).

Care: For continued use of the dynamometer for years to come be sure to take good care of it. Keep unit in protective carrying case when not in use to prevent water damage and accidental dropping.

Wipe off any water or liquid that might get on dynamometer to prevent the spring from rusting.

USING THE SMEDLEY HAND DYNAMOMETER

Grip Size To size grip, open latch on right side of dynamometer so grip could rotate freely.

Rotate grip clockwise for larger grip and counter-clockwise for smaller grip. As a grip size reference, measure the distance between grip base and base of frame. Record this measurement.

When appropriate grip size is reached, close latch and administer grip test.

Tests should not be administered while grip can rotate freely. Grip size should always be the same for a single individual. Procedure when administering grip strength tests, face gauge outward, facing away from patient.

Stabilize patient's arm and instruct patient to gradually increase effort until maximal effort is attained. This might take three seconds or so. Record reading and reset dial. Repeat test three times, allowing the patient a brief period to rest between each trial and use.

INDICATION OF USE & PRECAUTIONS

Those devices are designed for professional use by qualified operators with necessary skills. The use of the product involves direct contact with patient's skin; always wear protective gloves, do not use on dull skin, handle carefully to avoid accidental abrasion. Clean and sterilize before each use; use neutral non-aggressive disinfectant solution (PH7), do not use solutions containing chloride ions or high alcohol content. This product is intended to provide easy-to-read indicative measurements.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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