# GIMA

PROFESSIONAL MEDICAL PRODUCTS

## **OCCHIALI INGRADENTI BINOCULAR LOUPES LOUPES BINOCULAIRES LUPAS BINOCULARES LUPAS BINOCULARES ΜΕΓΕΘΥΝΤΙΚΕΣ ΔΙΟΠΤΡΕΣ LOUPES**

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**PRODUCT CARE**

- Do not expose your product to liquid or moisture or to humidity.
- Do not touch the surface of lens with hand or hard objects.
- Never clean the surface of lens with any hard and dirty materials.
- Do not look into the sun or into a powerful light source with this loupe.

***Purpose and features of the product*****Purpose**

SL series Binocular Loupes is designed for all possible medical needs. They are widely accepted by Stomatology, ENT, Gynecology, Cardiology, Brain surgery and Plastic surgery etc..

**Features**

- Multi-coated, high grade glass lenses produce crisp, clear images.
- The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- The angle of the loupes can be adjusted between 0°~40° with ease and fatigue is minimized even after prolonged use.
- 6mm up and down displacement. Pupil distance can be adjusted individually which could be suited for all kind of faces.
- Adapts to any Head band and Headlight of Zumax
- Smart design, convenient for mobile diagnosis.

***Operating Instruction***

First step, please fasten the strap to the frame of glasses.

Put on the loupe and check that the pupil distance is symmetrically set.

Now look at a well it three-dimensional object. Adjust one eyepiece after the other until you see a round image with no shadows and a focused object in each. Then screw in both set screws.

The position of the optics is adjusted over a wide range by means of the hinge mount. The ease of movement of the hinges can be adjusted by means of the hinge screws.

***Maintenance of product*****Maintenance of lens**

- The instrument not in use should be covered and put into the packing box.
- Do not disassemble the optical lens by unskilled or unqualified technicians.

**Cleaning**

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.

**Caution: Be careful to prevent the solvent into the instrument**

## Specifications

Magnification	Working distance (mm)	Field of view dia. (mm)
2.5x	340/420	80/90
3.5x	340/420	55/60

### Environment requirements

Operation	Environment Temperature	+10°C – +40°C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa

**Shipping** Common conveyance

Storage	Environment Temperature	+10°C – +40°C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa

	Product code		Lot number		Medical Device
	Keep in a cool, dry place		Keep away from sunlight		Manufacturer
	Caution: read instructions (warnings) carefully		Consult instructions for use		Medical Device compliant with Regulation (EU) 2017/745
	Date of manufacture		Imported by		Authorized representative in the European community
	Model number				

### GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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