GIMA

PROFESSIONAL MEDICAL PRODUCTS

MONOFILAMENTO MONOFILAMENT MONOFILAMENT EINZELFADEN MONOFILAMENTO MONO-FILAMENTO MONOMHMATIO

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ENGLISH

MONOFILAMENT

Show the filament to the patient and touch it to his/her hand or arm so that he/she knows it does not hurt.

Use the 10 gram filament to test sensation at the indicated sites on each foot as shown. Apply the filament along the perimeter of, and NOT directly on, an ulcer, callus, scar or necrotic tissue.

The monofilament is an instrument for a preliminary diagnosis of sensory neuropathy which can be caused, for example, by diabetes. The foot is one of the parts of the body which are affected by the damage caused by neuropathy as the peripheral sensory perception is compromised; the monofilament can detect a patient's ability to perceive moderate pressures.

The MONOFILAMENT bends at a pressure of 10 grams for at least 50/60 times before bending to a lower force; it is recommended to change the wire after the use on 6/15 patients and allowed to rest for at least 24 hours before being reused. It can work optimally at a temperature of 20/26°C

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GIMA

with humidity 30/50% without deforming. With temperature higher than 33°C and with humidity at 80%, the filament can be deformed and make the diagnosis less accurate. The filament was tested to determine the accuracy of bending to 10 grams and the deformation force has been recognized at +/- 10%. Substitute definitely the filament after 150/200 bendings.

Hold the filament by exercising a moderate pressure perpendicular to the skin and use a smooth motion when testing.

Hold the filament perpendicular to the skin and use a smooth motion when testing. Use a 3-step sequence that includes (1) touch the skin, (2) bend the filament, and (3) lift from the skin (See Figures 1- 3). Do not use rapid movement. The approach, skin contact, and departure of the filament should be approximately 1.5 seconds in duration.

Ask the patient to respond "yes" when the filament is felt. If the patient does not respond when you touch a given point on the foot, continue on to another site.

When you have completed the sequence, REPEAT on the area(s) where the patient did not indicate feeling the filament.

Use the filament in a random sequence.

LOSS OF PROTECTIVE SENSATION AT ANY ONE OF THE EIGHT SITES (shown on the diagram above) INDICATES A FOOT AT HIGH RISK.

It is recommended to disinfect the product after each use on a patient by using a non-alcoholic antibacterial disinfectant.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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