## **LUXAMED® LuxaScope Auris LED**

Deutsch .....	2
English.....	8
Français.....	14
Italiano .....	20
Español.....	26
Polski.....	32
Português .....	38
Suomi .....	44
EMC Guidelines.....	50

**[www.luxamed.de](http://www.luxamed.de)**

**Hersteller:** LUXAMED GmbH & Co. KG \*  
Daniel-Weil-Str.3 \* 89143 Blaubeuren \* Germany  
Telefon: +49(0)7344 92905 0 \* Fax: +49(0)7344 92905 10  
E-Mail: info@luxamed.de \* Internet: www.luxamed.de

- DEUTSCH -

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Recommended separation distances between portable and mobile RF communications equipment and the A1.416.114/A1.416.214/A1.416.314

The A1.416.114/A1.416.214/A1.416.314 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the A1.416.114/A1.416.214/A1.416.314 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters ) and the A1.416.114/A1.416.214/A1.416.314 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF IEC 61000-4-3	3 Vrms 80MHz to 2.5GHz	3V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the SYSTEM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p style="text-align: center;">80MHz to 800MHz</p> <p style="text-align: center;">800MHz to 2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p>

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RF emissions CISPR 11	Class B	The A1.416.114/A1.416.214/A1.416.314 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	

**Table 2:Guidance and manufacturer's declaration – electromagnetic immunity**

The A1.416.114/A1.416.214/A1.416.314 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1.416.114/A1.416.214/A1.416.314 should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Note:  $U_t$  is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the A1.426.114/A1.426.214/A1.426.314

The A1.426.114/A1.426.214/A1.426.314 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the A1.426.114/A1.426.214/A1.426.314 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the A1.426.114/A1.426.214/A1.426.314 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

#### Guidance and manufacturer' declaration - electromagnetic emissions - for all LuxaScope Auris LED A1.416.114, A1.416.214, A1.416.314

**Table 1:Guidance and manufacturer's declaration – electromagnetic emissions**

The A1.416.114/A1.416.214/A1.416.314 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1.416.114/A1.416.214/A1.416.314 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The A1.416.114/A1.416.214/A1.416.314 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause terference in nearby electronic equipment

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- ENGLISH -

**To ensure proper use of this LUXAMED® product, please read this user manual carefully and retain for future reference.**

## LUXAMED® LuxaScope Auris LED

### **Purpose:**

To check for signs of ear disease or abnormal development

### **Description:**

A systematic inspection of the external ear canal, surrounding tissue, ear canal, and tympanic membrane.

**Model:** A1.416.114, A1.416.214, A1.416.314, A1.416.914  
 A1.426.114, A1.426.214, A1.426.314, A1.426.914

### **Accessories:**

1. Otoscope ×1
2. Disposable ear tips, grey,  
     2.5 mm ×10  
     4.0 mm ×10
3. User manual ×1
4. Power Supply GTM41078-0605-USB ×1: only for A1.426.114/214/314
5. Plugs: EU/UK/US ×1: only for Power Supply GTM41078-0605-USB
6. lithium-ion rechargeable battery (10440 350mAh):  
     only for A1.426.114/214/314 (embedded in the handle)

### **Optional Accessories:**

1. USB Cable ×1: only for A1.426.114/214/314
2. Alkali cells(AAA/LR03) ×2: only for A1.416.114/214/314

**Table 3: Guidance and manufacturer's declaration – electromagnetic immunity**

The A1.426.114/A1.426.214/A1.426.314 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1.426.114/A1.426.214/A1.426.314 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the SYSTEM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance</b> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{GHz}$ Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 Vrms 80MHz to 2.5GHz	3V/m	

**Table 2:Guidance and manufacturer's declaration – electromagnetic immunity**

The A1.426.114/A1.426.214/A1.426.314 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1.426.114/A1.426.214/A1.426.314 should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV lines to lines ±2 kV lines to earth	±1 kV lines to lines ±2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) For 0,5 cycle 40% UT (60% dip in UT) For 5 cycle 70% UT (30% dip in UT) For 25 cycle <5% UT (>95% dip in UT) For 5 cycle	<5% UT (>95% dip in UT) For 0,5 cycle 40% UT (60% dip in UT) For 5 cycle 70% UT (30% dip in UT) For 25 cycle <5% UT (>95% dip in UT) For 5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of A1.426.114/A1.426.214/A1.426.314 requires continued operation during power mains interruptions, it is recommended that the A1.426.114/A1.426.214/A1.426.314 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

## 1. LUXAMED® instrument handle

The LUXAMED® instrument handle is only to be used with LUXAMED® instruments. These instruments are only to be used by authorized, trained personnel. Use the device in accordance with its designated use and for its intended purpose only.

### Power supply:

2 Alkali cells (AAA/LR03) or the LUXAMED® instrument handle can be operated by the LUXAMED® loading cover with integrated charging technology and a LUXAMED® rechargeable battery. Charging time:

approx. 4 h.

The LUXAMED® loading cover signals, via state of charge indicator, that the battery is completely charged (LED extinguishes).

NOTE: with rechargeable battery fully charge before first startup!

When using the charging cable to a computer it must meet the requirements of IEC 60601-1.

### Instrument connection:

Connecting port with groove (1).

### Switch on:

Press the push button (2) once.

### Adjust the brightness:

Adjust the brightness upwards: while the power is on, press and hold the push button (2) until the desired brightness is achieved.

Adjust the brightness downwards: press and hold the push button (2) again until the desired brightness is achieved.


### Switch off:

While the power is on, press the push button (2) once. In addition, this otoscope has an automatic switchoff function after 3 minutes.

### Changing batteries:

Screw off the battery cap (3) and remove the batteries. Insert the new batteries (AAA/LR03) with the positive pole (+) in direction of the instrument head. Re-screw the battery cap.

Never change batteries next to patient.

### Changing rechargeable battery:

Rechargeable battery is only allowed to be replaced by qualified and authorized

personnel.

In the event of a defective rechargeable battery, please return the handle to LUXAMED® or an authorized specialist. Replacement by inadequately trained personnel could result in a HAZARD (such as excessive temperatures, fire or explosion)

**Cleaning:**

The instrument handle can be cleaned with a damp cloth (use pH-neutral or slightly alkaline cleaning materials).

**Disinfection:**

Disinfection or disinfecting cleaning of the outer surfaces of the instrument handle simply by wipe disinfection. Use disinfectants, approved for use on medical devices made of plastic. Spray- and immersion-type disinfection is not approved.

**Sterilisation:**

Not approved for sterilisation.

**Disposal:**

This product must be disposed of at a separate collection of electric and electronic equipment.

**Storage and transport conditions:**

Ambient temperature between -10° C and +45° C.  
Rel. humidity between 30% and 90%.

**Technical Specifications:**

Nominal voltage 2.5V or 3.7V  
Nominal current 23mA  
Application part BF  
Protection class Class II  
Ingress protection IP20  
Environment usage: 5° C -40° C

10-75% humidity  
700-1070hPa

Power supply 2×AAA Alkaline batteries (2.5V otoscope)  
1×3.7V lithium secondary battery (3.7V otoscope)

**Warning:**

1. Please make sure that the otoscopes good contact with the power adapter while charging;
2. Strictly prohibit the use of this product when charging or unattended, users should unplug the USB plug before leaving.
3. When cleaning/disinfecting devices that have been used on patients, be sure to use appropriate
4. Personal protective equipment and follow all safety precautions/warnings.
5. Always remove batteries prior to cleaning or disinfecting any device (not rechargeable

**Table 1: Guidance and manufacturer's declaration – electromagnetic emissions**

The A1.426.114/A1.426.214/A1.426.314 is intended for use in the electromagnetic environment specified below. The customer or the user of the A1.426.114/A1.426.214/A1.426.314 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The A1.426.114/A1.426.214/A1.426.314 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The A1.426.114/A1.426.214/A1.426.314 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Compliance	

# EMC Guidelines

Below cables information are provided for EMC reference.

Cable	Max. cable length, Shielded/unshielded		Number	Cable classification
Charger USB cable to LuxaScope Auris LED (only for A1.426.114, A1.426.214, A1.426.314)	1m	unshielded	1 Set	DC Power

## Important information regarding Electro Magnetic Compatibility (EMC)

LuxaScope Auris LED needs special precautions regarding EMC and put into service according to the EMC information provided in the user manual and other documents; LuxaScope Auris LED conforms to this EN 60601-1-2:2007+AC:2010/IEC 60601-1-2:2007 standard for both immunity and emissions. Nevertheless, special precautions need to be observed:

- The LUXAMED® LuxaScope Auris LED is only to be used with a medical approved charger, plug and USB cable;
- The use of accessories and cable other than those specified, with the exception of accessories and cables sold by LUXAMED® of LuxaScope Auris LED as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY or decreased LIFESPAN of the LuxaScope Auris LED.
- LuxaScope Auris LED should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, The LuxaScope Auris LED should be observed to verify normal operation in the configuration in which it will be used.

## Guidance and manufacturer' declaration - electromagnetic emissions - for LuxaScope Auris LED A1.426.114, A1.426.214, A1.426.314

batteries). Make sure that no liquid penetrates into the housing.  
6. Use of any accessories or materials not indicated in the user's manual can degrade the minimum safety of the equipment.

7. During use metal components near instrument head can become warm. This is especially true if device is on for extended periods of time. Do not leave device on when not in use.

8. Oscopes are not serviced or maintained while in use with the patient.

9. Regarding any significant RISKS of reciprocal interference posed by the presence of the otoscopes during specific investigations or treatments.

10. Regarding potential electromagnetic or other interference between the otoscopes and other devices advice on ways to avoid or minimize such interference.

11. It is normal situation that the product state is from standby to working, which caused by the electrostatic interference in use. As long as manual off or auto off after 3 minutes, it can be back to standby.

12. In the diagnostic process, if patients have a burning sensation, user need to stop diagnosing temporarily and turn off the otoscope LED immediately.

13. Contraindication:

-Avoid direct contact with mucosa or injured shin.

-Do not shine into the eyes directly.

## 2. LUXAMED® LuxaScope Auris LED

### Purpose:

LUXAMED® LuxaScope Auris LED otoscopes are to be used for the examination of the outer ear canal.

Furthermore, the instrument can be used for other illumination purposes with non-invasive examinations.

### Warning:

Please make sure that, prior to any use, the ear tips are fixed correctly. NEVER use without ear tips.

During examination, ear tip should not remain more than 1 minute in the ear (skin contact). Inflation for pneumatic test must be done very carefully (risk of excess pressure and any resulting injuries).

Use the otoscopes only with LUXAMED® spare parts and accessories.

Never use for eye examination.

Please note that the slot for fixing the ear tips is sharp edged.

### Application:

Insert connector with guide balls (5) on the instrument handle.


### Fixing ear tips:

LUXAMED® LuxaScope Auris LED otoscopes are only to be inserted into the ear canal with fixed ear tips (grey).

For this, the ear tip can be fixed with the locking pin by slightly turning it clockwise

through the slot upon the tip of the otoscope (4).  
For any other illumination purpose no ear tip has to be used, as long as there is no direct contact with skin or mucosa.

The magnifying glass has an approx. 3-fold magnification.  
The glass can be pivoted in clockwise direction.

**Application of pneumatic test:**

The LUXAMED® LuxaScope Auris LED otoscope has a connection option for the insufflation bulb for pneumatic test (6).

**Cleaning:**

Cleaning with a soft, lint-free cloth, lightly moistened with alcohol.

Do not immerse in water or other liquids.

Ear tips are for single use only. Reprocessing disposable items may impair the functionality and safety of the product.

**Sterilisation:**

Not approved for sterilisation.

**Disposal:**

This product must be disposed of at a separate collection of electric and electronic equipment.

**Storage, transport and operating conditions:**

Ambient temperature between -10° C and +45° C.

Rel. humidity between 30% and 90%.

**To Obtain Warranty Service:**

Send item(s) postage paid to LUXAMED® , Attn: Repair Dept., Daniel-Weil-Str.3 \* 89143 Blaubeuren \* Germany

Please include your name and address, phone no., proof of purchase, and a brief note explaining the problem.

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LUXAMED GmbH & Co. KG - Daniel-Weil-Str.3 - 89143 Blaubeuren - Alemania