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Finally the blood pressure and pulse rate will be calculated and displayed on the LCD screen. The blood pressure classification indicator and irregular heartbeat symbol (if any) will blink on the screen. The result will be automatically stored in the monitor.

d. After measurement, the monitor will turn off automatically after 1 minute of no operation.

e. During measurement, you can press the "START/STOP" button to turn off the monitor manually.


Note: Please consult a health care professional for interpretation of pressure measurements.

7. Displaying stored results

- In StandBy Mode, press "MEM" button, the monitor will display sign of current group. The amount of results in current user memory zone will be displayed. See picture 7. Press "START/STOP" button to switch group, press "MEM" to confirm current group. Then LCD will display the average value of all results in the current user memory zone. See picture 7-1. If no result stored in the current user memory zone, LCD will display "0" for blood pressure and pulse rate. See picture 7-2.
- Press "MEM" button, LCD will display the average value of all the results which is measured from 5 o'clock to 9 o'clock in last 7 days in the current user memory zone. See picture 7-3. If no result stored from 5 o'clock to 9 o'clock in last 7 days, LCD will display "0" for blood pressure and pulse rate. See picture 7-4.


20

See picture 7-4.

21

- c. Press "MEM" button again, LCD will display the average value of all the results which is measured from 18 o'clock to 20 o'clock in last 7 days in the current user memory zone. See picture 7-5. If no result stored from 18 o'clock to 20 o'clock in last 7 days, LCD will display "0" for blood pressure and pulse rate. See picture 7-6.

- d. Press "MEM" button again, the most recent result will be displayed with date

and time stamp. See picture 7-7. Irregular heartbeat symbol (if any) and blood pressure classification indicator will blink at the same time. If the monitor has no result stored in the current user memory zone, the LCD will display "0" for blood pressure and pulse rate. See picture 7-8.

e. Press "MEM" button again to review the next result. See picture 7-9. In this way, repeatedly pressing the "MEM" button displays the respective results measured previously.

f. When reviewing the results, the monitor will turn off automatically after 1 minute of no operation. You can also press the "START/STOP" button to turn off the monitor manually.

Note: When the monitor displaying the measurement, the classification color indicator can be shown different color according to the systolic pressure and diastolic pressure. Refer to the "ASSESSING HIGH BLOOD PRESSURE FOR ADULTS" section.

8. Deleting measurements from the memory


When any result is displaying, keeping on pressing button "MEM" for three seconds, all results will be deleted. Press the button "START/STOP", the monitor will turn off.

9. Assessing high blood pressure for adults

The following guidelines for assessing high blood pressure (without regard to age or gender) have been established by the World Health Organization (WHO). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration. Consult with your physician for accurate assessment, and never change your treatment by yourself.

Classification of blood pressure for adults

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the color scheme and that the color scheme is meant only to discriminate between the different levels of blood pressure.

10. Technical alarm description

The monitor will show 'HI' or 'LO' as technical alarm on LCD with no delay if the determined blood pressure (systolic or diastolic) is outside the rated range specified in part SPECIFICACIONES. In this case, you should consult a physician or check if your operation violated the instructions.

The technical alarm condition (outside the rated range) is preset in the factory and cannot be adjusted or inactivated. This alarm condition is assigned as low priority according to IEC 60601-1-8.

The technical alarm is non-latching and need no reset. The signal displayed on LCD will disappear automatically after about 8 seconds.

11. Troubleshooting (1)

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD Display shows abnormal result	The cuff position was not correct or it was not properly tightened	Apply the cuff correctly and try again
	Body posture was not correct during testing	Review the "BODY POSTURE DURING MEASUREMENT" sections of the instructions and re-test
	Speaking, arm or body movement, angry, excited or nervous during testing	Re-test when calm and without speaking or moving during the test
	Irregular heartbeat (arrhythmia)	It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer

12. Troubleshooting (2)

PROBLEM	POSSIBLE CAUSE	SOLUTION
LCD shows battery symbol	Low Battery	Change the batteries
LCD shows "Er 0"	Pressure system is unstable before measurement	
LCD shows "Er 1"	Fail to detect systolic pressure	Don't move and try again
LCD shows "Er 2"	Fail to detect diastolic pressure	
LCD shows "Er 3"	Pneumatic system blocked or cuff is too tight during inflation	Apply the cuff correctly and try again. If the monitor is still abnormal, please contact the local distributor or the factory.
LCD shows "Er 4"	Pneumatic system leakage or cuff is too loose during inflation	Apply the cuff correctly and try again. If the monitor is still abnormal, please contact the local distributor or the factory.

25

Table 2
Enclosure Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

26

LCD shows "Er 5"	Cuff pressure above 300mmHg	
LCD shows "Er 6"	More than 3 minutes with cuff pressure above 15 mmHg	Measure again after five minutes. If the monitor is still abnormal, please contact the local distributor or the factory
LCD shows "Er 7"	Innes memory error	
LCD shows "Er 8"	Device parameter checking error	
LCD shows "Er A"	Pressure sensor parameter error	
No response when you press button or load battery	Incorrect operation or strong electromagnetic interference	Take out batteries for five minutes, and then reinstall all batteries

MAINTENANCE

- ⚠ Do not drop this monitor or subject it to strong impact.
- ⚠ Avoid high temperature and solarization. Do not immerse the monitor in water as this will result in damage to the monitor.

27

- If this monitor is stored near freezing, allow it to acclimate to room temperature before use.
- ⚠ Do not attempt to disassemble this monitor.
- It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.
- Clean the monitor with a dry, soft cloth or a soft cloth squeezed well after moistened with water, diluted disinfectant alcohol, or diluted detergent.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user's appropriately qualified technical personnel to repair those parts of equipment which are designated repairable can be supplied.
- The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years, and the cuff integrity is maintained after 1,000 open-close cycles of the closure.
- The monitor requires 6 hours to warm from the minimum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20°C.
- The monitor requires 6 hours to cool from the maximum storage temperature between uses until the monitor is ready for its INTENDED USE when the ambient temperature is 20°C.
- Not servicing/maintenance while the monitor is in use.

28

- 12 It is recommended the cuff should be disinfected 2 times every week if needed (For example, in hospital or clinic). Wipe the inner side (the side contacting the skin) of the cuff with a soft cloth squeezed after being moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.

EXPLANATION OF SYMBOLS ON UNIT

	Follow instructions for use
	Caution: read instructions (warnings) carefully
	Type BF applied part
	WEEE disposal
	Medical Device complies with Directive 93/42/EEC

29

	Manufacturer
	Date of manufacture
	Authorized representative in the European community
	Serial number
	Covering Protection rate
	Keep in a cool, dry place
	Keep away from sunlight
	Product code
	Lot number

30

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1
Emission

Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

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