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LEO BLOOD PRESSURE MONITOR

Use and maintenance book


ATTENTION: Operators must read and understand this manual completely before using the product.

REF CONTEC08A (GIMA 32902)

CONTEC MEDICAL SYSTEMS CO., LTD No.112 Qianhuang West Street, Economic & Technical Development Zone, Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA
Made in China

EC REP Proline GmbH, Bremerstr. 56, 40239 Dusseldorf Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaltaly.com - export@gimaltaly.com
www.gimaltaly.com

Chapter1 Functions and Purpose

1.1 Description of functions:

The Sphygmomanometer apply to measure the non-invasive blood pressure and SpO₂ of human (adult, children, neonate), uses three-user mode, each user can store 100 items records of measurement results at most. Each record includes detailed measuring time, systolic pressure, diastolic pressure, average pressure, pulse rate and record number, etc. With 2.8 inch color LCD screen, clear interface, the function of data review is complete. User can implement ON/OFF, manual measuring, system setup, parameters change and other operations with seven buttons which are located on the front panel of the device.


The Sphygmomanometer uses audible and visual prompt, when the battery power is low, the buzzer will intermittent buzzing and LCD screen displays "Low Power" to prompt user replace batteries. When the measurement data exceeds the set prompt limit, the font color of measurement results will change to red and the audible prompt will occur, user can open and close the prompt sound according to needs. With timing shutdown function, if there is no operation and SpO₂ measurement, the device will automatically turn off after 2 minutes. With USB interface, Users can send measurement results to PC. Refer to the help or explication of the related software for specific operation.

1.2 Purpose:

In Europe, The Electronic Sphygmomanometer is applied to Blood Pressure (BP) measure and monitor for adult, pediatric, and neonatal. It has an optional SpO₂ Measurement Function. The product is suitable for being used in family, hospital, oxygen bar, community healthcare, physical care in sports (It can be used before or after doing sports, and it is not recommended to use the device during the process of having sport) and etc. In other areas, The Electronic Sphygmomanometer is intended to measure the systolic, diastolic and mean blood pressure as well as pulse rate via non-invasive oscillometric technique in which an inflatable cuff is wrapped around the upper arm. It can be used on adult, pediatric and neonatal individuals.

Chapter2 Main Unit

The production is in the package. Open the package and confirm whether the production is whole.

Accessories:

Specification: limb circumference 22-32 cm (middle part of upper arm), please select suited cuff when measuring children or other limb circumference.

Optional Accessories:

AC adapter

Input: voltage: AC 100 V~240 V frequency: 50 Hz/60 Hz Rated current: AC 150 mA
Output: DC5.0 V±0.2 V 1.0 A

SpO₂ probe: Integrated SpO₂ probe (This part is only suitable for European Union market)

A. SpO₂ measurement

Range:0%~100% Error: 70~100%:±2%;Below 70%:unspecified Resolution: 1%

Note: because SpO₂ probe measurements are statistically distributed, only about two-thirds of SpO₂ probe measurements can be expected to fall within ±Arms of the value measured by a CO-OXIMETER.

B. Pulse rate measurement

Range:30 bpm~250 bpm Error: ±2 bpm or ±2% (select the larger) Resolution: 1bpm

C. Optical sensor: red light(wavelength: 660 nm, output power less than 6.65 mW) infrared light(wavelength: 880 nm, output power less than 6.75 mW). Optical sensors are light-emitting components that affect other medical devices that use this wavelength range. This information may be useful to clinicians who perform optical therapy.

D. Error in weak filling condition: SpO₂ and pulse rate can be shown correctly when pulse-filling ratio is 0.4 %. SpO₂ error is ±4%; when measuring range is 30 bpm~100 bpm, pulse rate error is ±2 bpm; when measuring range is 100 bpm~250 bpm, pulse rate error is ±2%.

Note:

◎ The optional probe of the Sphygmomanometer is an integrated SpO₂ probe, the measuring part is integrated with the probe;

◎ The service life of the integrated SpO₂ probe is three years.

Cuff:

Choose the right cuff based on the patient's upper arm circumference, there are several suitable cuffs(range of limb circumference, middle part of upper arm)

the range of limb circumference is 6-11 cm

the range of limb circumference is 10-19 cm

the range of limb circumference is 18-26 cm

the range of limb circumference is 22-30 cm

the range of limb circumference is 22-43 cm

the range of limb circumference is 32-43 cm

△ Note △

● The cuff is a consumable. Calculate by measuring 6 times a day(3 times each morning and evening), the service life of the cuff is about 1 year.(using our experimental conditions)

● In order to correctly measure blood pressure, please replace the cuff in time.

● If the cuff leaks, please contact our company to buy a new one. The cuff purchased separately does not include the airway tube plug. When replacing, please do not throw the airway tube plug away, install it on the new cuff.


△ Note △

When the product and accessories described in this manual are about to exceed the period of use, they must be disposed according to relevant product handling specification. If you want to know more information, please contact our company or representative organization.

Chapter3 External Interfaces

△ Note △

When removing NIBP cuff, please take plug at the front of the windpipe to pull out.

△ Note △

All analog and digital equipment connected to this device must be certified to IEC standards(such as IEC60950: Information technology equipment-Safety and IEC60601-1: Medical electrical equipment-Safety), and all equipment should be connected to in accordance with the requirement of the valid version of the IEC60601-1 system standard. The person connecting the additional equipment to the signal input and output port is responsible for whether the system complies with the IEC60601-1 standard.

Chapter4 Battery/AC Adapter Installation

The product can use battery or AC adapter as power source.

4.1 Battery Installation

- ① Demount the battery cover in the direction of the arrow.
- ② Install "AA" batteries according to polarities.
- ③ Slide to close the battery cover.

Icon : the battery power will exhaust. Replace with four new batteries (the same sort) at the same time. Test while low power may cause data deviation and other problems.

Turn the unit off before replacing the batteries.

△ Note △

When the battery reaches the end of its life, or if the battery is found to have odor, deformation, discoloration or distortion, stop using the battery and dispose of the used battery in accordance with local regulations, otherwise it will cause environmental pollution.

4.2 Usage of power adapter

1. Connect the sphygmomanometer and the power adapter. Insert the power adapter connector into the power adapter socket on the device.
2. Please insert the power plug of the adapter into the AC 100 V~240 V socket.

△ Note △

The device can be disconnected from the power supply network by unplugging the adapter plug.

When cut off the power supply, first cut off the connection of power socket and the regulated power supply, then cut off the connection of regulated power supply and the sphygmomanometer.

Please be sure to use dedicated medical grade power adapter.

△ Note △

When regulated power supply and batteries are both used at the same time, the battery power will not be consumed.

Switch regulated power supply and battery as power supply when the device is off, otherwise, the device may shutdown due to power failure.

The device can be used normally after it is turned on, without waiting for the device to be prepared.

Chapter5 Button Functions

All the operations to the Electronic Sphygmomanometer are through the buttons. The names of the buttons are above them. They are:

- [ON/OFF] ON/OFF button. Press this button to turn on/off the device.
- [START/STOP] Press it to inflate the cuff and start a blood pressure measurement. During measuring, press it to cancel the measurement and deflate the cuff.
- At all levels interface, the three buttons correspond respectively with the text prompts below the LCD screen, pressing any button will carry on the corresponding function, such as [UP][MENU][ENTER][DOWN] etc.

Chapter6 Setting the Date and Time

It is necessary to set date and time after power on.

The Electronic Sphygmomanometer can automatically stores measurement results with date and time.

If batteries power exhausts or removed, the time to stop.

At the moment, please reset date and time.

The Electronic Sphygmomanometer stores measure results of three users automatically, and up to 100 items for every user. If the date and time are set correctly, the date and time when measuring will be correct in the memory, otherwise it may not be correct. The results can be uploaded to PC via USB and processed with the PC software.

1. There are two modes of time setting:

- (1) When using the Sphygmomanometer for the first time or after the Sphygmomanometer has been placed without power supply for a certain time(more than 3 minutes), after power on, there is a prompt of time error on the main interface, set date and time with , and .
- (2) Press button on the main interface to enter system menu, then enter item, the current time will be displayed on the screen. Set date and time with , and .

2. After setting, select option and press button to confirm the setting value. If you do not want to change the time, select option and press button to return to the previous menu.

△ Note △


The range of year is from 2010 to 2099. When the year reaches 2099, pressing the button will return to 2010.

Chapter7 About Unit

There are two units: "mmHg" and "kPa".

The default is: "mmHg".


Enter the submenu in the , then select option to switch units between "mmHg" and "kPa".

Chapter8 User Switch

The Electronic Sphygmomanometer stores the measure results of three users automatically, and up to 100 items for every user.


Press button in main interface to switch users. Or press item in menu to switch users.

△ Note △

When the is set to , current user can be switched under main interface; when set to a certain user, it will not be able to switch under main interface.

User type can be set to adult, pediatric and neonatal three different kinds, setting method is as follows:

Chapter9 Over-limit Prompt Function

The Sphygmomanometer has two kinds of reminding methods:the technical parameter over-limit prompt and the physiological parameter over-limit prompt.

9.1 Physiological parameter over-limit prompt

The sphygmomanometer has the function of over-limit prompt, the user can press button to enter system menu, select option to enter its interface, then set the limit value of blood pressure. When the BP measurement result is higher than the high limit or lower than the low limit and the prompt is ON, the physiological prompt will occur; in interface, select option to enter its interface, when the SpO₂ measurement result is higher than the high limit or lower than the low limit and the prompt is ON, the physiological prompt will occur.

In the state of physiological prompt, press any button to cancel the prompt and it does not affect the next prompt; the prompt can be closed permanently with prompt switch of the prompt setup menu until the prompt switch be opened again.

9.2 Technical parameter over-limit prompt

When power is about to exhaust and prompt is ON, the prompt will occur. This prompt can not be canceled unless being closed or the power replaced.

△ Note △ Restore factory settings

If users want to restore the factory settings,select ,then the data will be restored to the initial state.

Chapter10 The Usage Method of Sphygmomanometer

10.1 Accurate Measurement Way

Measurement in quiet and relaxing state.


1. Adopt a comfortable sitting position, use back and arms to support the body.

2. Place your elbow on a table, the palm faces up and the body is relaxed.

3. The cuff is level with your heart.

4. Feet flat on the floor, and do not cross your legs.

△ Advice △

Try to measure your blood pressure at the same time each day with the same arm and the same pose for

consistency.

The high and low location of cuff will cause changes in measure results.

Do not touch the sphygmomanometer, cuff and windpipe during measure.

Measure

Chapter11 Memory Function

The sphygmomanometer is designed to store the blood pressure, pulse rate values and the date and time when measured, which are up to 100 groups. If there have been stored 100 groups, the earliest results will be deleted when saving the 101 group of measurement results.

11.1 Review Memory Values

1. In the main interface (interface when starting-up), press **【MEMORY】** button to review the latest measured values in big-font with the serial number from 1 to 100.

2. Press **【UP】 / 【DOWN】** button to circularly switch the former measurement values.

*The right figure shows that there is no measurement result can be displayed.

3. Press **【LIST】** button to display the data list interface.

End to display the measurement values:

Press **【EXIT】** button to return to the main interface or hold **【ON/OFF】** button to turn the device off.

11.2 Delete Memory Values

Users can delete all memory values of a user instead of separately delete one memory value.

1. Press **【MENU】** button to enter the system menu, select **【DELETE DATA】** option to enter its interface, select the user whose data to be deleted, after confirming again, all measurement results of the selected user will be deleted.

2. Finish Operation

Select **【EXIT】** to return to the previous menu, or hold **【ON/OFF】** button to turn the device off.

Chapter12 SpO₂: Measurement Function

(This chapter is only suitable for European Union market)

Precautions during SpO₂: Measurement:

⚠ Note ⚠

Make sure the nail covers the light. The probe cable should be on the backside of the hand. Improper probe placement or improper contact with the test site will influence the measurement.

SpO₂ value always displays in the fixed place.

©The test site shall not use external coloring agent (such as nail polish, dyestuff or color skin care products, etc.), otherwise it will affect the measurement.

©The fingers which are too cold or too thin may affect the measure accuracy, please insert the thicker finger such as thumb or middle finger deeply enough into the probe.

©The SpO₂ probe is suitable for children and adults (not suitable for infants and newborns). The device may not applicable for all patients. If you are unable to achieve stable readings, stop using it.

©Data averaging and signal processing will delay the SpO₂ displaying and data values transmitting. The update time of measurement data is less than 30 seconds, when signal attenuation, weak perfusion or other interference appears, it will result in time increasing of dynamic data averaging, which depends on the PR value.

©The PLETH waveforms are not normalized, which is used as the indicator for signal incompleteness. So the accuracy of the measured values may decrease when the waveform does not tend to smooth and stable. When the waveform tends to smooth and stable, the reading is optimal value, and the waveform at the moment is the most standard one.

©The temperature for the contact surface of the device with the body is less than 41°C, and this temperature value is measured by a temperature measuring device.

©The device does not provide over-limit alarm function, so it is inapplicable for using in places where need such function.

©The SpO₂ probe has been calibrated before leaving factory. It does not need to be calibrated during maintenance.

©SpO₂ probe is calibrated to show functional oxygen saturation.

©The SpO₂ probe and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between. Make sure the optical path is free from any optical obstacles like rubberized fabric, to avoid inaccurate measurement.

©As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.

©Percentage modulation of infrared signal as the indication of pulsating signal strength, patient simulator has been used to verify its accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions, compare them with the known SpO₂ and PR values of input signal.

©The claim for SpO₂ accuracy should be supported by clinical research measurements covering the entire spectrum. By artificially inducing to different stable oxygen levels, make it in the range of 70% ~ 100 % of SpO₂. Use secondary standard SpO₂ measuring equipment for comparison to collect SpO₂ values together with the tested product, compose paired data groups for accuracy analysis.

©Clinical report records data of 12 healthy volunteers, including 6 females and 6 males. Volunteers' age ranges from 21 to 29. Skin color is distributed from dark to light, including 3 dark black skins, 2 medium black skins, 5 light skins, 2 white skins.

©When using the device, please keep it away from the instruments that can generate strong electric or magnetic field. Use the device in an inappropriate environment may cause interference to the surrounding radio equipment or affect its working.

©If necessary, please log on our company's official website to download the list of SpO₂ probes and extension cords that can be used in conjunction with this device.

⚠ Warning ⚠

©Check if the cable of SpO₂ probe is in normal condition before measuring. After unplugging the SpO₂ probe cable from the socket, "SpO₂%" and "bpm" on the screen will disappear.

©Do not use the SpO₂ probe once the package or the probe is found damaged. Instead, you shall return it to the vendor.

©The supplied SpO₂ probe is only suitable for use with this device. This device can only use the SpO₂ probe described in this manual. It is the responsibility of the operator to check the compatibility of the device and the SpO₂ probe (and extension cable) before use. Incompatible accessories may result in device performance decreasing or cause injury to the patient.

©SpO₂ probe is a medical product that can be used repeatedly.

©The measured value may be normal seemingly for the testee who has anemia or dysfunctional hemoglobin (such as carboxyhaemoglobin (COHb), methaemoglobin (MetHb) and sulphaemoglobin (SuHb)), but the testee may appear hypoxia, it is recommended to perform further assessment according the clinical situations and symptoms.

©Pulse oxygen has only reference significance for anemia and toxic hypoxia, as some patients with severe anemia still show better pulse oxygen measurements.

©Measurement accuracy can be affected by the interference of electrosurgical equipment.

©Do not install the SpO₂ probe on an extremity with arterial catheter or receiving intravenous injection.

©Do not perform SpO₂ measuring and NIBP measuring on same limb at one time, because obstruction of blood flow during NIBP measuring may adversely affect the reading of SpO₂ value.

©Excessive movement (active or passive) of the subject or severe activity can affect the measurement accuracy.

©Excessive ambient light may affect the measuring results, such as surgical light (especially xenon light sources),

bilirubin lamp, fluorescent lamp, infrared heater and direct sunlight, etc. In order to prevent interference from ambient light, make sure to place the probe properly and cover the probe with opaque material.

©The measured value may be inaccurate during defibrillation and in a short period after defibrillation, as the SpO₂ probe not have defibrillation-proof function.

©The person who is allergic to silicone, PVC, TPU, TPE or ABS can not use this device.

©For some special patients, it should be a more prudent inspecting in the measurement part. The probe can not be clipped on the edema and tender tissue.

©Do not stare at the luminescent component directly when the device is turned on (infrared light is invisible), even if for maintenance purpose, or it may have bad influence to the eyes.

©Uncomfortable or painful feeling may appear if using the SpO₂ probe ceaselessly, especially for the microcirculation barrier patients. It is recommended that the measurement should not be taken at the same position for over 2 hours. Continuous, long measurements may increase the risk of unwanted changes in the skin characteristics, such as exceptionally sensitive, reddish, blistering or oppressive necrosis, especially for neonates or the patients with perfusion disorder and change or immature skin form. It should be paid special attention to check the placement position of the probe according to the skin quality change, correct optical alignment and attachment method. Check the attachment position periodically and change it when the skin quality decreases. A more frequent check may be required due to the difference of patient's state.

©Some models of functional tester or patient simulator can measure the accuracy of the device that reproduces the calibration curve, but it can not be used to evaluate the accuracy of this device.

©Please refer to related medical literature for detailed clinical restrictions and contraindications,

©This device is not used for treatment purpose.

©Do not use the SpO₂ probe during MRI and CT scanning, as the induced current may cause burns.

©When the device is ON, if power interrupt for more than 30s, the SpO₂ probe needs no operation after the power is restored, after the device is turned on, ensure that the SpO₂ probe can be used normally.

©The probe can be used before/after sport, but not recommended to use during exercising.

©Patient simulator has been used to verify the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator.

©Do not use the SpO₂ probe during MRI and CT scanning, as the induced current may cause burns.

©When the device is ON, if power interrupt for more than 30s, the SpO₂ probe needs no operation after the power is restored, after the device is turned on, ensure that the SpO₂ probe can be used normally.

©The probe can be used before/after sport, but not recommended to use during exercising.

©Patient simulator has been used to verify the pulse rate accuracy, it is stated as the root-mean-square difference between the PR measurement value and the value set by simulator.

Chapter13 SpO₂: Measurement Method

(This chapter is only suitable for European Union market)

1) Attach the SpO₂ probe to the appropriate site of the patient finger as the following figure.

Place SpO₂ probe

2) Plug the connector of the SpO₂ probe cable into the USB socket in the lower right of the device. The main interface will switch to SpO₂ interface. This operation brings no affection to other functions. If no finger inserted into the probe or finger placement is improper, the SpO₂ interface displays finger-out prompt; when the probe cable is unplugged from the USB socket, the SpO₂ interface will display probe-off prompt, a few seconds later, it automatically return to the main interface.

3) The pulse sound can be set to on or off during measuring the SpO₂. The operation steps are: Press **【MENU】** to enter the system menu, select **【POROPROMT SETUP】** to enter the setting interface, and select **【SpO₂ PROMPT】** to enter the SpO₂ prompt interface, in which you can set the pulse sound to on or off. Note: Please make sure the prompt switch is set to ON state before turning on the pulse sound. (set the prompt switch to on or off in the **【POROPROMT SETUP】** interface)

⚠ Note ⚠

SpO₂ display range: 0% ~ 100%, PR display range: 30 bpm (beats/min) ~ 250 bpm (beats/min)

A SpO₂ probe fault may exist in the following situations, please stop using it and contact our after-sales service:

©After reconnecting the probe cable with the USB port multiple times, the device main interface does not show any response;

©When probe cable is connected normally and the finger properly placed in the probe, the SpO₂ interface still displays finger-out or probe-off prompt;

©The SpO₂ interface prompts probe error;

©Percentage modulation of infrared signal as the indication of pulsating signal strength, patient simulator has been used to verify its accuracy under conditions of low perfusion. SpO₂ and PR values are different due to low signal conditions, compare them with the known SpO₂ and PR values of input signal.

©Measurement limitation

During operation, the accuracy of SpO₂ readings can be affected by:

- High-frequency electromagnetic interference such as interference from electrosurgical apparatus connected to the system.
- Intravenous dyestuff.
- Excessive patient movement.
- External light.
- Improper SpO₂ probe installation or incorrect contact position of the patient.
- Temperature of SpO₂ probe(optimal temperature range: 28°C ~ 40°C).
- Place the SpO₂ probe on an limb that has a blood pressure cuff, arterial catheter, or intravascular line.
- Concentrations of dysfunctional hemoglobin, such as carboxyhemoglobin (COHb) and methemoglobin (MetHb).
- SpO₂ is too low, Bad circular perfusion of the part being measured.
- Intravascular staining agents (such as indocyanine green or methylene blue), skin pigmentation.
- It is required to use SpO₂ probe which is provided by our company, contact with our sale department when necessary.

Chapter14 Error Message

Error message will be displayed in the screen if there is something wrong when measuring. The causes and solutions are shown as follows:

Error Message	Causes	Solutions
Self-test failure	Function abnormal	Please contact us
System failure	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 10)
Loose cuff	Cuff plug fall off	Make sure the cuff plug is securely inserted in the windpipe (refer to Chapter 10)
Air leakage	Air pressure error	Refer to the troubleshooting
Weak signal	The pulse signal is too weak or the cuff is loose.	Correctly connect cuff (refer to Chapter 10)
Overpressure	Cuff is blocked or squeezed	Correctly connect cuff (refer to Chapter 10)
Excessive movement	The signal extent is too big owing to the arm or body moving or other reasons when measuring	Keep arm, body still, measure again
Over range	Saturated signal	Keep arm, body still, measure again
Time out	It takes too much time	Keep arm, body still, measure again

Chapter15 Troubleshooting

Abnormal Phenomenons	Causes	Solutions
BP measurement values too high or	Cuff is not connected correctly. Talk or move arms when measuring	Correctly connect cuff (refer to Chapter 10) Keep quiet and restart a measurement

too low.	The turnup clothing presses the arm	Take off the clothing which presses the arm, and restart a measurement
No pressure	Cuff leakage	Buy a new cuff
	The cuff windpipe is not correctly connected with cuff	Correctly connect
	Cuff is not inflated	Stop using the device and contact us
Cuff deflates in short time	Loose cuff	Correctly apply cuff
	It can not carry on measurement when press the measurement button	Switch on the power once again and restart a measurement
	Power off suddenly when inflating	No use for a long time, the power of batteries can be exhausted owing to the changed temperature Replace all four batteries with new ones.
Hold the on/off button but can not start the device	Power of batteries can be exhausted	Replace all four batteries with new ones.
	The battery polarities is reversed	Check the battery installation for proper placement of the battery polarities.
Cuff inflation start before press the measurement button or never	Pull out the cuff to deflate. Stop using the device and contact us.	Pull out the cuff to deflate. Stop using the device and contact us.
Cuff never deflation	Pull out the cuff to deflate. Stop using the device and contact us.	Pull out the cuff to deflate. Stop using the device and contact us.
Air pressure error	No deflation or deflation error or inflation without stop	Pull out the cuff to deflate. Stop

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
REF CONTEC08A (GIMA 32902)

CONTEC MEDICAL SYSTEMS CO., LTD No.112
Qinhuang West Street, Economic & Technical
Development Zone, Qinhuangdao, Hebei Province,
PEOPLES REPUBLIC OF CHINA
Made in China

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