# GIMA

ARTICOLI & APPARECCHI PER MEDICINA

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

[gima@gimaitaly.com](mailto:gima@gimaitaly.com) - [export@gimaitaly.com](mailto:export@gimaitaly.com)

[www.gimaitaly.com](http://www.gimaitaly.com)

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## ECG ACQUISITION SYSTEM DMS 300-BTT02

### User manual

**ATTENTION:** The operators must carefully read and completely understand the present manual before using the product.

**REF** 33320 - 33321

Shenzhen Creative Industry Co., Ltd.  
2/F Block 3 Nanyou Tian'an Industry Town  
518054 Shenzhen, GD - P.R.C.  
Made in China

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestrasse 80, 20537 Hamburg - Germany

0123

## **CONTENTS**

<b>Preface</b> .....	<b>3</b>
General Information .....	3
EMC Declaration .....	4
<b>Chapter 1 Equipment Overview</b> .....	<b>7</b>
1.1 Intended Use .....	7
1.2 Specifications .....	7
1.3 Hardware Connection .....	8
1.4 System Requirements .....	8
1.5 USB security key driver installation .....	9
1.6 Resting 12-Lead ECG Software Installation .....	11
<b>Chapter 2 System Settings</b> .....	<b>14</b>
2.1 Launch the PC-ECG Software .....	14
2.2 Initial Screen .....	14
2.3 System Settings .....	15
<b>Chapter 3 Data Sampling</b> .....	<b>20</b>
3.1 Prepare Patient Skin .....	20
3.2 Start a New Patient Test .....	22
3.3 Set ECG Components .....	27
3.4 Quick ECG Collecting .....	28
3.5 Multi-Hour ECG Monitoring .....	28
<b>Chapter 4 Data Analysis</b> .....	<b>28</b>
4.1 12-lead Simultaneous ECG .....	28
4.2 FCG (Frequency ECG) Analysis .....	33
4.3 High Frequency ECG (HFECG) .....	37
4.4 QT Dispersion .....	38
4.5 Signal Averaged ECG (SAECG) .....	39
4.6 Vector ECG (VCG) .....	43
4.7 Timed Vector ECG (TVCG) .....	46
<b>Chapter 5 Other Functions of the Software</b> .....	<b>47</b>
5.1 Database Manager .....	47
5.2 Reanalysis .....	48
5.3 Cardio Card .....	49
5.4 Report .....	50
5.5 Exercise .....	51
<b>Chapter 6 Report Interpretation</b> .....	<b>55</b>
6.1 12-lead Simultaneous ECG .....	55
6.2 Frequency ECG (FCG) .....	58
6.3 High Frequency ECG .....	59
6.4 QT Dispersion .....	60
6.5 Hour ECG .....	61
6.6 Signal Averaged ECG (SAECG) .....	62
6.7 Vector ECG .....	63
6.8 Timed Vector ECG .....	64
<b>Appendix</b> .....	<b>65</b>

## PREFACE

### General Information

- Welcome to the DMS 300-BTT02 ECG Acquisition System. This system includes DMS 300-BTT02 ECG Acquisition Box and Resting 12-Lead ECG Software.
- In the following, we will introduce all the functions and uses of ECG Acquisition System specifically so that you can get familiar with it and grasp it quickly.
- FAQ are answered in this manual for you to solve all questions more easily.
- If you want to have a quick look at a certain function of ECG Acquisition System, please refer to the index page.
- The product bears the CE marking “CE 0197” indicating its conformity with the provision of the EC Directive 93/42/EEC concerning medical devices and fulfills the essential requirements of Annex V of this directive. It is class IIa (MDD) equipment.
- The product complies with the electromagnetic immunity requirements of standard IEC 60601-1-2/EN 60601-1-2 “Electromagnetic Compatibility- Medical Electrical Equipment”.
- DM Systems (Beijing) Co., Limited quality management system complies with the standard EN ISO 13485.

### Package List

A complete set of ECG Acquisition System should include the following items: (please check everything in the box on the first open and contact us if anything is missing.)

#### Unit

DMS 300-BTT02 ECG Acquisition Box	51201012
-----------------------------------	----------

#### Accessories

Patient Cable	20202520
Software CD (Resting 12-Lead ECG Software)	50101303
Security Key	11902019
Operator's Manual	50207203

### Safety and Effectiveness Considerations

- The following safety and effectiveness issues are to be considered prior to the usage of the ECG Acquisition System.

- Safety protection type of this ECG Acquisition System is DEFIBRILLATION-PROOF TYPE CF APPLIED PART.
- The reliability of the ECG Acquisition System is dependent upon conformance with the operation instructions, as detailed in this manual.
- The ECG Acquisition System will not cause abnormal operation of a patient's cardiac pacemaker or other electronic stimulator.
- Only a physician should give final interpretation.

### Warnings and Contraindications

- This equipment is not intended for use with high frequency surgical units. Disconnect the patient from the device before using the high frequency surgical unit.
- This equipment is not intended for use with electro surgery equipment.
- Do not remove unit covers. The ECG Acquisition System does not contain any user serviceable components.
- The ECG Acquisition System is not suitable for intracardiac application.
- DO NOT immerse the device in any fluid, place fluids on top, or attempt to clean the unit with any liquid detergents. This may cause an electrical hazard.
- The patient cable should not be kneaded or bent. Or the cable will be broken. The connector of the patient cable should not get in touch with the earth line or other bare metal parts.
- The ECG Acquisition System is defibrillation-proof. Therefore it is not necessary to remove the ECG electrodes during defibrillation.
- During defibrillation, do not touch the patient, the electrodes or lead wires. If defibrillation be carried out during ECG sampling, this would cause a signal delay 10 seconds.

- The electrodes should be in accordance with the requirements of ANSI/AAMI EC12 and Directive 93/42/EEC.
- The used electrodes should be disposed of in accordance with the requirements of environment protection.
- The equipment is used in accordance with the instructions given in this manual.
- Keep the recorder and patient cable clean, especially the components that touch patients.

**Equipment Symbols**

The following symbols appear on the equipment.

Table Equipment Symbols	
Symbols	Meaning
	Caution
	DEFIBRILLATION-PROOF TYPE CF APPLIED PART
	Direct Current
	Proper disposal of equipment
	CE Mark. (according to MDD93/42/EEC)
	Authorized Representative in the European Community
	Manufacturer

**EMC Declaration**

The ECG Acquisition System needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents; Portable and mobile RF communications equipment can affect the ECG Acquisition System. All cables and maximum length of cables, and other accessories with which the manufacturer of the ECG Acquisition System claims compliance with the requirements, Accessories that do not affect compliance with the requirements of these sub clauses need not be listed. Accessories and cables may be specified either generically or specifically.

**Warning:**

Transducers and cables sold by the manufacturer of the ECG Acquisition System as replacement parts for internal components need not be listed.

The use of accessories and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of The ECG Acquisition System as replacement parts for internal components, may result in increased emissions or decreased immunity of The ECG Acquisition System.

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The ECG Acquisition System is intended for use in the electromagnetic environment specified below. The customer or the user of the ECG Acquisition System should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The ECG Acquisition System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The ECG Acquisition System is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided the following warning is heeded. Warning: This ECG Acquisition System is intended for use by healthcare professionals only. This equipment/ system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the ECG Acquisition System or shielding the location.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The ECG Acquisition System is intended for use in the electromagnetic environment specified below. The customer or the user of the ECG Acquisition System should assure that it is used in such an environment.			
<b>IMMUNITY test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) for 0,5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ECG Acquisition System requires continued operation during power mains interruptions, it is recommended that the ECG Acquisition System be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Power frequency magnetic fields should be t levels characteristic of a typical commercial or hospital environment.

**Guidance and manufacturer’s declaration – electromagnetic immunity**

The ECG Acquisition System is intended for use in the electromagnetic environment specified below. The customer or the user of the ECG Acquisition System should assure that it is used in such an electromagnetic environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the ECG Acquisition System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	<p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{3} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{3} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**NOTE 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ECG Acquisition System is used exceeds the applicable RF compliance level above, the ECG Acquisition System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ECG Acquisition System.

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### Recommended separation distances between portable and mobile RF communications equipment and the ECG Acquisition System

The ECG Acquisition System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ECG Acquisition System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ECG Acquisition System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0,01	N/A	0,117	0,233
0,1	N/A	0,369	0,738
1	N/A	1,167	2,333
10	N/A	3,689	7,379
100	N/A	11,667	23,333

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Chapter 1 EQUIPMENT OVERVIEW

### 1.1 Intended Use

- The ECG Acquisition System is capable of sampling, recording patient's resting ECG. This system is applicable to the heart disease for the medical treatment institution.
- The ECG Acquisition System is intended for use under the direct supervision of a licensed health care practitioner.
- The device is to be used on an adult population, typically symptomatic.

### 1.2 Specifications

Voltage:	DC 5.0V, 200mA (USB Interface).
Sample rate:	10240Hz.
Sample resolution:	16 bit.
Lead:	12-lead standard ECG.
Input circuit current:	$\leq 0.1\mu\text{A}$ .
Input impedance:	$\geq 5\text{ M}\Omega$ .
Noise electric level:	$\leq 15\mu\text{Vp-p}$ .
Standard sensitivity:	10mm/mV, $\pm 5\%$ .
Anti-polarization voltage:	if add $\pm 300\text{mV}$ DC polarization voltage, the change of sensitivity is not more than $\pm 5\%$ .
CMRR:	$\geq 89\text{ dB}$ .
Frequency response:	0.5 - 150 Hz (+0.4 dB, -3.0 dB).
Lower frequency feature:	time constant is not less than 3.2s.
50Hz filter:	$\geq 20\text{dB}$ .
Equipment safety type:	Defibrillation-proof Type CF Applied Part.
Working environment:	Temperature: 5 - 40°C,
Relative humidity:	30 - 80%,
Atmospheric pressure:	860 - 1060hPa.
Storage and Transport environment:	Temperature: -40 - 50°C,
Relative humidity:	$\leq 93\%$ ,

Atmospheric pressure:  
Dimensions:  
Weight:  
Data Connector:

500 - 1060hPa.  
Approximately 120 × 85 × 28 mm  
Approximately 190g  
USB

### 1.3 Hardware Connection

Model: DMS 300-BTT02  
Patient Cable Model: IEC10/L-15P-M  
Electrode Model: 01-3633

Connect to computer USB port

Patient cable port

Connect to patient cable

### 1.4 System Requirements

CPU:	Pentium PIII or Celeron 800 or above
Memory:	at least 512M
Hard disk:	above 30G
VGA Card:	at least 4M memory
CD-ROM:	CD-ROM or DVD
Operation System:	Windows XP/Vista/7
Resolution:	1024x768 pixels

You can change these settings in Windows XP/Vista/7 in display property. If you have other questions in changing settings, please refer to Windows XP/Vista/7 help.

The computer should meet the requirements of EN60950-1.  
The printer should meet the requirements of EN60950-1.  
The computer should be placed outside the PATIENT ENVIRONMENT.

### 1.5 USB security key driver installation

1. Put the “DMS Software Installation CD” into CD-ROM. Click Sentinel Protection Installer File. The installation process starts. Please follow the wizard to do the installation..

Click “Upgrade” to continue.( Figure 1.5-1).

Figure 1.5-1

Click “Next” to continue

Figure 1.5-2

2. Make a checkmark in the option “I accept the terms in the license agreement”. Click “Next”. (Figure 1.5-3)


Figure 1.5-3

Click “Next” to continue

Figure 1.5-4

Click “Install” to continue

Figure 1.5-5

3. Click “Finish”. Plug the security key into the USB port of PC. The operation system will automatically drive the security key. It does not need to restart computer. (Figure 1.5-6).


Figure 1.5-6

### 1.6 Resting 12-Lead ECG Software Installation

1. The installations of Resting 12-Lead ECG Version 4.0 and Standard 4.0 are the same process.
2. Put the “DMS Software Installation CD” into CD-ROM.

Click “Resting12LeadECG.exe”.The installation process starts (Figure 1.6-1). Please follow the wizard to do the installation.

Figure 1.6-1

Input password “dmsdms”, click “OK”, Click “Next”. (Figure 1.6-2)


Figure 1.6-2

Input user name and organization. Click “Next”. (Figure 1.6-3).

Figure 1.6-3

3. If you want to make change of the patient data directory, click “Change” bar to select the desired directory. Click “Next”. (Figure 1.6-4).

The default installation directory of the software is C:\ Resting12Lead. It is recommended to select the drive with bigger hard disk. If the windows system is installed on the C:\ drive, it is better to save the patient data on other drive. This will prevent the data damage from the possible operation system reinstallation.


Figure 1.6-4

Click “Finish” when installation completes.


Figure 1.6-5

## Chapter 2 SYSTEM SETTINGS

### 2.1 Launch the PC-ECG Software

#### 2.1.1 Before You Begin

Before start Resting 12-Lead ECG, please check the following things:

- Check if the ECG acquisition devices have been correctly connected.
- Check if the ground line has been properly connected.
- Check if the USB security key has been plugged into the computer USB port.
- Check if the patient cable has been properly hookup on patient.
- Check if the windows system settings are correct. (Refer to section 1.5.1 “System requirements” for detailed information.).

#### 2.1.2 Launch the Resting 12-Lead ECG

Double click Resting 12-Lead ECG icon on the desktop with left mouse button to launch the software.

If it is the first time to use this software, the language is English. It allows customers to translate the software into their local language.

For Resting 12-Lead ECG, please refer to “Appendix B: How to translate the software into local language version” for detailed information.

#### NOTE

If it is the first time to use Resting 12-Lead ECG, Most of the function modules are grey color and inactive. Each function has a unique password. You must enter these passwords before start to use the corresponding function. Go to “Setting” on the main screen. Click “Password setting” and enter the passwords.

Please contact with your distributor for the passwords.

### 2.2 Initial Screen

#### 2.2.1 Resting 12-Lead ECG


Figure 2.2.1-1

**Table 2.2.1-1 Initial Screen Description**

Items	Description
<b>A</b>	Graph Button: include New Patient, Report, Tools and Exit. The new patient button include Standard ECG Collecting, Quick ECG Collecting, Multi-Hour ECG Monitoring, Simple Stress. Tools button include System setting, Leads Placement, Make Cardio Card, Import data from PDA, Output ECG data and Output VCG data.
<b>B</b>	Patient: Show the patient information
<b>C</b>	Main Function:Include 8 Main Function, ECG, QTD, FCG, HFECG, Hour, SAECCG, VCG, TVCG
<b>D</b>	Sync Data: Sync the patient database. When you reinstall the software, if you want to add the previous patient data to the present software, click “Sync Data” button, you can see the entire previous patient data come into the new software. The data directory must be same with the previous software.
<b>E</b>	Database: User can use this new type database list to search, compare, modify, delete, import, output patient.
<b>F</b>	Used Disk: it indicates the present hard disk the software occupied. If it is beyond 90%, it is better to delete some unuseful files or backup some files to release enough space for proper performance of the software.
<b>G</b>	Patient information: Show the patient ECG at the left bar, show the other information at the right bar.

**Table 2.2.1-2 Initial Screen Symbol Description**

Items	Description
<b>ECG</b>	Electrocardiography.
<b>FCG</b>	Frequency ECG.
<b>HFECG</b>	High Frequency ECG.
<b>QTD</b>	QT Dispersion.
<b>Hour</b>	Hour ECG
<b>SAECG</b>	Signal Averaged ECG, also called Ventricular Late Potential. (VLP).
<b>VCG</b>	Vector ECG.
<b>TVCG</b>	Timed Vector ECG.

**NOTE**

The Resting 12-Lead ECG 4.0 standard version only has the ECG and VCG functions. The initial screen is the same with Resting 12-Lead ECG 4.0 version. But all the function modules except for the ECG and VCG are grey and inactive.

**2.3 System Settings**

All the PC-ECG software require setup. Click “Setting” button on the main screen to display the system settings’ interface.

### 2.3.1 General Information


Figure 2.3.1-1 General information

It allows user to enter name, phone number, fax number and address of hospital. The hospital information will appear on the heading of the final report pages.

It allows user to enter physician's name, phone number and address. The information will appear on the heading of the report pages.

Set the directory to save patient data. It is recommended to select the drive with enough hard disk space.

### 2.3.2 Interface Setting


Figure 2.3.2-1 Physician

#### Custom Item

It provides four user-defined items. You may input any information such as E-mail address. This will appear on the heading of the final reports.

#### Language

Click the dropdown list to select the software language.

#### Data format

You can select the desired date format by clicking the dropdown list.

### 2.3.3 Signal Processing

In this window you are allowed to set software filter, sample time and COM port (Figure 2.3.3-1).


Figure 2.3.3-1 Signal Processing

#### Filter

A software filter is designed to help eliminate noises from ECG waves.

- Enable Main Filter: click the dropdown list and select 50Hz or 60 Hz.
- Enable ADS: It is strongly recommended to select this option.

This technology will guarantee stable baseline of ECG.

#### NOTE

If there are many noises in the ECG strips, first check if the ground line has been properly connected. Good grounding is the guarantee to get an interference-free ECG. If noises still exist, you need to turn on the software filter.

#### Sound

It is recommended to make a checkmark before “Enable Heart Rate Beep”. It will sound beep each time of the heart beat during data sampling.

#### Sample Time

Set the default sampling time of “12-lead simultaneous ECG” and “3-lead Simultaneous VCG”.

#### Device

Select “Bluetooth device” if you are using Bluetooth ECG transmitter and receiver or adapter.

#### COM Select

Select the corresponding COM port which connects to the ECG sampling devices.

#### NOTE

This is very important for ECG sampling. If you do not select a correct COM port, signal sampling will be failed.

### 2.3.4 Analysis Setting


Figure 2.3.4-1 Analysis Setting

#### Heart Rate Channel

Select a lead for heart rate calculation. For Wilson Lead, select any channel from 12 leads. For Frank Lead, select any channel from X, Y and Z. It is recommended to select the lead with higher amplitude of “R wave” for heart rate calculation (Figure 2.3.4-1).

#### QT dispersion

Set the normal rang of QT dispersion. The default value is 20-50 ms. User may change the value according to his own experiences.

#### Setting

Set the arrhythmias according to your experiences.

- **SVE Premature %**  
Set the SVE premature rate. The default value is 25%. It means that if the heart rate premature rate is equal or over 25%, the beat will be classified as SVE.
- **Pause Time**  
Set pause time. If the pause time is beyond the setting value, it will be considered as heart pause. The default value is 2.5s.
- **Bradycardia**  
The default value is 60bpm.
- **Tachycardia**  
The default values is 100bpm.

#### Auto Option

- **Statement & Condition:** With checkmarks before these two options, the system will give an auto conclusion and interpretation.
- **Set Measurement Point:** If the option is selected, the system will automatically set the measurement points when ECG sampling finishes. It is recommended not to make a checkmark of this option and it will allow you to manually set the measurement points.

#### HFECG Mode

You can select the display mode as “Lengner 6 Lead” or “12 Leads”.

### 2.3.5 Report and Password Setting


Figure 2.3.5-1 Report Setting

#### Printer Type

If the letters are reverse on the report, try to use another printer type.

#### Print Color

Select black or color to print corresponding report.

#### ECG Report

- Landscape (3\*4+2)
- Landscape (12)
- Landscape (3\*4+2) with Table
- Table page

If you make a checkmark in the box before table page, you will get a table page printout.

#### Orientation

- Landscape
- Portrait

#### Print Grid

If you make a checkmark in the box before "Without Grid", it will print only ECG and the grid will not be printed. This will increase the speed of printing, but you need to use grid paper.

If you want both the ECG and the grid printed, do not make a checkmark in the box.

#### Password Setting

If it is the first time to use the software, all of the function modules are inactive; it needs to enter the corresponding passwords to activate all these functions (Figure 2.3.9-1).

The passwords are in the envelope together with the products.

### 2.3.6 System Manager

The icon of “System Manager” is in the Low Right of the interface of “System Setting”.

Single click the icon with left mouse button.

Figure 2.3.6-1 System Manager

The password is “dms”

User can make changes of the hardware gain setting when you enter into “System Manager”, the default setting is “100%”.


Figure 2.3.6-2 Make change of gain

**NOTE:**

Only the authorized person can use this function. Maladjustment of the gain will result in faulty ECG amplitude calculation.

## Chapter 3 DATA SAMPLING

### 3.1 Prepare Patient Skin

#### Electrode

The quality of the electrodes should be in accordance with the requirements of ANSI/AAMI EC12 and Directive 93/42/EEC.

This model of ECG electrodes is recommended as follows: Manufacturer: ConMed Corporation  
Model : 01-3633

1. Require patient to rest, lying on his/her back before putting electrodes on patient. Patient must remain in the same position for resting without movement or talking.

2. Using standard 12-lead cables, follow the procedure below: For acquisition of the standard ECG leads 4 electrodes must be applied on the limbs and 6 on the chest. The limb electrodes should be placed above wrists and ankles. Each lead is marked as follows:

**Chest Electrode Locations**

<b>C1:</b>	Fourth intercostal space at the right border of the sternum
<b>C2:</b>	Fourth intercostal space at the left border of the sternum
<b>C3:</b>	On the 5th rib, midway between C2 and C4
<b>C4:</b>	Fifth intercostal space at the left mid-clavicular line
<b>C5:</b>	At the horizontal level of C4 at left anterior axillary line
<b>C6:</b>	At the horizontal level of C4 at the mid-axillary line
<b>C7*:</b>	At the left posterior axillary line in the 5th intercostal space
<b>C8*:</b>	At the left scapular line in the 5th intercostal space
<b>C3R*:</b>	Opposite C3, on the right side of the thorax
<b>C4R*:</b>	Opposite C4, on the right side of the thorax

\* means Additional standard leads.

**Limb Electrode Locations**

<b>R:</b>	Right Arm, on the medial side of wrist or lateral side between the shoulder and elbow.
<b>N:</b>	Right Leg, on the medial side above the ankle.
<b>L:</b>	Left Arm, on the medial side of wrist or the lateral side between the shoulder and elbow.
<b>F:</b>	Left Leg, on the medial side above the ankle.

**Lead Placement**


Figure 3.1-1 Lead Placement

Figure 3.1-1 shows the corresponding relationships between IEC and AAMI nomenclature and color coding standard. The patient cable

(Model: IEC10/QMD-15P-N) is in accord with IEC nomenclature and color coding standard..

**NOTE**

It is important that the electrodes are not tangled and the leads are properly connected. The patient's arms and legs are totally relaxed.

Also, check that the computer is plugged into a properly grounded socket. For extremely hairy chests, nestle the electrodes in the hair to maximize adhesion or shave the areas, or use adhesive tape over the electrodes.

**NOTE**

If the patient's skin is dirty, first cleanse with 70% Ethyl Alcohol (Rubbing Alcohol) and dry the skin completely.

**NOTE**

A thickness of the ECG baseline or flat line will usually mean the electrode is not conducting properly. Check the electrode on the right ankle first - it is the grounded electrode and is most important to having good conductivity.

**NOTE**

Ground wire should be properly connected. One end of the ground wire should be connected with the out cover of the main unit. The other end should be connected with the public ground wire or water pipe.

**NOTE**

The electrodes should be in accordance with the requirements of ANSI/AAMI EC12 and Directive 93/42/EEC.

**NOTE**

The used electrodes should be disposed of in accordance with the requirements of environment protection.

**NOTE**

The patient cable should not be kneaded or bent. Or the cable will be broken. The connector of the patient cable should not get in touch with the earth line or other bare metal parts.

**3.2 Start a New Patient Test**

Click "New Patient" icon and then "Standard ECG Collecting" on the main screen and you will see the following picture. (Figure 3.2-1). You have two different options which are "12-lead Simultaneous ECG", and "3-lead Simultaneous VCG". "12-lead Simultaneous ECG" use standard 12 leads hookups. "3-lead Simultaneous VCG" uses Frank X, Y, Z hookups. User can select to perform all of these functions or any one of them. The program will guide you to do the tests step by step.


Figure 3.2-1 Sample Style

**12-lead simultaneous ECG**

You can change the sampling time according to your needs. But if you want to sample data for FCG (Frequency ECG) analysis, it is strongly recommended sampling 90 seconds.

**3-lead Simultaneous VCG**

If you want to sample data for Vector ECG, SAECG (Signal Averaged ECG) and TVCG (Time Vector ECG) analysis, you need to apply Frank X, Y, Z hookups, If you want to test SAECG, the sampling time should not be less than 240 seconds.

User can select to perform all of these functions or any one of them.

If user select “12-lead Simultaneous ECG”, when sampling finishes, the ECG, FCG, HFECG, QTD, button will be lighted up.

If user select “3-lead Simultaneous VCG”, SAEKG, VCG and TVCG will be lighted up.

### 3.2.1 Input Patient Information

Click “Next” to Patient Information.


Figure 3.2.1-1 Patient Information

In this screen, it allows user to enter some basic information such as name, age, weight, date etc. “Last Name” and “Refer Physician” must be entered.

#### Type

Select the correct ECG acquisition device you are using. The default setting can be made in System Setting. (Click “Setting” button on the main screen of the software.) “Demo” option is only used for demonstration, not for real ECG sampling. If “Demo” is selected, it will simulate the real ECG sampling process.

#### MI

Middle name.

#### Directory

It is automatically generated by software.

#### Custom Item 1,2,3,4

User can enter any information in the custom item. The information will appear on the heading of the final report. You may also set these items at the “System Settings/Patient Info Custom”.

#### Copy From Selected Patient

Click “Copy From Selected Patient” button, it allows user to select a patient from previous data and start a new test. The patient’s information need not to be entered in this way. The old data will not be covered and the new data will get a new file name automatically. From the file name such as 20050512102 (May 12, 2005), we can see the date when the test was performed.

**Clear**

Click “Clear” to delete the previous information and type in the new one.

**Other Information Tab**

- Risk Factors
- Symptom Class

This information is only for doctor’s references and will not affect the analysis results.


Figure 3.2.1-2 Other Information

**3.2.2 12-lead Simultaneous ECG Sampling**

After filling in the patient information, single click “Next” with left mouse button to continue.


Figure 3.2.2-1 Wilson Lead System

Follow the wizard, connect patient cable to patient. Click “Next” to continue.


Figure 3.2.2-2  
ECG Sampling Screen

### ECG Sampling

At the beginning of ECG sampling, the ECG strips are yellow color and are not recorded by computer. Ensure the heart rate is identified by the program. You have two methods to check with this.

- The beep sound can be heard with each heart beat. It means that the heart rate has been identified by program. Make a checkmark before the box “Enable Beep”.
- Also you may enable beep sound in system setting. If you want to make “Enable Beep” in default, go to “Tools/System Setting/Signal Processing/Sound” to make changes of this setting.

When you see the ECG strips are stable and the heart rate has been identified, click “Start” button on the right bottom of the screen. At this time, the ECG strips become green color. Pay attention only the green color ECG can be recorded.

#### NOTE

To get an inference-free ECG strip, please check the following.

- Check if the heart rate has been identified by the program. You can check if the heart rate is displayed or the beep sound can be heard with each beat. If the heart rate is not identified, the probably reason is that QRS amplitude of the channel which is used for heart rate calculation is too small to be identified. Try to select another channel with higher QRS amplitude for heart rate calculation in field G.
- Make a checkmark before “Enable Main filter” and “Enable ADS Filter”. “Main filter” will help to filter the 50Hz AC interference. “ADS Filter” will help to get smooth and stable ECG strips. You may make these settings in “System Setting” and it will turn on the filter in default for each testing.

### Start

Click “Start” button on the right bottom of the screen to start ECG sampling. ECG waves turn into green. When it arrives the presetting time, ECG sampling will stop automatically. Click “Exit” to continue

**Pausa**

Click “Pause” to stop sampling. Click “Restart” to start a new sampling if you find some noises or other unexpected things happened such as electrodes loose.

If you only sample “12-lead Simultaneous ECG”, after sampling finishes, it will go to “Set ECG Components” screen (See section 3.3 “Set ECG Components”).

**3.2.3 Vector ECG (VCG) Sampling**

**Frank Lead System**

Frank X, Y, Z hookups should be applied on patient for vector ECG sampling.

<b>C1:</b>	Red X (-) Right mid-axillary line in the fifth intercostal space.
<b>C4:</b>	Brown X (+) Left mid-axillary line on the same horizontal line as C1.
<b>C3:</b>	Green Y (+) Mid-clavicle line on the same horizontal level as C1 and C4.
<b>C6:</b>	Purple Y (-) Any point of neck, avoid carotid artery and jugular vein.
<b>C2:</b>	Yellow Z (+) Mid-sternum, on the same horizontal level as C1 and C4.
<b>C5:</b>	Black Z (-) Center of the spine on the horizontal level as 1 and 2.

Click “Next” to continue (Figure 3.2.3-1).


Figure 3.2.3-1 Frank Lead System


Figure 3.2.3-2 VCG Data Sampling

### VCG Data Sampling

At the beginning of the sampling the ECG strips are yellow color and not recorded by computer.

Check if the heart rate has been identified by the program. You can check if the heart rate had been displayed or the beep sound can be heard with each beat. If the heart rate is not identified, it is probably that QRS amplitude of the channel which is used for heart rate calculation is too small to be identified. Try to select another channel with higher QRS amplitude for heart rate calculation.

When the ECG waves are stable and the heart rate has been identified, click “Start” button to sample data. This time the ECG waves become green color and can be recorded.

When it arrives the presetting time, it will stop sampling automatically. Click “Exit”. It will go to “Set ECG Components” screen. (Refer to section 3.3 “Set ECG Components”).

### 3.3 Set ECG Components


Figure 3.3-1 Set ECG Component

### Set ECG Components

After acquisition of an ECG in the automatic mode, the global measurement points for

- P1: P onset
- P2: P offset
- Q: QRS onset
- S: QRS offset
- T1: T onset
- T2: T offset

Can be adjusted manually.

Put the mouse on the cursor, hold and drag the cursor to the correct place. After you set all of the ECG components, click “Finish” to exit. It will go back to the PC-ECG software main screen. Now the function buttons become active. You can click each function button to perform corresponding data analysis.

### 3.4 Quick ECG Collecting

If you want to sample only 12-lead Simultaneous ECG data, you have a very quick method (Figure 3.3-1).

Click “Quick ECG” icon on the main screen, you will see the above picture. Fill in some necessary information (Refer to section 3.2.1 “Input Patient Information”).

### 3.5 Multi-Hour ECG Monitoring

If you want to sample multi hour 12-lead Simultaneous ECG data, you click Multi-Hour ECG Monitoring.

## Chapter 4 DATA ANALYSIS

### 4.1 12-lead Simultaneous ECG

#### ECG Review

Single lick “ECG” icon on the main screen with left mouse button, to enter into the 12-lead simultaneous ECG analysis (Figure 4.1-1).

The hart rate and RR interval are displayed on top of each beat with unit of BPM (Beat Per Minute) and ms (millisecond). The red marker indicates the active beat which can be edited by user.

Figure 4.1-1  
12-lead simultaneous ECG

### Next

Click “Right” Key to choose the beat on the right of the present one.

### Previous

Click “Left” Key to select the beat on the left of the present one.

### Edit

Click “Edit” button to label the beat as VE, SVE, Pause, Normal, Paced, Artifact, Failure to Capture or Failure to Sense. If VE is select, there will be a “V” displayed on the top of that beat. If SVE is selected, there will be an “S” displayed on the top of that beat. If Pause is select, there will be a “P” displayed on the top of that beat. If Normal is select, there is no label displayed on the top of that beat. If Artifact is select, the beat will be deleted and the triangle marker on the top of the beat will disappear. That means the beat has been deleted.

If Paced is select, there will be a “Paced” displayed on the top of that beat. If Failure to Capture is select, there will be a “FTC” displayed on the top of that beat. If Failure to Sense is select, there will be a “FTS” displayed on the top of that beat.

### Speed & Gain

It allows user to print the 8-second ECG strip in various speeds and gains. Click the drop down list of the speed and gain to select. You have three selections of speed which are 1.25cm/s, 2.5cm/s and 5cm/s. You have three selections of gain which are 2.5mm/mV, 5mm/mV and 10mm/mV. When 10mm/mV is selected, you can choose from 10mm/mV 2page, 10mm/mV 3page and 10mm/mV 4page. Click the 1, 2, 3 and 4 next to the speed and gain to see other leads.

### Preview & Print

Click “preview” icon to preview the report. Click “Print” to print the 8-second ECG strip.

### Ruler

Point and click on the “Ruler” icon and the mouse pointer turns into a small ruler with a little red cross. When take measurements, put the red cross at the start point, click and drag the mouse to the end point.

A box will mark the area you have measured. The results will be shown in the “Measurement Result” window.

### Labeling strip

To label a strip, click on the down arrow at the end of the label field. A series of label descriptions will appear. Point and click the desired description which will be printed out on the strip.

To enter your own label, simply mouse click in the label field. Using your keyboard, type in a description and press the ENTER key on the keyboard.

To add more labels to your label database, click on the “Edit” bar next to the “Clear” bar (Figure 4.1-1).


Figure 4.1-2  
Edit Data Dictionary

**Edit Data Dictionary**

Point and click on the blank field and type in the description. After you finish typing in the description, click on the “Add” bar. The new label will be added to the list. Click “OK”.

If you want to delete a label from list, simply click on the label you want to remove. It will be highlighted in blue. Click “Delete” bar. Click “OK”.

**NOTE**

1. Accuracy of amplitude measurements of waveform: 25µV (amplitudes ≤ 500µV) or 5% (amplitudes > 500µV).
2. Interval measurements on biological ECGs should meet the requirement of Table 4.1-1.
3. Stability of measurements against NOISE should meet the requirement of Table 4.1-2.

**Table 4.1-1 Acceptable mean differences and standard deviations for global durations and intervals for biological ECGs (ms)**

Global Measurement	Acceptable mean difference	Acceptable standard deviation
<b>P-duration</b>	±10	15
<b>PQ-duration</b>	±10	10
<b>QRS-duration</b>	±10	10
<b>QT-interval</b>	±25	30

**Table 4.1-2 Disclosed changes of measurements caused by NOISE on ECGS**

Global Measurement	Type of added NOISE	Disclosed differences	
		Mean ms	Standard deviation ms
<b>P-duration</b>	High frequency	±10	8
<b>P-duration</b>	Line frequency	±10	8
<b>P-duration</b>	Base-line	±10	8
<b>QRS-duration</b>	High frequency	±6	5
<b>QRS-duration</b>	Line frequency	±6	5
<b>QRS-duration</b>	Base-line	±6	5
<b>QT-interval</b>	High frequency	±12	10
<b>QT-interval</b>	Line frequency	±12	10
<b>QT-interval</b>	Base-line	±12	10

**Complex**

Click Complex tab, you can see 12 leads complex beat in one screen Use mouse left click to select which beat you want to adjust, use tab key to change the active.

Cursor, left and right key to move the position.

Figure 4.1-3 Complex

### Complex Big

Click complex big tab show you a single channel complex beat.

Figure 4.1-4 Complex

### Page Scan

Click Page Scan tab you can see a 90 seconds ECG data.


Figure 4.1-5 PageScan

### Interpretation

Click Interpretation tab you can review the ECG parameter, and change the conclusion, print report.


Figure 4.1-6 Interpretation

## 4.2 FCG (Frequency ECG) Analysis

### 4.2.1 CAD Locating Diagnosis

Single click “FCG” icon with left mouse button to enter into Frequency ECG analysis (Figure 4.2.1-1).


Figure 4.2.1-1 Frequency Domain Cardiogram

According to the distribution of the positive ECG power spectrum on 12 leads, CAD (Cardiac Artery Disease Locating Diagnosis) system will automatically give the location of CAD and myocardial ischemia.

If it is beyond the middle line in any lead, it means the abnormality in that lead. The combination of abnormality in different lead may indicate the different ischemia location of the heart. Please refer to the following suggestion for ischemia location.

Table 4.2.1-1 Ischemia Location Description

Leads	Location
V1+V2+V3+V4	Anteroseptal
V2+V3+V4+V5	Anterior
II+aVF+V1+V2	Inferior Posterior
I+aVL+V3+V4+V5+V6	Anterolateral
I+aVL+V5+V6	Lateral
I+aVR+aVL+V6	I
II+aVR+aVF	II
III+aVL+aVF	III
I+II+aVR+V5	aVR
I+III+aVL	aVL
II+III+aVF	aVF

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**NOTE:**

If more than one of above conditions is met at the same time, select the one that comes first from above.

If more than 10% heart beats are abnormalities, do not perform the FCG analysis. Because the abnormal beats will significantly affect the analysis results and the results are not reliable.

**Table 4.2.1-2 CAD Locating Diagnosis Screen Description**

Items	Description
<b>A</b>	Cardiac Electric Axis.
<b>B</b>	Automatic interpretation. You may also type in your own interpretation. Click “Save” bar in field C to save the changes.
<b>C</b>	Data dictionary, please refer to J item in the “Data Table Screen Description”. Click “Auto” to get the automatic analysis results. Click “Clear” to delete the results and enter the customized conclusion. Click “Load” to check the previous saved comments. To add more labels to your data dictionary, click on the “Edit Dictionary” bar.
<b>D</b>	Preview and Print report. Close and exit FCG.
<b>E</b>	The limb lead and chest lead heart map.

**4.2.2 Two Leads Compare**


Figure 4.2.2-1  
Two Leads Compare

**Frequency Spectrum Power Phase shift (PHS)**

The spectrum power of lead II and V5 are compared on the top of this screen. The X-axis is frequency and Y-axis is power spectrum.

Phase Shift is the time difference over the frequency axis and represents dys-synchronization of the dynamics of different parts of the heart.

PHS: Great deviation and fluctuation of phase shift. Poor conduction function of blood vessel. Change of blood dynamics.

**Vxy**

Cross-Correlation.

**Hxy**

Transfer function (Modulus).

**Coherence (CHR)**

These are the degree of mutual match of the electrical activity of different parts of the heart, the coherence of amplitude frequency and phase shift, and the characteristic quality or stability of the heart system.

- RSR: Insufficient blood to myocardium.
- CSR: Sometimes arrhythmia.

**Impulse Response**

Impulse Response is, essentially, the use of a unit impulse as input to excite a system. The resulting output is the Impulse Response. It also refers to the transferability of the response.

It can be affected by the physical parameters of the cardiac muscle, the damping characteristic between the cardiac muscle and blood flow, and viscosity of the blood. Delayed or abnormal response can suggest ischemia or latent infarction.

- PDN: The main peak of impulse response is upside down. Poor conduction function.
- MUP: Multiple main peaks. Poor conduction function or left ventricular malfunction.

**4.2.3 Precordial Lead Power**

Figure 4.2.3-1  
Precordial Lead Power

**Precordial Lead Power**

Power Spectrum is a frequency analysis of the ECG power, namely a power distribution of the ECG signal over each frequency. FCG graph is Power in mW (vertical axis) versus Frequency in Hz (horizontal axis). FCG is lined up from low frequency to high frequency as peaks, peak 1 to peak 4 is obvious and peaks after 5 are superimposed and hard to differentiate.

Power Spectrum is the most important in the Ultra-Phase analysis and can give significant pathological information about arrhythmia, fibrillation, poor myocardial function and insufficiency of main heart power.

To zoom out the picture, right click on any place of the screen.  
 To zoom in the picture, left click on any place of the screen.

### 4.2.4 Index Table


Figure 4.2.4-1 Index Table

#### Index Table

The indexes are independent of the conventional ECG and are unique because they are more sensitive. The patient with a normal conventional ECG may have some abnormal cardio spectrum index. The index with “+” sign means abnormality.

Table 4.2.4-1 Diagnosis Table Description

Items	Description
<b>R21</b>	The 2nd peak is greater than the 1st peak. Insufficient blood to myocardium.
<b>R43</b>	The 4th peak is greater than the 3rd peak. Insufficient blood to myocardium.
<b>R51</b>	The 5th peak is greater than the 1st peak. Insufficient blood to myocardium, compensation mechanism performed.
<b>R53</b>	The 5th peak is greater than the 3rd peak. Insufficient blood to myocardium, compensation mechanism performed.
<b>LO1</b>	1st peak is low. Historical injury of myocardium.
<b>LO3</b>	3rd peak is low. Historical injury of myocardium.
<b>LOA</b>	Average amplitude of the first 4 peaks is low. Myocardial ischemia. Insufficient blood to myocardium.
<b>HIA</b>	Average amplitude of the first 4 peaks is high. Ventricular hypertrophy.
<b>PDN</b>	The main peak of impulse response is upside down. Poor conduction function.
<b>MUP</b>	Multiple main peaks. Poor conduction function or left ventricular malfunction.
<b>RSR</b>	Insufficient blood to myocardium.
<b>CSR</b>	Sometimes arrhythmia.
<b>PHS</b>	Great deviation and fluctuation of phase shift. Poor conduction function of blood vessel. Change of blood dynamics.

### 4.3 High Frequency ECG (HFECG)

ECG signal visible on a conventional ECG have amplitudes of a millivolt and are characterized by frequencies less than 100Hz. The ordinarily invisible components have amplitudes in the microvolt range and are characterized by frequency from about 150Hz to about 250Hz. Deviation of these high frequency components from a normal pattern can be indicative of myocardial ischemia or myocardial infarction.

Single click “HFECG” icon with left mouse button to enter into high frequency ECG analysis screen (Figure 4.3-1).

HFECG is to detect high frequency components of QRS such as notches, slurs and beatings.


Figure 4.3-1 High Frequency ECG

Table 4.3-1 HFECG Screen Description

Items	Description
A	The ECG wave between the two broken lines is the present one. Click  button or using left and right arrow keys on keyboard to select other waves. Right mouse click on the ECG strip and select to show other leads.
B	The enlarged present ECG.
C	The statistic data of Notches, Slurs and Beadings. You can make changes of each value by double click on that value..
D	Select from the two display modes. If you select “12 lead diagnostic”, click “Change channel” to check other channels. There is no automatic result.
E	Type in the customized conclusion. Click “Save” to save the comments.

<b>F</b>	Data dictionary: click on the arrow down button and select the term from list. Click “Insert” bar to add the selected term to the “Suggested Interpretation”. Click “Save” button. Click “Compute” bar to get the automatic analysis results. Click “Clear” bar to delete the results and enter again. Click “Prior” bar to show the previous saved comments. To add more labels to your data dictionary, click on the “Edit” bar. Click on the label field and type in the description. Click “Add” bar. You will see the new label has been added to the list. Click “OK”. If you want to delete a certain label from list, simply click on the label you want to remove. It will be highlighted in blue. Click “Delete” bar and the description will be removed. Click “OK”.
<b>G</b>	Preview and print report. Click “Close” to exit HFECG analysis.
<b>H</b>	Mouse point and drag the bar to the desired location to change gain. Also you may use the up and down arrow key on key board to change gains.

### 4.4 QT Dispersion

Single click “HFECG” icon with left mouse button to enter into QT dispersion analysis (Figure 4.4-1).


Figure 4.4-1 QT Dispersion

Table 4.4-1 QT Dispersion Screen Description

Items	Description
<b>A</b>	Speed: 25mm/s and 50mm/s.
<b>B</b>	Gain: 2.5mm/mv, 5mm/mv, 10mm/mv and 20mm/mv.
<b>C</b>	T wave: To select the T end or T max as the measurement point.
<b>D</b>	Data table: automatic results.
<b>E</b>	QT chroma map: use different color to display the P-QRS-T amplitudes.

### Set QT Components

The results of QTD are the average value of the consecutive 5 waves. Use arrow up and down key on your keyboard to select the wave for QTD analysis.

Mouse point on the channel you want to make settings. The waves become green color. There are three vertical lines on the selected wave. The red vertical line is active. Use Tab key on keyboard to switch between these three lines. Use left and right arrow key on keyboard to move it to the desired location.

### Ruler

Point and click on the “Ruler” icon and the mouse pointer turns into a small ruler with a little red cross. Put the red cross at the start point, click and drag the mouse to the end point. A box will mark the area you have measured and the results will be shown in the “Measurement Result” window.

Click the “Ruler” icon again. The mouse pointer will change to normal.

### Preview

Click “Preview” to preview the report.

### Print

Click “Print” to print.

### Close

Click “Close” to exit QTD analysis.

**Table 4.4-2 Data Table Description**

Items	Description
<b>Max QT</b>	Maximum QT in the selected 5 waves.
<b>Min QT</b>	Minimum QT in the selected 5 waves.
<b>average QT</b>	Averaged QT in the selected waves.
<b>Max QTc</b>	Maximum QTc in the selected 5 waves.
<b>Min QTc</b>	Minimum QTc in the selected 5 waves.
<b>QTc</b>	Averaged QTc in the selected 5 waves.
<b>Max QTD</b>	Maximum QTd in the selected 5 waves.
<b>Max QTcd</b>	Maximum QTcd in the selected 5 waves.
<b>QTcd Averaged</b>	QTcd in the selected 5 waves.
<b>QTDR-QT</b>	Percentage of QT to RR interval.
<b>AdQTd</b>	Maximum QTd of border channels.
<b>Qtad</b>	The early part QTd.
<b>Qted</b>	The late part QTd.
<b>JTd</b>	Averaged JTd in the selected 5 waves.
<b>JTcd</b>	Averaged JTcd in the selected 5 waves.

## 4.5 Signal Averaged ECG (SAECG)

The SAECG is also called Ventricular Late Potential (VLP). The SAECG screen display consists of two parts. The top display is the R-R Interval graph. Each vertical line represents the heart rate of one beat. The longer the line the slower the heart rate is. Moving the cursor bar with a mouse click or with the arrow keys to the desired location, and using “Reject” button to reject the R-R interval can delete individual R-R intervals.

The large data display is one minute of ECG data. Each horizontal sweep represents 12 seconds of ECG data. Placing the mouse cursor on the desired location and clicking can eliminate undesired ECG beats. Reversing this process will cause rejected beats to be accepted. It is desirable to remove beats with artifact from the SAECG Late Potential file. You can access minutes 1 by pointing and clicking the mouse on the Page 1 tab, or you can change the minutes by pressing the PGUP and PGDN keys on the keyboard.

Single click “VE Chaos” icon with left mouse button to enter into Signal Averaged ECG analysis (Figure 4.5-1).


Figure 4.5-1  
Signal Averaged ECG

### Accept

Placing the cursor on the rejected beats and clicking “Accept” can cause the beats to be accepted

### Reject

Move the cursor to the desired vertical line on the R-R Interval graph, pointing and clicking “Reject” icon can delete individual R-R interval.

### Auto

Click “Auto”. Set the desired percentage of rejected beats whose R-R deviation is exceeding the percentage of the average heart rate. The default percentage is 10%. Place the cursor on the slide bar, click and drag to the desired percentage.

### Time

Click “Time” icon to activate the Time Domain screen display.

### Freq

Click “Freq” to activate the Spectral Frequency display.

### MIN R-R

Click “MIN R-R” and the screen cursor will be moved to the minimum R-R interval.

### MAX R-R

Click “MAX R-R” and the screen cursor will be moved to the maximum R-R interval.

### Left, Right, Up and Down Arrows

Click or press the ARROW keys on the keyboard to move the screen cursor. The Up and Down arrows move with larger increments than do the side to side arrows.

#### 4.5.1 Recommended Procedure to Perform SAECG Test

1. Frank X, Y, Z lead hookup should be applied.
2. Require patient to rest, lying on his/her back for 10 minutes before putting hookups on patient. Patient must remain in the same position for testing without moving or talking.
3. The data sampling time should be at least 5 minutes for SAECG analysis. Because it needs a total of 200 or more qualified beats for good SAECG analysis.
4. After sampling finishes, click SAECG button to perform analysis.  
Click on "Auto" button to set the desired percentage of rejected beats whose R-R deviation is exceeding the percentage of the average heart rate. The default value is 10%.
5. Placing the mouse cursor on the desired location and clicking can eliminate undesired ECG beats. Reversing this process will cause rejected beats to be accepted. It is desirable to remove beats with artifact and arrhythmia from the SAECG Late Potential file. You can access minutes 1 by pointing and clicking the mouse on the Page 1 tab, or you can change the minutes by pressing the PGUP and PGDN keys on the keyboard.
6. Click on "Time" button to perform the Time Domain Analysis set-up (See following: Time Domain Analysis).
7. Click on "Freq." button to perform the Frequency Domain Analysis. (See following Frequency Domain Analysis).

#### 4.5.2 Time Domain Analysis

There are four vertical lines and one horizontal fixed measurement on this screen. The active line will be highlighted in bright green. Use the Tab key on the keyboard to change the active screen line. Use the arrow keys on the keyboard to move the screen line to a proper position.

The top of the screen has 3 channels of a superimposed standard QRS. The first vertical line should be at the Q Point. If the line is not in this position, adjust it using arrow keys on the keyboard.

The second vertical line should be at the end of the QRS (J Point). If it is not, press the Tab key on the keyboard to activate the second line, use the arrow keys on keyboard to move the line to the end of the QRS.


Figure 4.5.2-1  
Time Domain Analysis

The bottom half of the screen shows the Filtered QRS. A computer algorithm locates the third vertical line to the right of the Filtered QRS. If you disagree with the location, press the Tab key on the keyboard to activate the line. Use the arrow keys on the keyboard for moving to the end of the Filtered QRS.

The vertical line at the beginning of the Filtered QRS moves simultaneously with the cursor for the beginning of the standard QRS.

Note that the Filtered QRS is usually a wider measurement than the standard QRS.

The calculation results are displayed in the data box on the right top of the screen.

To exit, single click the Close  with left mouse button at top right corner of the screen.

### 4.5.3 Frequency Domain Analysis

Click on “Freq.” icon and the SAECG Frequency 3-D Graph will come up on your screen.

View the Frequency Graph for Channel 1. You can move the active yellow vertical line at the end of the QRS with the arrow keys on your keyboard. Note that the Initial Window Width numbers will change in the data box directly above the 3D graph while 3D graph changes.

The information in the Data Box at the top right of the screen will also change.

Right mouse click and scroll down to show Ch. 2 or Ch3 then click.

The SAECG Frequency Graph will change with each channel. You will be able to view the changes in spectral power.


Figure 4.5.3-1  
Frequency Domain Analysis

Click on the Close icon  at the top right corner of the screen to back to the SAECG main screen.

Click on the Close icon  at the right corner of the screen to exit the SAECG analysis.

### 4.6 Vector ECG (VCG)

Single click “VCG” icon with left mouse button to enter into VCG analysis (Figure 4.6-1).


Figure 4.6-1 Vector ECG

Table 4.6-1 Vector ECG Screen Description

Items	Description
A	Vector loops of Frontal Plane.
B	Vector loops of Horizontal Plane.
C	Vector loops of right Saggital Plane.
D	Automatic analysis data. (Refer to “Interpretation of the parameters of field D” of this section).
E	Preview and print report.
F	Zoom in and print out any plane with different gains. Select the plane and loops you want to print out.
G	Select gain.
H	After selecting the desired gain then click on “Print Zoom Report” bar to print.
I	The vector ECG can be transformed to 18 leads. Make a checkmark on “Display 12 leads”, then click “Display” bar to see the standard 12 leads ECG. Make a checkmark on “Display 6 leads”, and then click on “Display” bar. You will see lead V3R, V4R, V5R, V7, V8, and V9 lead.

**Frontal Plane**


Figure 4.6-2 Frontal Plane

In this screen you can see the animation play of the vector loops. To set the play speed, mouse point on the little bar of the play speed and drag to the desired location. Click on the down arrow and select gain.

After you set the play speed and gain, click “Play” button to start. P loop is blue color, QRS loop is green color and the T loop is yellow color.

Click on “Right Saggital Plane” or “Horizontal Plane” tab and you can observe the corresponding animation play. The operation method is the same with the “Frontal Plane” above mentioned.

### 3D (Three Dimensional Vector Loops)

This function allows you to observe the vector loops in three dimensions. You have three display styles which are line, fill and transparent. If there is a checkmark before the “Line” style, the other two selections will become grey and can’t be selected. Another clicking on the checkmark of the “Line” will remove the checkmark and the other two selections will become available. You may select any of them. It is recommended to select both of the “Fill” and “Transparent” (Figure 4.6-3).


Figure 4.6-3 3D (Three Dimensional Vector Loops)

Make a checkmark in the box before “Rotation”, the whole picture will rotate and all directions of the vector loops can be seen. Clicking “Restore” bar can cause the vector loop to go back to the original location. Another clicking on the checkmark in the box before “Rotation” will stop rotation

#### Zoom

You may zoom in or out each plane by mouse point on the little bar and drag to the desired location.

#### Vista

You can select to view 3D, X-Y, Y-Z and X-Z plane.

#### Print

If you want to print screen, click “Save”. Then click “Preview”. Then click “Print”. With the same method you can print any of the vector loops displayed on screen.

#### Exit

Click on the Close icon  at the right corner of the screen to exit the Vector ECG analysis.

### 4.7 Timed Vector ECG (TVCG)


Figure 4.7-1  
Timed Vector ECG

With this function you can see the consecutive vector ECG on the time axis. The top half of the screen shows the X, Y, Z three channels ECG. The bottom half of the screen shows the “X-Y”, “X-Z” and “Z-Y” three planes vector ECG. The right of the screen shows the loops of the Frontal Plane, Horizontal Plane and Right Sagittal Plane.

#### Previous

Click “Previous” button to select the beat on the left side of the present one.

#### Next

Click “Next” button to choose the beat on the right side of the present one.

#### Ruler

Click on the “Ruler” icon and the mouse pointer turns into a small ruler with a little red cross. When take measurements, places the red cross at the start point and drag the mouse to the measurement’s end point. A box will mark the area you have measured and the results will be shown in the “Measurement result” window.

The functions of other operation icons please refer to section 4.1 “12-lead Simultaneous ECG”.

## Chapter 5 OTHER FUNCTIONS OF THE SOFTWARE

### 5.1 Database Manager

Figure 5.1-1  
Patient Database

#### Search

You may search for any patient with many different query conditions such as directory, ID, name, sex, age, telephone number, record date, refer physician and hardware type.

#### Clear

Click "Clear Condition" bar to delete all the existing query conditions.

#### Delete

Select a patient you want to delete from the patient list, click "Delete" bar. The patient will be removed from the patient database.

#### Modify

It allows user to modify patient's information. Select the patient you want to make modification, then click "Modify" bar. You can make any changes of the patient information. Click "Update" button to save the changes. Click "Close" to ignore the change. Click "Reset" to re-input the information.

#### Compare

This function allows you to make comparison of any two patients' data or any two times tests of the same patient. Select the desired patient by mouse click. Press "Ctrl" key on keyboard and hold, click another patient. Click "Compare" bar. You can see the two patient data are displayed on the same screen. You can type in your own suggestions, preview and print the comparison report.

#### Print

Click "Print" to print the patient list.

#### Export

This function allows you to save data to other directory or burn into CD. First, select the patients you want to remove from patient database. (Left click on one of the patient, then press and hold "Shift" key; left click on the last one you want to remove. In this way, you can select many patients one time.) Click "Export" button, select the directory you want to save data. The data will be moved to the indicated directory.

If you want to burn the data into CD, you must install special CD burning software. Burning the patient files to the CD.

**Import**

If you want to review the data you have exported to other directory or burned to CD, you must import these data back to the patient database.

Just click the “Import” button and locate the directory to save the patients. (See Figure 5.1-2). Click “Open” button and select the patients you want to import to the database.

Figure 5.1-2 Import Data

**NOTE**

In the process of importing patients back to the database, do not select any single patient file (See Figure 5.1-2), click “Open” bar directly and you will see all patient names are listed. Select the patient names you want to import, click “Import” button to continue.

**5.2 Reanalysis**


Figure 5.2-1 Reanalysis

Click on the icon of “Reanalysis” to enter into the reanalysis module.

### Set ECG Components

In the reanalysis interface, user should set ECG components. You can set ECG components from single lead or all leads.

- P1: P onset
- P2: P offset
- Q: QRS onset
- S: QRS offset
- T1: T onset
- T2: T offset

Click “Set VCG Measurement Point” button and you can set ECG components for X, Y or Z lead. Click on reanalysis bar to reanalyze the patient’s ECG data.

### 5.3 Cardio Card

Figure 5.3-1 My Cardio Card

It allows to a make a patient ECG card just like the ID card. Print the 12 lead ECG strips, fold in half and print it out. Then seal it with plastic seal machine.

#### Wave Location

Click on the arrow button to make adjustment of the ECG strip location.

#### Information

Enter patient information.

#### Interpretation

You can type in interpretation in the blank. If you make a checkmark in the little box before the “Interpretation” in the field of “Print Selection”, the interpretation will be printed out.

If you make modification of the “Interpretation”, the content will only appear on the ECG card but not be saved in the database.

#### Print Selection

Set the top and left margin of the printout. Make a check mark of the “Grid” if you want the grids be printed out. Make a checkmark of the “Interpretation” if you want the interpretation be printed out.

**Print**

Click “Print” to print the report.

**Preview**

Click “Preview” to preview the print.

**Save Setting**

Click “Save Setting” to save all of the settings you made.

**Exit**

Click “Exit” to close the function.

**5.4 Report**


Figure 5.4-1 Report

Click on the icon of “Report” on the main screen.

Make a checkmark in the report you want to print out. The grey color items are inactive and can not be selected. It is because that functions are not performed and there is no suitable data.

Click “Print All” if you want to print all the reports.

Click “Send E-mail” button. The selected reports can be sent to other doctors by e-mail via Microsoft Outlook.

GDT is only used for Germany customers.

You can create the ASCII file (ASCIIIECG.txt) by clicking “Save ASCIIIECG” button. You can view the ASCIIIECG.TXT file by entering the patient data directory and opening the patient file.

### 5.5 Exercise

This function is designed to work with the “Two-step Exercise Test”. Click on “Exercise” icon of the initial screen. The following picture will pop up.

Figure 5.5-1 Input Patient Information

#### Input patient Information

Input patient information such as Name, ID, Age, Sex, and ECG sampling device type etc. There are four “Custom Item” fields which allow user to enter customized information. (The setting of Custom Item please refer to section 2.3.7 “Patient info custom Item”). The customized information will be displayed on the heading of the printout report (Figure 5.5-1).

Click on “Clear” icon to delete all the information of current patient.

Click on “Sample” icon to start exercise ECG sampling.


Figure 5.5-2 Exercise ECG Sampling

**Table 5.5-1 Exercise ECG Sampling Screen Description**

Items	Description
<b>A</b>	Display the resting ECG. The ECG of beginning several seconds will be selected as the resting ECG which will be displayed in 5 seconds after sampling starts.
<b>B</b>	The present ECG waves. Make change of the ECG gain and speed by clicking on the “Gain” and “Speed” button on the right bottom corner of the screen.
<b>C</b>	Heart Rate (HR).
<b>D</b>	% Target HR: the percentage of the present HR and the target HR.
<b>E</b>	Target Heart Rate.
<b>F</b>	Heart rate of the resting ECG.
<b>G</b>	Delta Heart Rate: The difference of the present HR and the resting HR.
<b>H</b>	Time: ECG sampling time. The max sampling time is 15 minutes.
<b>I</b>	Left arrow and right arrow: select ECG as the resting ECG in case the resting ECG was not displayed automatically. Click on the right arrow bar to select the next beat. Click on the left arrow bar to select the previous beat.

**Resample**

If you are not satisfied with the sampled ECG, click “Stop” button to stop the present sampling. Click “Resample” button to start a new sampling.

**T-HR**

Click “T-HR” button. It will pop up the “Target Heart Rate” window. The formula is  $(220 - \text{Age}) * 85\% = \text{Target HR}$ . User can make change of the age and coefficient (See Figure 5.5-3).


Figure 5.5-3 Target Heart Rate

**Menu**

Click “Menu” button, it will pop up the picture shown by following picture (Figure 5.5-4). It allows user to set the print time. For example, make a checkmark in the little box after “HR 80”, it will make a time mark at the point when the HR is 80bpm for later print.

“End+1”: It will make a time mark at 1 minute after recovery stage starts.

“End+2”: It will make a time mark at 2 minutes after recovery stage starts.


Figure 5.5-4 Printing Setting

**Recovery**

Make a mark at the beginning of the recovery stage. Usually Click on this icon after the target HR is arrived.

**Save**

During ECG sampling, you may save the present ECG strips by clicking “Save” icon. These ECG strips can be printed after sampling finished. It allows 10 ECG strips to be saved. If the ECG strips are more than 10, some of them will be covered.

**Gain**

Make change of ECG gain.

**Speed**

Make change of ECG speed.

**Pgscan**

When ECG sampling finishes or the “Stop” icon is clicked on, click “Pgscan” to show the page scan window (Figure 5.5-5). The green line indicates the time when “Recovery” button is pressed. The red lines indicate the time when “save” button is pressed. The orange lines indicate the time set in the “Menu”.


Figure 5.5-5 Page scan

Move the mouse to the lead name on the left top screen, the mouse will change to the hand shape. Right click on the lead name to select any channel for display.

## Print

To preview and print ECG reports. Click “Print” button, it will pop up the following picture(Figure 5.5-6).


Figure 5.5-6 Print Report

The left half screen shows the settings made in the “Menu”. The right half screen shows the saved ECG strips.

Click the magnifying glass icon to preview ECG strips. Make a checkmark in the items you want to print. Click “Print” to print the selected ECG strips.

### NOTE

Print the patient reports before exit the program. The patient data will not be saved after exit the program.

# Chapter 6 REPORT INTERPRETATION

## 6.1 12-lead Simultaneous ECG


Figure 6.1-1 12-lead Simultaneous ECG

On the top right of the picture, it shows the speed and gain which is 25mm/s and 10mm/s.

On the far right of the picture, it shows PA, QA, SA, STM and TA which are described as follows:

**Table 6.1-1 Interpretation of Parameter**

Items	Description
<b>PA</b>	P wave amplitude of the averaged ECG.
<b>QA</b>	Q wave amplitude of the averaged ECG.
<b>SA</b>	S wave amplitude of the averaged ECG.
<b>STM</b>	The maximum ST segment.
<b>TA</b>	T wave amplitude of the averaged ECG.

On the left bottom corner shows the Frequency, Sample Time, Heart Rate, P Duration, QRS Duration, T Duration, PR Interval, QT interval, QTc Interval, P Axis, QRS Axis and T Axis.

On the right bottom corner shows the physician suggestion.


Figure 6.1-2 12-lead Simultaneous ECG

On the top of each beat shows the heart rate and RR interval with the unit of “Beat Per Minute” and “ms”.


Figure 6.1-3 12-lead Simultaneous ECG

Table 6.1-2 Interpretation of Data Table

Items	Description
P+(mV)	The positive going P wave amplitude
P-(mV)	The negative going P wave amplitude
Q (mV)	Q wave amplitude
R (mV)	R wave amplitude
S (mV)	S wave amplitude
R' (mV)	R' wave amplitude
S' (mV)	S' wave amplitude
STj (mV)	ST segment amplitude at J point
ST20	(mV) ST segment amplitude at post J point 20ms
ST40 (mV)	ST segment amplitude at post J point 40ms

Items	Description
ST60 (mV)	ST segment amplitude at post J point 60ms
ST80 (mV)	ST segment amplitude at post J point 80ms
T+ (mV)	The positive going T wave amplitude
T- (mV)	The negative going T wave amplitude
Q (ms)	Q wave time
R (ms)	R wave time
S (ms)	S wave time
R' (ms)	R' wave time
S' (ms)	S' wave time

On the left bottom corner shows the Frequency, Sample Time, Heart Rate, P Duration, QRS Duration, T Duration, PR Interval, QT interval, QTc Interval, P Axis, QRS Axis and T Axis.

On the right bottom corner shows the physician suggestion.

### 6.2 Frequency ECG (FCG)


Figure 6.2-1  
Frequency ECG

The above half of the report is Cardiac Artery Disease Locating Diagnosis system.

If it is beyond the middle line in any lead, it means the abnormality in that lead. The combination of abnormality in different lead may indicate the different ischemia location of the heart.

PC-ECG Software 2.0 version has no Index table.


Figure 6.2-2  
Frequency Analysis

The first and second lists of the report are Frequency Spectrum power of 12-lead ECG. The third list is II and V5 dual lead correlation

### 6.3 High Frequency ECG


Figure 6.3-1  
High Frequency ECG

On the left bottom corner shows the statistic table of the notches, slurs and beadings of each lead.

On the right bottom corner shows the physician suggestion.

### 6.4 QT Dispersion


Figure 6.4-1 QT Dispersion

On the far right side of the QTD report, the data highlighted with black rectangle is the maximum value of the QT/QTm. The data highlighted with the white rectangle is the minimum value of the QT/QTm.

QT: interval between Q point and T end.  
 QTM: interval between Q point and T max point.

Table 6.4-1 Interpretation of Data Table

Items	Description
Max QT	Maximum QT in the selected 5 waves.
Min QT	Minimum QT in the selected 5 waves.
Average QT	Averaged QT in the selected waves.
Max QTc	Maximum QTc in the selected 5 waves.
Min QTc	Minimum QTc in the selected 5 waves.
QTc	Averaged QTc in the selected 5 waves.
Max QTd	Maximum QTd in the selected 5 waves.
QTd	Averaged QTd in the selected waves.
Max QTcd	Maximum QTcd in the selected 5 waves.
QTcd	Averaged QTcd in the selected 5 waves.
QTdr	Percentage of QT to RR interval.
AdQTd	Maximum QTd of border channels.
Qtad	The early part QTd.
Qted	The late part QTd.
JTd	Averaged JTd in the selected 5 waves.
JTcd	Averaged JTcd in the selected 5 waves.

### 6.5 Hour ECG


Figure 6.5-1 Hour ECG

Table 6.5-1 Interpretation of Data Table

Items	Description
<b>SDNN</b>	The standard deviation of normal to normal intervals.
<b>rMSSD</b>	The root-mean square of the difference of successive R-R intervals.
<b>SDSD</b>	The standard deviation of the difference of successive R-R intervals.
<b>pNN10</b>	The percent of R-R intervals differing > 10ms from the preceding one.
<b>pNN20</b>	The percent of R-R intervals differing > 20ms from the preceding one.
<b>pNN30</b>	The percent of R-R intervals differing > 30ms from the preceding one.
<b>pNN40</b>	The percent of R-R intervals differing > 40ms from the preceding one.
<b>pNN50</b>	The percent of R-R intervals differing > 50ms from the preceding one.
<b>pNN60</b>	The percent of R-R intervals differing > 60ms from the preceding one.
<b>pNN70</b>	The percent of R-R intervals differing > 70ms from the preceding one.
<b>Total Power</b>	Total power.
<b>VLF</b>	Very Low Frequency.
<b>LF</b>	Low Frequency.
<b>HF</b>	High Frequency.
<b>LF/HF</b>	Ratio of Low Frequency to high Frequency.
<b>Triangular Index</b>	The integral of the density distribution divided by the maximum of the density distribution.

### 6.6 Signal Averaged ECG (SAECG)


Figure 6.6-1 Signal Averaged ECG

Table 6.6-1 Interpretation of Data Table

Items	Description
<b>Std QRS (ms)</b>	The standard deviation of normal to normal intervals.
<b>Total QRS (ms)</b>	The root-mean square of the difference of successive R-R intervals.
<b>Under 40 (ms)</b>	Duration of the terminal filtered QRS complex remains below 40 microvolt.
<b>Last 40 (mv)</b>	The voltage in the last 40ms of the filtered QRS complex.
<b>SD of noise</b>	Noise.
<b>Cycle</b>	Beat number.

Common criteria defining of a late potential includes the following:

- The filtered QRS complex is longer than 114ms,
- The terminal filtered QRS complex remains below 40 microvolt for more than 38ms, and
- There is less than 20microvolt of signal in the last 40ms of the filtered QRS complex.

### 6.7 Vector ECG


Figure 6.7-1 Vector ECG

Form the top to bottom shows the Frontal Plane, Horizontal Plane and Saggital Plane information. Take the Frontal Plane as an example, on the left shows the QRS loop and the area percentage in the four quadrants. In the middle shows the standard ECG as well as the gain and speed. On the right shows the P loop and T loop. On the far right shows the auto calculated VCG data which is described as below.

Table 6.7-1 Interpretation of Data Table

Items	Description
<b>QRS</b>	The maximum QRS loop amplitude (mV), angle (Degree) and the loop shape or rotation direction.
<b>P</b>	The maximum P loop amplitude (mV), angle (Degree) and the loop shape or rotation direction.
<b>T</b>	The maximum T loop amplitude (mV), angle (Degree) and the loop shape or rotation direction.
<b>CW</b>	Clockwise.
<b>CCW</b>	Counter clockwise.
<b>8</b>	"8" shape loop.
<b>Angle (QRS-T)</b>	The angle of the maximum QRS amplitude and the maximum T amplitude (Degree).
<b>MaxA (QRS/T)</b>	Ratio of the maximum QRS amplitude to the maximum T amplitude.
<b>0.01s</b>	The amplitude and angle of QRS loop at 0.01 second.
<b>0.02s</b>	The amplitude and angle of QRS loop at 0.02 second.
<b>0.03s</b>	The amplitude and angle of QRS loop at 0.03 second.
<b>0.04s</b>	The amplitude and angle of QRS loop at 0.04 second.

**6.8 Timed Vector ECG**


Figure 6.8-1  
Timed Vector ECG

The top screen shows the X, Y, Z three channels ECG. The middle screen shows consecutive vector ECG on the time axis. The bottom screen shows the loops of the Frontal Plane, Horizontal Plane and Right Sagittal Plane.

On the top half of the screen shows the HRT tachogram. Below the tachogram shows the Average, Minimum and Maximum R-R interval value and how many beats are calculated.

On the right of screen shows the risk indications which include:

- TO >= 1,0
- TS <= 2,5 ms/beat
- Age >= 65 years
- Precious MI
- Mean HR > 75 bpm
- HRV SDNN <= 20 units
- LVEF < 30%.

## Appendix

### Appendix A:

#### Recommended Basic Procedure of Using Resting 12-lead ECG Software.

When it is the first time to use RESTING 12-LEAD ECG software, it is recommended to do according to following steps which will help you easily perform a new patient test and analysis.

1. Select the ECG sampling devices you are using and set the corresponding COM port which connects to the ECG sampling devices.  
Select "Bluetooth device" if you are using DMS 300-BTT01/ DMS 300-BTR01 Model Bluetooth ECG transmitter and receiver"
2. Click on "System Settings" on the main screen. Click "Report and Password Setting" and type in the corresponding passwords for all the 7 functions. The password is in an envelope together with the products. After you input the entire passwords, you can see all the ten function modules are turning into active.
3. Make some system settings. enter general hospital information, the information will be displayed on the heading of final reports.  
Click "Signal Processing". Make a checkmark in the box before "Enable Main filter", "Enable ADS" and "Enable Beep". Select filter from 50Hz and 60Hz.

Click "Report and Password Setting". Select the print mode, color, Orientation. Select the report format of the printout from the Landscape (3\*4+2), Landscape (12) and Landscape (3\*4+1) with table. Make a checkmark if you want the data table be printed out. Click "General Information". Set the Raw Data Directory. Click the question sign to select the directory. Select Language you want to use.

4. Start a new patient test. Click "New Patient" icon on the main screen. Select the preferable items from 12-lead Simultaneous ECG, 3-lead Simultaneous VCG and set the sampling time for each item.  
  
Click "Next". You will be asked to type in some patient information. Pay attention in the "Type" field, select the correct ECG sampling device you are using. Follow the wizard to finish ECG sampling.
5. After ECG sampling finishes, you need to set the ECG components. The data will be analyzed according to your settings.
6. Now you can perform all the ten functions by click on the corresponding icons on the main screen. For example, clicking on the "ECG" icon on the main screen allows user to perform 12-lead simultaneous ECG analysis. Click on "FCG" icon, you can perform the Frequency ECG analysis. Please refer to the related chapter for the different functions. You may see some of the function modules are grey and inactive. It is because you did not sample data for that function. For example, if you only sample data for 12-lead simultaneous ECG, the SAEKG, VCG and TVCG icons will be grey color and inactive.

### Appendix B:

#### Maintenance and Cleaning

To keep your ECG Acquisition System in good working condition, please observe the following:

##### Cleaning components

1. The external housing of the ECG Acquisition Box should be cleaned and disinfected as required and at least semi-annually.
2. Remove the Data Input Line before cleaning.
3. Clean external surfaces with a damp cloth moistened with mild soap solution; wipe dry with a clean cloth.

##### Maintenance of the patient cable

1. The patient cable should be cleaned regularly. Use cotton cloth to wipe the dirt. Do not use the detergent with the content of alcohol, or the cable will become rigidified and easy to break.
2. The patient cable should be checked if it has physical break each time when starting a new patient test.
3. Do not submerge, autoclave, or steam clean the patient cable.
4. The patient cable should be disinfected by sodium hypochlorite (bleach) solution 10% in water.

##### Maintenance of ECG Acquisition Box

1. Keep the ECG Acquisition Box clean.

2. The devices should not work in humid environment.
3. Do not use acetone, ether, Freon, petroleum derivatives, or other solvents to clean the devices.
4. Do not let soap or water flow into the devices or come into contact with patient connector pins.
5. The Box should be disinfected by sodium hypochlorite (bleach) solution 10% in water.

### **Disposal of Waste Materials**

1. The used electrodes should be disposed of in accordance with the requirements of environment protection.

### **Discarding of the Equipment**

1. Discard the ECG Acquisition Box, patient cable, and accessories according to local laws.
2. Please follow the state's recycling laws or your facility's recycling policy to ensure proper disposal of the ECG Acquisition System patient cable and accessories.

## **Appendix C:**

### **Trouble Shooting:**

1. When I start the RESTING 12-LEAD ECG software, it pop up "Security key not found". How to solve the problem?  
A: First check if the security key that plugged into the computer USB port is loose. Then check if the security key driver has been installed and worked.
2. Why the sampling ECG is just line or artifacts?  
A: First, check if the ECG lead wire is loose from the ECG acquisition box. Second, check if the earth line worked. Third, check if the correct sampling device is selected.
3. The printer does not work even after several minutes after I click the print button in the software.  
A: Check the printer cable connected to the computer. Reinstall the printer driver.
4. The printer cable is good, but the printer does not work at all.  
A: Maybe there are too many programs are running at the same time. Close these programs and restart windows. If the printer still does not work, contact with your dealer.
5. In case of the program abnormally exit, will the patient data will be lost?  
A: The patient data will not be lost. Restart computer and RESTING 12-LEAD ECG software and you can check the data
6. How can I do if the printer speed is too slow?  
A: Go to printer properties. Click "Advance" tab. Click "Print Processor" button. Select "Raw" in the default data type.
7. Why the heart rate of the sampled ECG is two times of the actual heart rate?  
A: Please check if the T wave is too high in the heart rate channel. Change to other lead with lower T wave for heart rate calculation.
8. Why the sampled ECG is too tiny?  
A: Different person has different ECG waves. You may increase ECG gain in the software.
9. What is the meaning of "Custom Item1" on the report? How to delete this information?  
A: User may enter the customized information in the four "Custom Item" fields and this information will appear on the heading of the report. User can edit this in "System Settings/Patient Info Custom Item".
10. Why all of the function modules are grey color and inactive when I start RESTING 12-LEAD ECG software?  
A: If it is the first time to use the software, all of the function modules are grey color and inactive including FCG, HFECG, QTD, PMECG, SAECG, VCG, TVCG and Hour ECG, Each function has a unique password. You must enter these passwords before start to use the corresponding function. Go to "System Settings" on the main screen. Click "Report and Password setting" and enter the passwords. Please contact with your distributor for the passwords.
11. On the far right of the QTD report, there are some data highlighted with black rectangle and some with white rectangle. What does this mean?  
A: The data highlighted with black rectangle is the maximum value of the QT/QTm. The data highlighted with the white rectangle is the minimum value of the QT/QTm.

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

### **GIMA WARRANTY CONDITIONS**

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons.

Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.