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**SILVER/SILVER CHLORIDE REUSABLE
ECG ELECTRODES AND
ACCESSORIES**

DIRECTIONS FOR USE

INDICATIONS

Surface ECG recording.

PRELIMINARY NOTES

Normally a surface ECG (at rest or during stress) consists of 12 leads: 3 limb bipolar recordings, 3 limb unipolar recordings and 6 precordial unipolar recordings. Therefore, 10 electrodes, 4 for the unipolar and bipolar derivations and 6 for the precordial derivations have to be placed on the patient. The following FIAB electrodes are available:

MODELS F9010SSC – F9010PSSC

Silver silver chloride peripheral plates electrode with universal connection to be applied to the limbs. These plate electrodes are fixed to the limbs by means of the rubber belts models F9011PG – F9011P – PG910/15.

MODELS F9023SSC - F9024SSC – F9024OSSC

Silver silver chloride clamp electrodes kept in position by an elastic mechanism

MODELS F9008SSC – F9009SSC - F9015SSC - F9016SSC

Silver silver chloride suction chest electrodes with universal connection provided with rubber balls. After their positioning on the skin surface, rubber balls must be squeezed to improve adherence. These chest electrodes are used for precordial recording at rest

MODELS F9002SSC – F9003SSC

As the models above, with soft rubber ring. In case of rest ECG the electrodes are positioned as follows:

- Upper limbs: about 10 cm above the wrists on the internal forearm face.
- Lower limbs: about 10 cm above the malleolus on the internal leg face.
- Precordials: (V₁) 4th left intercostal space at the right parasternal line. (V₂) 4th left intercostal space at the left parasternal line. (V₃) halfway between V₂ and V₄. (V₄) 5th left intercostal space at the midclavicular line. (V₅) 5th left intercostal space at the anterior axillary line. (V₆) 5th left intercostal space at the midaxillary line.

In case of a stress test, to avoid ECG noises caused by limbs movements, limb electrodes must be positioned on the back of the patient. The precordial electrodes must be positioned as follows:

- Right arm electrode: on the right shoulder-blade.
- Left arm electrode: on the left shoulder-blade.
- Right leg electrode: lowest intercostal space on the right posterior axillary line.
- Left leg electrode: lowest intercostal space on the left posterior axillary line.

INSTRUCTIONS

Clean the skin surfaces with a cotton flock soaked with ether. Spread conductive paste (type FIAB G005) on the above described skin areas with an amount enough to ensure a complete adherence of the metallic plate. In precordial areas shaves hairs if necessary. Put electrodes on the so treated area.

To ensure a stable position of the plate electrode, insert its special support in a hole of rubber belt, then place the belt around the limb and close it on the support choosing the most effective tension.

To connect the electrodes insert the plug of the patient cable in the socket of the electrode connection cable.

POSITIONING

Suction chest electrode (precordial derivations): place the electrode on the wished point by squeezing the rubber ball. Squeeze the rubber ball and press the electrode until a good adherence of the chest on the skin is obtained, then release the rubber ball.

Clamp electrode (peripheral derivations): place the clamp so that the conductive plate adheres on the skin of the arm or of the leg on the points previously selected.

Plate electrode: insert the special connection of the electrode in one of the holes of the rubber belt, roll the belt around the limb and fasten it on the electrode connection using the hole which gives enough tension.

ATTENTION

European and American standards use the following colours to identify patient cable leads:

- **EUROPEAN STANDARD:**
right arm: **RED** – left arm: **YELLOW** – right leg: **BLACK** – left leg: **GREEN**.
- **USA STANDARD:**
right arm: **WHITE** – left arm: **GREEN** – right leg: **BLACK** – left leg: **RED**.

WARNING

- **A non sterile device is supplied.**
- ECG devices and connection cables which will be connected to FIAB electrodes must be in accordance with regulations in force.
- Devices must be connected and activated by qualified personnel.
- Do not spread the electroconductive paste on wounded or scarred skin.
- Do not use these electrodes for monitoring during cardioversion or defibrillation procedures because time elapsing from D.C. shock to signal recovery is longer than 10 seconds.
- Do not use these electrodes for short or medium monitoring (monitoring during surgical procedures, Holter recording).
- The electrodes are in accordance with the biocompatibility ISO10993-1.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

CLEANING / DISINFECTION

To clean the electrodes, cleaning and disinfecting detergents commonly suitable for surgical devices are suggested. These may be, for example, benzalkonium chloride-based products. Before using these cleaning products read carefully their instructions. After cleaning wash with water.

N.B. Do not use ultrasonic cleaners. Don't use removers, don't sterilize with steam. It possible to clean with ethyl alcohol with concentration less than 10%.

PRODUCT LIFE

The Silver Silver chloride (Ag/AgCl) coating of electrodes is obtained with a galvanic proceeding. The silver silver chloride could be removed by a not suitable cleaning (for example by abrasive systems) and this could cause the registration of a non suitable ECG signal. In this case eliminate the products.

- Plate electrode: endless.
- Clamp electrode: intensive and protracted use eventually deteriorates the flexible part (spring).

- Suction electrodes: in presence of chaps on the rubber ball or failure to adhere to the skin, replace the rubber ball.

STORAGE

The product should be stored in its original packaging in rooms with ambient conditions of temperature and relative humidity, as specified on the package label. The overlap of weights on the packaging may damage the product.

GUARANTEE - NOTICE

FIAB SpA guarantees that the product complies with Directive 93/42/EEC and has been manufactured according to the procedures of FIAB Quality System certified ISO 13485. No responsibility may be ascribed to the producer who shall not be held liable for medical costs, direct or indirect damage due to lacking function or malfunction of the above products, when used differently from the instructions for use. We recommend to report opportunely any malfunction or defect of the product to FIAB Quality Assurance Service.

WASTE DISPOSAL

Waste coming from hospitals must be disposed of in accordance with regulations in force.

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