
## VITA UTILE

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
## USEFUL LIFE


The useful life of the device is established in years 5 from the commissioning. For this period Titanox s.r.l. guarantees the availability of spare parts and safe operation as long as the environmental and use conditions defined in the instructions for use are respected by the user..

# INSTRUCTIONS MANUAL

## DRY HEAT STERILIZING UNIT MOD. PASTEUR ELECTRIC FORCED HOT AIR

**CE**  
0476

**TITANOX** s.r.l.

FABBRICA ARTICOLI MEDICO SANITARI  
MEDICAL SANITARY ITEM FACTORY

26038 Torre dè Picenardi (CR) – ITALY – Via Canove, 2/A – Canove dè Biazzi

Tel. (0039) 0375 394065 (r.a.) – Fax (0039) 0375 394067

Home page: <http://www.titanox.com> – E-Mail [info@titanox.com](mailto:info@titanox.com)

## INTENDED USE


The device is designed to be used in sanitary environments for the purpose of sterilizing non-heat sensitive medical and surgical instruments. The process of disinfection is based exclusively on the thermal destruction of the microorganisms present upon the instruments at a typical temperature of about 180° C.

Effectively reaching a sterile state for the instruments inserted into the sterilizers depends upon multiple factors:

- The level of initial contamination of the instruments (total bacterial load);
- The type of contaminating microorganisms;
- The permeability of all the contaminated instrument parts by the hot air produced by the sterilizer.

For more information see the paragraph **\*\*attention\*\*** ahead.



 **The device must be used only and exclusively by qualified personnel.**

## INDEX

Pos. 1 Thermoregulator internal temperature 50 ÷ 210° C	Pos. 6 Warning light “resistances”
Pos. 2 Programming Timer 2 hours or Manual	Pos. 7 Ventilation motor fan (inside)
Pos. 3 Thermometer indicating internal temperature 0 ÷ 200° C	Pos. 8 Closure with key
Pos. 4 Connection plug with fuses	Pos. 9 Silicone gasket
Pos. 5 Warning light “timer” or “manual” connected	Pos. 10 Security thermostat 260° C

## SYMBOLS AND SAFETY WARNINGS

## WARRANTY CONDITIONS

- 1) The device is warranted for a period of one year from date of purchase.
- 2) Warranty covers the substitution or repairs free-of-charge of components with manufacturing defects.
- 3) The device will be repaired only at our factory. Charges, risks arising from the transport of the device shall be on purchaser's account.
- 4) In the event of repairs at purchaser's home, purchaser shall be charged fixed call costs covering partial reimbursement of travel and professional visit by personnel.
- 5) Warranty coverage excludes: internal lighting, damages caused by carelessness of purchaser, incorrect and improper uses and installations not conforming to warnings, indicated in these booklet instructions or in any case results from phenomenon unrelated to the normal working of the device.
- 6) The warranty expires when the device is tampered with or repaired by unauthorised personnel.
- 7) It is excluded the substitution and the extension of the warranty following a breakdown.
- 8) It is excluded any compensation for damages direct or indirect of any nature to persons or objects arising from use or suspension of use of the device.
- 9) The warranty expires immediately if the relative certificate shows alterations, erasing, or it's not issued or convalidated by us. The certificate must accompany the device, or handed to maintenance personnel for home-repairs.

The manufacturing company Titanox S.r.l. is responsible for the safety, reliability and performance of the device if:

- the assembly, the additions, the re-setting, the modification or repairs are carried out by personnel of Titanox S.r.l.;
- the electrical system to which it is connected conforms to safety norms in country of installation;
- the device is used in conformity to instructions of use and maintenance.

This liability expires immediately when the device is tampered with or repaired by unauthorised personal.

**For any further requirements of spare parts, repairs or checks, it's necessary contact directly the manufacturer: TITANOX S.r.l. - Via Canove 2/A – Loc. Canove de' Biazzi – 26038 Torre de' Picenardi (CR) – Italia - Tel. (0039) 0375 394065 – Fax (0039) 0375 394067 communicating the registration number of the device to repair.**

## ORDINARY MAINTENANCE

Before initiating any maintenance operation, make sure that:

- the device is not connected to the power source.
- the device is at ambient temperature.

Keep the internal parts and the grid perfectly clean. Even though they may change colour and become brownish, never use abrasive or inflammable products for cleaning them.

Keep the external parts perfectly clean that they have to be always specular to prevent corrosion and dust.

## PERIODIC MONTHLY MAINTENANCE

After removing the plug from the power socket, check that the fuses are not oxidised specially when the device is not used for a long time or kept in a humid ambient.

The power socket should not change colour or oxidise. If that happens, replace it immediately.

The power cable must be integral and it should not show cuts, abrasions or bending.

The resistances and the internal electrical system do not require any maintenance.

## PERMITTED USES

The device must be used to sterilize metal materials whose point of melting is higher than 300°C (surgical instruments, plates, and metal screw).

Inside the sterilizer place only metal containers without plastic parts and without textile material.

## FORBIDDEN USES

It is forbidden to place in the device items whose temperature of melting is lower than 300°C or it's not known to the operator.

## ACCESSORY PARTS

2 Internal shelves

1 Mains cable

## SELLING OFF

The sterilizing unit is made of various materials with mechanical, electro-mechanical and electronic parts.

The selling off has to be made according with the regulations in force in the utilizing Nation.

\*\*\*ATTENTION\*\*\*

**Laboratory tests following a cycle of 120 minutes at 180°C on Bacillus subtilis var niger ATCC 9372 spores demonstrated the efficiency of the sterilizer.**

Titanox does not guarantee and cannot ensure the effective achievement of the sterile state of the instruments placed in the sterilizer, according to the definition of sterile medical instruments foreseen by the EN 556 Norms. The user of the sterilizer therefore has the responsibility to conduct all of the confirmation procedures of the sterilization process and the necessary verifications in order to ensure the effective completion of every single sterilization cycle.

## TECHNICAL SPECIFICATIONS

<b>Models:</b>	<b>A3-213-400V</b>	<b>A3-214-535V</b>	<b>A3-215-670V</b>
Loaded Chamber	3 kg	4 kg	5 kg
<b>External sizes:</b>			
Width mm	570	705	835
Height mm	400	475	555
Depth mm	345	450	565
<b>Internal Sizes:</b>			
Width mm	405	535	670
Height mm	210	345	415
Depth mm	255	320	420
<b>Weights:</b>			
Net weight kg.	13 kg	22 kg	35 kg
Gross weight kg.	15 kg	25 kg	41 kg
<b>Electrical Characteristics:</b>			
Nominal tension (Voltage)	230 V	230 V	230 V
Nominal power (Watt)	450 W	950 W	1900 W
Nominal frequency (Hz)	50/60 Hz	50/60 Hz	50/60 Hz
Net's fuses (mm 5x20)	F5A-250 V	F5A-250 V	F10A-250 V
Ventilation motor (fan)	23 W	23 W	23 W

The device is in conformity to electrical safety norms provided for by the normative institutes and supplied with bipolar plug which assures perfect electrical grounding.

**NON-COMPLIANCE WITH INSTRUCTIONS DESCRIBED IN THIS BOOKLET SHALL FREE COMPANY TITANOX S.R.L. FROM ANY LIABILITY.**

<b>SYMBOLS</b>	<b>MEANING</b>
	Warning !
	See the annexed documentation
	Caution! Hot surface (max. 210 °C)
	It points the Manufacturer's name
	This symbol, together with the year, points the production date
	RAEE symbol to handle electrical and electronic devices wastes

Available operation:

**TIME-PROGRAMMED:** timer programmable up to 2 hours

**CONTINUOUS CYCLE:** manual (without timer)

## ENVIRONMENTAL CONDITIONS

- Ambient temperature from 5 to 40° C.
- Relative humidity max. 80% for temperatures up to 31° C with linear decrease up to 50% at the temperature of 40° C condensing included.
- Atmospheric Pressure from 500 to 1060 hPa.
- Voltage supply variation not higher than  $\pm 10\%$ .
- Value of transitory over-voltage in conformity to the installation category (which provides for 2500V).


## INSTALLATION

The device has been calibrated and tested at factory and as such does not require any further calibration or adjustments before installation and start-up.

Unpack the device and follow the next advices:

1. Position the device on a levelled flat, hard and smooth surface, established with non-inflammable material.
2. Allow space not lower than 10 cm from walls and surrounding furniture.
3.  Do not install the device nearby water sinks or similar to avoid contacts with water or substances which could cause electrical short-circuit to the system.
4. Install the device in a well aerated location, not near windows or external doors which could cause an unnatural circulation of the air inside the device and therefore compromise its correctly running.
5.  Do not install the device nearby heat sources or near other electrical devices.
6. Install the device in such a way that the power cable is never twisted or bent, but it should connect to the socket free and unhindered. Avoid to positioning the power cable nearby heat sources or other devices that shall caused damage on it in a long run.
7. In the event the device is placed on a trolley, always check that the lower part of the device is not enclosed or hindered, as to allow continuous and sufficient ventilation.

Once the device is correctly installed and power cable connected, it's ready for use.

## SAFETY WARNINGS

- The device is designed for use in internal locations.
- The device is not designed for use in presence of gasses or explosive vapours.
- No water or other liquids should be poured into device neither on its base.
- Before any cleaning or maintenance actions, the power cable must always first be removed.
- Make sure that the electrical system has electrical grounding and that it's in conformity to the safety norms in the country of installation.
- Do not remove any label or plate, in case of need, ask for more.
- Always insist on original spare parts.
- Do not open the door of the device until the thermometer of the internal temperature do not indicate a temperature lower than 30°C.
- Medical devices, before being placed in the sterilizer, should be washed and dried. The residual water steam may create deposits on the probes and resistances by altering the precision.

## TIME – PROGRAMMED RUNNING

1. Insert the supplied connection socket of the mains cable in the plug of the device (Pos. 4) and the feeder- plug in the wall mains socket, after checking the voltage value.
2. Turn the thermoregulator knob (Pos. 1) and set the desired temperature. This operation shall always be carried out when device is cold.
3. Place materials to be sterilized inside the device.
4. Lock up the door (Pos. 8).
5. Set, by turning the timer knob (Pos. 2) anti-clockwise, the sterilization time adding necessary time in order to bring the device and the material inside to the desired temperature. **The necessary time to reach the temperature of 180°C is about 30/40 minutes in case of empty chamber sterilizer, and about 60 minutes with Loaded Chamber well arranged on the trays.** This operation switch on the warning light (Pos. 5) indicating that the device is working and the warning light “resistance” (Pos. 6) indicating that the resistances are turned on in the circuit. The graduated thermometer (Pos. 3) allows to check the progressive increase of the internal temperature. When the internal temperature has reached the programmed value, the warning light “resistance” (Pos. 6) is turned off for then, after brief time, again turned-on indicating the re-connection of the resistances in the circuit. This operation allows to keep inside the device a constant temperature up to the end of the sterilization cycle. After programmed time has elapsed, the device turns itself off automatically.

In order to maintain an even temperature inside of the heating chamber and to allow a better circulation of the air, the device is featured with a ventilation motor fan internally wired.

## CONTINUOUS -CYCLE RUNNING

The device has also a continuous-cycle running. For use this kind of running proceed up to point 4 of the procedure for time-programmed running described above. Turn, at this moment, the timer knob (Pos. 2) anti-clockwise. The indefinite-time manual continuous-cycle is inserted with a click. To turn off the device, turn back the timer knob to position 0.

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