DISPOSABLE SENSOR

REF GSE0004 (GIMA 33754) GSE0004 (GIMA 33757)
 GSE0002 (GIMA 33755) GSE0002 (GIMA 35109)
 GSE0005 (GIMA 33756)

Contec Medical Systems Co., Ltd.
 No.112 Qinhuang West Street, Economic & Technical Development Zone,
 Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA
 Made in China

EC REP Prolinx GmbH, Brehmstr. 56, 40239
 Duesseldorf Germany

Imported by:
Gima S.p.A.
 Via Marconi, 1 - 20060 Gessate (MI) Italy
 gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Pulse Oximeter Probe user manual

Product name: Pulse Oximeter Probe

Product model: GSE0004 (GIMA 33754) GSE0002 (GIMA 33755) GSE0005 (GIMA 33756) GSE0004 (GIMA 33757) GSE0002 (GIMA 35109)

Scope of application:
 This product is used to match CONTEC Patient Monitor, Electronic Sphygmomanometer, Pulse Oximeter collect and transmit the SpO₂ signal from patient with continuance and no trauma. It is inapplicable to monitor the weak perfusion moving state and monitor for long, so check the measuring position or change for another position per 4 hours.

Taboo disease:
 Don't fix the product on the position with tissue injury. It is inapplicable for the patient or users allergic to PVC, TPU, TPE, ABS plastic.

Product performance:

- 1) The range of SpO₂ measurement: 70%~100%; Accuracy: 70~100%:±2%; Below 70%: unspecified.
- 2) The range of pulse measurement: 30~250bpm; Accuracy: ±2bpm or ±2%(select larger).
- 3) Optical Sensor:
 Red light (wavelength is 650~670nm, 6.65mW)
 Infrared (wavelength is 880~910nm, 6.75mW)

M33754-EM-Rev.2 - 10/24

Main configuration: Consisting of plug, cable and probe.

Power supply requirement: The special power is supplied by from the equipments of CONTEC Patient Monitor, Electronic

Sphygmomanometer, Pulse Oximeter which are applicable to the requirements of IEC60601-1.

Directions for use:

Note: This product is type BF applied part.

Sketch map	Model explanation	Applied crowds	Placement
	GSE0003 (2.3.10.00006) Adult disposable sensor	Weight>30Kg adult	Recommendatory placement: forefinger
	GSE0004 (2.3.10.00007) Pediatric disposable sensor	Weight10~50kg Pediatric	Recommendatory placement: forefinger
	GSE0005 (2.3.10.00008) Infant disposable sensor	Weight3~20kg Infant	Recommendatory placement: big toe
	GSE0002 (2.3.10.00009 & 2.3.10.00009) Neonatal disposable sensor	Weight<3kg neonate /weight>40kg adult	Recommendatory placement: neonate's foot, adult's forefinger

Figure 1

- 1) As **Figure 1**, the pulse oximeter probe of different types is applied to different crowds.
 - 2) Select proper probe and put recommendatory placement according to **Figure 1**.
 - 3) Arrange the cable along the back of hand when place the pulse oximeter probe.
 - 4) Connect pulse oximeter probe with pulse oximeter, Electronic Sphygmomanometer, or patient monitor and check if the operating procedure accords with the procedure introduced in user manual.
 - 5) Pulse Oximeter Probe GSE0002 (2.3.10.00002) needs the help of the FST0001 pulse Oximeter probe extension cable to be connected into the jack of the Pulse Oximeter CMS60D, CMS70A for normal use. Connect the FST0001 pulse Oximeter probe extension cable into the jack of the Pulse Oximeter CMS60D, CMS70A, then connect pulse oximeter probe to the other end of the FST0001 pulse Oximeter probe extension cable. Plug the sensor connector firmly into the FST0001 blood oximeter extension cable.
- Pulse Oximeter Probe GSE0002 (2.3.10.00002) needs the help of the FST0004 pulse Oximeter probe extension cable to be connected into the jack of the Electronic Sphygmomanometer CONTEC08A for normal use. Connect the FST0004 pulse Oximeter probe extension cable into the jack of the Electronic Sphygmomanometer CONTEC08A, then connect pulse oximeter probe to the other end of the FST0004 pulse Oximeter probe extension cable. Plug the sensor connector firmly into the FST0004 pulse Oximeter extension cable.
- Pulse Oximeter Probe GSE0002 (2.3.10.00002) needs the help of the FST0014 pulse Oximeter probe extension cable to be connected into the jack of the new Patient Monitor CMS8000 for normal use. Connect the FST0014 pulse Oximeter probe extension cable into the jack of the Patient Monitor CMS8000, then connect pulse oximeter probe to the other end of the FST0014 pulse Oximeter probe extension cable. Plug the sensor connector firmly into the FST0014 pulse oximeter extension cable.

Notice items:

- 1) pulse oximeter probe placement, the position without ductus arteriosus, BP cuff and vein input pipe is top-priority.
- 2) If the pulse oximeter probe can't monitor the state of pulsation, it shows that the position of probe is improper, or the position is too thick, too thin or having too deep pigment to reach a proper translucidity effect. If above things has happened, place the probe again or select probe of other type.

Explanation about graphs and symbols used on the product

	Follow instructions for use		Date of manufacture		Disposable device, do not re-use		Keep in a cool, dry place
	Keep away from sunlight		Type BF Applied Part		WEEE disposal	IPX2	Covering Protection rate
LOT	Lot number	SpO2	The pulse oximeter saturation (%)	BPM	Pulse rate (bpm)	REF	Product code
	Manufacturer	EC REP	Authorized representative in the community European	CE 0123	Medical Device complies with Directive 93/42/EEC		
	Temperature limit		Atmospheric pressure limit		Humidity limit		
MD	Medical Device	UDI	Unique device identifier	LATEX	Not made with natural rubber latex		

Guidance and manufacturer's declaration – electromagnetic emissions-for pulse oximeter probe

Guidance and manufacturer's declaration – electromagnetic emission		
The pulse oximeter probe is intended for use in the electromagnetic environment specified below. The customer of the user of the pulse oximeter probe should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment – guidance
Conducted and radiated RF EMISSIONS CISPR 11	Group 1	The pulse oximeter probe uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted and radiated RF EMISSIONS CISPR 11	Class B	
Harmonic distortion IEC 61000-3-2	Not Applicable	The pulse oximeter probe is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations and flicker IEC 61000-3-3	Not Applicable	

Guidance and manufacturer's declaration – electromagnetic immunity – for pulse oximeter probe

Guidance and manufacturer's declaration – electromagnetic immunity			
The pulse oximeter probe is intended for use in the electromagnetic environment specified below. The customer or the user of pulse oximeter probe should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields Should be at levels characteristic of a typical location in a typical commercial or hospital environment.

3) This pulse oximeter probe should be applied to the special medical equipment. Operator is responsible to check the compatibility. Incompatible fittings or device will influence the measuring result.

4) The disposal of scrap instrument and its accessories and packing (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.

Maintenance/cleaning/disinfection:

- 1) Check if the product is undamaged and clean before using.
- 2) This product is not allow to use disinfection liquid for disinfection, this probe belong to one-off products.

Note: Don't immerse the product in the liquid, and don't expose it under the strong ultra-violet radiation

Service life: Suggest this product use only once, don't use again.

Environment requirements:

- Transport and storage
- 1) Temperature: -20°C~+55°C
 - 2) Humidity: less than 95%
 - 3) Pressure: 500hPa~1060hPa
- Operating
- 1) Temperature: +5°C~+40°C
 - 2) Humidity: less than 90%
 - 3) Pressure: 700hPa~1060Pa

Statement:

- 1) pulse oximeter probe needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in User Manual and test report.
- 2) Portable and mobile RF communications equipment can affect pulse oximeter probe.

Warning:

- 1) The use of cables other than those specified, with the exception of cables sold by CONTEC as replacement parts for internal components, may result in increased emissions or decreased immunity of pulse oximeter probe.
- 2) pulse oximeter probe should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the pulse oximeter probe should be observed to verify normal operation in the configuration in which it will be used.
- 3) Improper usage can result in inaccurate measurement.
- 4) Using it under too strong light will cause inaccurate measurement, in case of that, please set a opaque stuff around the probe to cut light off.
- 5) You should move the probe to other position per 4 hours at least. Because the state of local skin can influence the ability of skin to enduring probe, it is necessary to replace the position of probe according to the state of patient. Please do that when skin integrity changes.
- 6) The dyestuff in blood vessel cab cause the inaccurate measurement.
- 7) The performance of pulse oximeter probe is influenced by movement easily, so it is not suitable for active patient to use it.
- 8) Don't fix the probe with belt or bundle it tightly, because the vein pulsation can cause inaccurate SpO₂ measurement.
- 9) Same as other medical equipment, the cable should be set properly to avoid entangling or asphyxiating patient.
- 10) Don't use it in the process of MRI scan, because the conductor current may burn the skin of patient, moreover, the probe will influence MRI image and MRI set will also influence the accuracy of SpO₂ measurement.
- 11) Don't change the product at will, otherwise the capability or accuracy of product will be influenced.
- 12) The probe is not intended for use during patient transport outside the healthcare facility.
- 13) DO NOT use the probe while the patient is being scanned by MRI or CT.

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

ARTICOLI & APPARECCHI PER MEDICINA

SENSORE MONOUSO

REF GSE0004 (GIMA 33754) GSE0004 (GIMA 33757)
GSE0002 (GIMA 33755) GSE0002 (GIMA 35109)
GSE0005 (GIMA 33756)

Contec Medical Systems Co., Ltd.
No.112 Qinhuang West Street, Economic & Technical Development Zone,
Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA
Made in China

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