User's manual

## **Head immobilizers and ankle immobilizers**

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## Index

en

General information	page	1
Warnings	page	1
Description of product	page	1
Operating instructions	page	2
Maintenance and cleaning	page	3
Spare parts	page	3

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
## 1. GENERAL INFORMATION

### 1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

### 1.2 Conservation of the instruction manual

The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

### 1.3 Symbols used

 General or specific warning

 See instructions for use

 Lot number

 Product code

 Product compliant with specifications of the Directive 93/42/CEE

 Single use

### 1.4 Servicing request

For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on tel. 0039 0521 541111, fax 0039 0521 541222, e-mail info@spencer.it or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY.

### 1.5 Demolition

Follow the current regulations.

### 1.6 Labelling

The serial number as indicated below can be found on each appliance and must not be removed or covered. In order to facilitate assistance please indicate or communicate the lot number (LOT) on the label.

## 2. WARNINGS

### 2.1 General warnings

- Before carrying out any kind of operation on the appliance, the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Regularly check the appliance.
- In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and the safety, the appliance must be immediately removed from service.

- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with. In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Handle with care.

### 2.2 Specific warnings

- The product must be used by trained personnel only.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- The head immobiliser or ankle immobiliser should not be exposed nor get in contact with heat sources or flammable agents.
- The device must be applied to the patient by at least two operators.
- Always check the integrity of all the parts of the device before use.
- Replace the belts immediately if worn out or damaged (if present).
- Select the fixation points of the belts accurately (if present).
- Position and regulate the device so not to obstruct any operations the rescuers may need to carry out and the use of the rescue appliances.

### 2.3 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

## 3. PRODUCT DESCRIPTION

### 3.1 Intended use

The Spencer head immobiliser and ankle immobiliser have been designed to increase the stability of the cranium and the spine and likewise the distal areas of the lower limbs during movement and transport of patients on the spine board. Their function is to eliminate the results caused any movement that may occur and in any direction with a guarantee of stability even if subjected to stress by vertical components (perpendicular to the fixing board). Designed to immobilise the head or the ankles when moving a patient who has suspected spinal injuries. The Spencer head immobiliser and ankle immobiliser, if used in combination with the spine board, restraint belts and rigid cervical collar, offer efficient immobilising which permits extended transport or transport over uneven surfaces even of patients in a critical state.

### 3.2 Main components (fig. 1 - 2 page I)

 Nylon belts to optimise stability of device

- B** Side support cushions with ear holes
- C** Semi rigid base for application of cushions with Velcro® system
- D** Removable cushion for use as paediatric head immobiliser
- E** Restraint belts for use with spine board
- F** Shaped surface for restraint of the ankles

### 3.3 Models

- SH00300A  
Fix Head Advanced universal head immobiliser, orange  
SH00310A  
Fix Head Advanced universal head immobiliser, black  
SH00201A  
Spencer Contour Anatomic universal head immobiliser, yellow/black  
SH00200A  
Spencer Contour Anatomic universal head immobiliser, blue  
SH00111C  
Super Blue Compact universal head immobiliser, yellow  
SH00112C  
Super Blue Compact universal head immobiliser, orange  
SH00117C  
Super Blue Compact universal head immobiliser, green  
SH00110C  
Super Blue Compact universal head immobiliser, blue  
SH00104C  
Super Blue Compact universal head immobiliser, black  
SH00240A  
Tango Fix Integrated adult/paediatric head immobiliser  
SH00250A  
U-Fix Adult head immobiliser  
SH00260A  
Pedi Fix Paediatric head immobiliser  
SH00262A  
Pedi Go Paediatric head immobiliser, only cushions, for Baby Go and Pedi Loc spine boards  
ST02605A  
Pedi Roll Paediatric head immobiliser  
SH00160A  
751 Inflatable disposable head immobiliser  
SH00130C  
755 Universal disposable head immobiliser  
SH00120C  
756 Stabilized disposable head immobiliser  
SH00150A  
FXA Universal ankle immobiliser  
SH00151A  
FXA Pro Anatomic universal ankle immobiliser

### 3.4 Technical data (chart 3.4 page 4)

## 4. OPERATING INSTRUCTIONS

### 4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packaged ensuring also that there are no risks of shocks, bumps or falls during the transport itself. Keep the original packaging for use in case of any further transport.

Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client.

The appliance must be stored in a dry place free from humidity. During storage take care not to put heavy materials onto the device. In no way and under no circumstances should the device be considered as a work top.

### 4.2 Preparation

On receipt of the product:

- remove the packaging and display the material so that all components are visible
- check that all the components/pieces on the accompanying list are present

The appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage.

Before putting the device in service make the following careful checks:

- general integrity (absence of cuts, holes, abrasions)
- general functionality
- Velcro® and belts tightness (if present)

If everything is regular, the device is ready to be used.

### 4.3 Functioning

#### 4.3.1 Application of head immobiliser Fix Head, Spencer Contour, Super Blue, Tango Fix, U-Fix, Pedi Fix (fig. 3 - 4 - 5 - 6 page II)

The product must be used by trained personnel only and must be applied to the patient by at least two operators.

Attach the device to the spine board before use using the belts supplied (E fig. 1 page I), making sure that the side of the spine board with Velcro® is topside. Check the position and the correct fixing of the belts before use with the patient. Load the patient on the spine board maintaining the correct alignment. Block the patient's head by placing the special anatomical cushions symmetrically on each side and without applying any pressure. Fix the two belts which restrain the head (A fig. 1 page I) by passing them over the anatomical cushions and uniting them with the buckles on the opposite side.

Proceed and attach the belts, adjust the tension and block them using the Velcro® system. Proceed in the same way with both belts making sure to block the chin strap first and then the forehead strap. Before transporting the patient it is necessary to proceed with the application of the special restraints for the spine board which will block the patients shoulders so as to avoid compression of the vertebrae in the neck area.

#### 4.3.2 Double use adult/paediatric head immobiliser: Spencer Contour and Tango Fix

Spencer Contour and Tango Fix head immobilisers can be used both for adult and paediatric patients.

They have a removable cushion in the centre of the base applied with Velcro® (fig. 7 page II). If using the head immobiliser with an adult patient, the cushion must be in position in the centre between the two cushions. When in use for a paediatric patient the cushion must be removed. This will also permit the two anatomical cushions to be positioned closer and the Velcro® permits fixing in the correct position.

The Tango Fix head immobiliser is characterised by the special "U" shape of the cushions. The cushion for immobilising the paediatric patient is situated inside the main cushion. Therefore in a confined space you have both adult and paediatric head immobiliser which will permit the correct immobilising for both paediatric and adult patient (fig. 8 page II).

### **4.3.3 Application of the Pedi Go, Pedi Roll, 751, 755, 756 head immobiliser**

The device must be applied only by trained operators. At least two operators are essential for the application.

The Pedi Go, Pedi Roll and 751 head immobilisers are extremely simple to use and to apply. When the patient is in the correct position on the spine board, maintain the correct spinal alignment and proceed to immobilise the head by positioning the cushions at the sides of the head. The Pedi Go and Pedi Roll should then be fixed to the spine board and to the patient's head using the special belts (fig. 9 page II).

The 751 head immobiliser can be fixed using 50 mm surgical tape (fig. 10 page II).

To immobilise the patient with the 755 and 756 an immobiliser base must first be placed on the spine board. Once the patient has been loaded and correctly aligned, with the patient's head positioned on the head immobiliser base proceed to block in position by lifting the extension on the sides until they are positioned on the sides of the head but do not exert any pressure. Fixing of the device should be carried out using the purpose made straps which block the device and the head in position and then maintain them in position on the spine board (fig. 11 - 12 page II).

### **4.3.4 Application of FXA, FXA Pro ankle immobiliser (fig. 13 - 14 page II)**

The device must be applied only by trained operators. At least two operators are essential for the application. Assemble the device on the spine board before use using the purpose made belts supplied (E fig. 1 page I). Load the patient in the correct alignment. Continue to position the ankles in the hollows on the cushion. Block the ankles in a stable, correctly aligned position using the top fixing wrap (fig. 13 - 14 page II).

## **4.4 Troubleshooting (chart 4.4 page 4)**

## **5. MAINTENANCE AND CLEANING**

### **5.1 Cleaning**

Failure to carry out cleaning operations can cause the risk of cross infection.

Do not use aggressive substances or solvents of any type. The device can be washed using delicate detergents or water and soap. At the end always dry the product with a soft and clean cloth. Using high pressure water will raise the risk of lesions. Ensure you have dried the device perfectly so to avoid the presence of residual detergents, which could deteriorate or compromise the integrity and duration of the device.

### **5.2 Maintenance**

#### **5.2.1 Precautionary Maintenance**

The person responsible for every day maintenance can only substitute the spare parts indicated on paragraph 6 "Spare Parts". All other substitutions or repairs can be carried out only by the manufacturer or by a centre authorised by the manufacturer. For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device does not require programmed servicing.

en

#### **5.2.2 Special servicing**

**Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations.**

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The device, if used as indicated in the following instruction manual, has an average life span of 5 years. The life span can be expanded only following a general revision of the product that must be carried out by the manufacturer or by a centre authorised by the manufacturer.

## **6. SPARE PARTS**

SH00303A

Base for orange Fix Head

SH00302A

Cushion for orange Fix Head

SH00301A

Kit head/chin straps for Fix Head

SH00242A

Base for Spencer Contour, Tango Fix, U-Fix and Pedi Fix

SH00214A

Cushion for yellow/black Spencer Contour

SH00212A

Cushion for blue Spencer Contour

SH00244A

Kit head/chin straps for Spencer Contour, Super Blue, Tango Fix, U-Fix, Pedi Fix, Pedi Go and Pedi Roll

SH00243A

Base for Super Blue

SH00100C

Cushion for yellow Super Blue

SH00101C

Cushion for orange Super Blue

SH00102C

Cushion for green Super Blue

SH00108C

Cushion for blue Super Blue

SH00103C

Cushion for black Super Blue

SH00241A

Adult/paediatric cushion for Tango Fix

SH00251A

Adult cushion for U-Fix

SH00261A

Paediatric cushion for Pedi Go

### 3.4 Head immobilizers and ankle immobilizers

Code	Fix Head SH0300A	Spencer Contour SH00200A	Super Blue SH00110C	Tango Fix SH00240A	Pedi Fix SH00260A	Pedi Roll ST002605A	FYA SH00150A	Ankle immobilizers SH00151A
Total weight (kg)	2,3	1,5	1	2,5	1,1	0,092	0,4	0,376
X-ray compatible	●	●	●	●	●	●	●	●
Immobilization directions	►▼►△	►▼►△	►▼►△	►▼►△	►▼►△	►▼►△	►▼►△	►▼►△
Floatation	●	All	All	●	●	Pedi Loc / Baby Go	All	All
Compatibility with spine board	All	All	All	All	Pedi Loc / Baby Go	Pedi Loc / Baby Go	All	All
Compatibility with scoop stretcher	●	●	●	●	●	●	●	●
Fixation straps	Velcro® / straps	2	2	2	2	2	1	1
Type of fixation	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps	Velcro® / straps
Reusable	●	●	●	●	●	●	●	●
<b>Base</b>								
Dimensions (mm)	375 x 10 x h275	390 x 15 x h260	384 x 21 x h290	445 x 28 x h280	445 x 28 x h280	345 x 120 x h90	345 x 120 x h90	345 x 120 x h90
Weight (kg)	1,2	0,5	0,3	0,5	0,5	0,4	0,4	0,4
Type	Semi rigid collection	Removable layer	Stuffed	Removable layer	Removable layer			
Material	Plastic compound	Expanded PE / Spandex®	Expanded PE / Spandex®	Polyethylene / rubber	Polyethylene / rubber	Closed cell expanded PE	Plastic compound	Plastic compound
Manufacturing system	Injection	Injection	Injection	Blow moulding / injection	Blow moulding / injection	Hand made	Injection	Injection
Colour	Orange	Blue	Blue	Black / yellow	Black / yellow	Blue	Orange	Black
<b>Head blocks</b>								
Dimensions (mm)	225 x 90 x h150	235 x 130 x h160	245 x 80 x h160	230 x 125 x h160	190 x 70 x h123	160x80xh80	160x80xh80	160x80xh80
Weight (kg)	1,2	1	0,5	0,5	1,2	0,6	0,6	0,6
Structure	Stuffed pillows	Stuffed pillows	Stuffed pillows	Rigid with soft pillows	Rigid with soft pillows			
Material	Plastic compound	Plastic compound	Plastic compound	Polyethylene / rubber	Polyethylene / rubber	Vinyl coated expanded		
Manufacturing system	Injection	Injection	Injection	Blow moulding / injection	Blow moulding / injection	Hand made		
Auricular inspection	●	●	●	●	●	●		
Colour	Orange	Blue	Blue	Black / yellow	Black / yellow	Blue	Orange	Black
Functioning temperature (°C)	from -20 to +40	from -20 to +40	from -20 to +40	from 20 to +60	from 20 to +60	from -20 to +40	from -20 to +40	from -20 to +40
Storage temperature (°C)	from -20 to +40	from -20 to +40	from -20 to +40	from 20 to +60	from 20 to +60	from -20 to +40	from -20 to +40	from -20 to +40
<b>PROBLEM</b>								
Excessive mobility of the head/ankles	The basis presents an excessive mobility in relation to the spine board it is strapped to							
Colour	Loose belts							
Functioning temperature (°C)	Anatomic supports positioned wrong							
Storage temperature (°C)	Wear and tear, improper use							
<b>CAUSE</b>								
Lesions to the device								

### 4.4

### REMEDY

Tighten the belts at the basis around the spine board	
Check the correct closure and tension of the fastening belts	
Ensure that the hollow part of the anatomic support is in contact with the head/ankles	
Put immediately the device out of service and contact the service centre	

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**FXA**

**FXA Pro**

We declare that the appliances conform with the Directive 93/42/CEE "Medical Devices".

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ISO 9001 Quality management systems - Requirements

EN ISO 13485 Medical devices - Quality management systems - Requirements for regulatory purposes

EN 980 Graphical symbols for use in the labelling of medical devices

CEI EN ISO 14971 Medical devices - Application of risk management to medical devices

Guarantee of Quality System for the production and the final control of the products certified by the notifying body TÜV SÜD Product Service GmbH.

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## **Warning**

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l. Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, layout or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.

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