# GIMA

PROFESSIONAL MEDICAL PRODUCTS

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## **DESCRIPTION**

This is a vacuum product which has a valve. After exhausting the air with pump, it will become solid to achieve the effect of restriction. Open the valve when not using it. Then it will fill with air and the vacuum mattress kit will be soft enough to put into a bag.

## **STEPS**

1. Open the valve, twist the knob to the “OPEN” setting.
2. Unfold it to distribute particles evenly throughout the product.
3. Insert the pipe of the suction accessory into valve and exhaust a little air first based on no move of the particles in the product. Twist the knob to the “CLOSE” setting.
4. Move the patient who need restrict onto the product and wrap up the patient. Restrict the safety straps.
5. Open the valve (twist the knob to “OPEN”), continue to exhaust air until the product become completely solid. Close the valve (twist the knob to the “CLOSE” setting) and remove the pipe.
6. There are handles for transportation around the product, they can be used to transit patient to stretcher or ambulance.

### **Warning**

1. The product must be used by trained personnel only
2. If any soften phenomena occurred during the restriction process, air should be pumped in sequentially and check the valve whether it is closed.
3. The device must always be checked before use to detect working anomalies or damages caused by transport or storage.
4. Two or more operators are needed to correctly use the product.

## **STORAGE AND MAINTENANCE**

1. Clean the product after every use.
2. Do not use aggressive chemical to clean the product.
3. Storage and operating temperature: could not be lower than -30 °C.

## Symbols

	Manufacturer		Product code
	Date of manufacture		Lot number
	Caution: read instructions (warnings) carefully		Keep in a cool, dry place
	Medical Device compliant with Regulation (EU) 2017/745		Keep away from sunlight
	Medical Device		Consult instructions for use

## GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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