# GIMA

PROFESSIONAL MEDICAL PRODUCTS

## **BARELLA A CUCCHIAIO SCOOPING STRETCHER CIVIÈRE DE RELEVAGE CAMILLA A CUCCHARA SCHAUFELTRAGHE MACA TIPO CONCHA**

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The atraumatic stretcher is made almost entirely of anodized anticorrosal aluminum suitable for resisting corrosion and at the same time to ensure the necessary strength and lightness of the whole.

The load-bearing elements are built using extruded aluminum tubes, of suitable dimensions to withstand the stresses caused by use. The various parts are assembled by means of steel rivets and in the least stressed points by aluminum rivets. The closing hooks are in galvanized iron and have been designed to ensure maximum safety and prevent accidental opening.

It can be divided longitudinally into two mirror parts, it is equipped with an extension system with relative locking systems to adapt it to the patient's height and is supplied with a set of 3 safety belts.

**THE CONSTRUCTIVE FEATURES OF THE ATRAUMATIC STRETCHER MEET THE REQUIREMENTS OF THE TECHNICAL STANDARD UNI EN 1865-1: 2015 AND ALL TECHNICAL TESTS HAVE BEEN DONE BY A NOTIFIED ORGANISATION WHO ISSUE THE "CE" CERTIFICATES.**

## INSTRUCTIONS

To load the injured person, one or both of the cast aluminum joints at the ends of the stretcher must be released. This operation is carried out by fully pressing button A. Open the stretcher and slide it under the patient, when closing the two poles, clic the shutters release as shown in fig. 2 and pay attention that the shutters are well closed.

To regulate the poles lenght rotate the small lever B up to 90° as fig. n. 2 and pull the poles as much as needed, then secure the small lever B and be sure that it reaches the catch holes.

Before moving the patient, it is advisable to secure him/her with the three fixing straps, located respectively at the level of the chest, of the pelvis and a few centimetres above the knee.

The stretcher can be stored open when at the minimum extension or folded in two.

To folded it, should pull completely the poles and bend them as seen in fig. 1.

For a correct use when lifting a patient, the stretcher should be pulled by the obtained sideways handles and not by the poles which could be dangerous for the patient and damage the stretcher.

To clean the stretcher use soap and water and/or eventually a suitable disinfectant.

At irregular intervals it would be convenient to oil the sliding poles with a little vaseline grease.

**CAUTION:** the utilisation of the stretcher should be done by skilled staff. Do not allow untrained persons to help during the use of the device, as they could cause injury to the patient or to themselves.

- Always check the integrity of the device before each use. In case of anomalies or damages which could compromise the functionality and safety of the device, and therefore of the patient and the operator, it is necessary to immediately remove the device from service.
- Should the device be malfunctioning, it is necessary to immediately use a similar product, in order to guarantee the continuity of the rescue operations.
- It is forbidden to use the device for any use other than that described in this manual.
- Do not alter or modify the device arbitrarily; modifications could cause unpredictable operation and damage to the patient or rescuers.
- The device must not be tampered with in any way (modification, adjustment, addition, repair), otherwise we are no longer responsible for its correct operation or for any damage caused by the device itself; furthermore, the CE certification and the product warranty are invalidated

## GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

	Caution: read instructions (warnings) carefully		Keep in a cool, dry place		Consult instructions for use
	Medical Device compliant with Regulation (EU) 2017/745		Keep away from sunlight		Product code
	Manufacturer		Medical Device		Date of manufacture
	Authorized representative in the European community		Lot number		

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