GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

SATURIMETRO PEDIATRICO OXY-PED OXY-PED PEDIATRIC FINGERTIP OXIMETER OXYMÈTRE PÉDIATRIQUE OXY-PED PULSOXIMETER FÜR KINDER OXY-PED SATURÓMETRO PEDIÁTRICO OXY-PED MEDIDOR DE SATURAÇÃO PEDIÁTRICO OXY-PED ΠΑΙΔΙΑΤΡΙΚΟ ΚΟΡΕΣΤΟΜΕΤΡΟ OXY-PED مقياس التأكسج خاص بالاطفال

Manuale d'uso e manutenzione
Use and maintenance book
Instructions de fonctionnement et entretien
Betriebs- und wartungs anweisungen
Manual de uso y mantenimiento
Manual de uso e manutenção
Εγχειρίδιο χρήσης και συντήρησης
دليل الإستعمال والرعاية

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.
ATTENTION: The operators must carefully read and completely understand the present manual before using the product.
AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.
ACHTUNG: Diese Anleitung muss vor dem Einsatz des Produkts aufmerksam gelesen und vollständig verstanden werden.
ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.
ATENÇÃO: Os operadores devem ler e entender completamente este manual antes de usar o produto.
ΠΡΟΣΟΧΗ: Οι χειριστές αυτού του προϊόντος πρέπει να διαβάσουν και να καταλάβουν πλήρως τις οδηγίες του εγχειριδίου πριν από την χρήση του.
الحذر: على العمال قراءة وفهم هذا الدليل بكامله قبل البدء باستعمال المنتج.

M34266-M-Rev.4.05.19

REF 34266

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Instructions to User

Read these instructions carefully before using this equipment. These instructions describe the operating procedures to be followed strictly.

Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Please peruse the relative content about the clinical restrictions and caution.

1. SAFETY

1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the oximeter.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The oximeter cannot be used together with the devices not specified in User's Manual.

1.2 Warnings

- Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- DO NOT use the oximeter while the testee measured by MRI and CT.
- To dispose the device, the local law must be followed.

1.3 Attentions

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it. When it is carried from cold environment to warm and humid environment, please do not use it immediately.

- DO NOT press the keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection to the oximeter is not permitted. Refer to User's Manual for instructions of cleaning and disinfection.
- DO NOT have the oximeter immersed in liquid. When it needs cleaning, please wipe its surface with disinfect solution by soft material. Do not spray any liquid on the device directly.
- The device should be kept out of the reach of children

2. OVERVIEW

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb of the blood, so-called the O₂ concentration in the blood. It is an important bio-parameter to the respiration. Many of the respiration disease will cause hypoxemia, even endanger the patient's health. As a result, monitoring the SpO₂ is indispensable in the clinical rescuing. The traditional method to measure SpO₂ is to analyze the sample of the patient's blood, so can get the partial pressure of oxygen and calculate the SpO₂ by use the bloodgas analyzer. This method can not be used to monitor continuously. For the purpose of measuring the SpO₂ more easily and accurately, GIMA developed the Fingertip Oximeter. The device can measure the pulse rate and blood perfusion index simultaneously.

The Fingertip Oximeter is compact, convenient to use and carry and with low power consumption. You just need to put the fingertip into the sensor of the device, the SpO₂ value will appear on the screen immediately.

2.1 Features

- Small in size and lightweight;
- Color OLED, various display modes, display directions adjustable;
- The device can accurately measure SpO₂ value, pulse rate value and perfusion index;
- The device will automatically start measuring after putting finger into rubber cushions;
- The device will power off automatically without signal for about 8 seconds;
- Audible & visual alarm function;
- Low voltage indication;

2.2 Major Applications and Scope

This device is applicable to home, hospital (including internal medicine, surgery, anesthesia, pediatrics, emergency room etc.), oxygen bar, the community medical center, alpine area and it also can be used before or after sports, and the like. The Fingertip Oximeter can detect SpO₂, pulse rate and blood perfusion index.

This device is not appropriate to be used for continuous monitoring.

2.3 Environment Requirements

Operating Temperature: 5°C~40°C

Operating Humidity: 30%~80%

Atmospheric pressure: 70kPa~106kPa

2.4 SpO₂ Common Knowledge

1. Meaning of SpO₂

SpO₂ is the saturation percentage of oxygen in the blood, so called O₂ concentration in the blood; it is defined by the percentage of oxyhemoglobin (HbO₂) in the total hemoglobin of the arterial blood. SpO₂ is an important physiological parameter to reflect the respiration function; it is calculated by the following method:

$$\text{SpO}_2 = \text{HbO}_2 / (\text{HbO}_2 + \text{Hb}) \times 100\%$$

HbO₂ are the oxyhemoglobins (oxygenized hemoglobin), Hb are those hemoglobins which release oxygen

2. Principle of Measurement

Based on Lambert-Beer law, the light absorbance of a given substance is directly proportional with its density or concentration. When the light with certain wavelength emits on human tissue, the measured intensity of light after absorption, reflecting and attenuation in tissue can reflect the structure character of the tissue by which the light passes. Due to that oxygenated hemoglobin (HbO₂) and deoxygenated hemoglobin (Hb) have different absorption character in the spectrum range from red to infrared light (600nm~1000nm wavelength), by using these characteristics, SpO₂ can be determined. SpO₂ measured by this oximeter is the functional oxygen saturation - a percentage of the hemoglobin that

can transport oxygen.

In contrast, hemoximeters report fractional oxygen saturation - a percentage of all measured hemoglobin, including dysfunctional hemoglobin, such as carboxyhemoglobin or methaemoglobin.

Clinical application of pulse oximeters: SpO₂ is an important physiological parameter to reflect the respiration and ventilation function, so SpO₂ monitoring used in treatment has become more popular. (For example, such as monitoring patients with serious respiratory disease, patients under anesthesia during operation and premature and neonatal infants) The status of SpO₂ can be determined in timely manner by measurement and will allow finding the hypoxemia patient earlier, thereby preventing or reducing accidental death caused by hypoxia effectively.

3. Factors affecting SpO₂ measuring accuracy (interference reason)

- Intravascular dyes such as indocyanine green or methylene blue.
- Exposure to excessive illumination, such as surgical lamps, bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight.
- Vascular dyes or external used color-up product such as nail enamel or color skin care.
- Excessive patient movement.
- Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- Exposure to the chamber with High pressure oxygen.
- There is an arterial occlusion proximal to the sensor.
- Blood vessel contraction caused by peripheral vessel hyperkinesias or body temperature decreasing.

4. Factors causing low SpO₂ Measuring value (pathology reason)

- Hypoxemia disease, functional lack of HbO₂.
- Pigmentation or abnormal oxyhemoglobin level.
- Abnormal oxyhemoglobin variation.
- Methemoglobin disease.
- Sulfhemoglobinemia or arterial occlusion exists near sensor.
- Obvious venous pulsations.
- Peripheral arterial pulsation becomes weak.
- Peripheral blood supply is not enough.

2.5 Caution

- A. The uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- B. For the individual patients, there should be a more prudent inspecting in the placing process. The device can not be clipped on the edema and tender tissue.
- C. The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance man, can not stare at the light.
- D. Testee can not use enamel or other makeup on the finger.
- E. Testee's fingernail can not be too long.
- F. The finger should be placed properly, or else it may cause inaccurate measurement.
- G. Do not shake the finger and keep at ease during using.
- H. Do not put wet finger directly into oximeter.
- I. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- L. Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO₂.
- M. Make sure the optical path is free from any optical obstacles like rubberized fabric.
- N. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- O. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.

3. TECHNICAL SPECIFICATIONS

A. Display mode:

OLED Display

B. Power supply requirement:

3.6V LiR2450 rechargeable cell×1

Supply voltage: 3.3V~4.2V

Continues working time: 4hours

C. Operating current: $\leq 40\text{mA}$

D. SpO₂ Parameter Specifications:

Measuring range: 35% - 100%

Accuracy: $\leq 3\%$ (during 70% - 100%)

E. Pulse Rate Parameter Specifications:

Measuring range: 30bpm-240bpm

Accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$ (which ever is greater)

F. Blood Perfusion Parameter Specifications:

Measuring range: 0%~20%

G. The performance under low perfusion condition

The accuracy of SpO₂ and PR measurement still meets the specification described above when the modulation amplitude is as low as 0.6%.

H. Resistance to ambient light interference:

The accuracy of SpO₂ and PR measurement still meets the specification described above when the device is tested by SpO₂ simulator (Fluke Biomedical Index 2 series) while setting the emulating interference of sun light and 50Hz/60Hz fluorescent light.

I. Dimensions: 56 mm (L) x 47 mm (W) x 32 mm (H)

Net Weight: 45g (including batteries)

4. ACCESSORI

A. A hanging cord

C. A charger

E. A charger cable

G. A User Manual

B. Two batteries

D. A power adapter

F. A pouch

H. Quality Certificate

5. INSTALLATION

5.1. Front view


Figure 1

5.2. Battery

1. Insert the coin cell into the battery compartment properly in the right direction (as shown in Figure 2).


Figure 2

2. Replace the cover.

- Please take care when you insert the coin cell, improper insertion may make the oximeter not work.

6. OPERATION

1. Open the clip as shown in Figure 3.

Figure 3 Put Finger into the Oximeter

2. Put finger into the rubber cushions of the clip (make sure the finger is in the right position), and then clip the finger.
3. The device will power on automatically in 2 seconds, and start to display software version number.
4. Next enter into data display screen (as shown in Figure 4). User can read the values and view the waveform from display screen.
“%SpO₂”: SpO₂ symbol; “99”: SpO₂ value;
“PR”: Pulse rate icon; “65”: Pulse rate value;
“♥”: Pulse beat symbol;
“|”: Pulse intensity histogram.


Figure 4

5. When the display is shown as Figure 4, press Display Key to switch display screen.

- Press Display Key once, display screen (as shown in Figure 4) will be flipped 180°.
- Press Display Key twice, display screen will be changed as Figure 5.
- Press Display Key three times, display screen (as shown in Figure 5) will be flipped 180°.
- Press Display Key four times, display screen will back to the screen as shown in Figure 4.
- Press Display Key circularly, display screen will be switched between the two screens (as shown in Figure 4 and Figure 5), and four directions display alternately.


Figure 5

6. Longtime press Display Key (about 2 seconds), display screen will be shown as Figure 6. Differences between Figure 6 and Figure 4 are as follows
- In Figure 4, SpO₂ and pulse rate are being monitored and displayed on the screen.
 - In Figure 6, SpO₂ and perfusion index are being monitored and displayed on the screen.


Figure 6

- When the screen displays as shown in Figure 6, press Display Key circularly, the display screen will be switched between the two screens (as shown in Figure 6 and Figure 7), and four directions display alternately.
- Display screen (as shown in Figure 6 or Figure 7) will return to the screen as shown in Figure 4 or Figure 5 if without operation in 10 seconds


Figure 7

9. Alarm Indicator

When measuring, if SpO₂ value and pulse rate value exceeds the preset alarm limits, the device will alarm automatically and the value on the screen exceeding limit will blink; at this time press Display Key to shut down the alarm. Exceeding pulse rate alarm limit: sound twice as an interval.

Preset alarm range:


SpO₂ alarm: Low limit: 90%

Pulse Rate alarm: High limit 120bpm - Low limit: 50bpm


7. INSTRUCTIONS FOR CHARGING THE CELL

The rechargeable coin cell is used to supply power for the oximeter. When it is run out, you need to use the provided charger to charge it, details is as follows:

1. Insert the cell into charger; pay attention to its polarity (as shown in the following figure).

2. Connect the charger and power adapter by charger cable in the following way.

3. The power adapter should be plugged into the AC mains with 100~250VAC, 50/60Hz.
4. Status of charging indicator
 - Green: full
 - Orange: no cell
 - Red: charging

Note:

- The charger will stop charging automatically when the coin cell is full (the indicator becomes green). Moreover, it has over heat protection function, when the coin cell temperature is higher than 45°C, it will also stop charging automatically.
- The charger can be connected to the USB port of a computer by the charger cable for charging.

8. REPAIR AND MAINTENANCE

- A. Please charge the battery when the low-voltage indicator lightens.
 - B. Please clean the surface of the device before using. Wipe the device with alcohol first, and then let it dry in air or wipe it dry.
 - C. Please take out the battery if the oximeter will not be used for a long time.
 - D. The best storage environment of the device is -20°C to 60°C ambient temperature, 10% to 95% relative humidity and 50kPa to 107.4kPa atmospheric pressure. Necessary servicing must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves
- The expected useful life (not guaranteed) of this device is 5 years.

High-pressure sterilization cannot be used on the device.

Do not immerse the device in liquid.

It is recommended that the device should be kept in a dry environment. Humidity may reduce the useful life of the device, or even damage it.

8.1 Battery Maintenance

Keep the both sides of coin cell clean.

- Low temperature may decrease the performance of coin cell, and low battery indicator may appear early. In such case, please put coin cell into pocket for warm before use, thus bring it back to normal condition.
- Do not let any conductive metal (such as tweezers) contact both sides of coin cell simultaneously to avoid short circuit.
- Charge the coin cell for 8~10 hours each time; ambient temperature should be 5°C~40°C.
- If the coin cell is full after charging, but its performance decreases apparently, it means the coin cell is exhausted, please change a new one.

9. TROUBLESHOOTING

Trouble	Possible reason	Solution
The SpO ₂ and Pulse Rate display instable.	1. The finger is not placed inside enough. 2. The finger is shaking or the patient is moving.	1. Place the finger properly and try again. 2. Let the patient keep calm.
The device can not turn on.	1. The coin cell is drained or almost drained. 2. The coin cell is not inserted properly. 3. The device's malfunction.	1. Charge the coin cell. 2. Reinstall coin cell. 3. Please contact the local service center.
The indicator light is off suddenly.	1. The device will power off automatically when it gets no signal for 8 secs. 2. The coin cell is almost drained.	1. Normal. 2. Charge the coin cell.

Declaration of Conformity:

The manufacturer hereby declares that this device complies with the following standards:

- IEC 60601-1:2005+A1: 2012,
- IEC60601-1-2:2014,
- IEC60601-1-11:2010, ISO 80601-2-61:2011 and follows the provisions of the council directive MDD93/42/EEC.

10. KEY OF SYMBOLS

	Keep in a cool, dry place		Medical Device complies with Directive 93/42/EEC
	Follow instructions for use		Manufacturer
	Serial number		Date of manufacture
	Product code		Type BF applied part
	Lot number		Keep away from sunlight
	WEEE disposal		Caution: read instructions (warnings) carefully
	Pulse oxygen saturation		Perfusion Index
	Pulse rate (beats per minute)		Low battery voltage

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons.

Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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$$SpO_2 = \frac{HbO_2}{(HbO_2 + Hb)} \times 100\%$$

HbO₂ are the oxyhemoglobins (oxygenized hemoglobin), Hb are those hemoglobins which release oxygen.

2. Principle of Measurement

Based on Lambert-Beer law, the light absorbance of a given substance is directly proportional with its density or concentration. When the light with certain wavelength emits on human tissue, the measured intensity of light after absorption, reflecting and attenuation in tissue can reflect the structure character of the tissue by which the light passes. Due to that oxygenated hemoglobin (HbO₂) and deoxygenated hemoglobin (Hb) have different absorption character in the spectrum range from red to infrared light (600nm~1000nm wavelength), by using these characteristics, SpO₂ can be determined. SpO₂ measured by this oximeter is the functional oxygen saturation - a percentage of the hemoglobin that can transport oxygen.

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- Intravascular dyes such as indocyanine green or methylene blue.
- Exposure to excessive illumination, such as surgical lamps, bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight.
- Vascular dyes or external used color-up product such as nail enamel or color skin care.
- Excessive patient movement.
- Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.
- Exposure to the chamber with High pressure oxygen.

- There is an arterial occlusion proximal to the sensor.
- Blood vessel contraction caused by peripheral vessel hyperkinesias or body temperature decreasing.

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