PALLONE RIANIMATORE MONOUSO
DISPOSABLE RESUSCITATOR BAG
INSUFFLATEUR DE REANIMATION MONOUSAGE
BALON DE REANIMACION MONOUSO
REANIMADOR DESCARTÁVEL
EINWEG-BEATMUNGSBEUTEL
ΑΝΑΖΩΓΟΝΗΤΗΣ ΜΙΑΣ ΧΡΗΣΗΣ
RESUSCITATOR FÖR ENGÅNGSBRUK
RESUSCITATOR DE UNICĂ FOLOSINȚĂ
ELDOBHATÓ ÚJRAÉLESZTŐ
ΡΕΑΝΙΜΑΤΟΡ ΖΑ ΕΔΗΟΚΡΑΤΗΑ ΥΠΟΤΡΕΒΑ

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PRODUCT NAME:

Besmed Disposable Resuscitator Set Model:

- RE-22415 - Single Use Resuscitator with 60 cm H₂O POP-OFF, 1600 ml, Adult
- RE-22515 - Single Use Resuscitator with 40 cmH₂O POP-OFF, 500 ml, Child
- RE-22615 - Single Use Resuscitator with 40 cmH₂O POP-OFF, 280 ml, Infant

PRODUCT DESCRIPTION

Besmed Disposable resuscitator set are manually operated and meant as an adjunct to artificial respiration and cardiopulmonary resuscitation. The resuscitator set can be used to ventilate apnoeic patients, and augment ventilation in spontaneously breathing patients undergoing oxygen therapy. The resuscitator set come in different sizes and are designed according to different compression capacities for infants, children, and adult patients.

Technical Specification:

Product Name	Disposable Resuscitator Set		
Manufacturer	Besmed Health Business Corp.		
Representative Photo			
Device type	Disposable		
Patient Usage type	Single use		
Patient group	Adult	Child	Infant
Accessories	Manometer, mask, reservoir, oxygen tube, PEEP valve and connector		
Pop-Off	W/ Pop-Off		
Deliverable pressure	60 cmH ₂ O	40 cmH ₂ O	
Bag volume	1600 ml	500 ml	280 ml
Reservoir volume	2500 ml		1000 ml
Suitable patient body weight	>40 kg	11 - 40 kg	5 - 10 kg

Stroke volume range; One hand	770 ml	300 ml	160 ml
Stroke volume range; Two hands	900 ml	350 ml	190 ml
Minimum delivered volume	>600 ml	>150 ml	
Maximum cycle rate	20 breaths/min		40 breaths/min
Oxygen concentration	With reservoir: > 85% - Without reservoir: > 35%		
Dead Space	< 5.5 ml		

Material	Resuscitator body:	PVC
	Mask:	PP + TPR
	Reservoir bag:	PVC
	Oxygen tubing:	PVC
	non-rebreathing valve:	PC
	Oxygen reservoir / gas intake socket:	PC
	Duckbill / gas intake / oxygen intake valve:	PVC
	Pop-off valve:	PC, PVC, stainless steel
	Manometer:	Silicone
	PEEP valve:	PC, Silicone, Stainless Steel (for spring)
Connector:	Silicone	
Package	Non-sterile package	
Operation temperature	25±5 °C, Room Temperature	
Storage temperature	-40 to 60°C, 30-60% RH	
Shelf Life	5 years	
Specifications	Patient connection port and face mask connectors	Connection stayed intact
	Dismantling and reassembling	Proper ventilation occurred after dismantling and reassembling the device
	Patient valve function after contamination with vomitus	There is no change in the performance of the device after contamination with vomitus.

Specifications	Drop test	No function variation in the device after dropping	
	Immersion in water	There is no change in the performance of the device after immersion in water.	
	Expiratory resistance	2.07 cmH2O	
	Inspiratory resistance	3.03 cmH2O	
	Dead space and rebreathing	4.78 ml	
	Minimum delivered volume	629.23 ml	
	Pressure limitation	There is no performance deviation.	
		The resuscitator without and with reservoir bag is more than 35 % and 85 %	
	Patient valve malfunction	1.07 cmH2O	
Attestation	Test item	Guideline	Result
	Cytotoxicity	ISO 10993-5	PASS
	Skin irritation test	ISO 10993-10	PASS
	Sensitization	ISO 10993-10	PASS
	Emissions of particulate matter	ISO 18562-2	PASS
	Emissions of VOCs	ISO 18562-3	PASS

INTENDED USER:

Medical personnel received the training of respiratory care.

INTENDED USE:

Besmed Disposable resuscitator set are designed for use as an adjunct to artificial respiration and cardiopulmonary resuscitation. The resuscitator can be used to ventilate the apneic patients and to augment ventilation and oxygen delivery to the spontaneously breathing patient.

DIRECTION FOR USE:

1. Prior to use, read the Instructions, cautions, and warnings.
2. Connect the oxygen supply tubing to a regulated oxygen source.

3. Adjust the gas flow so that the reservoir expands completely during inspiration and nearly collapses as the squeeze bag refills during exhalation.
4. Prior to connection to a patient, check the function of the resuscitator, preferably attached to a test lung, by observing that the intake, reservoir, and patient valves are allowing all phases of the ventilator cycle to occur.
5. As required, connect the resuscitation mask to the patient connector.
6. Follow accepted Advance Cardiac Life Support (ACLS) or Institution-approved for ventilation.
7. Compress the squeeze bag to deliver a breath. Observe the chest rise to confirm inspiration.
8. Release pressure on the squeeze bag to allow exhalation. Observe the chest fall to confirm exhalation.
9. During ventilation, check for:
 - Signs of cyanosis
 - Adequacy of ventilation
 - Airway pressure
 - Proper function of all valves
 - Proper function of reservoir and oxygen tubing
10. Should the non-rebreathing valve become contaminated with vomiting, blood, or secretion during ventilation, disconnect the device from the patient and clear the non-rebreathing valve as follows:
 - Rapidly compress the squeeze bag to deliver several sharp breaths through the non-rebreathing valve to expel the contaminate. If the contaminate does not clear.
 - Rinse the non-rebreathing valve in water and then rapidly compress the squeeze bag to deliver several sharp breaths through the non-rebreathing valve to expel the contaminate. If the contamination still does not clear, discard the resuscitator.

Warning:

1. Do not use the resuscitator set in toxic atmospheres.
2. Do not administer supplemental oxygen in the presence of open flames.
3. This device is intended for use by qualified medical and emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques.
4. Proficiency in the assembly, disassembly, and use of this device should be demonstrated before use on a patient.
5. Always monitor airway pressure with a manometer when ventilating a patient.
6. Only qualified personnel trained in the use of Positive End Expiratory Pressure (PEEP) should level the function of the resuscitator before using it on a patient.
7. If overriding the pressure relief valve, great caution must be taken not to allow the pressure in the patient's airways to become too high.

Caution:

1. If overriding the pressure relief valve, great caution must be taken not to allow the pressure in the patient's airways to become too high.
2. Do not attempt to disassemble the pressure relief valve assembly. Disassembly will damage the component.
3. Dispose of in accordance with local regulations procedures or hospital protocol

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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