# **CANNULA DI GUEDEL BESMED BESMED GUEDEL AIRWAY CANULE DE GUEDEL BESMED VÍA AÉREA BESMED GUEDEL**

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## **INSTRUCTION FOR USE OF BESMED GUEDEL AIRWAY**

### **1. INTENDED USE**

Besmed Guedel airway intended to provide a clear airway through the oral cavity and the pharynx by preventing blockage by the tongue. It is used to maintain the patency of the respiratory passages through the oral cavity and pharynx while patient remain in unconscious condition.

### **2. PRODUCT DESCRIPTION**

Besmed Guedel airway is an oropharyngeal airway to maintain or open a patient's airway. It preventing the tongue from covering the epiglottis, which could prevent the person from breathing. Flexible and anatomical shape to minimize airway trauma. In endotracheal intubation, guedel airway is inserted orally into the trachea to allow gas exchange, often via mechanical ventilation. The tube can be placed under direct visualization.

#### **2.1 Indications**

Besmed Guedel airway is intended to provide a clear airway through the oral cavity and the pharynx by preventing blockage by the tongue. It is used to maintain the patency of the respiratory passages through the oral cavity and pharynx while patient remain in unconscious condition.

#### **2.2 Contraindications**

- Avoid inserting a guedel airway in conscious or semiconscious patients because it may stimulate gagging and vomiting.
- Improper sizing airway, it can cause bleeding or close the glottis and airway.

#### **2.3 Intended User**

Medical personnel received the training of respiratory care.

#### **2.4 Risks/ Side Effects**

A potential risk associated with using guedel airways is the possibility of vomiting, which will lead to aspiration. It may also worsen airway obstruction if the airway is too small or inadequately sized (i.e., too small). Inappropriately sized airways may also cause laryngospasm (e.g., too large). The insertion of a guedel airway can also cause damage to the oral structure or dentition.

### 3. SPECIFICATION/MATERIAL/APPLICATION

#### 3.1 Device Family and technical specification

Representative Photo	Ref. No.	Product Description
	AW-61700	Guedel Airway, 40-110mm
	AW-61140	Guedel Airway, 40mm, Pink
	AW-61150	Guedel Airway, 50mm, Blue
	AW-61160	Guedel Airway, 60mm, Black
	AW-61170	Guedel Airway, 70mm, White
	AW-61180	Guedel Airway, 80mm, Green
	AW-61190	Guedel Airway, 90mm, Yellow
	AW-61110	Guedel Airway, 100mm, Red
	AW-61111	Guedel Airway, 110mm, Orange
	AW-61112	Guedel Airway, 120mm, Purple

#### 3.2 BESMED GUEDEL AIRWAY TECHNICAL SPECIFICATIONS

Product Name	Guedel Airway
Manufacturer	Besmed Health Business Corp.
Reference No.	AW-61700, AW-61140, AW-61150, AW-61160, AW-61170, AW-61180, AW-61190, AW-61110, AW-61111, AW-61112
Representative Photo	
Material	Airway body: LDPE
	Bite block: PP

Identification	Color coded			
Color	AW-61150	Blue	AW-61180	Green
	AW-61160	Black	AW-61190	Yellow
	AW-61170	White	AW-61110	Red
	AW-61112	Purple	AW-61111	Orange

Size		40 mm – 120 mm	
Usage type		Single use	
Patient population		Infant to Adult	
Package		Sterile package	
Sterilization method		Ethylene oxide	
Duration of contact		Limited ( $\leq$ 24 h)	
Specifications	Resistance to collapse of the buccal portion	83,46%	84,44%
	Patency of lumen	The steel ball ( $\phi$ 3.6 mm and $\phi$ 4.3 mm), diameter greater than 75% of minimum inside dimension of guedel airway, was able to pass through the guedel airways.	
Operataion temperature		25 $\pm$ 5 °C, Room Temperature	
Storage temperature		-20°C to 60°C	
Relative humidity		30 - 80 % RH	
Shelf life		5 years	
		Test Item	Guideline
Attestation	Functional test	ISO 5364:2016	PASS
	Cytotoxicity	ISO 10993-5	PASS
	Skin Irritation	ISO 10993-10	PASS
	Sensitization	ISO 10993-10	PASS
	Usability	IEC 62366-1	PASS

### 3.3 Application

Choosing the right size of the guedel airway by measuring against the side of the patient's face.

## 4. DIRECTION OF USE

4.1 Read this instruction manual carefully before use.

4.2 Wash hands with warm soapy water. Dry thoroughly. Put on gloves and aprons.

4.3 Patient should be assisted to a comfortable position, which is normally sitting or supported at an angle of 45°.

4.4 Adjust the height of the bed so you can work comfortably.

4.5 Select the proper size of the guedel airway by measuring from the first incisors to the angle of the jaw and select the suitable size to it to avoid over insertion.

4.6 Remove the package and insert into the patient's mouth upside down.

4.7 Once contact is made with the back of the throat, the airway is rotated 180 degrees, allowing for easy insertion.

4.8 Assure the tongue is secured.

### 5. WARNING/ CAUTION

- Read this instruction before use.
- Replace a new one when the package is open or component damaged prior to use.
- Failure to do so may result in serious injury to the patient and/or the user.
- Never use the product outside the recommended technical specifications (intended use).
- The guedel airway is a single use device. Please discard it after use.
- It is used only on unconscious patients. Patients who are conscious or semiconscious are not allowed to use it.
- Any serious incident that has occurred in relation to the device should be reported to Besmed or its Authorized Representative/competent authorities and the Regulatory Authorities of the country.

### 6. DISPOSAL

Dispose of in accordance with local regulations.

### 7. LEGEND OF SYMBOLS

	Medical Device compliant with Regulation (EU) 2017/745		Product code
	Importated by		Manufacturer
	Please read instructions carefully		Date of manufacture

	Lot number		Disposable device, do not re-use
	Expiration date (see box / package)		Keep in a cool, dry place
	Keep away from sunlight		Authorized representative in the European community
	Temperature limit		Sterile
	Medical Device		Disposable device, do not re-use
	Do not resterilize		Single sterile barrier system
	Not made with natural rubber latex		No-DEHP formulation
	Unique device identifier		

#### GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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