# GIMA

PROFESSIONAL MEDICAL PRODUCTS

**TEST STERILIZZAZIONE DI TIPO 5**  
**TYPE 5 STEAM INTEGRATOR**  
**INTÉGRATEUR À VAPEUR DE TYPE 5**  
**INTEGRADOR DE VAPOR TIPO 5**  
**INTEGRADOR QUÍMICO A VAPOR TIPO 5**  
**DAMPFINTEGRATOR TYP 5**  
**PARNÍ INTEGRÁTOR TYPU 5**  
**ΣΥΣΚΕΥΗ ΑΠΟΣΤΕΙΡΩΣΗΣ ΑΤΜΟΥ**  
**ΤΥΠΟΣ 5**  
**TYPE 5 STEAM INTEGRATOR**  
**ТΥΡПІН 5 НӨҮРҮІNTEGRAАТТОRІ**

**STEAM INTEGRATOR TIP 5**  
**5 TÍPUSÚ GŐZ INTEGRÁTOR**  
**INTEGRATOR PAROWY TYPU 5**  
**INTEGRATOR PENTRU STERILIZARE CU**  
**ABURI TIP 5**  
**TYPE 5 STOOMINTEGRATOR**  
**TYP 5 PARNÝ INTEGRÁTOR**  
**TYPE 5 DAMPINTEGRATOR**  
**ИП 5 ΠΑΡΕΗ ΙΗΤΕΓΡΑΤΟΡ**  
**5 TIPO STERILIZAVIMO TESTAS**

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**ENGLISH****Warnings:**

Upon receipt of goods please check whether the package was undamaged and ensure you about the expiry date of the materials.

**Introduction:**

This Indicator is composed of a special non-cellulose strip and a special indicating steam ink which is turning into end color under Steam sterilization conditions. The Indicator was developed based upon the use of non-toxic, heavy metal free inks. For more details please refer to the technical data sheet.

The dimensions of the Indicator is 55 mm x 25 mm.

The Type 5 Steam Integrator is giving a conclusion on the sterility of the container or pouch and its content.

**Use recommendations:**

This Indicator can be used inside pouches, trays and baskets wrapped in sterile barrier systems, as well as inside instrument containers to check on local Steam sterilization conditions.

**Intended use:**

This Indicator is recommended for the monitoring of Steam sterilization cycles.

**Warnings:**

Do not use this Indicator in other type of sterilization cycles with a different SV. (See actual product for the SV)

Any serious incident that has occurred in relation to the Steam Integrator should be reported to the manufacturer.

**Instructions for use:**

Place the Indicator inside the pouch, pack, basket, tray or container together the materials which will be sterilized.

Keep the Indicator away from the bottom of the object to prevent condensate to effect the results.

In rigid systems, as containers, place the indicator in each corner, or at least in two diagonally opposite corners to assure that Steam has reached the most critical points at those places. Place pouches sideways / vertically in order to avoid the plastic film side to block the Paper back side which might influence the effectiveness of the sterilizing agent.

**Results interpretation:**

If the indicator changes to the reference end color printed on the product, it indicates that the process reached the condition and the sterilization process was successful. If one or more sterilization parameters were not accomplished the indicator is not reaching its full end color. In the latter case the materials / instruments cannot be considered sterile and must be processed again.

**Storage conditions:**

Please take into consideration the storage conditions provided for at the box label of the product. Only in case the storage conditions were met, SP Medikal guarantee products compliance to standard within the expiry date.

Do not expose the Steam Integrators to sterilizing agents or other chemical products of any kind.

**Compliance to standards:**

The Indicators are compliant to standard ISO 11140 - 1: 2014.

**Shelf life:**

The Indicators have a limited shelf life and needs to be stored in the recommended conditions. The manufacturing date as well as lot number are indicated on the packaging label.

**GIMA WARRANTY TERMS**

The Gima 12-month standard B2B warranty applies.

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