Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403 Export: tel. +39 02 953854209/221/225 fax +39 02 95380056 gima@gimaitaly.com - export@gimaitaly.com www.gimaitaly.com

PROFESSIONAL MEDICAL PRODUCTS

TEST PER LA VISIONE STEREOSCOPICA – TEST FOR STEREOCOPIC VISION – EXAMEN DE LA VISION STÉRÉOSCOPIQUE – TEST FÜR STEREOKOPISCHES SEHEN – TEST PARA LA VISIÓN ESTEREOCÓPICA – TESTE PARA A VISÃO ESTEREOCÓPICA – ΤΕΣΤ ΣΤΕΡΕΟΣΚΟΠΙΚΗΣ ΟΡΑΣΗΣ –

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ENGLISH

Evaluates the cooperation between both eyes in forming a single image, starting from the two images that are supplied separately to the brain by the two eyes.

Can be used for children over one year of age and adults.

METHOD AND HOW TO USE:

The test consists of a card showing 4 shapes:

a ball, a fish, a star and a tree.

Familiar images have been used, which can be recognised and named both by adults and by very young children.

The person being tested is asked to indicate what he/she sees in the test and his/her eye movements are observed.

The ball is the only element visible using just one eye. The other shapes (each of which is placed at different depth levels) are perceived only by someone with binocular vision and visual acuity over 2-3/10. For adults and cooperative children the card should be held at a distance of 30-40 cm, with a request to indicate the shapes seen in it. If only the ball is seen, the result is pathological. If all four shapes can be identified, then the result is normal.

Briefly:

Positive result:

All objects are identified and named

Unclear result:

Only one of the images is recognised

Negative result:

None of the three-dimensional images is recognised.

INDICATIONS FOR EVALUATION OF THE TEST:

Suspected strabismus, amblyopia or reduced sight in one or both eyes or as a screening test for stereoscopic vision.

If a child does not cooperate, the test can be shown and the tester can observe whether the child tries to grasp one of the objects apart from the ball (the fish, the star or the tree). If this occurs, the test is normal.

If it does not happen, the test cannot be considered pathological for certain and may be repeated some time later.

Important: in the event of a negative or unclear result, it is advisable to carry out a complete eye examination.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the

warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic elds, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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