# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy  
Italia: tel. 199 400 401 - fax 199 400 403  
Export: tel. +39 02 953854209/221/225 - fax +39 02 95380056  
gima@gimaitaly.com - export@gimaitaly.com  
[www.gimaitaly.com](http://www.gimaitaly.com)

## **VIDEO STETOSCOPIO LYNX MULTIFUNZIONE + SpO<sub>2</sub> LYNX MULTI-FUNCTIONAL VISUAL STETHOSCOPE + SpO<sub>2</sub> STÉTHOSCOPE MULTIFONCTIONS VISUEL LYNX + SpO<sub>2</sub> VIDEO ESTETOSCOPIO LYNX MULTIFUNCIÓN + SpO<sub>2</sub>**

**Manuale d'uso - User manual  
Notice d'utilisation - Manual del usuario**

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## **USER NOTICE**

### **Announce**

This user manual is written and compiled in accordance with the council directive MDD93/42/EEC for medical devices and harmonized standards. This user manual includes special documents, which are under protection of copyright law. All rights reserved. Without written announcement from our company, the user manual should not be transferred, copied or translated into other language.

Owing to the forthcoming renovation or error contained in the manual, the specific products you received may not be totally in accordance with the description of this User Manual. We would sincerely regret for that. Our company reserves the final elucidative right.

### **Manufacturer's responsibility**

Our company is responsible for safety, reliability and performance of this equipment only in the condition that:

All assembly operations, expansion, change, modification and repairs are conducted by persons authorized by our company. The device is operated under strict observance of this manual.

- ⚠ Warning: This device is not intended for treatment. If the result is distrustful, please use other methods to verify immediately.
- ⚠ Warning: This device can't be used as normal ECG machine, its ECG function only must be used to detect the ECG waveform.

### **Warranty**

The unit can not be repaired by users themselves. All services must be done by the engineers approved by manufacturer. We warrant that each product we sell you is free from defects in labour and materials and shall conform to its product specifications as defined in the user documentation. If the product doesn't function as warranted during the warranty period, we will repair or replace it without charge. Misuse, improper maintenance may void the warranty.

### **Label presentation**

- ⚠ Warning: the information that should be noted to avoid injury to the patient and the operator.
- ⓘ Attention: important information that you should know

# CONTENT

<b>Chapter 1 Safety Guidance .....</b>	<b>28</b>
1.1 Safety Operation Explaining .....	28
1.2 Warning .....	28
1.3 Attention .....	29
<b>Chapter 2 Brief Introduction.....</b>	<b>30</b>
2.1 Overview.....	30
2.2 Features.....	30
2.3 Main Performance .....	30
2.4 Accessories .....	31
<b>Chapter 3 Appearance and Structure .....</b>	<b>31</b>
3.1 Front Panel Introduction.....	32
3.2 Side Panels Introduction .....	33
3.3 Rear Panel Introduction.....	34
<b>Chapter 4 Operating Guide .....</b>	<b>34</b>
4.1 Monitor Operation .....	34
4.2 Menu Operations:.....	36
4.3 Alarm Operation .....	39
4.4 Charge Operation .....	40
<b>Chapter 5 Product Specification.....</b>	<b>40</b>
<b>Chapter 6 Cleaning/Maintenance and Disinfecting .....</b>	<b>41</b>
6.1 Cleaning .....	41
6.2 Disinfecting.....	42
6.3 Maintenance.....	42
<b>Chapter 7 Key of Symbols.....</b>	<b>43</b>
<b>Chapter 8 Trouble Shooting .....</b>	<b>44</b>
<b>Appendix .....</b>	<b>45</b>

# Chapter 1

## SAFETY GUIDANCE

### 1.1 Safety Operation Explaining

This unit is internally powered equipment, the degree of shock protection is type B applied part  .

Type B applied part protection means that these patient connections will comply with permitted leakage currents, dielectric strengths of IEC 60601-1.

To avoid any possible danger, please operate the device according to safety guidance as following:

### 1.2 Warning

- The infrared is harmful to the eyes, so the user and the maintenance man should not stare at the light part of the SpO<sub>2</sub> probe (the infrared is invisible).
- For the special patients, there should be a more prudent inspecting in the placing process. The SpO<sub>2</sub> probe can not be clipped on the edema and tender tissue.
- Uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- Please don't use this device when charging.
- Explosive hazard - DO NOT use the device in environment with inflammable gas such as some ignitable anesthetic.
- Please don't use this device when MRI or CT scanning.
- The person who is allergic to rubber can not use this device.
- The user must use the oximeter probe which is provided by our company, don't use the parts which are not manufactured by our company.
- The device don't suit all user, if you can't get steady value, please stop using it.
- This device is not intended for treatment.
- The device can not be transported mixed with toxic, harmful, corrosive material.

### 1.3 Attention

- 🔔 Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- 🔔 The device may be a little different from the picture of manual, please make actual device as standard.
- 🔔 Before using, check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and performance. It is recommended that the device should be inspected once a month at least. When there is obvious damage, stop using the device before changing the damaged parts.
- 🔔 Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- 🔔 The device is designed for continuous operation and is 'ordinary'. If the device gets wet, please stop operating it.
- 🔔 Keep the device away from dust, vibration, corrosive substances.
- 🔔 Electromagnetic Interference-Ensure that the environment in which the device is operated is not subject to any sources of strong electromagnetic interference, such as MRI scan, mobile telephones, etc. Keep them far away.
- 🔔 Since the friction of body or clothes will bring high voltage static, and interfere the veracity of the tested data by the apparatus, it is avoided to be used in the environment of high voltage static.
- 🔔 This product is calibrated before leaving factory. Please don't measure this device with functional tester for the device's related information.
- 🔔 The device has the function of limits alarming, when the measured data is beyond the highest or lowest limit, the device would start alarming automatically on the premise of the alarming function is on.
- 🔔 The device has the function of alarming, this function can either be paused, or closed (default setting) for good, this function could be turned on through menu operation if you need. Please check the chapter 4.3 as a reference.
- 🔔 When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- 🔔 DO NOT operate keys on front panel with sharp materials.
- 🔔 When cleaning and sterilization, please refer to user manual related chapters for details.
- 🔔 The disposal of scrap instrument and its accessories and packings (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- 🔔 The device has normal useful life of three years since the first electrified use.

- 🔔 After the service life, please return the products to the manufacture or disposal the products according to local regulations.
- 🔔 Please refer to text for else notice proceeding.

## **Chapter 2**

### **BRIEF INTRODUCTION**

#### **2.1 Overview**

The device is a kind of multi-functional visual stethoscope. It is added with ECG function on the basis of the traditional medicinal stethoscope. It can display a real-time ECG waveform and calculate the heart rate, if you connect SpO<sub>2</sub> probe to this instrument, the pulse oxygen Saturation and pulse rate can also be monitored just through one finger. Doctors can gain more information about the patients. This stethoscope used high screen size LCD (liquid crystal display), which can display the realtime waveform and date.

The product is applicable for family, hospital, private Clinique, community medical treatment, oxygen station, sport health care (it is advised to be used before and after sport, and it is not advised to be used in the process of sport) and so on. Its multi-function is effective on emergency treatment.

#### **2.2 Features**

- Operation of the product is simple and convenient.
- Light weight, small size, easy to carry.
- Low power consumption.
- Multi-function, know more physiological information.

#### **2.3 Main Performance**

- Auscultation sound output
- ECG waveform display
- Heart rate display
- SpO<sub>2</sub> monitor through connecting with an external oximeter probe
- SpO<sub>2</sub> value display
- Pulse waveform, pulse rate and bar graph display
- Pulse rate sound indication
- Battery power display
- Alarm function
- Background light setting
- Contrast value setting

## 2.4 Accessories

- A User Manual
- A power adapter
- A data line
- An oximeter probe (type BF applied part )

**Note:** Type BF applied part means that the device has specifically safeguard for electroshock. Type BF applied part protection means that these patient connections should comply with permitted leakage currents, dielectric strengths of IEC 60601-1.

## Chapter 3 APPEARANCE AND STRUCTURE

Fig. 3-1 Front panel view


Fig. 3-2 Side panels view


Fig. 3-3 Rear panel view

### 3.1 Front Panel Introduction

ECG monitor interface is as Fig. 3-4:



Fig. 3-4 ECG monitor interface

SpO<sub>2</sub> monitor interface is as Fig. 3-5:

Fig. 3-5 SpO<sub>2</sub> monitor interface

Monitor interface is as Fig. 3-6:

Fig. 3-6 Monitor interface

## 3.2 Side Panels Introduction

### 1) Power Switch

function: Press the “Power switch” to turn on/off the device. This function is ineffective in menu.

### 2) Wheel Key

function: Menu operation (Alarm value setting, Contrast value setting, alarm and pulse control), Background light setting.

Menu operation: On any interface, press the wheel key to enter the main menu (consult 4.2 for particular content).

Background light setting: On any interface, move the wheel key to the left to open or close background light.

### 3) USB Socket: %SpO<sub>2</sub> /

function: 1.The socket for connecting SpO<sub>2</sub> probe  
2.The socket for charge

### 4) Charging Indicator Light

function: The indicator light shows green in charging state. It goes out when the charge has finished.

### 5) Earpiece

function: monitor auscultation-sound.

### **3.3 Rear Panel Introduction**

#### **1) Stethoscope Probe**

function: collect heart-beat signal.

#### **2) Electrode**

function: collect electrical signal of the body.

## **Chapter 4 OPERATING GUIDE**

### **4.1 Monitor Operation**

#### **4.1.1 Auscultation and ECG Monitor**

- 1) Let the patient laying down.
- 2) Daub some medical conductive paste to the surface of the three electrodes well-proportioned.
- 3) Turn on the device, it will enter into the ECG measurement interface (do not connect the oximeter probe to the device)
- 4) Hold the device and stick to the chest of patient. All the electrodes should stick to the kin upon heart. It can display the real-time ECG and heart rate while monitoring each auscultation section by earpiece (the device has the ability of automatic magnification, ECG waveform can be magnified automatically).
- 5) When the device is testing,press the “Power switch” shortly to pause the ECG waveform , so the operator can observe it carefully. Press it again to show the waveform again.
- 6) Press the “Power switch” Long to turn off the device after measuring.

#### **Attention:**

- Please check the device before using, and confirm that it can work normally.
- When testing the ECG, keep the device in tranquil environment.
- The device should be in a proper position, or else it may result in inaccurate measure.
- When testing the ECG, the electrodes must be assured to connect with the body thoroughly and stably.
- The device should not be used at the same location or limb which has arterial catheter or blood pressure cuff or receiving intravenous injection.
- If only to use the auscultation function, medical conductive paste is not needed, you can measure directly through the sensor which connects with the body thoroughly and stably.

### 4.1.2 SpO<sub>2</sub> and Pulse Rate Measurement

- 1) Connect the oximeter probe to the device. The probe is limited to be produced by our company; never replace it with the similar ones by other manufacturers.
- 2) Put the finger into the probe.
- 3) Press the “Power switch” long to open the device, it will enter the SpO<sub>2</sub> measurement Interface.
- 4) Under measuring interface, the user could read the correlative data from display screen.
- 5) Press the “Power switch” Long to turn off the device after measuring.

Fig. 4-1 SpO<sub>2</sub> monitor interface

#### Attention:

- ⚠ Please check the device before using, and confirm that it can work normally.
- ⚠ The fingernails of the testee should not be too long.
- ⚠ Testee can not use enamel or other makeup.
- ⚠ As to the fingers which are too thin or too cold, it would probably affect the normal measure of the patients' SpO<sub>2</sub> and pulse rate, please clip the thick finger such as thumb and middle finger deeply enough into the probe.
- ⚠ Do not shake the finger and keep the patient in a stable state during using process.
- ⚠ Fingernails and the luminescent tube should be on the same side (as Fig. 4-1).
- ⚠ Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight, etc.
- ⚠ There shouldn't be rubber fabric and so light barrier on the way of light, or else it may result in venous pulsation and inaccurate measure of SpO<sub>2</sub> and pulse rate.
- ⚠ The SpO<sub>2</sub> sensor and photoelectric receiving tube should be arranged in a way with the subject' s arteriole in a position there between.
- ⚠ The SpO<sub>2</sub> sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.

- 🔔 Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
- 🔔 If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
- 🔔 The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
- 🔔 Please read the measured value when the waveform on screen is equally and steady-going, this measured value is optimal value. And the waveform at the moment is the standard one.
- 🔔 The problem of overrating would display when the patient is suffering from toxicosis which caused by carbon monoxide, the device is not recommended to be used under this circumstance.

#### **Clinical restrictions:**

- As the measure is taken on the basis of arteriole pulse, minimal pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO<sub>2</sub> waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO<sub>2</sub> determination by this monitor may be inaccurate.
- The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO<sub>2</sub> measure.
- As the SpO<sub>2</sub> value serves as a reference value for judgement of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO<sub>2</sub> measurement

#### **Note:**

- 🔔 For the particular content of clinical restrictions and taboo disease, please read interrelated medicine literature as reference.

## **4.2 Menu Operations:**

After turning on the device, press the wheel key down under any test condition to enter the main menu as Fig.4-2. On the main menu interface, by moving the upper wheel key to the left or right can move the menu choice bar. Then press the wheel key down to choose corresponding option:


Fig.4-2 Main Menu

### 1) Contrast value setting

On the main menu interface, by moving the upper wheel key to the left or right can move the menu choice bar to “Alarm”, Then press the wheel key down, a small triangle will display on the right of “Alarm” (see figure 4-2), then move the wheel key to the left of right again to enter the alarm setting menu of figure. Choose “contrast” option, the parameter shows “025” (see figure 4-3), then move the wheel key to the left or right,each move of wheel key, the data will raise or descend for one time accordingly(setting range: 15~35, default is 25, the value will turn to 25 when the setting value is beyond the limits. After setting,press the wheel key down to save and exit the option.




Fig. 4-3 Contrast value setting

### 2) Alarm setting

Choose “Alarm” option, then press the wheel key down, a small triangle will display on the right of “Alarm” (see figure 4-4a), here move the wheel key to the left or right again to enter the alarm setting menu of figure 4-4b.

a

b

Fig. 4-4 Alarm setting

Choose the option which need change for setting, each move of wheel key, the data will raise or descend for one time accordingly. After setting, press the wheel key down to save and exit the option.

Number	Function	Explaining
1	SpO <sub>2</sub> alarm high-limit	Setting range: 0~100, default 99 (Unit: %)
2	SpO <sub>2</sub> alarm low-limit	Setting range: 0~100, default 85 (Unit: %)
3	PR alarm high-limit	Setting range: 0~254, default 120 (Unit: bpm)
4	PR alarm low-limit	Setting range: 0~254, default 40( Unit: bpm)
5	Exit	Exit alarm setting menu

### 3) The alarm state setting

Choose “Alarm sound” option, and change the state:choose “On” to open the alarm and choose “Off” to close the alarm. After setting,press the wheel key down to save and exit the option.

### 4) Pulse sound setting

Choose “Pulse sound” option, and change the state:choose “On” to open the alarm and choose “Off” to close the alarm. After setting,press the wheel key down to save and exit the option.

### 5) Enter monitor interface

Choose “Monitor” option,change the state:choose “On”, enter monitor interface after exiting menu as Fig. 4-5; choose “Off”, exit monitor interface after exiting menu.


Fig. 4-5 Monitor interface

## 6) Exit

Choose “Exit” option to exit menu.

## 4.3 Alarm Operation

Alarm includes the alarm of measure data's going beyond the limits, the alarm of low-voltage, the alarm of finger's out of position, the alarm of probe's out of position. Alarm PRI: the alarm of low-voltage > the alarm of finger's out of position = the alarm of measure data's going beyond the limits = the alarm of probe's out of position.

### Alarm function:

The alarm of low-voltage: When there is alarm of low-voltage, the battery mark will flicker; when the voltage is so low that the device can't work in gear, the device will turn off automatically after displaying “Low Power”.

The alarm of finger's out of position: when there is sound-alarm of finger's out of position, the screen will show “Finger Out”.

The alarm of probe's out of position: when there is sound-alarm of probe's out of position, the screen will show “Probe Out”.

If the sound alarm occurs, and the screen shows none of the three, it means the measure data's going beyond the limits.

### Alarm control:

When alarm is on, or during the period of alarming, you can pause it by moving the wheel key to the right, moving the wheel key to the right again, alarm resumes.

If you use the shortcut key to make alarm pause, the alarm will be resumed in about 60 seconds automatically. If you move the upper wheel key to the right again in 60s, it can renew alarm as well (under test condition).

If you want to close alarm for a long time, you must enter the menu and set the alarm state “off”.

## 4.4 Charge Operation

There are two kinds of charging methods:

- 1) Connect the device with computer by USB, then the device should be under charging state.
- 2) Connect the device with power supply by power adapter, then the device should be under charging state.

**Charging indicator light:** When the device is in the state of battery charging, the indicator light is on, when the battery power is full, the indication light would be off accordingly.

 For measuring exactly, when the battery energy show shortage, please recharge battery immediately.

## Chapter 5 PRODUCT SPECIFICATION

**Product Name:** Multi-functional visual stethoscope

**Safety:** Complies with IEC60601-1: 2005

**Classification:**

EMC: Group I Class B

Anti-electroshock Type: Internally powered equipment.

Anti-electroshock Degree: type B applied part 

Harmful Liquid Proof Degree: IPX1

Degree of Safety in Presence of Flammable Gases: Equipment not suitable for use in presence of flammable gases.

Working System: Continuous running equipment.

**Physics characteristic:**

Size: 87mm× 63mm×19 mm

Weight: 180g (including battery)

**Environment requirements**

- 1) Storage and Transportation Environment:
  - a) Temperature:- 10°C~+50°C
  - b) Relative Humidity: 10%~95%
  - c) Atmospheric pressure range: 50kPa~106kPa
- 2) Operating Environment:
  - a) Temperature: 5°C~40°C
  - b) Relative Humidity: 20%~80%
  - c) Atmospheric pressure range: 70kPa~106 kPa

**Battery type:** charging lithium battery

**Working-voltage:** 3.6 V DC ~ 4.2V DC

**Power:** < 0.3VA

The continuance of working hour: It can last more than 20 hours when the power is full.

LCD display: Resolution 160×96, display area 60×35mm

### **Main Parameters**

#### **1 Measurement of Heart Rate**

Measurement range: 30bpm~300bpm

Resolution: ±1 bpm

Accuracy: ±2 bpm

#### **2 Measurement of SpO<sub>2</sub>**

Measurement range: 0%~100%

Resolution: 1%

Accuracy: When the SpO<sub>2</sub> measurement range is 70%~100%, the permission of absolute error is ±2% below 70% unspecified

#### **3 Measurement of pulse rate**

Measurement range: 30bpm~250bpm

Resolution: 1 bpm

Accuracy: ±2 bpm or ±2% (select larger)

#### **4 Measurement Performance in Weak Filling Condition:**

SpO<sub>2</sub> and pulse rate can be shown correctly when pulse-filling ratio is 0.4%.

SpO<sub>2</sub> error is ± 4%, pulse rate error is ±2 bpm or ±2% (select larger).

#### **5 Resistance to surrounding light:**

The deviation between the value measured in the condition of man-made light and indoor natural light and that of darkroom is less than ±1%.

## **Chapter 6 CLEANING/MAINTENANCE AND DISINFECTING**

### **6.1 Cleaning**

Keep the outside surface of the device clean and free of dust and dirt, and clean exterior surface of the device(including LCD screen)with a dry, soft cloth.

If necessary, clean the chassis with a soft cloth soaked in a solution of soap or water, and wipe dry with a clean cloth immediately.

#### **Attention:**

-  Before cleaning, switch off the power.
-  Don't use strong solvent, such as acetone.
-  Never use an abrasive such as steel wool or metal polish.

-  The temperature of the water which is used for cleaning the device should be below 60°C.
-  Do not allow any liquid to enter the product, and do not immerse any parts of the device into any liquids.
-  Avoid pouring liquids on the device while cleaning.
-  Don't remain any cleaning solution on the surface of the device.

## 6.2 Disinfecting

After cleaning the device, wipe the surface of device with ethanol, self-air dry (or clean with a clean, dry cloth).

### Attention:

-  Never try to sterilize the equipment by high temperatures or high pressure steam sterilizing process.
-  Do not use E-beam or gamma radiation sterilization or other methods.

## 6.3 Maintenance

The user must check that the equipment does not have visible evidence of damage that may affect patient safety or device capability before using.

The recommended inspection interval is once per month or less. If damage is evident, replacement is recommended before using.

The device is frangible and must be handled with care.

These methods will prolong the life of the unit.

- 1 Please recharge the battery when the screen shows  .
- 2 Recharge the battery soon after the over-discharge. The device should be recharged every six months when it is no regularly used. It can extend the battery life following this guidance.
- 3 Please clean and disinfect the device after using to avoid infection.
- 4 The device needs to be calibrated once a year (or according to the calibrating program of hospital). It also can be performed at the state-appointed agent or just contact us for calibration.

The accuracy of the device is controlled by the equipment and can not be adjusted by user. If the result is distrustful, please use other methods to verify immediately or contact local distributor or manufacture to get help.

### Attention:

-  The device should be maintained by the eligible engineers who have been accredited by our company.
-  If the device is not used for a long time ,please put it in the casing after cleaning and disinfecting.

## Chapter 7

### KEY OF SYMBOLS

Signal	Description
SpO <sub>2</sub>	Pulse oxygen saturation(%)
PR	Pulse rate (bpm)
HR	Heart rate (bpm)
	Full-voltage
	Low-voltage
	Alarm-sound "Off"
	Alarm-sound "On"
	Close the pulse sound indication
	Open the pulse sound indication
SN	Serial number
-- and ---	Waiting for measuring data
PX1	Ingress of liquids rank
	Attention consult ACCOMPANYING DOCUMENTS
	WEEE (2002/96/EC)
	This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community

## Chapter 8

# TROUBLE SHOOTING

Trouble	Possible Reason	Solution
1. The SpO <sub>2</sub> and Pulse Rate can not be displayed normally.	1.The finger is not properly positioned. 2.The patient's SpO <sub>2</sub> value's going beyond the limits.	1. Place the finger properly and try again. 2.Try again. Go to a hospital for a diagnosis if you are sure the device works all right.
2. The SpO <sub>2</sub> and Pulse Rate are not displayed stably.	1. The finger is not placed inside deep enough. 2. The finger is shaking or the patient is moving.	1. Place the finger properly and try again. 2. Let the patient keep still.
3. ECG display astatically.	1. Electrodes do not connect with skin very well. 2.The device is shaking. 3. The patient is shaking or moving.	1. Recoat gel and press electrodes. 2.Keep the device on the body still. 3. keep the patient in a stable state.
4. The device can not be turned on.	1. The battery is drained away or almost drained away. 2. The malfunction of the device.	1. Recharge battery. 2. Please contact the local service center.
5. Noise.	Here is disturbance sources.	Remove the disturbance sources and other electrical instruments.
6. The device can not be used for full time after charge or the battery can not be full charged even after 10 hours charging time.	The battery is broken.	Please contact the local service center.

### Note:

-  If the trouble is still not got rid of, please contact warrant service agent as soon as possible.
-  For the damage aroused by client dismantling and mending, our company does not assume any responsibility.

**Disposal:** *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

## **GIMA WARRANTY CONDITIONS**

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production.

The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons.

Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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