
HEEL PROTECTOR

Ref. 680 011 Αναφ. 680 011
45*10*15 cm 150gr 45*10*15 εκ. 150 γρ.

Ref. 680 008 Αναφ. 680 008
22*15.12 cm 90gr 22 * 15,12 εκ. 90 γρ

Ref. 680 010 Αναφ. 680 010
22*15.12 cm 90gr 22 * 15,12 εκ. 90 γρ

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EN- Please read these safety instructions carefully before use, and keep them for future question.

1. Safety indications

- ♥ ALWAYS consult with your physician or supplier to determine proper adjustment and use of the device.
- ♥ Don't store heel support outside and don't let under the sun (Beware: don't let the product in a car in case of high outside temperature !!), store in a dry room in a flat position to prevent deforming the cushion.
- ♥ Store the heel support away from sharp objects and pointed edges to prevent tears and punctures of the cushion and cover.
- ♥ Use and store cushion in a flat unfolded position.
- ♥ Look for cracks before use.
- ♥ Don't use in extreme temperature conditions (above 38°C or 100°F or below 0°C or 32°F).
- ♥ No modification must be made to this device.
- ♥ ALWAYS check if the protector is properly and securely locked in place before using.
- ♥ In case you didn't use your cushion for a long time, let it be checked by the technical staff, before using.
- ♥ ALWAYS observe the weight limit on the labelling of your device.
- ♥ Use only as a protector and in a proper way: don't knee or step on the product.
- ♥ If your skin develops redness, discontinue the use of the cushion immediately and consult your doctor or therapist.
- ♥ Follow the advice of the professional who prescribed this medical device, respect duration and frequency of use.
- ♥ In case of discomfort consult your professional.
- ♥ Avoid contact with any sharp objects or pointed edges that may tear or puncture the cushion or the cover.
- ♥ Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

2. Intended use, indications and contraindications

It is easier to prevent pressure ulcers than to treat them. Bedsores can be avoided by changing frequently lying or seating positions to avoid stress on

the skin and to minimize the risk of pressure on vulnerable bony areas. Other strategies include skin care, regular skin inspections and good nutrition. Using a heel protector can help relieve pressure and make sure that the heel is correctly positioned in the chair or in the bed. It is one of the best ways to avoid bedsores, specifically reducing pressure, friction and shear from the heels, „floating“ the heels off the surface of the bed or chair.

Which patients are at risk of pressure ulcers?

Patients with limited or no mobility due to their age, an accident or illness, and who sit for long periods of time in the day. Patients who cannot reposition themselves independently when they get uncomfortable. These patients might slump to one side, forward or slide down the chair, needing help from carers or family members to regain a mid-line posture. Patients who have thin or fragile skin and who are at a high pressure injury risk. Patients who are hoisted for all transfers.

Indications:

- Patients with or at risk of developing minor pressure ulcers
- Pain relief, pressure redistribution, reduction of friction and shear forces
- Patients with significant bony prominences
- Post surgery and bedsores prevention.

Contraindications:

- Don't use this product in case of :
- Excessive sweating causing a risk of skin maceration
 - Severe skin defect (ulcers)
 - Inadequate supporting surface
 - Avoid layering incontinence products

3. Description and Fitting instructions

Ref 680010

The heel protector is a contoured moulded visco-elastic foam support. The foam increases the contact surface and maximizes pressure relief. The heel protector's shape ensures a good posture of the leg and completely offloads heels. This heel protector is delivered with a polyurethane cover.

Hook and loop closure.

Ref 680011

Coated cylindrical hollow core fibres provide soft comfortable support, while

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