AUDIOMETER

TRIANGLE

MULTILANGUAGE USER MANUAL

Document title: AU1S-Triangle User Manual PT-TR-NL-LT-NO
Code: AU1-MA011 _B
Revision: Rev. 02
Date: 2025.02.19

LYDMÅLER

TRIANGLE

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AVAILABLE TESTS
Manual pure tone audiometry, Automatic pure tone audiometry (Fixed intensity, Screening, Houghson Westlake)

AUDIOMETRY	
Stimulus	Pure tone, Warble
Masking	Narrow-band noise (NBN), White noise (WN)
Attenuator step	5 dB
Method of presentation	Continuous, Pulsed (cadence: 0.5 Hz, 1 Hz, 2 Hz)
Frequency accuracy	0.1%
Intensity accuracy	±3 dB between 125Hz and 4kHz; ±5 dB over 4kHz

Total harmonic distortion (THD)	AC: less than 2.5% BC: less than 5.5%	
Warble tone	Frequency of the modulating signal: 5 Hz Modulation waveform: sine wave Modulation range: ±12%	
NBN	Band: ½ octave, i.e.: - lower cut-off frequency $f_l = f/1.1892$ - upper cut-off frequency $f_u = f * 1.1892$ where f is the center frequency	
WN	Lower cut-off frequency: 100 Hz Upper cut-off frequency: 24kHz	
Compatible transducers		
Type	Manufacturer	Model
Supra-aural headphones	Radioear Corp.	DD45
Circum-aural headphones	Radioear Corp.	DD65
Insert earphones	Etymotic Research Inc.	ER-3C
Insert earphones	Radioear Corp.	IP30
Bone vibrator	Radioear Corp.	B71 ¹

PURE TONE AND WARBLE TONE AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES				
Freq [Hz]	AC DD45 [dB HL]	AC DD65 [dB HL]	AC ER-3C / IP30 [dB HL]	BC B71 [dB HL]
125	65	65	80	-
250	85	80	90	35
500	100	95	100	50
750	100	95	100	55
1.000	100	95	100	60
1.500	100	95	100	60
2.000	100	95	100	60
3.000	100	95	100	60
4.000	100	95	100	60
6.000	95	80	90	40
8.000	85	80	75	35

¹ Available if the optional “Bone conduction” license is activated

NARROW BAND NOISE (NBN) AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES			
Freq [Hz]	AC DD45 [dB EM]	AC DD65 [dB EM]	AC ER-3C / IP30 [dB EM]
125	40	40	60
250	65	60	75
500	80	75	85
750	85	75	90
1.000	85	80	90
1.500	85	80	90
2.000	85	80	90
3.000	85	80	90
4.000	85	80	90
6.000	85	70	85
8.000	85	70	80

WHITE NOISE (WN) MAXIMUM INTENSITIES		
AC DD45 [dB SPL]	AC DD65 [dB SPL]	AC ER-3C / IP30 [dB SPL]
100	85	95

REFERENCE EQUIVALENT THRESHOLD LEVEL FOR PURE TONE				
	DD45	DD65	ER-3C IP30	B71(*)
Ref. std.	ISO 389-1 (ANSI S3.6)	Vendor Tech. Specificat.	ISO 389-2 (ANSI S3.6)	ISO 389-3 (ANSI S3.6)
	Coupler: IEC 60318-3	Coupler: IEC 60318-1	Coupler: IEC 60318-5	Mastoid: IEC 60318-6
Freq [Hz]	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 1 µN]
125	47.5	30.5	26.0	-
250	27.0	17.0	14.0	67.0
500	13.0	8.0	5.5	58.0
750	6.5	5.5	2.0	48.5
1000	6.0	4.5	0.0	42.5
1500	8.0	2.5	2.0	36.5

2000	8.0	2.5	3.0	31.0
3000	8.0	2.0	3.5	30.0
4000	9.0	9.5	5.5	35.5
6000	20.5	21.0	2.0	40.0
8000	12.0	21.	0.0	40.0
(*) Calibration of bone vibrator (B71) refers to mastoid placement.				

SOUND ATTENUATION VALUES			
Freq. [Hz]	DD45(*) [dB]	DD65 [dB]	ER-3C [dB]
125	3.0	8.3	33.5
250	5.0	15.5	34.5
500	7.0	26.1	34.5
750	-	-	-
1000	15.0	32.4	35.0
1500	-	-	-
2000	26.0	43.6	33.0
3000	-	-	-
4000	32.0	43.8	39.5
6000	-	-	-
8000	24.0	45.6	43.5

(*) With MX41\AR or PN 51 cushions

PATIENT – OPERATOR COMMUNICATION
Talk-over through embedded microphone
Patient response trigger

PATIENT MANAGEMENT	
Max number of patients	100
Details stored	Patient details (first name, last name, date of birth, gender), date and time of test, left and right audiogram (AC and BC)

PHYSICAL SPECIFICATIONS	
Dimensions	(WxDxH) 160 × 217 × 48 mm / 6.3 × 8.5 × 1.9"
Weight	420 g / 14.8 oz
Display	LCD TFT 2.8" RGB, 240 × 320 pixels Viewing area 43.2 mm × 57.6 mm

Touchscreen	Capacitive
-------------	------------

SOCKETS ON THE REAR PANEL			
<i>Description</i>	<i>Connector</i>	<i>Pins</i>	<i>Specifications</i>
Power supply / Computer communication	USB type B		1 5 V DC 2 Data – 3 Data + 4 GND
L and R headphones	Audio jack, 1/4" mono		1 GND 2 Signal (8 V on 10Ω load)
Bone vibrator			
Patient response trigger			1 GND 2 Input Switch

INTERFACE WITH COMPUTER ²	
Connection	USB (no driver needed)
Minimum PC requirements	USB 2.0 port
Compatible software products	Inventis Maestro suite

On request Inventis will make available circuit diagrams, component part lists, descriptions, calibration instructions or other information that will help service personnel to repair those parts of the device that are designated by Inventis as repairable by service personnel.

² Available if the optional “Computer connectivity” license is activated

Electromagnetic compatibility

The instrument has been tested and respects the limits for electro-medical devices specified by IEC 60601-1-2 standards. These limits ensure reasonable protection against hazardous interference in typical medical installations.

The device generates, uses and radiates radio frequency energy. If not installed and used according to the instructions in this manual, it may interfere with other nearby devices. No guarantee is given that interference will not occur under certain conditions.

This device is suitable for use in professional healthcare facility environments, i.e. in hospitals, except for near active HF surgical equipment and RF-shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

The instrument should not be used adjacent to or stacked on other equipment. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

The existence of electromagnetic interference can be verified easily by turning the device off and back on again. If it is found that the device is indeed interfering with other equipment, try to solve the problem by adopting one of the following solutions:

- Change the orientation and/or position of the affected device.
- Move the two devices further away from each other.
- Contact the manufacturer or authorised service center for further assistance.

List of cables, transducers and accessories

Cables, transducers and accessories for which Inventis claims compliance with the IEC 60601-1-2 standard are those supplied with the device, in particular:

- Medical grade USB power adapter
- USB cable: 2m length, screened.
- AC transducers: 2m length, screened.
- BC transducers: 2m length, screened.
- Patient button: 2m length, screened.

The use of accessories, transducers and cables other than those specified, except for transducers and cables sold by the manufacturer as spare parts for internal components, may result in increased emissions or decreased electromagnetic immunity of the device and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

The instrument does not have any ESSENTIAL PERFORMANCE as per the IEC 60601-1 standard.

Note: All instructions necessary for maintaining compliance with regard to electromagnetic compatibility can be found in the maintenance section in this manual. No further steps are required.

Guidance and manufacturer's declaration – electromagnetic emissions		
Triangle is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment Guidance
RF emissions CISPR11	Group 1	Triangle uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	Triangle is suitable for use in professional healthcare facility environments and directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flickers emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
Triangle is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such an environment.			
Immunity Test	IEC 60601 test level	Compl. Level	Electromagnetic environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ⁽¹⁾ ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air ⁽¹⁾	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a professional healthcare facility environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a professional healthcare facility environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. 40% U_T (60% dip in U_T) for 5 cycles. 70% U_T (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. 40% U_T (60% dip in U_T) for 5 cycles. 70% U_T (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	Mains power quality should be that of a professional healthcare facility environment. If the user of Triangle requires continued operation during a power mains interruption, it is recommended that Triangle be powered with an uninterruptable power supply or its battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic field should be at levels characteristic of a professional healthcare facility environment.
Note: ⁽¹⁾ A lock or reboot of the device with no permanent damage is acceptable ⁽²⁾ U_T is the a.c. main voltage prior to application of the test level.			

Manufacturer's indications and declaration - electromagnetic immunity			
Triangle is intended for use in an electromagnetic environment as specified below. The customer or the user of Triangle should assure that it is used in such an environment.			
Immunity tests	Test level to IEC 60601	Level of Compliance	Electromagnetic environment – guidelines
Conducted RF IEC 61000-4-6	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	3 Vrms 0.15 – 80 MHz 6 Vrms in ISM bands between 0.15 - 80 MHz	Portable and mobile RF communications equipment should be used no closer than 30 cm (12 inches) to any part of Triangle, including cables specified by the manufacturer.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz	3 V/m 80 MHz - 2.7 GHz	
Radiated RF (from RF Wireless communication equipment) IEC 61000-4-3	9 V/m	9 V/m	
	704-787 MHz 5100 – 5800 MHz	704-787 MHz 5100 – 5800 MHz	
	27 V/m	27 V/m	
	380 - 390 Mhz	380 - 390 Mhz	
Proximity Fields IEC 61000-4-39	28 V/m	28 V/m	
	430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	430 - 470 Mhz 800 - 960 Mhz 1700 - 1990 Mhz 2400 - 2570 Mhz	
	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz	8A/m 30KHz CW 65 A/m 134.2KHz Pulse Mod 2.1KHz	
	7.5 A/m 13.56MHz Pulse Mod 50KHz	7.5 A/m 13.56MHz Pulse Mod 50KHz	Interference may occur in the vicinity of equipment marked with the following symbol: 
	<i>Note:</i> At 80 MHz and at 800 MHz, the higher frequency range is applied.		
<i>Note:</i> These guidelines may not be valid in all situations. Electromagnetic propagation is affected by absorption and reflection on encountering structures, objects and persons.			
<i>Note:</i> (1) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If measured field strength in the location Triangle is used in exceeds the applicable RF compliance level above, Triangle should be observed to verify its normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Triangle. (2) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			