Instruction Manual

Manuale d'istruzione

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Skin Stapler Remover

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WARNING
Please read the *instruction manual* carefully before using it!

1. Main Structure

The metal part of the remover is made of 1Cr18Ni9Ti, and the handle of the remover is made of ABS plastic.

2. Working Principle

After the handle was pushed by outside force, the upper cutter of the stapler remover force on the middle part of the stapler which is on the nether cutter of the remover. Under the mechanical transmission of the each part, two ends of the stapler are to be turned up and out of the skin surface and that achieve the aim of removing the staplers.

3. Capability Index

- The metal parts of the skin stapler removers have the good function of anti-eroded.
- The skin stapler remover is sterilized with ethylene oxide. The leftover of the ethylene oxide should be less than 10 mg/kg before the goods are released from factory.

4. Applicative area

Skin stapler remover is used to remove the staplers.

5. How-to-use

a. Insert the front part of the nether cutter of the remover to the underside of the stapler

b. Clamp the upper and nether handle tightly until they are fully closed. Stapler turns to be out of shape.

c. Take out of the completely unshaped stapler vertically from the skin surface operated.

6. Important Notice:

- The remover must be operated by the professional person who has the medical certificates.
- The remover could not be polluted when it was taken out of the package. Forbidden to use if the package broken.
- Please use the skin stapler remover within the period of validity.

7. Reserve and Transportation

The product should be kept in rooms with relative humidity $\leq 80\%RH$, with no erosive gases and good ventilation. Not to be stowed under other cargo, and handle with care.

8. Period of Validity


5 years.

9. Service after Sold

If you have any complains for this product, please contact with suppliers or manufacturers.

10. GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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Ningbo Advan Electrical Co., Ltd.
Corp. GmbH (Europe)
Gima S.p.A. - Via Marconi, 1