
HEINE Dermatoscopes

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49 (0) 81 52 / 38 - 0
Fax +49 (0) 81 52 / 38 - 202
E-Mail: info@heine.com · www.heine.com
med 0415 2016-05-23

HEINE Dermatoscopes

| | |
|------------|----|
| DEUTSCH | 3 |
| ENGLISH | 5 |
| FRANÇAIS | 7 |
| ESPAÑOL | 9 |
| ITALIANO | 11 |
| SVENSKA | 13 |
| NEDERLANDS | 15 |
| DANSK | 17 |
| NORSK | 19 |
| SUOMI | 21 |
| PORTUGUÊS | 23 |

HEINE Dermatoskopes

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Guidance and manufacturer's declaration – electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such environment.

| Emission test | Compliance | Electromagnetic environment – Guidelines |
|---|------------|---|
| RF emissions CISPR11 | Group 1 | The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic Emissions IEC 61000-3-2 | Class A | Symmetrical three-phase devices and other devices. |
| Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3 | Passed | |

Guidance and manufacturer declaration - Electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – Guideline |
|--|---|---|---|
| Electrostatic discharge (ESD) acc. to IEC 61000-4-2 | ± 6 kV contact discharge ± 8 kV air discharge | ± 6 kV contact discharge ± 8 kV air discharge | Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for mains cables ± 1 kV for input and output lines | ± 2 kV for mains cables ± 1 kV for input and output lines | The supply voltage quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth | ± 1 kV voltage phase – phase ± 2 kV voltage phase – earth | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | < 5% U _T , (>95% dip in U _T) for 1/2 period 40% U _T , (60% dip in U _T) for 5 periods 70% U _T , (30% dip in U _T) for 25 periods <5% U _T , (>95% dip in U _T) for 5 seconds | < 5% U _T , (>95% dip in U _T) for 1/2 period 40% U _T , (60% dip in U _T) for 5 periods 70% U _T , (30% dip in U _T) for 25 periods <5% U _T , (>95% dip in U _T) for 5 seconds | Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

Comment: U_T is the a.c. supply voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment Guidelines |
|-------------------------------|---|--------------------|---|
| Conducted RF IEC 61000-4-6 | 3 V _{eff} 150 kHz to 80 MHz | 3 V _{eff} | Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. |
| Radiated HF IEC 61000-4-3 | 3 V/m 80MHz to 2,5GHz | 3 V/m | <p>Recommended separation distance: $d = 3,5/3 \cdot \sqrt{P/W}$ $d = 3,5/3 \cdot \sqrt{P/W}$ 80 MHz to 800 MHz $d = 7/3 \cdot \sqrt{P/W}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

Note 1: At 80Hz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.
To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation.
If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

Recommended separation distances for portable and mobile RF communication equipment and the device

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|---|---|---|--|
| | 150 kHz to 80 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$ | 80 MHz to 800 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$ | 800 MHz to 2,5 GHz $d = 7/3 \cdot \text{SQRT}(P)$ |
| 0.01 | 0.1 | 0.1 | 0.2 |
| 0.1 | 0.4 | 0.4 | 0.7 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.7 | 3.7 | 7.4 |
| 100 | 11.7 | 11.7 | 23.3 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Technical specification

| | |
|---|---|
| Environmental conditions for operation | +10 °C to +35 °C 30 % to 75 % rel. humidity 700 hPa to 1060 hPa |
| Environmental conditions for storage | +5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa |
| Environmental conditions for transport | -20 °C to +50 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa |
| Nominal voltage | 3.0 V – 3.7 V |
| Nominal current | 440 – 760 mA |
| Protection class | internal power supply |
| IP-Code | IP20 |
| Device classification according to IEC 6247 | Group 2 |
| Applied part | Type BF (for contact plate) |
| HEINE mini3000® Dermatoscope | #109 (2,5 V) |

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Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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