Tegaderm™ Alginate Ag

CARATTERISTICHE

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BIBLIOGRAFIA

CHRONIC WOUNDS

| | | | | | |
|---|---------------------------|---|---|-------------------------------|--|
| Tegaderm™ Ag Mesh Dressing with Silver — vs — Silver sulfadiazine cream | Randomized Clinical Trial | The Efficacy of Silver Mesh Dressing Compared with Silver Sulfadiazine Cream for the Treatment of Pressure Ulcers | Wound reduction was greater for patients managed with Tegaderm™ Ag Mesh Dressing with Silver. The cost of treatment using Tegaderm™ Ag Mesh Dressing with Silver was significantly less than using silver sulfadiazine cream. | Aprig, C. et al. | J Med Assoc Thal. 2011;94(5): 559-565 |
| Tegaderm™ Ag Mesh Dressing with Silver | Case Series | Evaluation of 3M™ Tegaderm™ Ag Mesh Dressing with Silver in the Management of Wounds Where Sub-infectious Microbial Colonization is Suspected | Case studies of slow to non-healing wounds (a venous stasis ulcer, a dehisced surgical wound, and a Stage IV ischial pressure ulcer) were managed using Tegaderm™ Ag Mesh Dressing with Silver. All wounds were suspected of critical colonization. Management of these wounds with Tegaderm™ Ag Mesh Dressing with Silver over 3-4 weeks led to reduced signs of critical colonization with resumption of normal wound healing. | Brown-Etris, M. Punchello, M. | White Paper based on poster presented at 2007 WOCN 70-2009-7189-6 |
| Tegaderm™ Ag Mesh Dressing with Silver | Case Series | Use of 3M™ Tegaderm™ Ag Mesh Dressing with Silver in Patients with Recurrent Lower Extremity Ulcers | Tegaderm™ Ag Mesh Dressing with Silver was used on two patients, each with venous leg ulcers. Tegaderm™ Ag Mesh Dressing with Silver demonstrated progressive healing in a 6-8 week period versus wound deterioration with treatment of hydrogel, contact layer and zinc paste bandage. | Peltier, G. et al. | White Paper (2007) 70-2009-7525-1 |
| Tegaderm™ Ag Mesh Dressing with Silver | Case Series | Management of Recalcitrant Ulcers with 3M™ Tegaderm™ Ag Mesh Dressing with Silver | Two patients, one with a venous stasis ulcer and another with a recurrent ulcer on an index finger, had failed to improve with multiple treatment modalities, including silver therapies. Rapid closure occurred using Tegaderm™ Ag Mesh Dressing with Silver on these recalcitrant wounds. | Agbim, S. | White Paper (2006) 70-2009-7535-0 |

ANTIMICROBIAL EFFICACY

| | | | | | |
|--|-------------------------|--|---|--|--|
| Tegaderm™ Ag Mesh Dressing with Silver | <i>in vitro</i> Studies | The Power of Silver | High level summary of efficacy of Tegaderm™ Ag Mesh Dressing with Silver. Includes the efficacy of Tegaderm™ Ag Mesh Dressing with Silver against several pathogens including MRSA, MRSE, VRE, MDR, Pseudomonas aeruginosa, and Candida albicans. Includes information on how ionic silver works. | 3M | Product Monograph (2006) 70-2009-7137-5 |
| Tegaderm™ Ag Mesh Dressing with Silver — vs — ConvaTec AQUACEL® Ag Hydrofiber® Dressing with Ionic Silver | <i>in vitro</i> Studies | <i>In Vitro</i> Investigation of the Long-term Time-Kill Kinetics of 3M™ Tegaderm™ Ag Mesh Dressing with Silver Compared to Aquacel® Ag Hydrofiber® Dressing with Ionic Silver | The sustained activity of Tegaderm™ Ag Mesh Dressing with Silver against wound microorganisms demonstrate 14-day efficacy against commonly encountered wound pathogens including MRSA, VRE, anaerobic bacteria and yeast, and fungi. Tegaderm™ Ag Mesh Dressing with Silver compared favorably to ConvaTec AQUACEL® Ag Hydrofiber® Dressing with Ionic Silver. | Thomhill, G. Stahl, J. Opp, C. | White Paper (2006) 70-2009-7151-6 |
| Tegaderm™ Ag Mesh Dressing with Silver — vs — Smith & Nephew ACTICOAT® Dressing — vs — ConvaTec AQUACEL® Ag Hydrofiber® Dressing with Ionic Silver | <i>in vitro</i> Studies | <i>In Vitro</i> Assessment of the Microbicidal Activity of 3M™ Tegaderm™ Ag Mesh Dressing with Silver Using a Short-term Time-Kill Procedure | Tegaderm™ Ag Mesh Dressing with Silver performed better than Smith & Nephew ACTICOAT® Dressing when challenged with: Candida albicans, Enterococcus faecalis (MDR), Enterococcus faecium (VRE), Staphylococcus aureus (MRSA), and Staphylococcus epidermidis (MRSE). Tegaderm™ Ag Mesh Dressing with Silver was found to be significantly faster than ConvaTec AQUACEL® Ag Dressing in reducing the same microbial load. | Baxter, K. Stahl, J. Morse, D. Opp, C. | White Paper (2006) 70-2009-7152-4 |
| Tegaderm™ Ag Mesh Dressing with Silver | <i>in vitro</i> Studies | <i>In Vitro</i> Assessment of the Broad-Spectrum Activity of Ionic Silver in 3M™ Tegaderm™ Ag Mesh Dressing with Silver | Zone of inhibition studies demonstrating Tegaderm™ Ag Mesh Dressing with Silver inhibits the growth of pure cultures of multiple microorganisms. | Thomhill, G. Stahl, J. Opp, C. | White Paper (2006) 70-2009-7153-2 |