

Test Marker Cardiaco Combo a dischetto per sangue Intero, plasma o siero
Step Cardiac Combo Marker test disk for whole blood, plasma or serum
Test Marqueur Cardiaque Combo sur cassette sur sang total, plasma ou sérum
Marker Test Cardiaco Combo de disco para sangre entera, plasma o suero

MANUALE D'USO
OPERATOR'S MANUAL
MANUEL D'UTILIZATION
MANUAL DE USO

ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.

ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.

ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.



A creatine phosphokinase/creatinase or isoenzymes test system is a device intended to measure the activity of the enzyme creatine phosphokinase or its isoenzymes (a group of enzymes with similar biological activity) in plasma, serum or whole blood. Measurements of creatine phosphokinase and its isoenzymes are used in the diagnosis and treatment of myocardial infarction and muscle diseases such as progressive, Duchenne-type muscular dystrophy.

The One Step Cardiac Combo Marker Whole Blood, Plasma or Serum Test is a rapid assay test. The analytical sensitivity of Troponin I, CK-MB and Myoglobin is 1 ng/ml, 5 ng/ml and 50 ng/ml respectively. For professional use only.

PRECAUTIONS

The One Step Cardiac Combo Marker test kit should be stored at room temperature. The test device is sensitive to humidity as well as to heat. Perform the test immediately after removing the test device from the foil pouch. Do not use it beyond the expiration date.

SPECIMEN COLLECTION AND STORAGE

Whole Blood specimen collection: Collect an anticoagulated blood sample by using heparin as the anti-coagulant. Troponin I, CK-MB or Myoglobin is very unstable in serum or whole blood specimen. Whole blood or serum specimen must be tested as soon as possible.

Plasma/Serum specimen collection

1. Centrifuge whole blood to get plasma/serum specimen.
2. If specimens are not immediately tested, they should be refrigerated at 2-8 degrees C. Specimens should be at room temperature before running a test.
3. Specimens containing precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

Warnings

1. For in vitro diagnostic use only.
2. Do not eat or smoke while handling specimens.
3. Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
4. Avoid splashing or aerosol formation.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
7. Do not use the test kit if the pouch is damaged or the seal is broken.

PROCEDURE

1. Remove the test disk from the foil pouch, and place it on a flat, dry surface.
2. Holding the sample dropper above the test disk (Figure 1) and add 1 hanging drop into the Sample Well. After the drop is absorbed into the Sample Well, add another hanging drop, repeat the procedure until a total of 3 hanging drops have been added to the Sample Well. If specimen drops are added too quickly, specially for blood specimen, it may cause clogging of the Sample Well.
3. As the test begins to work, you will see purple color move across the Result Window in the center of the test disk. Note: If purple color dye does not begin to flow through the "Result Window" within 2 minutes, add one more drop of sample.
4. Interpret test results at 10 to 15 minutes. Do not interpret test results after 20 minutes.

Caution: The above interpretation time is based on reading the test results at room temperature of 15 to 30 °C. If your room temperature is significantly lower than 15 °C, then the interpretation time should be properly increased.

INTERPRETATION OF THE TEST

1. A purple band in the “C” area is the Control Line.
2. There are three Test Lines. cTnI as Troponin I, CKMB as CK-MB and MYO as Myoglobin.

Positive Result: The presence of purple color “C” control line together with any purple color of “T1”, “CKMB” or “MYO” test line, regardless of which band appears first, indicates a positive result (Figure 2). Note: Generally, the higher the analyte level in the specimen, the stronger the Test Line band color will be. When the specimen analyte level is close to but still within the sensitivity limit of the test, the color of the Test Line band will be very faint.

Note: Specimens containing very low levels of Myoglobin may develop color “MYO” band over 20 minutes.


Negative Result: The presence of only one purple color “C” control line indicates a negative result (Figure 3).

Invalid Result: If after performing the test no purple color control line is visible within the Result Window, this result is considered invalid (Figure 4). Not following the procedures correctly or using a test kit that has deteriorated can cause invalid results. It is recommended that the specimen be re-tested.

Note: A positive result will not change once you have established your answer at 20 minutes. However, in order to prevent any incorrect results, the test result should not be interpreted after 20 minutes. Some specimens with a high rheumatoid factor concentration may yield a nonspecific positive result.

LIMITATIONS OF THE TEST




Although the One Step Cardiac Combo Marker Test is accurate in detecting cTnI,CK-MB and MB, a low incidence of false results can occur. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.






REFERENCE


- 1, Andersen PT, Miller-Petersen J, Henneberg EW, et al, "IpermiogliHypermyoglobinemia After Successful Arterial Embolectomy," *Surgery*, 1987, 102:25-31.
- 2, Cummins B, Russell GJ, Chandler ST, et al, "Uptake of Radioiodinated Cardiac Specific Troponin I Antibodies in Myocardial Infarction," *Cardiovasc Res*, 1990, 24(4):317-27.
- 3, Delanghe J, Chapelle JP, el Allaf M, et al, "Quantitative Turbidimetric Assay for Determining Myoglobin Evaluated," *Ann Clin Biochem*, 1991, 28(Pt 5):474-9.
- 4, Donnelly R and Hillis WS, "Cardiac Troponin T," *Lancet*, 1993, 341(8842): 410-1.
- 5, Gibler WB, Gibler CD, Weinschenker E, et al, "Myoglobin as an Early Indicator of Acute Myocardial Infarction," *Ann Emerg Med*, 1987, 16:851-6.
- 6, Mair J, Smidt J, Artner-Dworzak E, et al, "Rapid Diagnosis of Myocardial Infarction by Immunoturbidimetric Myoglobin Measurement," *Lancet*, 1991, 337(8753):1343-4.
- 7, Parmacek MS and Leiden JM, "Structure, Function, and Regulation of Troponin C," *Circulation*, 1991, 84(3):991-1003.
- 8, "Troponin T and Myocardial Damage," *Lancet*, 1991, 338(8758):23-4, (editorial).
- 9, Serrano S, Chueca P, Carrasco E, et al, "Predictive Value of Myoglobin in Early Diagnosis of Acute Mtocardial Infarction," *Ann Emerg Med*, 1990, 19(8):953.
- 10, Silva DP Jr, Landt Y, Porter SE, et al, "Development and Application of Monoclonal Antibodies to Human Cardiac Myoglobin in A Rapid Fluorescence Immunoassay," *Clin Chem*, 1991, 37(8):1356-64.
- 11, Vrenna L, Castaldo AM, Castaldo P, et al, "Comparison Between Nephelometric and RIA Methods for Serum Myoglobin, and Efficiency of Myoglobin Assay for Early Diagnosis of Myocardial Infarction," *Clin Chem*, 1992, 38(5):789-90.

Simbologia / Index of symbols - TEST

	<p>Leggere e seguire attentamente le istruzioni per l'uso <i>Please read and follow instructions carefully</i></p>
	<p>Solo per uso diagnostico <i>in vitro</i> <i>For in vitro diagnostic use only</i></p>
	<p>Conservare tra 4° e 30° C <i>Store between 4° and 30° C</i></p>

	<p>Dispositivo monouso, non riutilizzare <i>Disposable device, do not re-use</i></p>
	<p>Data di scadenza (vedi scatola / bustina) <i>Expiration date (see box / package)</i></p>
	<p>Numero di lotto (vedi scatola / bustina) <i>Lot number (see box / package)</i></p>

	<p>Prodotto conforme alla Direttiva Europea n. 98/79/CE sui dispositivi diagnostici in vitro <i>Product complies with European Directive no. 98/79/EC on In Vitro diagnostic devices</i></p>
	<p>24525 Codice prodotto <i>Product code</i></p>

FABBRICANTE / MANUFACTURER:

GIMA S.p.A.
 Via Marconi, 1
 20060 Gessate (MI), ITALY