

COVID-19 ORAL FLUID ANTIGEN RAPID TESTS

- 24588 SALIVA COVID TEST - self test - box of 1
- 24590 SALIVA COVID TEST - professional - box of 20

The COVID-19 Antigen Rapid Test (Oral Fluid) is a qualitative rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2 Nucleocapsid (N) protein antigens in oral fluid specimens from individuals with suspected infection. It's a single-use test kit intended to detect the coronavirus SARS-CoV-2 through human oral fluid. This test is designed for self-testing (24588) or professional use (24590). Research shows that the SARS-CoV-2 virus is stable in oral fluid for long periods of time at warm temperatures.

- fast results: 15 minutes
- high accuracy
- convenient procedure
- easy interpretation
- simple operation
- relative specificity: 99,3%

Storage temperature: 2-30°C
Shelf life*: 12 months (24588) 24 months (24590).

*Gima stock may vary from 4 to 10 months (24588) 6 to 22 months (24590)



24588

STANDARD ACCESSORIES

- | | |
|--|--|
| Test devices | Biosafety bag |
| Collection devices (funnel, tube and tube tip) | Procedure card (24590 only). |
| Buffer | Multilingual box: GB, FR, IT, ES, DE, PT |
| Package insert | |

<p>Step 1: Specimen collection Remove the funnel and plastic tube; fit the funnel onto the tube. Deeply cough 3-5 times. Gently spit oral fluid into the funnel. The oral fluid (non-bubble) should just reach the height of scale line.</p>	<p>Step 2: Specimen preparation Add entire buffer to the tube with oral fluid. Fit the tube tip onto the tube. Gently squeeze the tube 10-15 times to mix well.</p>	<p>Step 3: Testing Place the test cassette on a flat and level surface. Invert the tube and add 2 drops of solution to the specimen well(S) of the test device and then start the timer.</p>	<p>Step 4: Read the results in 15 minutes. Do not interpret the result after 20 minutes. After test is completed, dispose all the components</p>
---	--	---	---

SARS-CoV-2 ANTIGEN RAPID TEST KIT - PROFESSIONAL

• 24585 SARS-CoV-2 ANTIGEN TEST KIT - box of 25
Rapid test kit (Colloidal gold immunochromatography) intended for the qualitative detection of SARS-CoV-2 antigen and variants in clinical samples (nasal swab). It includes detection of major variants, recognizing the N protein on the surface of the virus.

- non invasive, simple to use
- convenient, no device required
- rapid, gets results in 15 min.
- stable, with high accuracy
- inexpensive, cost-efficiency
- storage condition: 4~30°C
- 18 months shelf life since manufacturing date

(average shelf life of Gima stock: 4 to 16 months)
Kit includes 25 test cards, 25 nasal swabs, both individually packed and 25 single-dose dropper of diluent solution (3 ml).
Multilingual instructions: GB, IT, ES, PT, DE, RU.



- Hospital
- Test site
- Mass Screening
- Station
- Airport
- Hotel
- Corporation



<p>Operation Steps Step 1: use swab to collect sample</p>	<p>Step 2: peel off adhesive</p>	<p>Step 3: insert swab from well B to well A</p>	<p>Step 4: A) add 6 drops B) rotate the shaft, two rounds each direction.</p>	<p>Step 5: fold left side over, and wait for 15 minutes</p>
--	----------------------------------	--	---	---

SARS-CoV-2 ANTIGEN RAPID TEST KIT - SELF TEST

• 24583 SARS-CoV-2 ANTIGEN TEST KIT - box of 1
Rapid, lateral flow immunoassay test intended for self qualitative detection of SARS-CoV-2 nucleocapsid antigens with nasal swab. Self-collection must be always carried out by an adult. This test is intended for use in individuals with symptoms or other epidemiological reasons to suspect a COVID-19 infection. The diagnostic sensitivity and diagnostic specificity of the product is 95.9% (90.8-98.2%) and 100% (96.3-100.0%) respectively. Kit includes: 1 test card, 1 nasal swab, 1 single-dose dropper; operation card and instructions GB, IT.

- non-invasive: 2-3 cm insertion
- rapid, gets result in 15 mins.
- stable with high accuracy
- inexpensive

Testing results

	Positive (+)		Negative (-)	Invalid	
C: control line	C	C	C	C	C
T: test line	T	T	T	T	T
	Solid line	Faint line			



24583