

HISTOFREEZER® CRYOSURGICAL SYSTEM

- **30591 HISTOFREEZER® - kit with applicator 2 mm**
Two 80 ml aerosol bottles and 60 applicators 2 mm.
- **30592 HISTOFREEZER® - kit with applicator 5 mm**
Two 80 ml aerosol bottles and 52 applicators 5 mm.
- **30589 HISTOFREEZER® - large kit mixed applicators**
Two 80 ml aerosol bottles, 24 applicators 2 mm, 36 applicators 5 mm.
- **30585 HISTOFREEZER® - mini kit mixed applicators**
One 80 ml aerosol bottle, 16 applicators 2 mm, 16 applicators 5 mm.

The Histofreezer® portable cryosurgical system is a low cost, highly effective cryosurgical agent that provides safe and effective treatment for warts and other common benign skin lesions. This therapeutic treatment makes use of local freezing for the controlled destruction or removal of living tissue. Suitable for adults and children.

Multilingual box and manual in: GB, FR, IT, DE, ES, NL, PL. From web you can download: GR, PT, TH, Arabic. For professional use only.

9 types of benign lesions	% of resolution*
Verruca Vulgaris (Common Warts)	94%
Verruca Plantaris (Plantar Warts)	86%
Molluscum Contagiosum	100%
Lentigo (Age Spots)	100% - facial 99% - non facial
Condyloma (Genital Warts)	89%
Seborrheic Keratosis	95%
Acrochordon (Skin Tags)	100%
Verruca Plana (Flat Warts)	93%
Actinic keratosis (Facial and Non-Facial)	91% - facial 87% - non facial

*Summary of resolution rates by type of lesion from a number of clinical trials obtained with histofreezer + post-treatment of 27% salicylic acid / 17% lactic acid

Saturated 2 mm applicators will stay cold for 90 seconds.



30591

Saturated 5 mm applicators will stay cold for 120 seconds.



30592



30589



30585



2 mm
5 mm

2 mm disposable applicator for small lesions
5 mm disposable applicator for large lesions

CHARGE
custom applicator press valve until droplets start to form.

INVERT
canister for 15 seconds for optimal freezing temperature.

APPLY
a frozen applicator to lesion for 15-40 seconds.

CRYOMEGA® AND CRYOCLEAR® MULTI-USE CRYOSURGICAL DEVICES

- **30586 CRYOMEGA® with 16 g cartridge - non-refillable**
CryoOmega® dispenses liquid nitrous oxide (-89°C) the coldest portable cryogen and it is cleared for 16 common lesions including verrucae, skin tags, and seborrheic keratoses. Ideal choice in portable cryosurgery for physician offices looking for convenience and great first-time outcomes. It provides clinically effective freezing with pinpoint precision. No contact device for use on different patients at 1 cm distance. Up to 300 seconds treatment time per device (approx 30 treatments of 10 seconds). For medical healthcare professionals.

MADE IN THE USA



Type of Lesion	Freeze time range (sec)*
Verruca Vulgaris	10-20
Verruca Plantaris	10-20
Molluscum Contagiosum	3-10
Skin Tags	5-10
Lentigo	2-5
Actinic Keratosis	5
Verruca Plana	5-15
Seborrheic Keratosis	5-10
Genital Lesions	5-12

*Suggested freezing time as reported in literature

N₂O
-89°C


30586



Multilingual pack and instructions (GB, FR, IT, ES, DE, NL, PL, FI).

- **30587 CRYOCLEAR® with 16 g cartridge - non-refillable**
CryoClear® uses Carbon Dioxide which freezes at -79°C to a depth of no more than 2 mm into the skin, delivering a safe and clinically effective treatment to superficial skin indications. Over 95% high efficacy combined with less destruction of healthy tissue makes it perfect for treating age spots, sun spots, and skin tags. It offers a better and faster cryotherapy solution, even as combination treatment to chemical peels, facials, micro-needling or dermaplaining. No contact device for use on different patients at 0.5 cm distance. Up to 300 seconds treatment time per device (approx 60 treatments of 5 seconds). Suggested for use at aesthetics and cosmetic sites.

MADE IN THE USA



Indication	Freeze time range (sec)*
Actinic Keratosis (sun spots)	5-10
Lentigo (Age Spots)	5-10
Skin Tags	5-10

*Suggested freezing time as reported in literature

CO₂
-79°C

30587



Multilingual pack and instructions (GB, FR, IT, ES, DE, NL, PL, FI).