

GIMA 23916
GIMA 23918
GIMA 23921

CE 0197



Shandong Lianfa Medical Plastic Products Co., Ltd
No.1 Shuangshan Sanjian Road, Zhangqiu, Jinan City,
250200, Shandong P. R. China
Made in China




01-1328-II, 01-1330-II, 01-1332-II



Linkfar Healthcare GmbH
Niederrheinstraße 71, 40474 Düsseldorf, Germany



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



M23916-M-Rev.2.05.26

ENGLISH

BLOOD LANCETS

Model:

I, II, III

Specification:

14G-33G

Intended Use:

The Blood lancets is intended for capillary blood sampling in order to obtain a small blood sample for various tests. It is designed for use by both healthcare professionals and lay users, providing a safe and convenient method for blood specimen collection in clinical and home healthcare settings.

Contraindication:

Use with caution in individuals with coagulation mechanism disorders.


Information about the matching lancet device:

LINKFAR Lancing Devices/ Lianfa Lancing Devices/ OneTouch® UltraSoft®/ One-Touch® UltraVue/ FreeStyle Lite®/ Bayer Microlet®/ Bayer Microlet® Vaculance/ Glucolet Lancing Device/ BD Lancet Device/ Easy Touch® lancing device/ Mylife Unio (Ypsomed)/ Mylife Pura (Ypsomed)/ Beurer Lancing Device Autolet Lancing Device/ Glucolet Lancing Device/ SD Biosensor Lancing Device/ TaiDoc Lancing Device/ CareSoft Lancing Device/ Tempo Lancing Device/ Autolet/ MICROLET® NEXT/ Medifun/ SteriLance/ AIMSCO Lancing device

Product main structure:

The Blood Lancets primarily consists of lancet body, needle and cap.

The specific structure of the product is as follows:



1. Cap - 2. Needle - 3. Lancet body


Instructions for use

1. Wash hands with soap and warm water. Dry thoroughly.
2. Insert a new blood lancet into the lancet holder of lancet device.
3. Remove lancet cover (steps 1 & 2 above).
4. Consult the instructions for use of your lancet device.
5. Select a fingertip site, and proceed to puncture the skin.
6. After use, discard the used lancet in appropriate container.
7. Wash hands with soap and warm water again.

Wear gloves if use in healthcare facilities.

Precautions and cautions:

- Please read the instructions before use.
- Check the package before use.
- Do not use it if the package is damaged.
- Do not reuse lancets.
- Do not store in Lancing Devices.
- For use only on a single patient.
- It is forbidden to use after the expiration date
- Discard the entire device to biohazard container after use.
- Used Lancets should be treated as sharps waste and/or medical waste and disposed of properly following

**Clinical benefit:**

Help medical staff and patients complete capillary blood sampling in order to obtain a small blood sample for various tests.

Warning:

- Not intended for more than one use. Do not use on more than one patient.
- Improper use of blood lancets can increase the risk of inadvertent transmission of blood borne pathogens, particularly in settings where multiple patients are tested.
- The expiration time is 5 years. Using the product beyond the expiration date is strictly prohibited.
- For single use only. Do not use more than one patient. The device is deactivated after a single use, cannot be used more than once.
- The product has no therapeutic or diagnostic function.
- If the lancet's protective cap is damaged or lost, please do not use it.
- Use with caution for those with coagulation mechanism disorders.
- For obese users or individuals who require larger blood volumes, it is recommended to choose products with a smaller "G" value or large puncture depth for blood collection.

Shelf life:

Valid period of product is 5 years.

Storage and Transportation:

1. The device shall be stored at an environment where the temperature is between -20 °C and 40 °C and the relative humidity of less than 80%. The storage environment shall be free of corrosive gases, dry, away from sunlight, well-ventilated, and clean.
2. The device shall be handled with care during transportation and handling.

Disposal:

- The product must be disposed of after use.
- Products contaminated with blood must be disposed of in accordance with the valid regulations for hazardous materials.
- The product must not be used after the expiry date!
- Safe disposal must be carried out in accordance with current national regulations.

Adverse events:

We conducted a search for adverse events and state of the art literature on equivalent and similar devices in the market. No adverse events found for equivalent and similar products, indicating compliance with the current technical standards. If any adverse events occur during the use of the product, please contact us and we will resolve the relevant issues as soon as possible. Additionally, please note that any serious incidents related to the device should be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

Indice dei simboli - Symbol index - Index des symboles - Índice de símbolos- índice de símbolo - Symbolindex - Ευρετήριο συμβόλων - Indeks symboli - Index symbolů - Symbol index - Symboli-indeksi - Indeks simbolov - Index symbolov - Index de simbol - Symbool index - Indeks simbola - Szimbólum index - Symbolindeks - Indeks over symboler - Индекс на символа - Simbolių rodyklė - Simbolu rādītājs - Sūmbolite indeks - Указатель символов - تالوييوس تارشؤم

	IT - Fabbrikante GB - Manufacturer FR - Fabricant ES - Fabricante PT - Fabricante DE - Hersteller GR - Παραγωγός PL - Producent CZ - Výrobce SE - Tillverkare FI - Valmistaja SI - Proizvajalec SK - Výrobca RO - Producător NL - Fabrikant HR - Proizvođač HU - Gyártó DK - Fabrikant NO - Produsent BG - Производител LT - Gamintojas LV - Ražotājs EE - Tootja RU - Производитель AR - الشركة المصنعة
	IT - Data di fabbricazione GB - Date of manufacture FR - Date de fabrication ES - Fecha de fabricación PT - Data de fabrico DE - Herstellungsdatum GR - Ημερομηνία παραγωγής PL - Data produkcji CZ - Datum výroby SE - Tillverkningsdatum FI - Valmistuspäivämäärä SI - Datum proizvodnje SK - Dátum výroby RO - Data fabricației NL - Productiedatum HR - Datum proizvodnje HU - Gyártás dátuma DK - Fabrikationsdato NO - Fabrikasjonsdato BG - Дата на производство LT - Pagaminimo data LV - Izgatavošanas datums EE - Valmistamise kuupäev RU - Дата выпуска AR - تاريخ التصنيع
	IT - Dispositivo monouso, non riutilizzare GB - Disposable device, do not re-use FR - Dispositif pour usage unique, ne pas réutiliser ES - Dispositivo monouso, no reutilizable PT - Dispositivo descartável, não reutilizar DE - Für einmaligen Gebrauch, nicht wiederverwenden CZ - Jednorázový prostředek, nepoužívejte opakovaně DK - Engangsenhed, må ikke genbruges GR - Προϊόν μιας χρήσεως. Μην το χρησιμοποιείται εκ νέου. EE - Ühekordne, ärge kasutage seda mitu korda FI - Kertakäyttöinen laite, ei saa käyttää uudelleen HR - Uredaj za jednokratnu upotrebu, nemojte ponovo koristiti HU - Eldobható eszköz, ne használja újra LT - Vienkartinis prietaisas, nenaudokite pakartotinai LV - Vienreiz lietojama ierīce, nelietojiet to atkārtoti PL - Jedno urządzenie, nie używaj ponownie NL - Voor eenmalig gebruik, niet hergebruiken RO - Dispozitiv de unică folosință, a nu se refolosi SK - Zariadenie na jedno použitie, nepoužívajte opakovane SI - Za enkratno uporabo, ne uporabiti ponovno SE - Engångsanordning, får ej återanvändas NO - Engangsutstyr, ikke gjenbruk BG - Изделие за еднократна употреба, да не се използва повторно RU - Одноразовое изделие. Повторное использование запрещено AR - أداة أحادية الاستخدام، لا تستخدم من جديد
	IT - Sterilizzato per irradiazione GB - Sterilized using irradiation FR - Stérilisé par irradiation ES - Esterilizado por irradiación PT - Esterilizado por irradiação DE - Durch Bestrahlung sterilisiert GR - στείρεινεται με ακτινοβολήση PL - Sterylizowane napromienianiem CZ - Sterilizováno ozařováním SE - Steriliserad genom bestrålning FI - Steriloitu säteilyttämällä SI - Sterilizirano z obsevanjem SK - Sterilizované ožiaréním RO - Sterilizat prin iradiere NL - Gesteriliseerd door bestraling HR - Sterilizirano zračenjem HU - Besugárzással sterilizált DK - Steriliseret ved bestråling NO - Sterilisert ved bestråling BG - Стерилизиран чрез облъчване LT - Sterilizuojamas švitinant LV - Sterilizēts ar apstarošanu EE - Steriliseeritud - kiiritusega RU - Стерилизовано с использованием облучения AR - تعقيمها عن طريق التشعيع
	IT - Conservare in luogo fresco ed asciutto GB - Keep in a cool, dry place FR - À conserver dans un endroit frais et sec ES - Conservar en un lugar fresco y seco PT - Armazenar em local fresco e seco DE - An einem kühlen und trockenen Ort lagern GR - Διατηρείται σε δροσερό και στεγνό περιβάλλον PL - Przechowywać w suchym miejscu CZ - Skladujte na větraném a suchém místě SE - Förvara på svalt och torrt ställe FI - Säilytä kuivassa ja viileässä SI - Hraniti na suhem in hladnem mestu SK - Skladujte na chladnom a suchom mieste RO - A se păstra într-un loc răcoros și uscat NL - Koel en droog opslaan HR - Čuvati na hladnom i suhom mjestu HU - Száraz, hűvös helyen tárolandó DK - Opbevares køligt og tørt NO - Må oppbevares på et tørt og kaldt sted BG - Да се съхранява на хладно и сухо място LT - Laikyti vėsioje ir sausoje vietoje LV - Uzglabāt vēsā, sausā vietā EE - Hoida jahedas ja kuivas kohas RU - Хранить в сухом прохладном месте AR - يحفظ في مكان بارد وجاف
	IT - Limite di umidità GB - Humidity limit FR - Limite d'humidité ES - Límite de humedad PT - Limite de humidade DE - Feuchtigkeitsgrenzwert GR - Όριο υγρασίας PL - Granica wilgotności CZ - Limit vlhkosti SE - Fuktighetsgräns FI - Kosteusraja SI - Omejitev vlažnosti SK - Hranica vlhkosti RO - Limită de umiditate NL - Drempelwaarde vochtigheid HR - Granica vlažnosti HU - Páratartalom határérték DK - Grænse for luftfugtighed NO - Fuktighetsgrense BG - Ограничение на влажността LT - Drėgmės riba LV - Mitruma robeža EE - Niiskuse piirang RU - предел влажности AR - حد نسبة الرطوبة
	IT - Numero di lotto GB - Lot number FR - Numéro de lot ES - Número de lote PT - Número de lote DE - Chargennummer GR - Αριθμός παρτίδας PL - Kod partii CZ - Číslo šarže SE - Satsnummer FI - Eränumero SI - Številka partije SK - Číslo šarže RO - Număr de lot NL - Partijnummer HR - Broj serije HU - Tételszám DK - Batchnummer NO - Produksjonsserienummer BG - Номер на партида LT - Partijos numeris LV - Partijas numurs EE - Partii number RU - Номер партии AR - رقم الدفعة
	IT - Data di scadenza GB - Expiration date FR - Date d'échéance ES - Fecha de caducidad PT - Data de validade DE - Ablaufdatum GR - Ημερομηνία λήξεως PL - Data ważności CZ - Datum ukončení platnosti SE - Utgångsdatum FI - Viimeinen voimassaolopäivä SI - Rok uporabnosti SK - Dátum expirácie RO - Valabil până la data de NL - Vervaldatum HR - Datum isteka HU - Lejárati dátum DK - Udløbsdato NO - Utløpsdato BG - Срок на годност LT - Galiojimo laikas LV - Derīguma termiņš EE - Aegumiskuupäev RU - Срок годности AR - تاريخ انتهاء الصلاحية

	IT - Limite di temperatura GB - Temperature limit FR - Limite de température ES - Límite de temperatur PT - Limite de temperatura DE - Temperaturgrenzwert GR - Διατηρείται μεταξύ -10 και 49°C PL - Przechowyj pomiędzy 1 °C CZ - Uchovávejte při teplotě mezi a °C SE - Lagras mellan och °C FI - Säilytyslämpötila -10/49°C SI - Hranite pri temperaturi med in °C SK - Uchovávaťe pri teplote od do °C RO - A se păstra la temperaturi cuprinse între și °C NL - Bewaren tussen en °C HR - Čuvati između i °C HU - és °C között tárolandó DK - Grænse for luftfugtighed NO - Oppbevares mellom og °C BG - Да се съхранява между и °C LT - Temperatūros riba LV - Uzglabāt temperatūrā līdz °C EE - Temperatuuripiirang RU - Хранить при температуре °C AR - يحفظ بين و درجة مئوية
	IT - Attenzione - Leggere e seguire attentamente le istruzioni (avvertenze) per l'uso GB - Caution - read instructions (warnings) carefully FR - Attention - lisez attentivement les instructions (avertissements) ES - Precaución - lea las instrucciones (advertencias) cuidadosamente PT - Cuidado - leia as instruções (avisos) cuidadosamente DE - Achtung - Anweisungen (Warnungen) sorgfältig lesen GR - Προσοχή - διαβάστε προσεκτικά τις οδηγίες (εισατάξεις) PL - Ostrzeżenie - Zobacz instrukcję obsługi CZ - Pozor - Pečlivě si přečtěte a dodržujte pokyny (varování) k použití SE - Varsamhet - läs anvisningarna (varningar) noga FI - Huomio - Lue käyttöohjeet (varoitukset) ja noudata niitä huolellisesti SI - Pozor - Preberite in skrbno sledite navodilom (opozorilom) za uporabo SK - Pozor - Pozorne si prečítajte a dodržiavajte pokyny na použitie (výstrahy) RO - Atenție - Citiți și respectați cu atenție instrucțiunile (avertismentele) de utilizare NL - Opgelet - Lees en volg aandachtig de gebruiksaanwijzing (waarschuwingen) HR - Pozor - Pročitajte i pažljivo slijedite upute (upozorenja) za upotrebu HU - Figyelem - Figyelmesen olvassa el és kövesse a használati utasításokat (figyelmeztetéseket) DK - Forsigtig - Læs instruktioner (advvarsler) omhyggeligt NO - OBS! Les og følg anvisningene (advarslene) svært nøye BG - Внимание: Прочетете внимателно и спазвайте стриктно инструкциите (предупрежденията) за употреба LT - Dėmesio - perskaitykite ir atidžiai laikykites naudojimo instrukcijų (spėjimų). LV - Uzmanību - Izlasiet un uzmanīgi ievērojiet lietošanas instrukcijas (brīdinājumus) EE - Tähelepanu! Lugege kasutusjuhend (hoiatused) läbi ja järgige neid hoolikalt RU - Внимательно читайте инструкции ((предупреждения AR - الحذر - قراءة التعليمات (التحذيرات) بعناية
	IT - Identificatore univoco del dispositivo GB - Unique device identifier FR - Identifiant unique de l'appareil ES - Identificador de dispositivo único PT - Identificador exclusivo do dispositivo DE - Unique Device Identifier (Eindeutige Kennung des Geräts) GR - Μοναδικό αναγνωριστικό συσκευής PL - Unikalny identyfikator urządzenia CZ - Jediněčný identifikátor zařízení SE - Unik identifierare för enheten FI - Laitteen yksilöllinen tunniste SI - Enolični identifikator naprave SK - Jediněčný identifikátor zariadenia RO - Identificatorul unic al dispozitivului NL - Unieke identificatie van het apparaat HR - Jedinstveni identifikator uređaja HU - Az eszköz egyedi azonosítója DK - Unik identifikator for enheden NO - Unik identifikator for enheten BG - Уникален идентификатор на устройството LT - Unikalus įrenginio identifikatorius LV - Unikāls ierīces identifikators EE - Seadme kordumatu identifikaator RU - Уникальный идентификатор устройства AR - معرف فريد للجهاز
	IT - Codice prodotto GB - Product code FR - Code produit ES - Código producto PT - Código produto DE - Erzeugniscode GR - Κωδικός προϊόντος PL - Numer katalogowy CZ - Kód výrobku SE - Produktkod FI - Tuotekoodi SI - Koda izdelka SK - Kód výrobku RO - Cod produs NL - Productcode HR - Šifra proizvoda HU - Termékkód DK - Produktkode NO - Produktkode BG - Код на продукта LT - Prekės kodas LV - Produkta kods EE - Toote kood RU - Код изделия AR - كود المنتج
	IT - Dispositivo medico GB - Medical Device FR - Dispositif médical ES - Producto sanitario PT - Dispositivo médico DE - Medizinprodukt GR - Ιατροτεχνολογικό προϊόν PL - Wyrób medyczny CZ - Zdravotnický prostředek SE - Medicinteknisk produkt FI - Lääkinnällinen laite SI - Medicinski pripomoček SK - Zdravotnícka pomôcka RO - Dispozitiv medical NL - Medisch hulpmiddel HR - Medicinski uređaj HU - Orvostechnikai eszköz DK - Medicinsk udstyr NO - Medisinsk utstyr BG - Медицинско изделие LT - Medicininis prietaisas LV - Mediciniskā ierīce EE - Meditsiiniline seade RU - Медицинское изделие AR - جهاز طبي
	IT - Conservare al riparo dalla luce solare GB - Keep away from sunlight FR - À conserver à l'abri de la lumière du soleil ES - Conservar al amparo de la luz solar PT - Guardar ao abrigo da luz solar DE - Vor Sonneneinstrahlung geschützt lagern GR - Κρατήστε το μακριά από ηλιακή ακτινοβολία PL - Przechowywać z dala od światła słonecznego CZ - Skladujte mimo sluneční světlo SE - Skyddas från solljus FI - Säilytä auringonvalolta suojassa SI - Hraniti zaščiteno pred sončno svetlobo SK - Skladujte mimo slnečného svetla RO - A se păstra ferit de razele soarelui NL - Afgeschermd van zonlicht opslaan HR - Čuvati zaštićeno od sunčeve svjetlosti HU - Napfénytől védve tárolandó DK - Må ikke udsættes for sollys NO - Må oppbevares på et sted uten direkte sollys BG - Да се съхранява на място, защитено от слънчева светлина LT - Saugoti nuo saulės spindulių LV - Uzglabāt prom no saules gaismas EE - Hoida eemal päikesevalgusest RU - Беречь от попадания солнечных лучей AR - يحفظ بعيدًا عن ضوء الشمس
	IT - Sistema di barriera sterile singolo GB - Single sterile barrier system FR - Système de barrière stérile unique ES - Sistema de barrera estéril única PT - Sistema de barreira estéril único DE - Einzelbarriere-Sterilsystem GR - Μονό αποστειρωμένο σύστημα φραγμού PL - System pojedynczej bariery sterylnej CZ - Sterilní systém s jednou bariérou SE - Sterilt system med en barriär FI - Yhden esteen steriili järjestelmä SI - Sterilni sistem z eno pregrado SK - Sterilný systém s jednou bariérou RO - Sistem steril cu o singură barieră NL - Steriel systeem met enkele barrière HR - Sterilni sustav s jednom barijerom HU - Single barrier steril rendszer DK - Enkeltbarriere steril system NO - Enkelt steril barriersystem BG - Единична барьерна стерилна система LT - Vieno barjero sterili sistema LV - Viena barjera sterila sistēma EE - Ühe barjääriga steriilne süsteem RU - Не содержит фталата ДЭФ AR - نظام حاجز معقم واحد

	IT - Importato da GB - Imported by FR - Importé par ES - Importado por PT - Importado por DE - Eingeführt von GR - Εισαγώγη από PL - Importowane przez CZ - Dovezeno uživatelem SE - Importerad av FI - Tuojaa SI - Uvozil SK - Dovážal RO - Importat de NL - Geïmporteerd door HR - Uvezeno od strane HU - Importálta DK - Importeret af NO - Importert av BG - Внесено от LT - Importavo LV - Importēja EE - Importitja RU - Ввозимые AR - مستورد عن طريق
	IT - Leggere le istruzioni per l'uso GB - Consult instructions for use FR - Consulter les instructions d'utilisation ES - Consultar las instrucciones de uso PT - Consultar as instruções de uso DE - Gebrauchsanweisung beachten GR - Διαβάστε προσεκτικά τις οδηγίες χρήσης PL - Przeczytaj instrukcję użytkowania CZ - Přečtěte si návod k použití SE - Läs bruksanvisningen FI - Lue käyttöohjeet SI - Preberite navodila za uporabo SK - Prečítajte si návod na použitie RO - Citiți instrucțiunile de utilizare NL - Lees de gebruiksaanwijzing HR - Pročitajte upute za uporabu HU - Olvassa el a használati utasításokat DK - Se brugsvejledningen NO - Les bruksinstruksjonene BG - Прочетете инструкциите за употреба LT - Perskaitykite naudojimo instrukcijas LV - Izlasiet lietošanas instrukcijas EE - Lugege kasutusjuhendit RU - См. инструкцию AR - اقرأ بدقة وحرص تعليمات الاستخدام
	IT - Dispositivo medico conforme al regolamento (UE) 2017/745 GB - Medical Device compliant with Regulation (EU) 2017/745 FR - Dispositif médical conforme au règlement (UE) 2017/745 ES - Producto sanitario conforme con el reglamento (UE) 2017/745 PT - Dispositivo médico em conformidade com a regulamentação (UE) 2017/745 DE - Medizinprodukt im Sinne der Verordnung (EU) 2017/745 GR - Ιατρική συσκευή σύμφωνα με την ΚΑΝΟΝΙΣΜΟΣ (ΕΕ) 2017/745 PL - Wyrób medyczny zgodny z Rozporządzenie (UE) 2017/745 CZ - Zdravotnický prostředek v souladu s nařízením (EU) č. 2017/745 SE - Den medicintekniska produkten överensstämmer med förordning (EU) 2017/745 FI - Lääkinnällinen laite, joka vastaa asetusta (EU) 2017/745 SI - Medicinski pripomoček, skladen z uredbo (EU) 2017/745 SK - Zdravotnícka pomôcka v súlade s nariadením (EU) 2017/745 RO - Dispozitiv medical conform regulamentului (UE) 2017/745 NL - Medisch hulpmiddel in overeenstemming met verordening (EU) 2017/745 HR - Medicinski proizvod sukladan propisu (EU) 2017/745 HU - A 2017/745/EU rendeletnek megfelelő orvostechnikai eszköz DK - Medicinsk udstyr i overensstemmelse med forordning (EU) 2017/745 NO - Medisinsk utstyr som samsvarer med gjeldende regelverk (EU) 2017/745 BG - Медицинско изделие в съответствие с Регламент (ЕС) 2017/745 LT - Medicinos prietaisas, atitinkantis reglamentą (ES) 2017/745 LV - Medicīniska ierīce, kas atbilst Regulai (ES) 2017/745 EE - Määrusele (EL) 2017/745 vastav meditsiiniseade (UE) 2017/745 RU - Медицинское изделие отвечает положениям Директивы (UE) 2017/745 AR - جهاز طبي يتوافق مع التوجيه
	IT - Non utilizzare se l'imballaggio è danneggiato GB - Don't use if package is damaged FR - Ne pas utiliser si le colis est endommagé ES - No usar si el paquete está dañado PT - Não use se o pacote estiver danificado DE - Nicht verwenden, wenn das Paket beschädigt ist CZ - Nepoužívejte, pokud je obal poškozen DK - Må ikke bruges, hvis pakken er beskadiget GR - Μην το χρησιμοποιείτε αν η συσκευασία είναι κατεστραμμένη FI - Ei saa käyttää, jos pakkaus on vaurioitunut HR - Ne koristiti ako je pakiranje oštećeno HU - Ne használja, ha a csomagolás sérült NL - Nauwendokite, jei pakuoite pažeista LV - Nelietot, ja iepakojums ir bojāts PL - Nie używać, jeśli opakowanie jest uszkodzone NL - Niet gebruiken als de verpakking beschadigd is RO - A nu se utiliza dacă ambalajul este deteriorat SK - Nepoužívaťe, ak je obal poškodený SI - Ne uporabljajte, če je embalaža poškodovana SE - Använd inte en förpackning som är skadad NO - Ikke bruk hvis emballasjen er skadet BG - Да не се използва, ако опаковката е с нарушена цялост EE - Ärge kasutage, kui pakend on kahjustatud RU - В случае повреждения упаковки не использовать AR - لا تستخدم في حالة تلف الحزمة
	IT - Non-pirogeno GB - Non-pyrogenic FR - Apyrogène ES - Apirógeno PT - Não pirogénico DE - Nicht pyrogen GR - Μη πυρογόνο PL - Niepirogenny CZ - Niepyrogenní SE - Icke-pyrogen FI - Ei-pyrogeneinen SI - Nepirogen SK - Nepyrogénny RO - Apirogen NL - Niet-pyrogeen HR - Nepirogeno HU - Nem pirogén DK - Ikke-pyrogen NO - Ikke-pyrogen BG - Непириогенен LT - Nepirogeninis LV - Nepirogēns EE - Mittepirogeenne RU - Апириогенный AR - غير يبروجيني
	IT - Rappresentante autorizzato nell'Unione Europea GB - Authorized representative in the European Union FR - Représentant autorisé dans l'Union européenne ES - Representante autorizado en la Unión Europea PT - Representante autorizado na União Europeia DE - Autorisierter Vertreter in der Europäischen Union GR - Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση PL - Autoryzowany przedstawiciel w Unii Europejskiej CZ - Autorizovaný zástupce v Evropské unii SE - Auktoriserad representant i Europeiska unionen FI - Valtuutettu edustaja Euroopan unionissa SI - Pooblaščen zastopnik v Evropski uniji SK - Autorizovaný zástupca v Európskej únii RO - Reprezentant autorizat în Uniunea Europeană NL - Gemachtigde vertegenwoordiger in de Europese Unie HR - Ovlašteni predstavnik u Europskoj uniji HU - Hivatalos képviselő az Európai Unióban DK - Autoriseret repræsentant i Den Europæiske Union NO - Autorisert representant i EU BG - Упълномощен представител в Европейския съюз LT - Įgaliotasis atstovas Europos Sąjungoje LV - Pilnvarotais pārstāvis Eiropas Savienībā EE - Volitatud esindaja Euroopa Liidus RU - Уполномоченный представитель в Европейском Союзе AR - الممثل المعتمد في الاتحاد الأوروبي
	IT - Rappresentante autorizzato nel Regno Unito GB - Authorized Representative in the UK FR - Représentant autorisé au Royaume-Uni ES - Representante autorizado en el Reino Unido PT - Representante autorizado no Reino Unido DE - Autorisierter Vertreter im Vereinigten Königreich GR - Εξουσιοδοτημένος Αντιπρόσωπος στο Η.Β PL - Autoryzowany Przedstawiciel w Wielkiej Brytanii CZ - Autorizovaný zástupce ve Spojeném království SE - Auktoriserad representant i Storbritannien FI - Valtuutettu edustaja Isonsa-Britanniassa SI - Pooblaščen zastopnik v Združenem kraljestvu SK - Autorizovaný zástupca v Spojenom kráľovstve RO - Reprezentant autorizat în Marea Britanie NL - Geautoriseerde vertegenwoordiger in het VK HR - Ovlašteni predstavnik u UK HU - Meghatalmazott képviselő az Egyesült Királyságban DK - Autoriseret repræsentant i Storbritannien NO - Autorisert representant i Storbritannia BG - Оторизиран представител в Обединеното кралство LT - JK įgaliotasis atstovas LV - Pilnvarotais pārstāvis Apvienotajā Karalistē EE - Volitatud esindaja Ühendkuningriigis RU - Уполномоченный представитель в Великобритании AR - الممثل المعتمد في المملكة المتحدة