GIMA

PROFESSIONAL MEDICAL PRODUCTS

I.V. INFUSION STAND ON 5 WHEELS TROLLEY WITH HANDLE AND SHELF - S/S - 4 HOOKS - LOAD 20 KG

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
Jiangsu Saikang Medical Equipment Co., Ltd
No. 35 Lehong Road, Modern Agriculture Demonstration
Park, Zhangjiagang City, Jiangsu Province, China
Made in China

REF SKH041(16)

Gima 27816

EC REP EVERBIZ GmbH
LANDSBERGER STR. 155
80687 MUNCHEN GERMANY

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

**ENGLISH
INTENDED USE**

Intravenous (IV) stands or poles are medical devices that are designed to hang bags containing intravenous fluids or medicines which need to be administered to the patient.

1. Installing guide

1. **Install the lower rod to the base.**
 - Connect the lower rod to the 5-foot base by screws and fix.
2. **Install the plastic platform to the lower rod.**
 - Connect the platform to the lower by screws and fix.
3. **Install the upper rod to the platform.**
 - Connect the upper rod to the plastic platform by screws and fix.
4. **Install the 4 hooks on the top.**
 - Insert the 4 hooks on the top of the upper rod.

2. Use method

2.1 Use method of height adjustment

Adjust the knob, adjust the infusion stand to the appropriate height, check whether the screws are tight, ready to use.

3. Maintenance

1. In order to use the IV pole safely, it is necessary to perform regular safety inspections on the IV 1. In order to use the IV stand safely, it is necessary to perform regular safety inspections on the IV pole. It is recommended to check it every six months to ensure that the connections are not loose and all operations are normal.
2. When the IV stand frame ages and it reaches a certain servicing time, which the metal part of the IV pole and plastic parts of the stainless steel basin frame can be recycle.
3. Avoid scratching the panel with sharp-angle appliances or knives during use, and clean it frequently to keep it clean and dry.
4. If the panel is accidentally stained with stains, it is recommended to clean it. Do not use alkaline or corrosive chemicals to clean the IV pole, which will cause the stainless steel surface to rust.
5. Regularly check the casters to avoid hitting or overloading, which will cause damage to the casters.
6. It is forbidden to force the casters to brake when braking, which will cause damage to the brake system.
7. The wear of the tire surface of the wheel can be detected by visual inspection. Some thread spun yarn and other accumulations may be wrapped around the wheel, remove the bolts and nuts of the wheel, clean up the debris, and check whether the bearing of the wheel is damaged. If the parts are not damaged, can be reassembled and used.
8. Judging whether the caster is replaced according to the time used, wear and tear, you can replace the wheels and make sure that the axles bolt and nuts are well put together and connected. If possible, use a lock washers or locknuts to help you with this process.
9. For the casters equipped with brakes, it is necessary to regularly check whether the brakes are working properly. Check the brakes once a day or before each use. For a IV stand equipped with multiple brake casters, only one brake caster can be locked at a time which will allow you to try and push the IV pole also check whether one of the brake performance of each caster is good. If the function of the brake fails due to wear or damage to the wheels, immediately replace the wheels and then retest the brake system.
10. If the mechanism of the brake system of the caster is damaged, and the brake needs to be repaired or replaced, please contact our after-sales service department or authorized dealer provide by our company. Every time the brakes are replaced, the brake performance of the casters must be retested.

4. Packaging, Transport, Storage

1. The packaging of stainless steel I.V. stand is carried out according to the contract or product standard.
2. Avoid rushing out, violent vibration, and prevent from sunlight and rain during the transportation of the stainless steel I.V. stand.
3. The IV stand should be stored:
 - a) Ambient temperature: -40°C ~ +70°C.
 - b) Relative humidity: ≤ 95%.
 - c) Atmospheric pressure: 500 hPa ~ 1060 hPa.

5. Max Load

Hooks: max load for each hook 3kg
 Pump Tray: max load 10kg
 Handle: max load 10kg
 Overall Structure: max load 32kg

6. After-sales service

1. Please keep the files that comes with the machine and invoices of this product properly, and you need to present these files when the company performs warranty and maintenance for the product.
2. If there is any problem in the process of using, please contact our company in time, so that our company can provide you with accurate and fast technical support and maintenance services in a timely manner.
3. From the date of sale, if the product is broken or damaged due to the correct installation and use according to the specification, the product will enjoy one-year free warranty and lifelong maintenance service with the "certificate" or invoice.
4. Do not disassemble the internal device of this product by yourself to avoid unnecessary damage. If you find that the use is affected due to quality problems, please contact our after-sales service department or our authorized dealer.
5. From the date of purchase, if it is indeed damaged or does not work normally due to quality problems within one year, the company will provide to the users a free repairing of the product.
6. Lifetime service from manufacturer: JIANGSU SAIKANG MEDICAL EQUIPMENT CO LTD

Notice! Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

| | | | |
|--|---|--|--------------------|
| | Lot number | | Manufacturer |
| | Authorized Representative in the European Community | | Imported by |
| | Unique Device Identifier | | Manufacturing Date |

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| | | | |
|--|----------------------------|--|--|
| | Atmospheric pressure limit | | Keep in a cool, dry place |
| | Humidity limit | | Caution: read instructions (warnings) carefully |
| | Temperature limit | | Consult instructions for use |
| | Medical Device | | Medical device compliant with Regulation (EU) 2017/745 |
| | Keep away from sunlight | | Product code |