

ENGLISH

SILICONE HOLLOW FIBER CUSHIONS

WARNINGS

Before performing any operation, visually check the integrity of the device. If the product shows any anomalies, immediately contact the point of sale where the purchase was made. Erregi sas assumes no responsibility for the improper use of the device. In the event of a serious incident related to the use of the device, please report the incident to the manufacturer and the competent authority of the member state where the user and/or patient is located.

DESCRIPTION

Cushions for the prevention of pressure ulcers, made of cotton with hollow siliconized fiber padding. Some devices are equipped with fastening systems using Velcro closures. High resilience and excellent breathability allow for effective prevention of pressure ulcer formation.

INTENDED USE

Pressure ulcer prevention for adult patients aimed at reducing pressure peaks, minimizing shear and friction forces, and promoting air circulation.

INDICATIONS AND PRECAUTIONS FOR USE

Medical device indicated for the prevention of pressure ulcers, to be used following clinical evaluation by a specialist physician. Place the cushion on the wheelchair seat, ensuring that the Velcro fastenings, where present, are positioned toward the backrest. Then secure the cushion to a rigid part of the wheelchair frame using the Velcro fastenings. Although the device is intended for the prevention of pressure ulcers, it cannot control all the factors that contribute to their development.

It is recommended to regularly monitor the affected area, frequently reposition the patient, and pay particular attention to hygiene. For this purpose, it is advisable to consult protocols and guidelines for the prevention of pressure ulcers. Always use the device while wearing a protective garment. This device is not fire-resistant.

CONTRAINDICATIONS

No contraindications have been detected, except in cases of proven sensitivity to the component materials.

COMPOSITION

COVER: Cotton

FILLING: Siliconized hollow fiber

WASHING INSTRUCTIONS

			
WASH AT MAX 40°C	DO NOT BLEACH	DO NOT IRON	DO NOT TUMBLEDRY
	USE ANY SOLVENT EXCEPT TETRACHLORETHYLENE		

The manufacturer is not responsible for any changes in the technical specifications of the device in case of non-compliance with the washing instructions.

DISINFECTION

Use products for medical-surgical use.

STORAGE

Store in a place protected from light, heat, and moisture.

DISPOSAL

Dispose of the device in accordance with the applicable environmental protection regulations and separate waste collection.

	Product code		Medical Device
	Lot number		Manufacturer
	Medical Device compliant with Regulation (EU) 2017/745		Consult instructions for use
	Keep away from sunlight		Keep in a cool, dry place
	Date of manufacture		Not made with natural rubber latex
	Authorized Representative in Switzerland		Distributed by
	Unique Device Identifier		

DO NOT remove the label from the Medical Device, as it allows tracing the supplied product (Manufacturer, UDI, batch number, product code).

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Gima 28549 - 28550 - 28551 - 28552



Erregi S.a.s.

Via Adua 40
21040 Vedano Olona (Va)
commerciale@erregiausilimedici.it
www.ausilimedici.it
MADE IN ITALY



03/04 - 01/04 - 05/04 - 04/04



Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in Italy



Élite Med Group Sagl

Via Luigi Lavizzari 4
6850 Mendrisio
Switzerland

M28549-Rev.5-07.26

