PROFESSIONAL MEDICAL PRODUCTS

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**CLEANING AND STERILIZATION INSTRUCTION FOR HOSPITAL
PLASTIC HOLLOWARES**

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All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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Devices used for a variety of medical purposes, such as fluids container, body waste collector, holder/carrier for instruments, intended for medical use as support of the medical and/or surgical procedure in which they are used.

1. Cleaning

Clean before and after use. (The first and any time after use, the device must be cleaned with warm water, disinfected and rinsed if possible with distilled water and carefully dried).

2. Disinfection, rinsing and drying

Pay attention that only neutral cleaning agents and disinfectants without chlorine, ammonia, phenol, aldehydes that are suitable for cleaning and disinfecting plastic instruments are used.

Disinfectants must conform to the applicable national regulations for disinfectants (such as FDA approval and CE mark).

Dry all the parts with lint-free disposable cloths (do not use recycled cloths). Use oil-free compressed air to dry holes.

3. Sterilization

Non sterile, sterilize if required (European standard ISO 17665 i.e. steam sterilization for medical devices to be followed).

Sterilization is achieved by exposing products to saturated steam at high temperatures (121°C to 134°C). Product(s) are placed in device called the autoclave and heated through pressurized steam to kill all microorganisms including spores. The device's exposure time to steam would be anywhere between 3 to 15 minutes.

Examples of minimum temperatures and times established for adequate level of microbial lethality in sterilization processes (refer table – 1 of ISO/TS 17665-2).

Temperature (°C)	Time Minutes
121	15
126	10
134	03

4. The device showing damage or crack must be immediately replaced and not to be used

Recommended storage condition: device to be stored in clean and dry condition.

	Medical Device compliant with Regulation (EU) 2017/745		Product code		Non-sterile
	Caution: read instructions (warnings) carefully		Lot number		Keep away from sunlight
	Consult instructions for use		Manufacturer		Keep in a cool, dry place
	Date of manufacture		Temperature limit		Humidity limit
	Medical Device		Unique device identifier		Authorized Representative in the UK
	Authorized Representative in Swiss				

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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Gima S.p.A.