



IT

HANDHELD TYMPANOMETER SCREENING AUDIOMETER

EN

TIMPANI


FR

MULTILANGUAGE USER MANUAL

DE

Document title: IM1P-Timpani User Manual IT-EN-FR-DE-ES
Code: IM1-MA008_A
Revision: Rev. 02
Date: 2025.02.19

ES



TYMPANOMETER

TIMPANI

USER MANUAL



Read this manual thoroughly before using the device. Pay particular attention to Chapter 1 (“Safety: warnings and information”) and Chapter 2 (“Installation, power-up and power down”).



Internal inspections and repairs must only be performed by authorized personnel.

Copyright: INVENTIS S.r.l. has the copyright on this manual. It may not be copied, reproduced or altered, in total or in any part, without the written specific authorization from INVENTIS S.r.l..

Inventis ® is a registered trademark owned by INVENTIS S.r.l.



Summary

<i>Foreword</i>	v
<i>Chapter 1 Safety: warnings and information</i>	1
1.1 Operator Manual	1
1.2 Operator responsibilities	1
1.3 Intended purpose.....	2
1.4 Indications for use and end users.....	2
1.5 Medical conditions	2
1.6 Main features	2
1.7 Use cases.....	3
1.8 Precautions.....	3
1.9 Disposal	6
1.10 Conformity.....	6
1.11 Symbols on labels	7
<i>Chapter 2 Installation, power-up and power down</i>	9
2.1 Opening the pack and inspecting the contents.....	9
2.2 Precautions.....	9
2.3 Connections.....	10
2.4 Power-up, power-down and main screen	10
2.5 Connection to PC.....	11
<i>Chapter 3 Controls and exams</i>	13
3.1 Controls for direct access to exams.....	13
3.2 Tympanometry	14
3.3 Acoustic reflex	15
3.4 Audiometry	17
3.5 Patient Management	19
3.6 Settings	21
<i>Chapter 4 Docking station</i>	25
4.1 Installation and general precautions	25
4.2 Conformity.....	26

4.3	Symbols on labels	26
4.4	Connections	27
Chapter 5 Maintenance		29
5.1	Periodic checks.....	29
5.2	Maintenance of Transducers.....	30
5.3	Cleaning the probe.....	31
5.4	Cleaning the instrument	33
5.5	Replacing the battery.....	33
5.6	Repairs and technical assistance.....	33
Chapter 6 Troubleshooting		35

Foreword

Thank you for purchasing an Inventis audiology device.

Advantageously compact and lightweight, the Timpani tympanometer is a powerful and versatile screening device, ideal for fast and accurate examinations of the middle ear.

EN

The Inventis company has always considered the use of its devices in conjunction with computers to be a factor of key importance. Installing the Maestro software suite, available with or without proprietary database or as a Noah module, any Inventis audiology device can be connected to a computer, and all examinations conducted then archived in the user's own database.

Bear in mind also that Inventis has developed a complete line of audiology devices: in addition to middle ear analyzers, the company's product line includes a range of audiometers, REM and HIT hearing aid fitting devices, a wireless video otoscope, and much more.

For further information, and to report any problems of whatever description that may be encountered, contact the company at:



INVENTIS S.r.l.
Corso Stati Uniti, 1/3
35127 Padova, Italia
Tel.: 049.8962844 – Fax: 049.8966343
www.inventis.it info@inventis.it

Chapter 1

Safety: warnings and information

1.1 OPERATOR MANUAL

Be sure to read this manual through completely, so that all of the features offered by the instrument can be used to their full potential. Pay particular attention to this chapter as it contains warnings that are of fundamental importance in ensuring safe and correct use of the device.

The safety warning symbol illustrated below is used in this manual to draw the attention of the reader to information of particular importance in matters of safety, and to guard against incorrect use.



1.2 OPERATOR RESPONSIBILITIES

The Timpani tympanometer guarantees consistent and dependable performance only when used in accordance with the instructions and procedures described in this manual.

Should the device need to undergo repair or maintenance work, it must be disconnected from the electrical power supply and not used again until after the repair work has been completed. Defective or faulty parts must only be replaced with original spare parts supplied by Inventis and all repairs must be carried out exclusively by Inventis or by staff authorized by Inventis. No parts of the device must be modified or replaced without authorization from Inventis.

The user assumes full responsibility for any malfunction resulting from improper use or operation, likewise from maintenance or repair work performed by third parties other than Inventis or Service Centers approved by Inventis. Inventis and approved Service Centers will answer for the performance and reliability of the equipment only if:

- adjustments, modifications or repairs have been entrusted to persons authorized by Inventis;

- the electrical system and earthing of the installation are in compliance with the specifications of standards for electro-medical appliances.

1.3 INTENDED PURPOSE

The Timpani tympanometer is a medical device intended to measure the biomechanical characteristics of the patient's middle ear, to help the operator assess his or her functional conditions for screening purposes.

Timpani is also a pure tone audiometer: by generating and offering the patient sound stimuli of different types and intensities, it helps the operator to assess the patient's auditory sensitivity for the purpose of screening.

1.4 INDICATIONS FOR USE AND END USERS

Timpani is intended for use by healthcare ENT professionals in hospitals, ENT clinics and audiology clinics as a tool for hearing screening programs and as an aid to the diagnosis of possible hearing disorders.

There is no limitation of the patient population in using the device. Always be sure to perform an otoscopy before using the device.

These tests must be conducted in a quiet environment to avoid artifacts.

1.5 MEDICAL CONDITIONS

Conditions of impaired sensitivity of the auditory system or any conditions in which the auditory system is thought to play a role in diagnosis.

1.6 MAIN FEATURES

The Timpani is a portable device that can be used to conduct middle ear screening tests simply, swiftly and accurately. Available with a range of optional licenses, the device is able to meet the needs of private medical practices, clinics and hospitals alike.

The device features:

- backlit color display with touchscreen interface, allowing graphic illustration of test results;
- compact and ergonomic design, lightweight construction;
- long endurance, with built-in rechargeable lithium battery;
- interaction with computer using Maestro software.

Depending on the licenses activated, the main functionalities available are:

- 226 Hz tympanometry test;
- 1000 Hz tympanometry test (with *1 kHz Probe Tone* license);
- ipsilateral acoustic reflex test with 226 Hz probe tone and stimuli:
 - o 1000 Hz with *Reflex – Basic (Stimuli @ 1 kHz)* license
 - o 500 Hz, 1000 Hz, 2000 Hz and 4000 Hz with *Reflex – Plus (Stimuli @ 0.5, 1, 2, 4 kHz)* license.
- Pure-tone audiometry test (with *Tonal Audiometry* license).

Other accessories available are a dedicated docking station and a portable thermal printer.

EN

1.7 USE CASES

Timpani can be used to conduct automatic tympanometry tests at low frequency (226 Hz) and at high frequency (1000 Hz, only with *1 kHz Probe Tone* license) and ipsilateral acoustic reflex tests (only with *Reflex – Basic* or *Reflex – Plus* license). Moreover, it is also possible to perform the audiometry test using pure tones by activating the *Screening Pure Tone Audiometry* license.

These tests must be conducted in a particularly quiet environment in order to avoid artifacts.

The Timpani tympanometer is intended for use by persons who have a detailed knowledge of the procedures associated with the tests supported by the instrument; the operator must therefore be either an audiometrist (or a technician having the requisite levels of audiological knowledge) or a medical practitioner in possession of specific skills (ENT or audiology specialist).

1.8 PRECAUTIONS



Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

To ensure correct and safe use of the device, the following precautions must be observed.

1.8.1 Installation and general precautions



Make certain that the required ambient conditions are met (during transport, storage, and operation):

<i>Operation</i>	<i>Temperature between +15°C and +35°C Relative humidity: between 30% and 90% (no condensation) Pressure: between 700 hPa and 1060 hPa</i>
<i>Transport and storage</i>	<i>Temperature between -10°C and 50°C Relative humidity: 90% max, no condensation Pressure: between 500 hPa and 1060 hPa</i>
<i>Warm-up time</i>	<i>1 minute</i>



The device will not be protected if exposed during use to flammable anesthetic gases or similar products. Risk of explosion.



Avoid installing and using the device near sources of strong electromagnetic fields: these could interfere with the operation of the equipment.



Use only original accessories supplied by Inventis, unless specifically indicated otherwise.



Use only medical grade power adapters, certified to IEC 60601-1, with the following specifications:

<i>Internal battery</i>	<i>Rechargeable Li-Ion, 18650 standard, 3.7V 2.6Ah</i>
	<i>Endurance: Minimum 4h continuous use</i>
	<i>Auto-off time: 5 minutes</i>
	<i>Stand-by time: 1 minute</i>
<i>External power adapter</i>	<i>Recharge time: from PC, standard USB port: 10h max; from dedicated power adapter: 3h max</i>
	<i>Input 100-240Vac 50/60Hz, 0.3-0.15A</i>
	<i>Output 5Vdc 1.4A</i>
	<i>Compliant with the IEC 60601-1 standard.</i>



The Timpani is a medical device: if connected to a computer (or any external device) located within the "patient area" (as defined in IEC 60601-1), this likewise must be a medical device, or protected by an isolating transformer, in order to ensure that the

combination of computer (external device) + tympanometer is in compliance with IEC 60601-1.



The Timpani must be installed and operated taking account of the information regarding electromagnetic compatibility (EMC) provided at the end of this manual.



The proximity of portable and mobile appliances used for RF communications can affect the operational efficiency of the instrument box. Refer to the information regarding electromagnetic compatibility (EMC) provided at the end of this manual.

1.8.2 Calibration



The calibration should be performed at least once every 12 months, and whenever a transducer is replaced.



The calibration of the instrument is valid only for the transducers supplied. If a transducer is replaced, the instrument must be recalibrated.



The calibration is valid for transducers supplied with the equipment, if connected directly to the instrument, without any interposition of extension leads. If the transducers are not connected directly to the device, a new calibration procedure will be required before the instrument is used.



If the transducer selected is not calibrated, an alert will appear in the test screens. It will not be possible to present any stimulus to the patient using transducers that have not been calibrated.



Take note of the calibration interval indicated. Use of the instrument beyond the calibration interval expiry date can lead to unreliable diagnoses.

1.8.3 Hygiene



The earpieces of the tympanometer probes are disposable, likewise that of insert earphones; do not use the same earpiece for different patients. Dispose of earpieces after use.



Disinfect the headphone pads between one patient and the next, following the procedure described in Chapter 5: Maintenance

1.8.4 Use



The instrument can generate tones at an intensity potentially damaging to the patient. Take particular care to set the intensity of the tone correctly before it is presented.



Never insert the probe tip into the ear canal without affixing an ear tip as this may damage the patient's ear canal.



Be sure to insert the probe tip in a way, which will assure an air tight fit without causing any harm to the patient. Using a proper and clean ear tip is mandatory.



When conducting audiometry using insert earphones, do not insert or in any way try to conduct measurements without proper foam tip in place.



Do not perform service or maintenance while using the device on a patient.

1.9 DISPOSAL

Like any other electronic device, the Timpani tympanometer contains certain extremely hazardous substances, albeit in extremely small proportions. If released into the normal refuse collection system without suitable preliminary treatment, these substances cause serious damage to the environment and to health. At the end of its service life, accordingly, each component of the appliance must go through a sorted collection process: this means that the user should deliver (or dispatch) waste items to the sorted collection centers set up by local authorities, or alternatively, hand them back to the dealer when purchasing a new appliance of the same or similar type.

Thanks to the sorted collection of waste items and the subsequent processing, recovery and disposal operations they undergo, appliances can be made from recycled materials, and any negative impact of improper waste management on the environment and on health can be suitably limited.

1.10 CONFORMITY

The Timpani tympanometer is a class IIa medical device, according to Annex VIII of Medical Device Regulation (MDR) 2017/745/EU.

The Inventis Quality Management System has been certified by leading assessment body TÜV as compliant with ISO 13485 standard.

1.11 SYMBOLS ON LABELS



Warning: the use of this device requires certain precautions.

To ensure safe use, consult the accompanying documentation.



Follow the instructions for use.



Device serial number. The number is made up of 13 alphanumeric characters indicating the model, series, year of manufacture and serial number:

- characters 1-5: Inventis product code
- characters 6-7: year of manufacture ("20" denotes 2020)
- characters 8-13: pro serial number



Catalog code



Name and address of manufacturer.



Medical Device



Device with applied parts, Type B (IEC 60601-1)



The device emits radio frequency



Product conforms to European Community Medical Device Regulation (MDR) 2017/745/EU. Class IIa device; number of notified body: 0123 (TÜV SÜD Product Service GmbH).

Rx Only

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.



The product is subject to the requirements of Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). In the event of this product being sold and/or scrapped, it must not be disposed of as ordinary household or industrial waste, but collected separately.



Do not reuse. Components bearing this mark can be used once only, and must not be reused thereafter.

(01)80541873807472
(21)IM1PA18200595



UDI code

Chapter 2

Installation, power-up and power down

2.1 OPENING THE PACK AND INSPECTING THE CONTENTS

Upon receiving the pack, check that the box is not damaged and that the parts contained are neither damaged nor defective.

Having made the various connections, carry out a further visual inspection before switching on, to check for possible damage.

Should the instrument or any of its parts or accessories appear to be damaged or defective, contact the dealer or Inventis service.



Keep the packaging materials, in case the instrument should need to be sent to the dealer or to Inventis for any reason.

2.2 PRECAUTIONS

Installation of the Timpani tympanometer is easy but needs to be done carefully: incorrect installation could lead to safety issues while using the system.

Like any other electrical or electronic device, the tympanometer will emit electromagnetic waves. Whilst the level of emissions is guaranteed to be within statutory limits, other electronic devices operating in the immediate vicinity could be affected, if particularly sensitive to electromagnetic interference. If this should occur (interference is verifiable by switching off the instrument and then switching on again) it may be possible to remedy the problem by adopting one or more of the following solutions:

- change the orientation and/or the position of the device affected by interference
- distance the device from the tympanometer
- plug the affected device into a power socket on a circuit other than the circuit to which the tympanometer is connected
- consult the manufacturer or a service center for assistance

2.3 CONNECTIONS

The Timpani can be connected either to a PC for recharging and transferring test data, or to the power adapter supplied. Use only the USB cable supplied with the product. Moreover, it is possible to recharge Timpani and transfer tests through the optional Docking station (see Chapter 4).

As long as it is connected to a power source, the appliance will remain active in recharge or trickle charge mode.



Use only the medical grade power adapter supplied with Timpani, certified according to the IEC 60601-1 standard



Make certain that the electrical power supply and ground connections comply with the applicable standards for electro-medical devices. Risk of electric shock

Plug all transducers and parts into the respective sockets as indicated in the following table:

Connector	Detachable part
PHONES	Air conduction headphones
PATIENT RESP.	Patient response button
	USB cable for connection to personal computer


2.4 POWER-UP, POWER-DOWN AND MAIN SCREEN

Turn on the instrument by pressing and holding the on/off button; it can be switched off by pressing and holding the same button.



At the startup of the instrument, a pressure initialization is performed: to complete initialization successfully, hold the tympanometer still with the probe positioned in free space.

A few seconds after power-up, the display of the instrument will show the main screen, illustrated below:



EN

Swipe your finger to the left on the display for the settings and to manage the patient memory.



2.5 CONNECTION TO PC

The Timpani tympanometer can be interfaced with a personal computer after installing the Inventis Maestro software.¹ Connect the Timpani tympanometer to a USB port of the computer using the cable provided (USB cable with A / mini B connectors) or plug it into the docking station (which in turn is connected to the computer via the USB cable provided).

¹ Maestro Spring 2019 version (1.09.0) or later



Use the supplied cable to connect the Timpani tympanometer to one of the USB ports of the computer

After a few seconds, the connected device will be recognized by the operating system.

Refer to Maestro user manual for more details about the software.

Chapter 3

Controls and exams

EN


3.1 CONTROLS FOR DIRECT ACCESS TO EXAMS

Button	Function
	Perform Tympanometry
	Perform Tympanometry and then the Acoustic Reflex test
	Perform Audiometry
	Perform Automatic Audiometry
	Save the current session to the patient memory
	Delete the current session
	Print the current session to the thermal printer (if configured and available)

3.1.1 Function buttons

Button	Function
	Return to the main screen
	Select the ear to test (in the example the right ear was chosen)
	Delete the current test

3.2 TYMPANOMETRY




Button	Function
226 Hz ☐ 1000 Hz	Probe tone selection
START TYMPANOMETRY	Force the tympanometry to start.

If the probe is detected as being correctly inserted in the ear of the patient, with a compliance measurement that is stable and within the specified range, the test will begin automatically; alternatively, the start of the test can be forced.

In the event of pressure losses being detected, the instrument will repeat the pressurized scan up to three times before generating an alert to indicate that there is a problem. Should it prove impossible to conduct the test as the result of pressure loss, the problem may be overcome by replacing the eartube with another of different size and/or adjusting the position and angle of the probe in the auditory canal.


3.2.1 Results

Indication	Information
	Tympanometry. The unit of measurement of the vertical axis (acoustic admittance) is expressed in: - ml (equivalent volume of air), with 226 Hz probe tone - mmho, with 1 kHz probe tone The horizontal axis indicates meatus pressure in relation to ambient pressure, expressed in daPa.
Press.: -25 daPa	Pressure at the tympanogram peak.
ECV: 0.71 ml	Ear Canal Volume: compliance in ml measured at the maximum value of the pressure range selected for the scan. This value is also called the "equivalent volume".
Comp.: 0.52 ml Comp.: 0.38 mmho	Compliance: Amplitude of the tympanogram peak measured against ECV. The unit of measurement reflects that of the tympanogram.
Grad.: 96 daPa	Gradient of the tympanogram: width of the tympanogram at 50% of the compliance value (only for the 226 Hz probe tone)

If the test has not been able to determine one or more of the above values, a double dash "--" will appear in place of the number.

3.3 ACOUSTIC REFLEX


At the end of the Tympanometry test, the instrument automatically performs the reflex test at the peak tympanometry pressure. If the Tympanometry test has not been performed, the Acoustic Reflex test is performed at atmospheric pressure.



Button	Function
	Start the Acoustic Reflex test.

3.3.1 Results

Button	Function
	Display the reflex plot at the specific frequency. Button information: ✓: reflex threshold detected, ✗: reflex threshold not detected

Indication	Information
	Reflex performance. The green segment indicates the duration of the stimulus.
0.5 kHz 75 dBHL 140 daPa	Reflex frequency, reflex level, pressure at which the reflex search was performed.

3.4 AUDIOMETRY

3.4.1 Common indicators

The following indicators are common to both manual and automatic tone audiometry tests.

Indication	Information
	Patient response button not pressed
	Patient response button pressed
	Headphones
	Headphones with active stimulus
	Insert earphones

EN


3.4.2 Manual audiometry




Button	Function
	Save the current point
   	Increase/decrease Level and Frequency
	Send the stimulus

Indication	Information
40 dB HL	Stimulus level
1 kHz	Stimulus frequency. If valid data are associated with the current frequency, this is highlighted in the color representing the side.

3.4.3 Auto threshold (automatic audiometry)




Button	Function
	Start the test
	Suspend the test
	Stop the test




The automatic test includes an initial familiarization phase, to train the patient on the threshold determination procedure, and this is followed by the actual test at all enabled frequencies. Stimulation takes place for a duration of 1.7 seconds and is then followed by a pause of random duration lasting between 1.7 and 2.5 seconds.

3.4.4 Results




3.5 PATIENT MANAGEMENT

The Patient management screen allows adding (or modifying) patients and reviewing stored exams. The first time the Patient Management screen is accessed, Timpani asks for a PIN to prevent data access from unwanted accesses. You can choose either to enter the PIN or disable data protection.



Message prompt at the first Patient Management screen access




Patient Management screen



It is recommended to enable the PIN for patient data protection

Button	Function
	Return to the main screen
	Create new patient
	Delete all saved patients
	Return to the list of patients
	Side of the test saved
	Delete the current patient
	Print the tests of the current patient

3.6 SETTINGS



EN

Button	Function
	Go back to the main screen
	Access the info screen, with serial number of the device, calibrated transducers, firmware version and other information for service

3.6.1 User-settable parameters

General configuration parameters for the device are listed below. The availability of some configuration parameters depends on the licenses active on the device.

- **Language:** Interface language. Default value: English (may vary based on destination)
- **Date and time:** Access the menu to adjust date and time and its format.
- **Automatic test repetition:** Enables/disables the possibility of repeating the test, re-inserting the probe in the ear (without having to manually delete the previous test). Default value: disabled.
- **Tympanometry:** Access the menu for tympanometry settings. Allows selecting the pressure range used to conduct the test: Standard [-400; +200] daPa or Reduced [-300; +100] daPa. Default value: Standard.
- **Acoustic reflex:** Access the menu for the Acoustic Reflex test settings:

- Frequency selection: the available stimulus frequencies can be selected individually: 0.5 kHz, 1 kHz, 2 kHz, 4 kHz. Default value: all frequencies enabled.
 - Test mode: sets the mode in which the test is conducted, namely fixed intensity or threshold search. Default value: threshold search.
 - Test configuration:
 - Select the level of the stimulus in dB HL (in fixed intensity mode). Default value: 90 dB HL
 - Selection of the initial and final level, selection of the variation rate in steps of 5 or 10 dB. Default value: 75-95 dB, 5 dB steps.
 - Reflex sensitivity: sensitivity with which identification of the reflex (variation of compliance) occurs, normal (0.04 ml) or strong (0.06 ml). Default value: normal (0.04 ml).
 - Data polarity: sets the method for representing data in the chart: negative polarity (decrease in compliance caused by the reflex is represented with a deflection of the reflex curve), positive polarity (decrease in compliance caused by the reflex is represented with an increase in the curve). Default value: negative.
- **Pure Tone Audiometry**: Access the menu for the Audiometric test settings:
- Frequency selection: selection of stimulus frequencies in the range 125 Hz – 8 kHz. The value 1kHz cannot be unselected. Default value: all frequencies enabled.
 - Stimulus mode: sets the stimulation mode, allowing the user to choose between continuous stimulation or pulsed 1 Hz stimulation. Default value: continuous.
 - Default intensity: sets the intensity of the stimulus from which to start the frequency change for the manual examination. Default value: 40 dB HL.
 - Automatic frequency skip function: enables / disables automatic switching to the next frequency after a value is saved. Default value: disabled.
 - Switch mode: allows the interrupter button to be used as a button (stimulation active as long as the key is pressed) or switch (stimulation is activated when the key is pressed and deactivated when it is pressed again). Default value: pulsating.

- AC transducer: sets the type of transducer for air conduction, allowing the user to choose between supra-aural and insert earphones. Default value: supra-aural headphones.
- **Data security**: Access the menu to modify the PIN and enable/disable it
- **Display brightness**: Set the display brightness between 20% and 100%. Default: 80%.
- **Printer**: Access to the print options menu:
 - Print patient details: allows enabling the printing of the patient's personal details. Default value: enabled.
 - Print reflex graphs: allows printing the acoustic reflexes in graph form. Default value: disabled.

3.6.2 License menu

Access the menu to enable additional licenses.

Chapter 4

Docking station

The docking station, available upon request, allows the user to store away the Timpani after use, recharging the device and transferring the data to the computer.¹

EN



Place the Timpani stably on the docking station to ensure correct communication.

4.1 INSTALLATION AND GENERAL PRECAUTIONS

For a complete list of applicable precautions, please refer to the *Installation and general precautions* paragraph in Chapter 1. In addition to the listed warnings, please be aware that the following applies for Timpani Docking station:

¹ Requires Maestro Summer 2020 version (1.10.0) or later



The Timpani Docking station is an accessory of a medical device: if connected to a computer (or any external device) located within the "patient area" (as defined in IEC 60601-1), this likewise must be a medical device, or protected by an isolating transformer, in order to ensure that the combination of computer (external device) + tympanometer is in compliance with IEC 60601-1.



The Timpani Docking station must be installed and operated taking account of the information regarding electromagnetic compatibility (EMC) provided at the end of this manual.

4.2 CONFORMITY

The Timpani Docking station is a class I accessory, according to Annex VIII of Medical Device Regulation (MDR) 2017/745/EU.

4.3 SYMBOLS ON LABELS



Warning: the use of this device requires certain precautions.

To ensure safe use, consult the accompanying documentation.



Follow the instructions for use.

Device serial number. The number is made up of 13 alphanumeric characters indicating the model, series, year of manufacture and serial number:

- characters 1-5: Inventis product code
- characters 6-7: year of manufacture ("20" denotes 2020)
- characters 8-13: pro serial number



Catalog code



Name and address of manufacturer.



Product conforms to European Community Medical Device Regulation (MDR) 2017/745/EU. Class I accessory.

Rx Only

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.



The product is subject to the requirements of Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). In the event of this product being sold and/or scrapped, it must not be disposed of as ordinary household or industrial waste, but collected separately.

(01)08054187380761
(21)IM1PB21000000



UDI code

4.4 CONNECTIONS

EN

Plug the docking station into a power socket using the power adapter supplied by Inventis and to the computer using the cable supplied (USB cable with A / mini B connectors). The two USB ports on the back of the instrument are interchangeable; both can communicate with the computer and also power the device. It is not necessary for both ports to be connected.



Use only the medical grade power adapter supplied by Inventis, certified according to the IEC 60601-1 standard



Make certain that the electrical power supply and ground connections comply with the applicable standards for electro-medical devices. Risk of electric shock

Chapter 5

Maintenance

EN

The Timpani tympanometer does not require any special periodic maintenance other than calibration and normal cleaning, both of which are described in this chapter. The instrument must be switched off before commencing any kind of cleaning operation.

The performance and safety of the instrument will be assured as long as the recommendations for care and maintenance indicated here are correctly observed.



Apart from replacing the battery, the inspection and servicing of internal components must be left entirely to technicians approved by Inventis.



Transducers are manufactured utilizing ultra-fragile diaphragms that could be damaged in the event of impact. Handle with care during maintenance operations.

5.1 PERIODIC CHECKS



The procedure described under this heading must be carried out when the instrument is used for the first time each day.



The tests must be conducted with the instrument positioned for normal use.

Before switching on the instrument, make certain that there is no sign of damage visible on the equipment, including the accessories and the external power adapter; make a visual inspection of the power cable and connectors to verify the integrity of the insulation, and ensure that they are not subject to any kind of mechanical loading or stress that could cause damage; check that all parts and cables are properly connected.

Verify correct operation of the probe and check the pressure. To this end, carry out the following sequence of steps:

- Fit a new earpiece to the probe;
- Select the tympanometry test;

- Check that the probe is identified as being open;
- Start the test in manual mode and check that the internal pump carries out pressurization cycles until the point, after a few seconds, that a pressure loss alert is generated, then press OK;
- Occlude the probe by placing a finger over it;
- Check that the probe is identified as being closed;
- Start the test manually and make certain it is carried out in a few seconds, showing a tympanometry graph that appears blank, with $ECV < 0.2$ ml;
- If calibration cavities of 0.5 ml, 2.0 ml and 5.0 ml capacity are available, run a tympanometry test on each one and check that the ECV value obtained is compatible in each case;
- If the optional acoustic reflex license is installed:
 - o Select the reflex test, holding the probe open;
 - o Check that the probe is identified as being open;
 - o Start the test manually and check that the cycle is performed as expected, given the reflex configuration selected; bringing the tip of the probe closer to the ear, in noiseless conditions, the stimuli should be audible.



If any part or transducer has any malfunction, consult Chapter 6 “Troubleshooting”.

Also, check that the calibration interval has not expired: the date is shown on the info screen accessible from the setup menu.



Calibration must be entrusted to technicians approved by Inventis. The operation should be performed at least once every 12 months, and whenever a transducer is replaced.

5.2 MAINTENANCE OF TRANSDUCERS



Do not use liquids or sprays to clean the tympanometer.

Do not allow dust to collect on the transducers. Also:

- The cushions of the headphones are made of biocompatible material but are not sterile. To prevent the spread of infections and to ensure their biocompatibility, they must be sanitized before being used on a new patient using wipes dampened with denatured alcohol or a microfiber cloth dampened with denatured alcohol.

- The ear tips of the probe and of the insert earphones are intended to be inserted in the patient's ear canal. They are made of biocompatible material and disposable: use once only and dispose of in accordance with current health and safety regulations.



Earpieces are not sterile. Reusing unsterilized earpieces can cause ear infections.



The headphone cushions can be repeatedly cleaned as described in the “Maintenance of Transducers” paragraph. In the event of any malfunctioning after any cleaning operation, contact an Inventis service technician.



Though the headphone cushions can be repeatedly cleaned, always check that their characteristics and integrity are maintained. To do so, it is sufficient to perform the tests described in the “Periodic checks” paragraph. As soon as any failure is encountered, contact an Inventis service technician to verify whether your transducer needs to be replaced.




To avoid damaging the headphones, do not press them against a flat surface. This can create a vacuum and damage the transducer (suction cup effect).

EN

5.3 CLEANING THE PROBE

To guarantee accurate compliance measurements, the three channels incorporated into the probe must be kept properly clean. In effect, these three channels carry the compliance measurement system, the stimulus speaker, and the pressurization system.


As illustrated below, the probe comprises a main body rigidly attached to the instrument, a tip (to which the earpiece is fitted) and a screw collar, which keeps the tip of the probe firmly secured to the body.




The procedure for cleaning the probe will now be described.

Remove the earpiece then loosen the screw collar so that the tip of the probe can be separated from the main body.

Use thin nylon threads to clean the three channels in the tip of the probe. Insert the thread into each of the channels in turn, from the base, and push through until it can be pulled out from the top.



Having cleaned the channels thoroughly, the tip of the probe must be refitted. Offer the tip of the probe to the main body, being careful to align the guide pin A presented by the body with the hole B located in the tip, as illustrated in the figure below. Retighten the screw collar.



Clean the outer surface of the probe using a lint-free soft cloth moistened with water and mild detergent; if the probe is to be sanitized, moisten the cloth with a 3% solution of hydrogen peroxide.



Do not immerse the probe or part of it in any kind of liquid

In the event of the probe being broken or malfunctioning, contact an Inventis service technician. The replacement of the probe must be entrusted exclusively to Inventis or to a service technician authorized by Inventis. If the probe is replaced, it must be calibrated before being used with the instrument.

5.4 CLEANING THE INSTRUMENT

Clean the instrument using a lint-free soft cloth moistened with water and mild detergent; if it is to be sanitized, moisten the cloth with a 3% solution of hydrogen peroxide.

5.5 REPLACING THE BATTERY

If the device fails to operate as long as expected despite a full charge, the battery may be damaged or dead.



The incorrect battery replacement can potentially be dangerous. Take particular care to replace it.

EN

Purchase a new battery from an Inventis-approved dealer; proceed to replace the existing battery from the instrument as described below:

- Switch off the instrument and disconnect it from the USB cable;
- Position it face down (display directed downwards) on a soft surface;
- Undo the screw retaining the flap of the battery compartment;
- Remove the battery. Separate the connectors without tugging; ease apart using tweezers;
- Fit the connector of the new battery;
- Position the lead inside the compartment below the screw and locate the new battery in its housing, then close the flap and secure with the retaining screw.

Recharge the instrument completely before use.



All parts mentioned in the manual are designed especially for use with this instrument. Only parts supplied by Inventis should be connected to the device.

5.6 REPAIRS AND TECHNICAL ASSISTANCE

Before contacting the service department, make certain that all the possible solutions in Chapter 6 *Troubleshooting* have been tried.

Parts that are to be returned to the manufacturer must be cleaned and sanitized, following the directions in this manual. Transducers must be shipped in a closed and sealed transparent bag.

Should the instrument need to be sent to the service department or returned to the dealer, it is important that the original packing be used, and that all accessories and transducers are enclosed.

Chapter 6

Troubleshooting

EN

Problem	Possible cause	Solution
No probe tone	Holes of probe tip blocked	Unscrew the probe tip and clean inside
No pressurization message: “Pressure loss”	Probe not correctly positioned and tightened	Check that the tip of the probe is properly tightened
	Probe not inserted tightly in ear / unsuitable earpiece	Change the earpiece and reinsert the probe Change the angle of the probe in the ear
Compliance measurements affected by noise	Probe not properly positioned	Reposition the probe so as to minimize vibrations
	Holes of probe tip blocked	Unscrew the probe tip and clean inside
No signal from a transducer	Transducer not connected properly	Check that the transducer is connected properly
	Transducer damaged	Contact Inventis service department or dealer
Unable to establish a direct connection between the PC and the Timpani or docking station	Problems with USB connection	Check the USB connection between instrument and computer
	USB cable damaged	Change the USB cable (USB A – mini B standard)
Unable to transfer the data to the PC via the docking station	Instrument not positioned correctly in the docking station	Check the positioning of the instrument Check that the contacts are clean Check connections

Problem	Possible cause	Solution
The instrument does not switch on	Battery flat	Connect the instrument to a power source and switch on the device
The display remains blank (LED on)	Instrument in stand-by	Touch the screen or press the power button
	Display damaged	Contact Inventis service department or dealer
Battery does not recharge	USB cable damaged	Change the USB cable (USB A – mini B standard)
	Adapter damaged	Contact Inventis service department or dealer
	Instrument not positioned correctly in the docking station	Check the positioning of the instrument Check that the contacts are clean Check connections
	Battery damaged	Replace the battery - Contact Inventis service department or dealer
A test cannot be accessed	Optional test not activated	Contact the Inventis technical service department to obtain the license, supplying the serial number of the device
<i>message:</i> “Hardware error”	Non-fatal internal error	Press OK to continue; if the problem persists, contact the Inventis service department
<i>message:</i> “Serious error”	Fatal internal error	Restart the instrument; if the problem persists, contact the Inventis service department