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Blood Glucose Test Strips

Use only with GIMACARE Multi-Functional Monitoring System.

Warnings

- For in vitro diagnostic use (for use outside of the body only).
- For single use only.
- Healthcare professionals and other users should handle everything that comes into contact with human blood carefully to prevent transmitting infectious diseases, including sanitized objects.
- Please read this sheet and your GIMACARE Multi-Functional Monitoring System Owner's Manual before you use this test strip. Use only GIMACARE Blood Glucose Test Strips with GIMACARE Multi-Functional Monitoring System to obtain accurate results, and be covered by the manufacturer's warranty.
- Results may be inaccurate when testing on patients with abnormally low blood pressure or those who are in shock.
- For patients with impaired peripheral circulation, collection of capillary blood from the approved sample sites is not advised as the results might not be a true reflection of the physiological blood glucose level. It may apply under the following circumstances: severe dehydration as a result of diabetic ketoacidosis or due to stress hyperglycemia, hyperosmolar non-ketotic coma, shock, decompensated heart failure NYHA Class IV or peripheral arterial occlusive disease.
- Keep test strips and lancets away from small children. If swallowed, consult a doctor immediately for advice.

Intended Use

- GIMACARE Blood Glucose Test Strips, when used together with GIMACARE Multi-Functional Monitoring System, allow your blood glucose levels to be measured by yourself at home or by healthcare professionals. It is intended to be used for the quantitative assessment of glucose (sugar) in fresh capillary whole blood samples from fingertips, and from venous, arterial and neonatal whole blood sample. Professionals may use test strips to test capillary, venous, arterial and neonatal blood sample; home use is limited to capillary whole blood testing. This system is not intended for use in the diagnosis or screening of diabetes mellitus.

Test Principle

- Your system measures the amount of glucose in whole blood. The test is based on the measurement of electrical current generated by the reaction of the blood glucose with the reagent of the test strip. The meter measures the current, calculates the glucose in the blood, and displays the result on the screen. The strength of the current produced by the reaction depends on the amount of blood glucose in the blood sample.

Limitations

- Lipemic Effects:** Blood triglycerides up to 3000 mg/dL (33.9 mmol/L) do not affect the results significantly, but may affect results at higher levels.
- Xylose:** Do not test blood glucose during or soon after a xylose absorption test. Xylose in the blood can give false elevated results.
- Metabolites:** Dopamine, L-Dopa, methyldopa, uric acid, ascorbic acid, and acetaminophen at normal blood concentration do not significantly affect blood glucose test results.
- There are no significant interference in the presence of galactose, maltose, or fructose observed in blood glucose tests.
- Altitude Effects:** Altitudes up to 10,742 feet (3,275m) do not affect test results.
- Use only heparin for anticoagulation of fresh capillary or venous whole blood.
- Hematocrit:** The hematocrit level is limited to between 0% and 70%. Please ask your healthcare professional if you do not know your hematocrit level.
- Neonatal Use:** This test strip can be used for the testing of newborns.

The following compounds when determined to be in excess of their limitation and tested with the GIMACARE Blood Glucose Monitoring System may produce elevated glucose results: Summary of substances and concentrations in excess of limitation with interference.

Substance	Limiting Concentration (mg/dL)	Therapeutic / Physiologic Concentration Range (or Upper Limit) (mg/dL)
Acetaminophen (Paracetamol)	> 20	0.45 - 3
Ascorbic Acid	> 5	> 5
Pralidoxime iodide	> 5	~ 10 (IV Dose 500 mg)
Uric Acid	> 10	2 - 8

Storage and Handling

- Do **NOT** use the test strips if they have expired.
- Store the test strips in a cool, dry place between 2°C and 30°C (35.6°F and 86°F) and between 10% and 90% relative humidity.
- Store the test strips in their original vial **ONLY**. Do **NOT** transfer them to a new vial or any other containers. (For strip vial only)
- Use each test strip immediately after taking it out of the vial or individual foil packet. Close the vial immediately after taking out a strip. (For strip vial only)
- Keep the vial closed at all times. (For strip vial only)
- Keep the test strips away from direct sunlight. Do **NOT** store the test strips in high humidity.
- Do **NOT** touch the test strips with wet hands.
- Keep the vial clean, dry, and alter the test strip

Chemical Components

- Glucose dehydrogenase (E. coli) 8%
- Electron shuttle 55%
- Enzyme protector 8%
- Stabilizing agents 29%
- Other inactive ingredients 29%

Strip Appearance

- Absorbent Hole**
Apply a drop of blood here. The blood will automatically be absorbed.
- Confirmation Window**
This is where you confirm if enough blood has been drawn into the absorbent hole of the strip.
- Test Strip Handle**
Use this to hold the strip to insert the test strip into the slot.
- Contact Bars**
Insert this end of the test strip into the meter. Push it in firmly until it will go no further.

- ATTENTION:** The front side of the test strip should face up when inserting the test strip. Test results might be wrong if the contact bar is not fully inserted into the test slot.

PLEASE WASH AWAY YOUR HANDS BEFORE PERFORMING ANY TESTING

- STEP 1**
Insert the test strip fully into the slot of the meter until it will go no further. When the strip is fully inserted, the meter will do several self-checks.

- Collect a blood sample with the test strip. A sufficient quantity of blood is required for the test to provide accurate results. Touch the blood drop with the absorbent hole of the test strip, and wait until the blood is completely absorbed. The meter will start counting down. Do **NOT** apply a smeared blood sample.

- STEP 3**
After a few seconds, the meter will display your test result. The last reading will be automatically saved in the meter. Turn it off by removing the test strip and throw away the used test strip.

- NOTE:** Please refer to your Owner's Manual for more information. The used lancet and test strips are potentially biohazardous. Please dispose of them carefully according to your local regulations.

- Reading Your Result**
The blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining blood glucose levels close to normal can reduce the risk of diabetes complications by up to 50%. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes. Your blood glucose readings deliver plasma equivalent results and are displayed either in milligrams of glucose per deciliter of blood (mg/dL) or in millimoles of glucose per liter of blood (mmol/L).

- Reference values**

Time of day	Normal plasma glucose range for people without diabetes (mg/dL or mmol/L)
Fasting ^{1,2} and before meal	Less than 100 mg/dL (5,6 mmol/L)
2 hours after meal	Less than 140 mg/dL (7,8 mmol/L)

- ¹: American Diabetes Association. Classification and Diagnosis of Diabetes: Standards of Medical Care in Diabetes-2022 Jan; 45(Supplement 1): S17-S38. <https://doi.org/10.2337/dc22-5002>
- ²: Fasting is defined as no caloric intake for at least 8 hours.

- Please consult your doctor to determine a target range that works best for you. Questionable or inconsistent results**

- If your test results are unusual or inconsistent with how you are feeling:
 - Make sure you confirm the window of the test strip is completely filled with blood.
 - Check the expiration date of the test strips.
 - Check the performance of your meter and test strips with the control solutions.
- Unusually high or low blood glucose levels may be symptoms of a serious medical condition. If most of your test results are unusually high or low, please contact your healthcare professional.

Quality Control Testing

- Our control solutions contain a known amount of glucose that can react with test strips. If you are concerned about the meter or test strips are not working properly, you can check the performance of the meter, test strips and your technique by comparing the control solution results with the range printed on the label of test strip vial or on the test strip package. Please refer to the Owner's Manual for complete testing instructions.
- Δ The reference range of the control solutions may vary with each new vial or package of test strips. Make sure you check the range on the label of your current vial or on the current package.

Additional Information

- Always wear gloves and follow your facility's biohazard control policy and procedures when performing tests involving patient blood samples. Use fresh whole blood samples only. Professionals may use test strips to test capillary, venous, arterial and neonatal blood sample. Sample Size: 0.5 µL Reaction Time: 5 seconds System Measurement Range: 10 to 800 mg/dL (0.56 to 44.4 mmol/L) Hematocrit Range: 0% to 70%

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