CATETERE ESTERNO PER UOMO + FASCIA HYDROCOLLOID
KIT MEN'S EXTERNAL CATHETER + HYDROCOLLOID STRIP
KIT CATHÉTER EXTÉRIEUR POUR HOMME + RUBAN HYDROCOLLOÏDAL
KIT CATÉTER EXTERNO HOMBRE + TIRA HIDROCOLOIDAL
КОМПЛЕКТ МЪЖКИ ВЪНШЕН КАТЕТЪР + ХИДРОКОЛОИДНА ЛЕНТА
SOUPIRAVA PÁNSKÝ EXTERNÍ KATETR + HYDROKOLOIDNÍ PÁSEK
KIT EXTERNER KATHETER FÜR HERREN + HYDROKOLLOID-STREIFEN
KIT ΑΝΔΡΙΚΟΣ ΕΞΩΤΕΡΙΚΟΣ ΚΑΘΕΤΗΡΑΣ + ΥΔΡΟΚΟΛΛΟΕΙΔΗ ΤΑΙΝΙΑ
KIT MIESTEN ULKOINEN KATETRI + HYDROKOLLOIDILUSKA
KIT FÉRFI KÜLSŐ KATÉTER + HIDROKOLLOID CSÍK
KIT EXTERNE KATHETER VOOR HEREN + HYDROCOLLOÏDE STRIP
ZESTAW CEWNIK ZEWNĘTRZNY MĘSKI + PASEK HYDROKOLOIDOWY
KIT CATETER EXTERNO MASCULINO + FAIXA HIDROCOLOIDE
KIT CATETER EXTERN BÄRBAȚI + BANĂ HIDROCOLOIDĂ
KIT EXTERN KATETER FÖR HERR + HYDROCOLLOID STRIP

Gima 28724 - Gima 28725

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ENGLISH

INSTRUCTIONS FOR USE

External male catheter with Hydrocolloid Strips is intended for collection of urine in males suffering from incontinence and provides secure and gentle fixation for external catheters. MEC and related accessories are compliant with Regulation (EU) 2017/745.

 Please read the instruction before using the product

Description:

The non-invasive non-sterile, Male External Catheter, are intended for use to divert urine in male patients with urinary incontinence. This catheter is only for external and single use. This product contains Latex.

Intended Use:

This catheter is only for external and single use. Diverts urine from penis to extension tubing and urinary collection bags.

Indications for use:

Male external Catheters (MEC) is intended to be designed to safely conduct urine through funnel tube, keeping fluid pathway open.

Hydrocolloid Strips are ideal for replacing regular adhesive tapes for long-term attachment, without causing skin irritation- latex free.

Conforms to the skin, absorbs body fluids, Provides secure adherence, Easy to apply and remove

How To apply the Male External Catheter:

1. Clean the Area before applying the Hydrocolloid strip.
2. The hair around the penis shaft may be trimmed before application.

3. Peel off one of the two paper cover strips.
4. Wind the sticky side of the strip around the penis shaft in a spiral motion starting from the back of the base and working upwards.
5. Do not make the Hydrocolloid strip too tight to prevent circulation problems.
6. The strip must overlap.
7. Before applying the MEC remove the other paper protection strip.
8. Apply the MEC leaving about 2cm between the top of the glans and the end of the MEC.
9. Attach the urine bag to the catheter.

Removal of the catheter:

1. Disconnect the catheter from the urine bag.
2. Gently remove the adhesive strip
3. Roll the catheter along length of the penis and lift gently off.

Warning:

- Male External Catheter is made of medical grade natural latex. Do not use the product if you are sensitive to latex. In case of suspected latex sensitivity stop using the catheter and refer to your doctor.
- Do not reuse the catheter.
- The Catheter must be Changed every 12 to 24 hours.

Use limitations:

Check Whether the package is in good condition before use, and confirm the package mark, production batch number or production date, and use within the validity period.

Recommendation:

- Do not apply the adhesive strip too tightly, to avoid circulation restriction.

Storage conditions

• Before use, please store the device in the original packaging away from contaminated areas, dust, sunlight, excessive moisture and temperature and damaging chemicals. The product shall be stored in a ventilated, dry room away from heat source and free of corrosive gas. Always inspect product and conduct a user seal check before use as specified in the user instructions.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

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