GIMA

PROFESSIONAL MEDICAL PRODUCTS

ULTRA

PREGNANCY TEST

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TEST DE GROSSESSE

TEST DE EMBARAZO

TESTE DE GRAVIDEZ

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GIMA 29103

Hangzhou AllTest Biotech Co., Ltd
#550, Yin Hai Street, Hangzhou Economic &
Technological Development Area, Hangzhou,
310018 P.R. China - Made in China

FHC-U103H

MedNet EC-REP GmbH Borkstrasse 10,
48163 Muenster, Germany

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China

INTENDED USE

The hCG Pregnancy Enhanced Sensitivity Rapid Test Midstream is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin in urine to aid in the early detection of pregnancy.

PRINCIPLE

HCG Pregnancy Enhanced Sensitivity Rapid Test Midstream is a rapid, one-step lateral flow immunoassay in midstream format for the qualitative detection of human chorionic gonadotropin (hCG) in urine to aid in the early detection of pregnancy. The test utilizes a combination of antibodies including a monoclonal hCG antibody to selectively detect elevated levels of hCG. The assay is conducted by adding urine to the hydrophilic stick and obtaining the result from the colored lines.

REAGENTS

The test contains anti-hCG particles and anti-hCG coated on the membrane.

PRECAUTIONS

Please read all the information in this package insert before performing the test.

- Do not use after the expiration date printed on the foil pouch.
- Store in a dry place at 2-30°C or 35.6-86°F. Do not freeze.
- Do not use if pouch is torn or damaged.
- Keep out of the reach of children.
- For in vitro diagnostic use. Not to be taken internally.
- Do not open the test midstream foil pouch until you are ready to start the test.
- The used test midstream should be discarded according to local regulations.

STORAGE AND STABILITY

The kit can be stored at room temperature or refrigerated (2-30 °C). Do not open the pouch until ready for use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

The urine specimen should be collected in a clean and dry container. A first morning urine specimen is preferred since it generally contains the highest concentration of hCG; however, urine specimens collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

If the urine specimen cannot be detected immediately, it should be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS PROVIDED	
• Test Midstream	• Package Insert
MATERIALS REQUIRED BUT NOT PROVIDED	
• Timer	• Specimen Collection Container


INSTRUCTIONS

Allow the test, urine specimen to reach room temperature (15-30°C) prior to testing.

1. Remove the midstream from the foil pouch and test it immediately in one hour.
2. Take down the cap of the midstream, hold the midstream so as to place the absorbent tip in the urine stream or place the absorbent tip (≥2/3) in urine in a clean cup for **at least 15 seconds**.

NOTE: Do not urinate on the Result Window.

3. Cover the cap on the testing midstream, then lay down the midstream on a clean and stable desk, start the timer immediately.
4. **Read the result at 3 minutes;** don't interpret the result after 10 minutes.

READING THE RESULTS

POSITIVE: Two colored lines appear. One line should be in the control line region (C) and another line should be in the test line region (T). One line may be lighter than the other; they do not have to match. This means that you are probably pregnant.

NEGATIVE: One colored line appears in the control line region (C). No line appears in the test line region (T). This means that you are probably not pregnant.

INVALID: The result is invalid if no colored line appears in the control line region (C), even if a line appears in the test line region (T). You should repeat the test with a new test midstream.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

LIMITATIONS

This test might give wrong results. Talk to your doctor before making medical decisions.

1. Some drugs with hCG (such as Pregnyl, Profasi, Pergonal, APL) can cause false positives. But alcohol, birth control pills, painkillers, antibiotics, and non-hCG hormone therapies shouldn't affect it.
2. If your urine is very dilute, it might not have enough hCG. If you still think you're pregnant, test again with first morning urine after 48 hours.
3. Low hCG levels (less than 50mIU/mL) can happen right after implantation. If the test is weakly positive, retest with first morning urine after 48 hours.

4. False positives can happen due to conditions like trophoblastic disease or cancers (testicular, prostate, breast, lung), cause elevated levels of hCG.^{1,2} Don't use hCG to diagnose pregnancy unless other causes are ruled out.
5. False negatives can happen if hCG levels are too low. If you still suspect pregnancy, test again with first morning urine after 48 hours. If it's still negative, see a doctor.
6. This test only suggests pregnancy. A doctor must confirm it after checking all findings.

EXTRA INFORMATIONS

How it works: It detects hCG, a hormone produced during pregnancy, in your urine.

When to test: The test can detect hCG 6 days before your missed period (5 days before day of the expected period). Use first morning urine for the best results, as it has the highest hCG concentration.

Accuracy: Over 99% accurate compared to another urine hCG test in a clinical evaluation.

Sensitivity: Detects hCG at 10 mIU/mL or higher.

If pregnant: See your doctor to confirm and discuss next steps.

If not pregnant: If your period doesn't start within a week, retest. If still negative and no period, see a doctor.

Disposal: Dispose of the used test according to local regulations.

BIBLIOGRAPHY

1. Dawood MY, BB Saxena, R Landesman Human chorionic gonadotropin and its subunits inhydatidiform mole and choriocarcinoma, Obstet. Gynecol. 1977; 50(2): 172-181

2. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross "Ectopicproduction of human chorionic gonadotropin by neoplasms", Ann. Intern Med. 1973; 78(1): 39-45

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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