PROFESSIONAL MEDICAL PRODUCTS

SPECULUM VAGINALE STERILE MONOUSO
STERILE VAGINAL SPECULUM FOR SINGLE USE
SPÉCULUM VAGINAL STÉRILE À USAGE UNIQUE
ESPÉCULO VAGINAL ESTÉRIL DE UN SOLO USO
ESPÉCULO VAGINAL ESTERILIZADO PARA UTILIZAÇÃO ÚNICA
VAGINALSPEKULUM ZUM EINMALIGEN GEBRAUCH
ΑΠΟΣΤΕΙΡΩΜΕΝΟΥ ΚΟΛΠΟΔΙΑΣΤΟΛΕΑ ΜΙΑΣ ΧΡΗΣΗΣ
STERYLNEGO WZIERNIKA DOPOCHWOWEGO
VAGINÁLNÍ SPEKULUM NA JEDNORÁZOVÉ POUŽITÍ
VAGINALSPEKULUM FÖR ENGÅNGSBRUK
STERIILIN KERTAKÄYTTÖISEN EMÄTTIMEN SPEKULAN
STERILNI VAGINALNI SPEKULUM ZA ENKRATNO UPORABO
STERILNÉHO JEDNORAZOVÉHO VAGINÁLNEHO ZRKADLA
SPECULUM VAGINAL STERIL DE UNICĂ FOLOSINȚĂ
STERIEL VAGINAAL SPECULUM VOOR EENMALIG GEBRUIK
STERILINI JEDNOKRATNI VAGINALNI SPEKULUM
STERIL ELDOBHATÓ HÜVELYI TÜKÖR
STERILT VAGINALT SPEKULUM TIL ENGANGSBRUG
СТЕРИЛЕН ВАГИНАЛЕН СПЕКУЛУМ ЗА ЕДНОКРАТНА УПОТРЕБА
VIENKARTINĖ STERILAUS MAKŠTIES SKĖTIKLIO
STERILU MAKSTS LĀPSTIŅU VIENREIZĒJAI LIETOŠANAI
ÜHEKORDSE STERIILSE SPEEKULUMI

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ENGLISH

Product Name: Sterile Vaginal Speculum For single use

Indications: This product is suitable for diseases that require gy-

neological examination.

Intend purpose: The product is used in gynecological examination and surgery.

Contraindications: None.

Caution:

A. When using the product, it should strictly implement the requirements of aseptic operation specifications and relevant regulations, and can be used by trained doctors, nursing staff. In the course of use, the operator or patient should be treated in a timely manner according to medical regulations.

B. If you find any possible quality doubt during use, please stop using it immediately and notify our company quickly.

C. Sterilized by EO gas. Product is sterile.

D. For single use. Discard after use.

E. Please checking the individual package before using it.

F. Don't use if package is opened or damaged.

G. Sterile and validity for 5 year. Don't use if overdue.

H. Discard after use. After use, the products shall be collected, transported, stored and disposed to ensure that the products are harmless to human health and the environment. Or the product can be disposed in accordance with the local medical regulations for waste disposal.

Patient group: It is suitable for people who require gynecological examination and surgery.

Intended users: It can be used by trained doctors, nursing staff.

Clinical benefits: The product plays an auxiliary role in the treatment of diseases.

Any residual risk and any undesirable side-effects: none

Install information: none

Performance characteristics: The product is provided in a sterile state and sterilized with ethylene oxide. Single-use

Physical performance:

a. The entire surface of the plastic part should be smooth and transparent, and the inner and outer sides of the upper and lower leaves should be smooth without obvious burrs.

b. The upper and lower leaves should match well. When the expander is in a closed state, the upper and lower leaves should basically match, and the deviation of the upper and lower leaves should not exceed 2 mm.

c. When the expander is expanded to its maximum position, a static load of 30 N is applied at a distance of 20 mm from the front end of the upper blade, and after 1 minute, no component shall break.

d. When the expander is expanded to its maximum position, a vertical downward load of 20 N is applied at a distance of 20mm from the front end of the upper blade, and the allowable deformation is not more than 30 mm.

Structural composition: The product are composed of upper bill and lower bill etc.

Method:

Open the primary package, close on upper bill and lower bill, insert into the vagina, press the handle, expand the field of vision, and lock the dilator, and then observe the patient's condition.

Storage: The product should be stored in a well-ventilated room temperature.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

All serious accidents concerning the medical device supplied by us must be reported to the manufacturer and competent authority of the member state where your registered office is located.

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