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): Proper cuff wearing icon. If the cuff is worn properly at the start of the measurement, this icon will be displayed.

): Improper cuff wearing icon. If the cuff is not worn correctly at the start of the measurement, this icon will be displayed and marked in the measurement result.

The measurement value can be stored automatically, using [memory function] (refer to Chapter11).

*Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your physician.

 **Note** 

Wait at least 4-5 minutes between measurements

When repeatedly measuring, the accurate blood pressure value may not be measured due to congestion in the arm. Please measure after the blood flow is smooth.

- When some factors affect the measurement results in measurement process, error messages will appear on the screen, you can obviate the malfunction and restart a measurement.
- The minimum value of the patient's physiological signal is the minimum limit that the device can measure. The device may obtain inaccurate measurement results when it is operated below the minimum amplitude or minimum value of the patient's physiological signal.

③In the state of the physiological parameter over-limit prompt is not triggered, press any button to carry on the corresponding button function; in audio prompt state, press any button (except 【ON/OFF】 button) to clear up the audio prompt.

④Take off the cuff, press 【ON/OFF】 button to turn the device off.

*The device will automatically turn off after 1 minute in which there is no operation to the device, even if you forget to turn the power off.

Chapter11 Memory Function

The sphygmomanometer is designed to store the blood pressure, pulse rate values and the date and time when measured, which are up to 100 groups. If there have been stored 100 groups, the earliest results will be deleted when saving the 101 group of measurement results.And the “OVERFLOW”will display on the bottom right of the screen,prompting that data record is full. You can delete the data.

11.1 Review Memory Values

- In the Main interface (the interface entered when startup), press the MEMORY button to display the average of the three latest measurements in large font. At this time, the icon  will appear on the screen.
 - Continue to press the DOWN button to view the latest measurement result. In the Memory interface, the measurement records are numbered from 1 to 100. Press UP and DOWN buttons to view other measurement records.
 - Press 【LIST】 button to display the data list interface.
 - Press 【TREND】 button to display trend interface.
- End to display the measurement values:

Press 【EXIT】 button to return to the main interface or hold 【ON/OFF】 button to turn the device off.

11.2 Delete Memory Values

Users can delete all memory values of a user instead of separately delete one memory value.

- Press 【MENU】 button to enter the system menu, select【DELETE DATA】option to enter its interface, select the user whose data to be deleted, after confirming again, all measurement results of the selected user will be deleted.

2. Finish Operation

Select 【EXIT】 to return to the previous menu, or hold 【ON/OFF】 button to turn the device off.

Chapter12 Bluetooth Transmission

The device can send the stored blood pressure values to a master device via Bluetooth (see the software help for details).

The user can search the device on the master device and connect them via the Bluetooth, then data can be transmitted. The data that have been transmitted to the master device will be marked with green dots in the review interfaces. In the Memory interface, the marker is located in the upper right corner of the right side of the serial number; in the BP List interface, the marker is located in the left side of the serial number.

Note: If the system time is incorrect, the device processes the measurement data but does not upload them, and the data will be marked with red dots.

About communication:

This product is a wireless device for low-power data communication systems in the 2.4GHz frequency band.

Due to the use of radio signals, it may be intentionally or unintentionally intercepted by third parties.

Do not use for confidential matters or life-critical applications.

It is strictly prohibited to connect other devices arbitrarily during device operation to prevent potential safety risks and data leakage.

Strictly protect all network confidential information to prevent unauthorized access

 **Note** 

When using this product near equipment operating in shared frequency bands, such as wireless routers, electronic microwaves, and wireless instruments, electromagnetic interference may occur. If so, please turn off other equipment or move this product to an area free from electromagnetic interference.

This section only applies to products with Bluetooth function.

Chapter13 Error Message

Error message will be displayed in the screen if there is something wrong when measuring. The causes and solutions are shown as follows:

Error Message	Causes	Solutions
Self-test failure System failure	Function abnormal	Please contact us
Loose cuff	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 10)
Air leakage	Cuff plug fall off	Make sure the cuff plug is securely inserted in the windpipe (refer to Chapter 10)
Air pressure error	Air pressure error	Measure again and do not move during measuring.
Weak signal	The pulse signal is too weak or the cuff is loose.	Correctly connect cuff (refer to Chapter 10)
Overpressure	Cuff is blocked or squeezed	Correctly connect cuff (refer to Chapter

		10)
Excessive movement Over range Saturated signal	The signal extent is too big owing to the arm or body moving or other reasons when measuring	Keep arm, body still, measure again
Time out	It takes too much time	

Chapter14 Troubleshooting

Abnormal Phenomenons	Causes	Solutions
BP measurement values too high or too low.	Cuff is not connected correctly.	Correctly connect cuff (refer to Chapter 10)
	Talk or move arms when measuring	Keep quiet and restart a measurement
	The turnout clothing presses the arm	Take off the clothing which presses the arm, and restart a measurement
No pressure	Cuff leakage	Buy a new cuff
	The cuff windpipe is not correctly connected with cuff	Correctly connect
Cuff deflates in short time	Loose cuff	Correctly apply cuff
It can not carry on measurement when press the measurement button		Switch on the power once again and restart a measurement
Power off suddenly when inflating	No use for a long time, the power of batteries can be exhausted owing to the changed temperature	Replace all four batteries with new ones.
Hold the on/off button but can not start the device	Power of batteries can be exhausted	Replace all four batteries with new ones.
	The battery polarities is reversed	Check the battery assembly for proper placement of the battery polarities.
Other phenomenon		Switch on the power once again and restart an operation. Replace the batteries. If no, please contact us.

Chapter15 Keys and Symbols

Your device may not contain all the following symbols.

Signal	Description	Signal	Description
	Caution: read instructions (warnings) carefully		Follow instructions for use
SIS	Systolic pressure	DIA	Diastolic pressure
INFO	Information	PR	Pulse rate (bpm)
IP20	Enclosure protection grade	CEM	Electromagnetic compatibility
	Manufacture date	P/N	Material code of manufacturer
	Type BF applied parts		Recyclable
	Open the Prompt sound indication		Close the Prompt sound indication
	Lot number		Expiration date
	This way up		Fragile, handle with care
	Keep in a cool, dry place		Atmospheric pressure limit
	Temperature limit		Humidity limit
	Manufacturer		Unique device identifier
	Low battery		Batteries Power
	1.No NIBP data to review 2.An indicator of signal inadequacy		1.No Pulse rate 2.An indicator of signal inadequacy
	WEEE disposal		Medical Device compliant with Directive 93/42/EEC
	Serial number		Class II applied
	Authorized representative in the European Union		Interface for connecting cuff
	MR Unsafe,can not be used in MRI		Socket for power adapter
	Medical device		Artery indicator label
	Keep away from sunlight		Imported by
	Cuff tied properly		Product code
	Irregular pulse		Average value icon
	Improper cuff wearing		Large movement during measurement

Chapter16 Maintenance, Cleaning and Keeping

***Please do obey the precautions and correct operating methods in this user manual. Otherwise, we will not responsible for any fault.**

 **Warning** 

Remove the batteries before cleaning. The accessories and main unit must be separated for cleaning.

Maintenance is not allowed during device using.

Do not squeeze the rubber tube on the cuff.

 **Caution** 

- High pressure disinfection to the device and accessories is not allowed.
- Do not let water or cleaning agent flow into the socket to avoid device damage.
- Do not soak the device and accessories in liquid.
- If any damage or deterioration of the device and accessories is found, please do not use it.

Maintenance:

- Clean the device and accessories regularly. It is recommended to clean them every one month.
- Before cleaning the device, remove batteries and disconnect it from the AC power. The accessories and main unit must be separated for cleaning. Do not maintain or repair the device during use.
- When cleaning the device, dip a clean cloth in isopropyl alcohol (70%), wring it out fully, and wipe the main unit, cuff and cuff windpipe separately for about 3 minutes, then use the other clean cloth moistened with distilled water, wring it out fully and respectively wipe the main unit, cuff and cuff windpipe for about 2 minutes. Repeat above 5 times until there is no obvious residual cleaning agent. Avoid isopropyl alcohol or water entering the main unit during cleaning. After cleaning, place the product in a dry and ventilated place to dry.
- Visual inspection to make sure the product is thoroughly cleaned. If any residue exists, please repeat the entire process described above.
- Clean the accessories before applying to a different patient.
- The device should be inspected and calibrated periodically (or obey the requirements of the hospital). It is available to inspect in the state specified inspection institution or by professional personal, or you can contact our company.Long press "USER" button in main interface for 5s to enter the calibration interface.

 **Advice** 

- Do not use gasoline, volatile oil, diluent, etc. to wipe the device.
- Do not clean or wet the cuff.

Storage:

 **Advice** 

- Do not expose the device in direct sunlight for long time, otherwise the display screen maybe damaged.
- The basic performance and safety of the device are not affected by the dust or cotton wool in home environment, while the device shall not be placed where with high temperature, humidity or dusty.
- Aged cuff may result in inaccurate measurement, please replace the cuff periodically according to the user manual.
- To avoid device damage, keep the device out the reach of children and pets.
- Avoid the device close to extreme high temperature such as fireplace, otherwise the device performance may be affected.
- Do not store the device with chemical medicine or corrosive gas.
- Do not place the device where there is water.
- Do not place the device where with slope, vibration or impact
- Take the batteries out if the device is not to be used for three months or longer.

Chapter17 NIBP Specification


Name	Electronic Sphygmomanometer	
Display mode	2.8" color LCD Display	
The degree of protection against ingress of water	IP20	
Measurement method	Oscillometric method	
Working mode	Automatic	
Operation mode	Continuous operation	
Cuff pressure range	0-300 mmHg(0-40 kPa)	
Accuracy	Pressure: ±3 mmHg(±0.4 kPa)	
Maximum pressure	Not exceed 300 mmHg (40 kPa)	
Measurement Range	Blood Pressure	SYS: 30~270 mmHg(4~36 kPa) DIA: 10~220 mmHg(1.3~29.3 kPa)
	Pulse: 40~240 bpm	
Error	Blood Pressure	The BP value measured by the device is equivalent with the measurement value of Stethoscopy, perform clinical verification in accordance with the requirements in ISO 81060-2, whose error meets the followings: Maximum mean error: ±5 mmHg Maximum Standard deviation: 8 mmHg
	Pulse	±5 bpm or ±5% select larger
Resolution	Pressure: 1mmHg(0.1 kPa)	
	Pulse: 1 bpm	
Operating temperature/humidity	+5 ℃~+40 ℃ 15 %RH~85 %RH(Non-condensing)	
Transport	Transport by general vehicle or according to the order contract, avoid pounded, shake and splash by rain and snow in transportation.	
Requirements for transport and storage environment(between two transports)	Temperature: -20 ℃~+55 ℃; Relative humidity: ≤95 %((no condensation)); No corrosive gas and drafty.	
Atmospheric pressure	700 hPa~1060 hPa	
Power supply	4 "AA" alkaline batteries, AC Adapter(AC, 100 V-240 V, optional)	
Rated current	600 mA	

Battery life	When the temperature is 23 ℃C, limb circumference is 270 mm, the measured blood pressure is normal, 4 "AA" alkaline batteries cab be used about 300 times.		
Dimensions	130(L)*110(W)*80(H) mm		
Unit Weight	300 gram(without batteries)		
Safety classification	Class II equipment (power supplied by power adapter)/Internally powered equipment (power supplied by batteries) Type BF applied part		
Service life	5 years from the date of manufacture		
Date of manufacturer	See the label		
Accessories	Standard Configuration: Cuff the range of limb circumference is 22-32 cm (middle part of upper arm) Model:IGN0066 Manufacturer:Contec Medical Systems Co., Ltd. Optional Configuration: AC Adapter Model:CMS0105 Manufacturer:Contec Medical Systems Co., Ltd. Cuff		
	Model	Limb circumference(cm)	Manufacturer
	IGN0009	10-19	Contec Medical Systems Co., Ltd.
	IGN0048	18-26	
	IGN0083	22-30	
	IGN0074	22-43	
	IGN0008	32-43	
	Basic performance	The error limit under environmental conditions is less than or equal to ± 0.4 kPa (± 3 mmHg)	
NIBP repeatability should be less than or equal to 0.4 kPa (3 mmHg)			
Bluetooth Specification	Working frequency: 2402 MHz ~ 2480 MHz Bandwidth: 2MHz Modulation mode: GFSK Transmitting power: 0 dBm, +4 dBm Receiving sensitivity: -93 dBm		

 **Disposal:** The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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- © Para medir corretamente a tensão arterial, substitua a braçadeira atempadamente.