AUDIOMETERS

BELL & HARP

MULTILANGUAGE USER

MANUAL

Document title: AU1D-Bell Harp User Manual IT-EN-FR-DE-ES
Code: AU1-MA003_A
Revision: Rev. 01
Date: 2022.11.24

BELL & HARP

AUDIOMETRI

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(*) Type A-E if the free field equivalent filter for AC output is enabled

CE CERTIFICATE	
MDR 2017/745/EU Classification	Class IIa
Classification rule (Annex VIII, 2017/745)	10
Notified body	TÜV SÜD Product Service GmbH Ridlerstrasse 65 D-80339 München, Germany
Notified body number	0123

AVAILABLE TESTS				
	<i>Bell Basic</i>	<i>Bell Plus</i>	<i>Harp Basic</i>	<i>Harp Plus</i>
Pure tone audiometry	•	•	•	•

Autothreshold	•	•	•	•
Speech audiometry	-	-	•	•
ABLB	-	-	•	•
SISI	-	-	•	•
DLI	-	-	•	•
Tone decay	-	-	•	•
Master hearing aid	-	-	•	•
Stenger	-	-	•	•
TEN test	-	-	-	opt.
QuickSIN	-	-	-	opt.

AVAILABLE SIGNALS				
Type	<i>Bell Basic</i>	<i>Bell Plus</i>	<i>Harp Basic</i>	<i>Harp Plus</i>
Pure tone	•	•	•	•
Warble tone	•	•	•	•
2 external inputs for speech audiometry (EXT1 and EXT2)	-	-	•	•
MIC input for live speech audiometry	-	-	•	•
Internal input (flash memory) for speech audiometry	-	-	-	•
Narrow-band noise (NBN)	•	•	•	•
White noise (WN)	-	•	•	•
Speech noise (SN)	-	-	•	•
Threshold equalizing noise (TEN)	-	-	-	•

SIGNALS SPECIFICATIONS	
Attenuators step	Bell: 5 dB Harp: 1 and 5 dB
Presentation mode	Continuous Pulsed, with rate of 0.5, 1 or 2 Hz
Frequency accuracy	0,1 %
Intensity accuracy	±3 dB between 125 Hz and 4 kHz ±5 dB above 4 kHz
Total Harmonic Distortion (THD)	AC: less than 2,5 % BC: less than 5,5 %
Warble tone	Frequency of the modulating signal: 5 Hz Modulation waveform: sine wave Modulation range: ±12%
NBN	Band: ½ octave, i.e.: - lower cut-off frequency $f_l = f / 1.1892$ - upper cut-off frequency $f_u = f \cdot 1.1892$

	where f is the centre frequency
WN	Lower cut-off frequency: 100 Hz Upper cut-off frequency: 24 kHz
SN	As specified in IEC 60645-1 §13
External signals	EXT1 and EXT2 input: max 3 Vrms
TEN	Lower cut-off frequency: 350 Hz Upper cut-off frequency: 6500 Hz

AVAILABLE OUTPUTS			
Output	Bell Basic	Bell Plus	Harp Basic/Plus
Air conduction			
TDH-39 / DD45 headphones or ER-3 / ER-5 / IP30 insert earphones	●	●	●
HDA-280 headphones	●	●	-
Bone conduction (B-71 bone vibrator)	-	●	●
Free field	-	-	●
Insert masking earphone	-	●	●

PURE AND WARBLE TONE AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES						
Freq. (Hz)	AC TDH39 DD45 (dB HL)	AC HDA280 (dB HL)	AC ER-3 IP30 (dB HL)	AC ER-5 (dB HL)	BC B-71 (dB HL)	FF (*) (dB HL)
125	80	90	90	90	-	75
250	100	100	105	100	45	85
500	110	115	110	110	65	95
750	115	120	115	120	70	95
1,000	120	120	120	120	75	95
1,500	120	120	120	120	80	95
2,000	120	115	120	115	80	95
3,000	120	115	120	115	75	95
4,000	120	115	110	110	75	95
6,000	105	100	95	100	55	90
8,000	95	90	90	90	50	85
12,500	-	80	-	-	-	-

(*) The maximum levels refer to the “normal” intensity range for free field. In case of extended range, add 10 dB to each value.

SPEECH AUDIOMETRY MAXIMUM INTENSITIES				
	AC TDH39 DD45 (dB HL)	AC ER-3 ER-5 IP30	BC B-71 (dB HL)	FF (dB HL)

		(dB HL)		
Type A	100	100	60	Norm.: 75 Ext.: 85
Type AE	80	100	60	Norm.: 75 Ext.: 85

EXTERNAL SIGNALS LEVEL INDICATOR (only for Harp Basic/Plus)	
Type of indicator	VU-meter
Dynamic range	+3...-20dB
Input Voltage at 0 dB	1.5 Vrms
Speed of the volume following	Increase: 60 dB/s Decrease: 60 dB/s

NARROW BAND NOISE (NBN) AVAILABLE FREQUENCIES AND MAXIMUM INTENSITIES							
Freq. (Hz)	AC TDH39 DD45 (dB EM)	AC HDA280 (dB EM)	AC ER 3 IP30 (dB EM)	AC ER 5 (dB EM)	BC (dB FL)	Insert mask (dB EM)	FF *) (dB HL)
125	60	65	70	65	-	80	75
250	80	85	85	85	105	95	85
500	95	100	95	95	110	105	95
750	100	105	100	100	110	110	95
1.000	105	120	105	100	110	115	95
1.500	105	110	105	100	110	115	95
2.000	105	105	105	100	110	115	95
3.000	105	105	105	100	100	115	95
4.000	105	105	100	100	100	110	95
6.000	100	95	90	95	90	95	90
8.000	90	95	80	80	90	95	85
12.500	-	70	-	-	-	-	-

(*) The maximum levels refer to the “normal” intensity range for free field. In case of extended range, add 10 dB to each value.

WHITE NOISE (WN) MAXIMUM INTENSITIES					
	AC TDH39 DD45	AC ER-3 ER-5 IP30	BC	Insert Mask	FF
WN	110 dB SPL	100 dB SPL	100 dB FL	100 dB SPL	Norm: 95 dB SPL Ext: 105 dB HL
WN Type A (speech)	90 dB EM	80 dB EM	45 dB EM	80 dB EM	Norm: 75 dB HL Ext: 85 dB HL

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ACOUSTICAL SAFETY OF THE DEVICE	
Alert condition	Tone intensity higher than 100 dB HL (IEC 60645-1, §5.2)
Safety measures in alert condition	1) The operator has to press “Higher dB” button to increase the intensity over 100 dB HL 2) Warning on the display 3) “Normally on” function disabled

COMPATIBLE TRANSDUCERS		
Type	Manufacturer	Model
Supra-aural headphones	Telephonics Corp.	TDH39
Supra-aural headphones	Radioear Corp.	DD45
(only for Bell) Supra-aural headphones	Sennheiser elec. GmbH	HDA-280 ¹
Insert earphones	Etymotic Research Inc.	EARTone 3
Insert earphones	Etymotic Research Inc.	EARTone 5
Insert earphones	Radioear Corp.	IP30
Bone vibrator	Radioear Corp.	B71
Insert earphone for masking	INVENTIS S.r.l.	IME-100

PATIENT – OPERATOR COMMUNICATION				
	Bell Basic	Bell Plus	Harp Basic	Harp Plus
Talk-over through:				
- Embedded microphone	•	•	•	•
- External microphone (not included).				
Talk back through:				
- Embedded speaker	-	-	•	•

¹ HDA-280 and ER-5 supported but not available anymore.

- External speakers (not included)				
- Patient microphone				
Patient response trigger	●	●	●	●

DATA MANAGEMENT		
	<i>Bell</i>	<i>Harp</i>
Internal memory capacity	100 patients	100 patients (*)
Integrated thermal printer (paper size: 112 mm)	-	Optional
Communication to the computer	USB	
Management software	Maestro, installed with or without proprietary database or as Noah module	

(*) Only pure tone audiometry exams can be stored in the internal archive.

DISPLAY	
Type	Graphical colour TFT LCD
Size	diagonal 4,3", 95 mm x 54 mm
Resolution	480 x 272
Available viewing modes	Normal / High contrast

MECHANICS	
<i>Models without the integrated printer</i>	
Size (LxDxH)	32 x 32 x 9 cm / 12.6 x 12.6 x 3.5 in
Weight (device only)	1.8 Kg / 4.0 lbs
Weight (including transducers and package)	4.2 Kg / 9.3 lbs
<i>Model with the integrated printer</i>	
Size (LxDxH)	32 x 39 x 9 cm / 12.6 x 15.4 x 3.5 in
Weight (device only)	2.3 Kg / 5 lbs
Weight (including transducers and package)	4.7 Kg / 10.4 lbs

SOCKETS ON THE REAR PANEL		
<i>Description</i>	<i>Type</i>	<i>Connector</i>
Power supply	In	DC plug
L and R headphones	Out	2 audio jack, 1/4" mono
Bone vibrator (not in Bell Basic)	Out	Audio jack, 1/4" mono
Insert masking earphone (not in Bell Basic)	Out	Audio jack, 1/4" mono
Patient response trigger	In	Audio jack, 1/4" mono
External microphone for talk-over or live speech tests (Harp)	In	Audio jack, 3.5 mm mono

USB	In - Out	USB type B
<i>Only on Harp audiometers</i>		
L and R Free Field output	Out	2 RCA
Monitor headset	Out	Audio jack, 3.5 mm stereo
Microphone for live speech tests	In	Audio jack, 1/4" mono
Patient microphone for talk-back	In	Audio jack, 1/4" mono
EXT1 and EXT2	In	2 RCA

SPECS OF INPUT FACILITIES	
<i>Input</i>	<i>Electrical property</i>
Power supply	Internal pin +6V, external pin 0V
Patient response	Switches 3V to logical input (switch current: 10mA) Impedance: 10K Ω
External input for speech aud.	Sensitivity: 3mV at max volume and 0Vu Impedance: 10K Ω . Freq. response: 75-12000Hz +/- 3dB
Microphone	Electret or 200 Ω dynamic microph. Impedance: 47K Ω Freq. response: 100-12KHz +/- 3dB Electret Bias: 2.2V through 2.2K Ω
Talk over	
Talk back	

SPECS OF OUTPUT FACILITIES		
<i>Output</i>	<i>Available Voltage</i>	<i>Nominal Impedance</i>
TDH 39	8Vpp	10 Ω
DD45	8Vpp	10 Ω
ER-3, ER-5, IP30	8Vpp	10 Ω
B 71	8Vpp	10 Ω
HDA 280	8Vpp	37 Ω
Free Field (SPKF)	8Vpp	100 Ω
Insert mask	8Vpp	32 Ω
Monitor	1Vrms by 32 Ω load	

SOUND ATTENUATION VALUES		
<i>Freq</i>	<i>TDH 39 / DD45^(*)</i>	<i>ER 3 / ER 5 / IP30</i>
<i>[Hz]</i>	<i>[dB]</i>	<i>[dB]</i>
125	3.0	33.5
250	5.0	34.5
500	7.0	34.5
750	-	-
1000	15.0	35.0
1500	-	-
2000	26.0	33.0
3000	-	-
4000	32.0	39.5
6000	-	-

8000	24.0	43.5
(*) With MX41\AR or PN 51 cushions		

EQUIVALENT FREE FIELD SPEECH AUDIOMETER		
	TDH39	DD45
	IEC 60645-1 / ANSI S3.6	IEC 60645-1 / ANSI S3.6
Freq.	G _F -G _C	G _F -G _C
[Hz]	[dB]	[dB]
125	-17,5	-21,5
160	-14,5	-17,5
200	-12,0	-14,5
250	-9,5	-12,0
315	-6,5	-9,5
400	-3,5	-7,0
500	-5,0	-7,0
630	0,0	-6,5
800	-0,5	-4,0
1000	-0,5	-3,5
1250	-1,0	-3,5
1600	-4,0	-7,0
2000	-6,0	-7,0
2500	-7,0	-9,5
3150	-10,5	-12,0
4000	-10,5	-8,0
5000	-11,0	-8,5
6300	-10,5	-9,0
8000	1,5	-1,5
Coupler: ANSI S3.7 (NBS-9A) IEC 60318-3 (6ccm)		

REFERENCE EQUIVALENT THRESHOLD LEVEL FOR PURE TONE						
	TDH39	DD45	HDA280	ER3 / ER5 / IP30	B71(*)	Free Field
Freq.	ISO 389-1 (ANSI S3.6)	ISO 389-1 (ANSI S3.6)	ISO 389-1 (ANSI S3.6***)	ISO 389-2 (ANSI S3.6)	ISO 389-3 (ANSI S3.6)	ISO 389-7 (ANSI S3.6)
	Coupler: ANSI S3.7-1995 (NBS-9A) IEC 60318-3 1998 (6ccm)f			IEC 60318-5	IEC 60318-6	refers to 0 deg incidence
[Hz]	dB [re 20 μPa]	dB [re 20 μPa]	dB [re 20 μPa]	dB [re 20 μPa]	dB [re 1 μN]	dB [re 20 μPa]
125	45.0	47.5	38.5	26.0	-	22.1
250	25.5	27.0	25.0	14.0	67.0	11.4
500	11.5	13.0	13.0	5.5	58.0	4.4
750	7.5 (8.0**)	6.5	9.0	2.0	48.5	2.4
1000	7.0	6.0	7.5	0.0	42.5	2.4
1500	6.5	8.0	9.5	2.0	36.5	2.4
2000	9.0	8.0	8.0	3.0	31.0	-1.3

3000	10.0	8.0	6.5	3.5	30.0	-5.8
4000	9.5	9.0	9.5	5.5	35.5	-5.4
6000	15.5	20.5	19.0	2.0	40.0	4.3
8000	13.0	12.0	18.0	0.0	40.0	12.6
10000	-	-	18.5	-	-	13.9
12500	-	-	30.5	-	-	12.3
(*) Calibration of bone vibrator (B71) refers to mastoid placement. (**) Value to apply for ANSI S3.6 (***) Value to apply for ANSI S3.6						
Reference Standard: ISO 389-1, ISO 389-2, ISO 389-3, ISO 389-5, ISO 389-8, ISO 389-7/A1 e ANSI S3.6						

REFERENCE EQUIVALENT THRESHOLD LEVEL FOR SPEECH IEC 60645-1					
	TDH39	DD45	ER3 / ER5 / IP30	B71	Free Field
	Coupler: ANSI S3.7 (NBS-9A) / IEC 60318-3 (6ccm)		IEC 60318-5	IEC 60318-6	
	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 1 µN]	dB [re 20 µPa]
Linear Type A	20.0	20.0	20.0	55.0	20.0
Equ. FF Type AE	20.0	20.0	-	-	-
For stimulus: LINE, MIC, DIGITAL, SPEECH NOISE, WHITE NOISE					

REFERENCE EQUIVALENT THRESHOLD LEVEL FOR SPEECH ANSI S3.6					
	TDH39	DD45	ER3 / ER5 / IP30	B71	Free Field
	Coupler: ANSI S3.7 (NBS-9A) / IEC 60318-3 (6ccm)		IEC 60318-5	IEC 60318-6	refers to 0 deg incidence
	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 20 µPa]	dB [re 1 µN]	dB [re 20 µPa]
Linear Type A	19.5	18.5	12.5	55.0	14.9
Equ. FF Type AE	19.5	18.5	-	-	-
For stimulus: LINE, MIC, DIGITAL, SPEECH NOISE, WHITE NOISE					

Annex B:

Specifications of the amplifier / speakers for free field

In order to reach the values of the intensities for free filed output specified in the Appendix A (both for “normal” and “extended” range), the amplifier and the loudspeakers (or the active loudspeakers) must meet the following specs.

Characteristics of the dual channel amplifier:

	Normal range	Extended range
Minimum output power	2 x 60 W rms	2 x 200 W rms
Passband	90 Hz – 20 kHz	

Characteristics of the loudspeakers:

	Normal range	Extended range
Minimum output power	2 x 60 W rms	2 x 200 W rms
Sensitivity (at 1 m for 1 W, 1 kHz)	> 88 dB SPL	> 98 dB SPL
Passband	90 Hz – 20 kHz	

Annex C:

Audiometric symbols

Air Conduction

	Right	Left	Bilateral
HL threshold without masking	○	×	⊗
HL threshold with masking	△	□	
UCL threshold	E m ┐ or ┐	≡ m ┐ or ┐	

Bone Conduction

	Right	Left
HL threshold without masking	< / >	> / <
HL threshold with masking	┐ / ┐ ◁ / ▷ or ▷ / ▷	┐ / ┐ ▷ / ▷ or ▷ / ▷

The direction of the symbols used for BC and the symbol used for masked BC can be selected in the device set-up window.

Free Field

	Right	Left	Bilateral
HL threshold			

In case of absence of response, a diagonal arrow pointing towards bottom is added in the lower part of the symbol.

Annex D:

Electromagnetic Compatibility

Bell & Harp have been thoroughly tested and respects the limits for electro-medical devices specified by IEC 60601-1-2 standards. These limits ensure reasonable protection against hazardous interference in typical medical installations.

The instrument generates, uses and radiates radio frequency energy. If not installed and used according to the instructions in this manual, it may interfere with other nearby devices. No guarantee is given that interference will not occur under certain conditions.

This instrument is suitable for use in professional healthcare facility environment, i.e. in hospital environments, except for near active HF surgical equipment and RF shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

Bell & Harp should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Bell & Harp should be observed to verify normal operation in the configuration in which it will be used.

The existence of electromagnetic interference can be verified easily by switching the instrument off and back on again. If it is proven that the device is indeed interfering with other devices, try to solve the problem by adopting one of the following solutions:

- change the orientation and/or position of the affected device;
- move the two devices further away from each other;
- contact the manufacturer or authorised service organisation for further assistance.

List of cables, transducers and accessories

Cables, transducers and accessories with which Inventis claims the compliance with the IEC 60601-1-2 standard are those ones supplied with the device, in particular the followings:

- 6Vdc Medical Grade Power
- Power supply cable (maximum length: 1.8 m)
- TDH39 or DD-45 transducers, cable 2 m, double signal shielded
- Insert earphones, ER-3 (Etymotic) or IP30 (RadioEar)
- Bone vibrator B-71 transducer with twisted not shielded double signal cable 2m
- IME-100 Insert masking earphone with shielded cable 2m

- Monitor headset with boom microphone, manufacturer: Sennheiser, model: PC3
- Talk back microphone with 2m shielded cable
- Patient response switch (manufactured by Inventis) with 2 m shielded cable
- USB cable, shielded, maximum length: 2 m

The use of accessories, transducers and cables other than those specified, except for transducers and cables sold by the manufacturer as spare parts for internal components, may result in increased emissions or decreased electromagnetic immunity of the device and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of Bell & Harp, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

The instrument has no ESSENTIAL PERFORMANCE as related to IEC 60601-1

Note: All necessary instruction for maintaining compliance with regard to electromagnetic compatibility can be found in the maintenance section in this manual. No further steps are required.

Guidance and manufacturer's declaration – electromagnetic emissions		
Bell & Harp are intended for use in the electromagnetic environment specified below. The customer or the user of Bell & Harp should assure that they are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR11	Group 1	Bell & Harp use RF energy only for their internal function. Therefore, their RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	Bell & Harp are suitable for use in professional healthcare facility environment and directly connected to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flickers emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity			
Bell & Harp are intended for use in the electromagnetic environment specified below. The customer or the user of Bell & Harp should assure that they are used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance

Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ⁽¹⁾ ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air ⁽¹⁾	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a professional healthcare facility environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a professional healthcare facility environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. $40\% U_T$ (60% dip in U_T) for 5 cycles. $70\% U_T$ (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	$< 5\% U_T$ ⁽²⁾ ($> 95\%$ dip in U_T) for 0,5 cycle. $40\% U_T$ (60% dip in U_T) for 5 cycles. $70\% U_T$ (30% dip in U_T) for 25 cycles. $< 5\% U_T$ ($> 95\%$ dip in U_T) for 5 s.	Mains power quality should be that of a professional healthcare facility environment. If the user of Bell & Harp requires continued operation during power mains interruptions, it is recommended that Bell & Harp be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic field should be at levels characteristic of a professional healthcare facility environment.
Note: ⁽¹⁾ A lock or reboot of the device with no permanent damage is acceptable ⁽²⁾ U_T is the a.c. main voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
Bell & Harp are intended for use in the electromagnetic environment specified below. The customer or the user of Bell & Harp should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compl. Level	Electromagnetic environment – Guidance

Conducted RF IEC 61000-4-6	3 Vrms 0.15 MHz to 80 Mhz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz	3 Vrms 0.15 MHz to 80 Mhz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer than 30cm (12 inches) to any part of Bell & Harp, including cables specified by the manufacturer. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 Mhz to 2,7 Ghz	3 V/m 80 Mhz to 2,7 Ghz	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular / Cordless) telephones and land mobile, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If measured field strength in the location in which Bell & Harp are used exceeds the applicable RF compliance level above, Bell & Harp should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating Bell & Harp.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity	
Function to verify for freedom from unacceptable risk	Acceptance pass/fail criteria
Sound generation operating correctly	No sound from trasnducers exceeding 80dB; a lock or reboot of the device is acceptable

AUDIOMETERS

BELL & HARP

MULTILANGUAGE USER MANUAL

Document title: AU1D-Bell Harp User Manual BG-RO-HU
Code: AU1-MA003_C
Revision: Rev. 01
Date: 2022.11.24

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AUDIOMETERS

USER MANUAL

Read this manual thoroughly before using the device. Pay particular attention to Chapter 1 (“Safety: warnings and information”) and Chapter 2 (“Installation, power up and power down”).

Internal inspections and repairs must only be performed by authorised personnel.

Copyright: INVENTIS S.r.l. has the copyright on this manual. It may not be copied, reproduced or altered, in total or in any part, without the written specific authorization from INVENTIS S.r.l..

Inventis ® is a trademark of INVENTIS S.r.l..

QuickSIN™ is copyrighted by Etymotic Research Inc. and licensed to INVENTIS S.r.l., 2013.

Summary

<i>Foreword</i>	v
CHAPTER 1: Safety: warnings and information	1
Operator Manual	1
Operator responsibilities	1
Intended purpose	2
Indication for use and end users of the device	2
Medical conditions	2
Precautions	2
Disposal	5
Conformity	6
Symbols on labels	6
CHAPTER 2: Installation, power up and power down	9
Precautions	9
Connections	9
Power-up and main screen	11
Power-down	11
Connection to PC	11
CHAPTER 3: Controls	13
Function buttons	13
Control of channels	13
Common controls	14
Gain control and playback of word lists (Harp model only)	15
Patient-operator communication	15
CHAPTER 4: Maintenance	17
Periodic checks	17
Maintenance of Transducers	18
Cleaning the instrument	19
Replaceable parts	19

EN

Loading the print paper 20
Repairs and technical assistance..... 20
CHAPTER 5: Troubleshooting 21

Foreword

Thank you for purchasing an Inventis audiometer.

Bell and Harp audiometers offer a combination of innovative features placing them at the top of their class: a generously proportioned colour graphics display, internal flash memory, elegant design, and notable ease of use.

EN

The Inventis company has always considered the use of its devices in conjunction with computers to be a factor of key importance. Maestro software, available *with or without proprietary database* or as *Noah module*, can be used to connect any Inventis audiology device to a computer, so that examinations can be archived in a database, and the examination currently in progress can be displayed on the computer screen. We would also remind you that Inventis has developed a complete line of audiology devices: in addition to these middle ear analyzers, the company's product line includes a range of audiometers and a wireless video otoscope.

For further information, and to report any problems of whatever description that may be encountered, contact the company at:

INVENTIS S.r.l.

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35127 Padova, Italia

Tel.: +39 049 8962844 – Fax: +39.049.8966343

www.inventis.it

info@inventis.it

Safety: warnings and information

OPERATOR MANUAL

We recommend that you read this manual thoroughly to obtain maximum efficiency and safety from your audiometer. It is particularly important that you read this chapter in its entirety. It contains essential warnings and important information on the safe and correct use of the instrument.

In this manual, the safety symbol shown below draws your attention to information that is particularly important for safe and correct use.

OPERATOR RESPONSIBILITIES

The Bell and Harp audiometers guarantee consistent and dependable performance only when used in accordance with the instructions and procedures described in this manual.

Should the device need to undergo repair or maintenance work, it must be disconnected from the electrical power supply, and not used again until after this work has been completed. When defective or failed parts must be renewed, accept only original replacement parts supplied by INVENTIS S.r.l. Any repairs must be entrusted exclusively to Inventis or to service technicians approved by Inventis.

No parts of the device must be modified or replaced without authorization from Inventis.

The user of the device assumes full responsibility for any malfunction resulting from improper use or operation, likewise from maintenance or repair work performed by third parties other than INVENTIS S.r.l. or Service Centers approved by Inventis. INVENTIS S.r.l. and approved Service Centers will answer for the performance and reliability of the equipment only if:

1. adjustments, modifications or repairs have been entrusted to persons authorized by Inventis;

2. the electrical power supply and ground connections of the system comply with applicable standards for electromedical devices.

INTENDED PURPOSE

Bell/Harp medical device is an audiometer. An audiometer is a device that helps the operator in defining the patient's auditory sensitivity by generating and delivering to the patient sound stimuli of different types and intensities for diagnostic purposes.

INDICATION FOR USE AND END USERS OF THE DEVICE

Bell/Harp is intended for use by healthcare ENT professionals in hospitals, ENT clinics and audiology offices in conducting hearing evaluations and assisting in diagnosis of possible otologic disorders. There is no patient population restriction in the use of the device; always be sure to perform an otoscopy before using the device.

These tests must be conducted in a quiet environment to avoid artifacts.

MEDICAL CONDITIONS

Conditions of impaired sensitivity of the auditory system or any conditions in which the auditory system is thought to play a role in diagnosis.

PRECAUTIONS

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

To ensure correct and safe use of the audiometer, the following precautions must be observed.

Installation and general precautions

Make certain that the required ambient conditions are met (during transport, storage, and operation):

Operation	Temperature: between 15°C (59°F) and 35°C (95°F)
	Relative humidity: between 30% and 90% (non-condensing)
	Pressure: between 700 hPa and 1060 hPa
Transport and storage	Temperature: between -10°C (14°F) and 50°C (122°F)
	Relative humidity: max. 90% non-condensing
	Pressure: between 500 hPa and 1060 hPa
Warm-up time	1 minute

EN

Bell and Harp audiometers will not be protected if exposed during use to flammable anaesthetic gases or similar products. Risk of explosion.

Avoid installing and using Bell and Harp audiometers close to any sources of strong electromagnetic field: these could interfere with the operation of the appliance.

Use only original detachable parts supplied by INVENTIS S.r.l., unless specifically indicated otherwise.

Use only power adapters intended for medical equipment provided with the device, certified to IEC 60601-1, with the following specifications:

- model without printer:

Main unit: 6V, 1.67A d.c.

External adapter: SL POWER MENB1010A0603F02

100-240Vac 50/60 Hz 0.3-0.1A (included) responding to IEC 60601-1 standard

- model with printer:

Main unit: 6V, 3.2A d.c.

External adapter: SINPRO MPU31-103

100-240Vac 50/60 Hz 0.9-0.34A (included) responding to IEC 60601-1 standard

Bell and Harp audiometers are medical devices: if connected to a computer (or any external device, such as a CD player) located within the “patient area” (as defined in IEC 60601-1), this likewise must be a medical device, or protected by an isolating transformer, in order to ensure that the computer (external device) + audiometer combination is in compliance with EN standard 60601-1.

Bell and Harp audiometers can be used in conjunction with a soundproof booth to conduct tests under optimum acoustic

conditions. Before connecting the audiometer to a soundproof booth, check that the sockets are compatible with the specifications prescribed for each connector.

Bell & Harp need special precautions regarding EMC and need to be installed and put into service according to the EMC information provided at the end of this manual.

Use of portable and mobile RF communications equipment can affect the correct operation of Bell & Harp devices. Make reference to the EMC information at the end of this manual.

The power adapter cable is considered the means to disconnect the device from the mains supply.

Do not position the device so that it is difficult to disconnect the device from the mains supply.

Calibration

The calibration should be performed at least once every 12 months, and whenever a transducer is replaced.

The calibration of the audiometer is valid only for the transducers supplied with the appliance. If a transducer is replaced, the audiometer must be recalibrated.

The calibration of the audiometer is valid for transducers supplied with the audiometer, if connected directly to the instrument, without any interposition of extension leads and without the passage from connectors to panel (as habitually occurs in soundproof booth installations). If the transducers are not connected directly to the audiometer, a new calibration procedure will be required before the instrument is used.

In each test window, when you select a not-calibrated transducer, the background of the 'output' area will be displayed in red colour. Moreover, you will not be able to send any stimulus through not-calibrated transducers.

Take note of the calibration interval indicated for the audiometer. Use of the instrument beyond the calibration interval expiry date can lead to unreliable diagnoses.

Hygiene

The eartips of insert earphones are disposable, likewise that of the masking earphone; do not use the same eartip for different patients. Dispose of eartips after use.

Disinfect the cushions of audio cups between one patient and the next, following the procedure described in CHAPTER 4: Maintenance.

Use

The audiometer can generate tones at an intensity potentially damaging to the patient. Take particular care to set the intensity of the tone correctly before it is presented.

When conducting audiometry using insert earphones, do not insert or in any way try to conduct measurements without proper foam tip in place

The desktop microphone shall be placed at a distance not greater than 15 cm from the mouth to reduce the background noise

Maintaining the previous intensity of the stimulus when changing frequency, transducer or stimulation side can result in potentially harmful signals being presented to the patient.

To present a stimulus signal stronger than 100 dB HL, the operator must first press the F1 function button (HIGHER dB), which is active only when the intensity of the stimulus reaches 100 dB HL.

EN

DISPOSAL

Like all electronic devices, your audiometer contains extremely small quantities of certain hazardous substances such as cadmium or mercury. If such substances are allowed to enter the normal waste disposal cycle without suitable preliminary treatment, they can cause damage to the environment and to health. All parts of the audiometer must therefore be disposed of separately.

At the end of its life, take (or have taken) the disused instrument to a civic waste disposal and recycling facility, or return it to the reseller against the purchase of an equivalent new instrument.

Separate waste collection and the subsequent operations of treatment, recycling and disposal facilitate the manufacture of new appliances from

recycled materials, limiting any negative impact on the environment and public health that might otherwise derive from improper disposal.

CONFORMITY

Bell and Harp audiometers are class IIa medical device, according to Annex VIII of Medical Device Regulation (MDR) 2017/745/EU.

The Inventis Quality Management System has been certified by leading assessment body TÜV as compliant with ISO 13485 standard.

SYMBOLS ON LABELS

Name and address of the manufacturer

This symbol means that, for a safety use of the device, it is important that the user refers to the warnings contained in this manual.

Refer to instruction manual for use

This symbol means this products is covered by the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). It is required not to dispose this product as unsorted municipal waste, but to collect it separately.

Device with applied parts of type BF (IEC 60601-1).

DC power supply

Product conforms to European Community Medical Device Regulation (MDR) 2017/745/EU. Class IIa device; number of notified body: 0123 (TÜV SÜD Product Service GmbH).

Rx only

Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner

Medical device

IP20 *IP (Ingress Protection) Code: this device is protected against the access of objects sized > 12.5 mm, not protected against liquid.*

REF *Catalogue number*

Serial number of the device. The number is made up of 13 alphanumeric characters indicating the model, series, year of manufacture and serial number. In particular, the number comprises these segments:

SN

- *first 5 characters: Inventis product code*
- *characters 6 and 7: year of manufacture ("10" denotes 2010)*
- *characters 8... 13: incremental number*

EN

UDI code

(01)08054187380068(21)AU1DC16200570

CHAPTER 2:

Installation, power up and power down

Whilst the installation of a Bell or Harp audiometer is a relatively simple procedure, it should be entrusted to a person with the requisite skills. If the installation is not performed correctly, the system could be affected by safety problems when in use.

This chapter describes the procedure for installing the system.

Keep the packaging materials, in case the audiometer should need to be sent to the dealer or to Inventis for any reason.

EN

PRECAUTIONS


Like any other electric or electronic device, the Bell and Harp audiometers will emit electromagnetic waves. Even though its emissions are within the limits of standards, other electronic devices close to the audiometer might be affected, if particularly sensitive to electromagnetic interference.

Should this occur, check just by switching the audiometer OFF and ON and try to eliminate the interference using one or more of the following solutions:

- change the orientation and/or the position of the device affected by interference;
- distance the affected device from the audiometer;
- plug the affected device into a power socket on a circuit other than the circuit to which the audiometer is connected;
- consult the manufacturer or a service centre for assistance.

CONNECTIONS

All connection points for detachable parts are located on the rear panel, as also is the power switch. This section refers to the Harp model. In the case of the Bell audiometer, certain connectors are not included.

Plug all detachable parts into the respective sockets as indicated in the following table:

Connector	Detachable part
 INSERT MASK	Insert earphone for masking
 BONE VIBR.	Bone vibrator
HEADPHONES	AC headphones: Right R and Left L
INSERT	AC Insert earphones: Right R and Left L
 SPKR	Free field speaker: Right R and Left L
MIC	Operator microphone for live speech tests
TALK BACK	Patient microphone
TALK OVER	Microphone for operator to patient communications or live speech test
MONITOR	Monitor headphones
EXT	External inputs 1 and 2 for speech audiometry with external audio source
 PATIENT RESP.	Patient response switch
USB	USB cable for the connection to the PC

6V---
1.67A
⊖⊕⊖⊕

Connect the power adapter to the input and to the relative power cord, which should be plugged into an a.c. mains socket supplying the voltage indicated on the adapter label.

All connections must be made with the appliance switched off, that is to say with the power switch set to position 0.

Use only the power adapters included in the package, certified to IEC 60601-1.

Make certain that the electrical power supply and ground connections comply with the applicable standards for electro-medical devices. Risk of electric shock

EN

The rear panel also lists the available transducers (*COMPATIBLE TRANSDUCERS*), depending on the device model.

POWER-UP AND MAIN SCREEN

Once all cables and leads are connected, the instrument can be switched on from the rear panel.

A few seconds after power-up, the display of the audiometer will show the main screen.

POWER-DOWN

Press the switch on the rear panel to safely power down the device.

CONNECTION TO PC

Bell and Harp audiometers can be interfaced with a computer with Maestro software, by way of a USB connection. The connection is plug-and-play, with no special drivers required for installation purposes: a few seconds after plugging in, the operating system will recognize the devices and install the drivers automatically.

Use the supplied cable (standard USB A/B cable) to connect the Bell/Harp audiometer to one of the USB ports of the computer

Refer to Maestro user manual for more details about the software.

CHAPTER 3:

Controls

FUNCTION BUTTONS

There are four buttons associated with configurable soft keys, whose functions are indicated at the bottom of the display window, immediately above the buttons. Moreover, there are six dedicated function buttons:

EN

Control	Operation
BACK 	Go back to the previous window
NEW PATIENT 	Clear all examinations conducted in the current session
SAVE EXAM 	Save the current test to the patient memory of the instrument
ERASE DATA 	Delete data pertinent to the current graph of the examination
TO PC 	Send the current examination to a computer, on which Maestro software must be installed
HELP 	Provide general information on the functions of the current window

CONTROL OF CHANNELS

Control	Operation
RIGHT LEFT  	Select stimulation on either side

	Select the output type
	Select the input signal
	Switch on the channel
	Enable the “Normally ON” option (the output will be permanently active, unless the channel is switched off by touching the interrupter)
	Attenuate the output (a knob for each channel)

COMMON CONTROLS

Control	Operation
	Select the pulsed mode for presentation of the stimulus signal
	Select the pulsed mode for presentation of the stimulus signal
	Enable tracking (keep the same difference in dB between stimulus and masking)
	Enable lock (the masking tone will be switched on / off together with the stimulus tone)
	Change to the previous/next frequency. In speech audiometry, count the words recognized and not recognized by the patient
	Store the hearing level

GAIN CONTROL AND PLAYBACK OF WORD LISTS (HARP MODEL ONLY)

Control	Operation
	Show the volume control window for communications between patient and operator
	Show the gain control window for the speech material
	Show the volume control window for monitor signals used by the operator
	Harp Plus: Start playing the current list, or pause the list if already playing. The Led remains permanently alight during Play, and blinks as long as Pause is selected
	Harp Plus: Terminates playback of the current list
	Harp Plus: - Bring up the list preceding the current list, if no list is playing - Jump back either to the previous word (if the list is indexed), or by 5 seconds, if a list is playing
	Harp Plus: - Bring up the list that follows the current list, if no list is playing - Jump forward either to the next word (if the list is indexed), or by 5 seconds, if a list is playing

EN

PATIENT-OPERATOR COMMUNICATION

Control	Operation
	Enable the operator to talk to the patient (by way of the transducer currently selected). Communication is enabled by pressing and holding the button disabled on release

	Enable the patient to communicate with the operator. The voice of the patient can be heard: 1) through the speaker incorporated into the chassis; 2) through headphones or external speakers connected to the “Monitor” socket on the rear panel
	Light up whenever the patient presses the response button

CHAPTER 4:

Maintenance

Bell and Harp audiometers do not require any special periodic maintenance other than calibration, checks and normal cleaning, all of which are described in this chapter.

The performance and safety of the instrument will be maintained if the recommendations for care and maintenance given in this chapter are observed.

The instrument must be switched off and disconnected from the power supply before commencing any kind of cleaning operation.

EN

The inspection and servicing of internal components must be left entirely to technicians approved by INVENTIS S.r.l., included the internal coin cell replacement.

Transducers are manufactured utilizing ultra-fragile diaphragms that could be damaged in the event of impact. Handle with care during maintenance operations.

Do not perform any service or maintenance operation when the device is used on a patient.

PERIODIC CHECKS

The procedure described under this heading must be carried out when the instrument is used for the first time each day.

The tests must be done with the audiometer in its installation position.

- Before switching on the instrument, check that no sign of damage is visible in any part of the device, including detachable parts and the external power supply; double check visual integrity of insulation of the mains cable and the connectors, and verify that they are not exposed to any kind of mechanical load that could involve damage; verify that all the parts and cables are properly connected.
- Check subjectively that the air conduction and bone conduction output is equal on both channels and all frequencies, in order to do

this 10 or 15 dB, just enough to hear is applied. The person who carries out this check should have good hearing

- Check at a level of 60dB in AC and 40dB in BC that there is no distortion, noise or parasitic signals in any of the frequencies
- Check that the interrupter key(s), the patient response switch(es) and the keyboard indicators function correctly
- Check that the attenuator knobs function correctly without noises or interference between channels
- Check the speech audiometry inputs by doing a speech test with every speech input
- Check the headband strain of headset and of bone vibrator
- Check the communication with the patient

If any part or transducer has any malfunction, consult the “Troubleshooting” Chapter.

Always check that the calibration interval has not elapsed: the expiry of the interval is indicated at the bottom of the main presentation screen.

Calibration must be entrusted to technicians approved by INVENTIS S.r.l.. The operation should be performed at least once every 12 months, and whenever a transducer is replaced.

MAINTENANCE OF TRANSDUCERS

Do not use liquids or sprays to clean the audiometer.

Do not allow dust to collect on the transducers. In addition:

- The headphone cushions are made of biocompatible material but are not sterile: to prevent the spread of infection and guarantee the biocompatibility of the material, whenever the headphones are to be worn by a new patient, the DD45/TDH-39 cushions must be wiped with denatured alcohol wipe or denatured alcohol with a microfiber cloth
- The ear tips of insert earphones and the masking insert earphones (IME-100) are intended to be inserted in the patient’s ear canal. They are made of biocompatible material and disposable: use once only and dispose of in accordance with current health and safety regulations

- The bone vibrator comes into contact with the patient's uninjured skin. It is made of biocompatible material, but it is not sterile: whenever this transducer has to be worn by a new patient, it must be wiped with denatured alcohol wipe or denatured alcohol with a microfiber cloth.

The eartips of insert earphones are not sterile. The use of unsterilized earpieces can cause ear infections.

The bone vibrator and the headphone cushions can be repeatedly cleaned as described in the "Maintenance of Transducers" paragraph. In the event of any malfunctioning after any cleaning operation, contact an Inventis service technician.

Though the bone vibrator and headphone cushions can be repeatedly cleaned, always check that their characteristics and integrity are maintained. To do so, it is sufficient to perform the tests described in the "Periodic checks" paragraph. As soon as any failure is encountered, contact an Inventis service technician to verify whether your transducer needs to be replaced.

To avoid damaging the DD45/TDH39 headphones, do not push it against a flat straight surface as this can create vacuum and cause a damage to the transducer (suction cup effect).

EN

CLEANING THE INSTRUMENT

To prevent the accumulation of dust on the audiometer, always cover the instrument with the protective cloth when not in use. Dust the audiometer regularly to ensure it is kept clean.

All parts not mentioned specifically in the previous section can be cleaned using a lint-free soft cloth moistened with a solution of water and mild detergent; in case of sanitization, moisten the cloth with hydrogen peroxide at a 3% concentration. The device allows multiple cleanings without degradation of basic safety or performances; always check that the device characteristics and integrity are maintained. To do so, it is sufficient to perform the tests described in the "Periodic checks" paragraph. As soon as any failure is encountered, contact an Inventis service technician to verify whether any parts need to be replaced.

REPLACEABLE PARTS

The transducers and detachable parts can be disconnected from the appliance. Should a fault develop in any one of these devices, the audiometer must be

switched off and isolated from the power supply, and the defective item then disconnected from the appliance.

All detachable parts of the audiometer are designed specifically for use with the appliance. Only parts supplied by Inventis should be connected to the audiometer.

LOADING THE PRINT PAPER

To load a new roll of paper in the printer, the rear cover must first be opened and the empty spool removed.

The heat-sensitive side of the paper is located externally of the winding mechanism: the paper must be positioned so that the heat-sensitive face does not come into contact with the rubber pinch roller.

Raise the green lever and offer the paper to the bottom part of the roller. Feed the paper through until it emerges from the top part at the required height. Lower the green lever.

Make certain that the paper is correctly aligned. If this is not the case, follow the manual loading instructions to correct the alignment.

Slip the paper through the slot in the cover and close the cover.

REPAIRS AND TECHNICAL ASSISTANCE

Before contacting the service department, make certain that all the possible solutions in the “Troubleshooting” Chapter have been tried.

All the parts being returned to the manufacturer for repair and service shall be cleaned and sanitized. Transducers should be sealed in a transparent bag.

Important: should the instrument need to be sent to the Inventis service department or returned to the dealer, make certain that the original packing is used, and that all detachable parts and transducers are enclosed.

CHAPTER 5:

Troubleshooting

Problem	Possible cause	Solution
No signal from a transducer	Transducer not connected to the correct output	Connect the transducer to the correct output
	Transducer damaged	Contact Inventis service department or dealer
No signal from patient response button when pressed	Wrong connection	Connect the patient response button to the correct socket
	Patient response button damaged	Contact Inventis service department or dealer
No signal from monitor headphones	Wrong connection	Connect the monitor headphones to the correct socket
	Headphone volume too low	Press the MONITOR button to increase the volume
Patient's voice cannot be heard	Problems with TALKBACK input connection	Check the connection to the TALKBACK input
	TALKBACK volume too low	Adjust the TALKBACK volume control
Connection between PC and audiometer	Problems with USB connection	Check the USB connection between instrument and computer

EN

cannot be established	USB cable damaged	Change the USB cable (standard USB A/B cable)
INT1 or INT2 input cannot be selected in speech audiometry mode	Internal memory not enabled	Verify that the audiometer model has an internal memory
	No speech material stored in the internal memory	Use <i>ATIT</i> software to check that the internal memory contains speech material
Date and time prompt appears when instrument is powered up	Internal battery failure	Contact Inventis service department or dealer to have the battery replaced
Display remains blank	Instrument not powered up correctly	Check the connection to the power supply and make certain the switch is in position “1”
	Adapter damaged	Contact Inventis service department or dealer
You cannot access to a test	Optional test not enabled	Contact your reference technical service to obtain the licence, communicating the device serial number

Where the audiometer is used in conjunction with a soundproof booth, check that the connections both inside the booth and between the booth and the instrument are correct and secure.

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