PALLONE RIANIMATORE MONOUSO
DISPOSABLE RESUSCITATOR BAG
INSUFFLATEUR DE REANIMATION MONOUSAGE
BALON DE REANIMACION MONOUSO
REANIMADOR DESCARTÁVEL
EINWEG-BEATMUNGSBEUTEL
ΑΝΑΖΩΟΓΟΝΗΤΗΣ ΜΙΑΣ ΧΡΗΣΗΣ
RESUSCITATOR FÖR ENGÅNGSBRUK
RESUSCITATOR DE UNICĂ FOLOSINȚĂ
ELDOVHATÓ ÚJRAÉLESZTŐ
ΡΕΑΝΙΜΑΤΟΡ ЗА ЕДНОКРАТНА УПОТРЕБА

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REF

GIMA 34277 RE-22415
GIMA 34248 RE-22515
GIMA 34249 RE-22615


Besmed Health Business Corp.
No. 5, Lane 116, Wu-Kong 2nd Road, Wu-Ku District,
New Taipei City 24888, Taiwan
Made in Taiwan

Casus Europe B.V.
Lange Viestraat 2b, 3511 BK Utrecht
The Netherlands

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

ENGLISH**Instruction for Use****Trade Name:**

Disposable Resuscitator set

1 Product Description

- Manual Resuscitator sets (Disposable) are manually operated and meant as an adjunct to artificial respiration and cardiopulmonary resuscitation. Manual Resuscitator sets (Disposable) is meant for artificial inspiration and expiration by simple manual compression. It is composed of 4 components: A. Non-rebreathing (duckbill) valve, B. Resuscitator bag, C. Reservoir valve, and D. Oxygen reservoir, the last two components being optional and to be removed if additional oxygen is not supplied.
- When compressing down the resuscitator, a positive pressure is created, which forces air into a patient's respiratory tract, causing inspiration. During this time, the intake valve is closed,

2. Intended Purpose

Manual Resuscitator sets (Disposable) are designed for use as an adjunct to artificial respiration and cardiopulmonary resuscitation. The resuscitator can be used to ventilate the apneic patients and to augment ventilation and oxygen delivery to the spontaneously breathing patient.

2.1 Indications

Device Models	Reference number	Product Description	
	RE-22415	Single Use Resuscitator 1600 ml W/O POP-OFF, Oxygen tubing, Mask, Oxygen reservoir, Adult	
	RE-22515	Single Use Resuscitator 500 ml W/ 40 cmH2O POP-OFF, Oxygen tubing, Mask, Oxygen reservoir, Child	
	RE-22615	Single Use Resuscitator 280 ml W/ 40 cmH2O POP-OFF, Oxygen tubing, Mask, Oxygen reservoir, Infant Single Use Resuscitator with Handstrap 280 ml W/ 40	
Product Group Name		Manual Resuscitator Sets	
Manufacturer		Health Business Corp.	
Device Name		Disposable Resuscitator Set	
Representative Photo			
Device type	Disposable		
Patient Usage type	Usò per un singolo paziente		
Accessories	Manometer, Mask, Oxygen reservoir, oxygen tube, Disposable PEEP valve, Intake valve, and connector		
Connector	15F/22mm		
Oxygen reservoir volume	2.5L		0.6L
Patient group	Adult	Child	InfantBambino (1-3 anni)
Resuscitator bag volume	1600 mL	500ml	280ml
POP-OFF	W/ POP-OFF		
Deliverable pressure	60 cmH ₂ O	40 cmH ₂ O	
Reservoir volume	2500ml		1000ml
Suitable patient body weight	>40 kg	11 - 40 kg	<10 kg
Stroke volume range; One hand	770ml	300ml	160ml
Stroke volume range; Two hands	900ml	350ml	190ml
Minimum delivered volume	>600 ml	>150 ml	
Maximum cycle rate	20 breaths/min		40 breaths/min
Oxygen concentration	With reservoir: > 85% Without reservoir: > 35%		
Dead Space	< 5,5 ml		
Material	Resuscitator body:	PVC	
	Resuscitator body:	Natural rubber	
	Mask:	PVC, and PP + TPR	
	Reservoir bag:	PVC	
	Oxygen tubing:	PVC	
	Disposable Pressure Manometer L Type:	PC, TPR	
	Non-rebreathing valve:	PC	
	Oxygen reservoir:	PE	

		Duckbill / gas intake valve:	PVC
		Pop-off valve:	PC, PVC, stainless steel
		PEEP valve:	PC, Silicone, (Stainless Steel for spring)
Package		Non-sterile package	
Latex/ DEHP content		No	
Operation temperature		25±5 °C, Room Temperature	
Storage temperature		da -40°C to 60°C	
Relative humidity		Da 30 - 70% UR	
Shelf Life		5 years	
Specification	Patient connection port and face mask connectors	Connection stayed intact	
	Dismantling and reassembling	Proper ventilation occurred after dismantling and reassembling the device	
	Patient valve function after contamination with vomitus	here is no change in the performance of the device after contamination with vomitus	
	Drop test	No function variation in the device after dropping	
	Immersion in water	There is no change in the performance of the device after immersion in water.	
	Expiratory resistance	2,07 cmH2O	
	Inspiratory resistance	3,03 cmH2O	
	Dead space and rebreathing	4,78ml	
	Minimum delivered volume	629,23ml	
	Pressure limitation	There is no performance deviation.	
	Supplementary and delivered oxygen concentration	The resuscitator without and with reservoir bag is more than 35 % and 85 %	
	Patient valve malfunction	1,07 cmH2O	
Attestation		Test item	Guideline
		Functional test	EN ISO 10651-4:2023
		Conical connector	EN ISO 5356-1:2015
		Cytotoxicity	EN ISO 10993-5:2009
		Sensitization	EN ISO 10993-10:2023
		Skin Irritation	EN ISO 10993-23:2021
		Emissions of particulate matter	EN ISO 18562-2
		Emissions of VOCs	EN ISO 18562-3
		Usability	IEC 62366-1
			PASS

2.2 Indications

Manual Resuscitator sets (Disposable) is a device intended to use as an adjunct to artificial respiration and cardiopulmonary resuscitation for emergent patient ventilation.

2.3 Intended User:

Medical personnel received the training of respiratory care.

2.4 Contraindications

This product does not apply to patients who do not suffer from illness or rescue need of emergent respiratory support.

2.5 Risk/Side Effects

An inappropriate inspiratory gas pressure can cause hazardous situation too. A high pressure may put too much stress to skin which can harm alveoli.

The amount of air/oxygen delivered should be very well controlled. The air entering the patient also has access to the stomach via the esophagus, which can inflate if the resuscitator is squeezed too hard (causing air flow that is too rapid for the lungs to absorb alone) or too much (causing excess air to divert to the stomach). This can lead to vomiting and subsequent aspiration of stomach contents into the lungs.

Holding the mask properly on the patient's face is very important while providing ventilation via resuscitators as pressure needed to force-inflate the lungs would be released to the environment. This is difficult when a single user attempts to maintain a face mask seal with one hand while squeezing the bag with other.

3. Instruction for Use:

1. Prior to use, read the Instructions, cautions, and warnings.
2. Check all the connection parts of the product, and tighten all the connections. Avoid leakage or loss when the product is used.
3. Choose the right size of the mask as per patient's face.
4. Connect the oxygen supply tubing to a regulated oxygen source.
5. Adjust the gas flow so that the reservoir expands completely during inspiration and nearly collapses as the squeeze bag refills during exhalation.
6. Prior to connection to a patient, check the function of the resuscitator, preferably attached to a test lung, by observing that the intake, reservoir, and patient valves are allowing all phases of the ventilator cycle to occur.
7. As required, connect the resuscitation mask to the patient connector.
8. Follow accepted Advance Cardiac Life Support (ACLS) or Institution-approved for ventilation.
9. Compress the squeeze bag to deliver a breath. Observe the chest rise to confirm inspiration.

10. Release pressure on the squeeze bag to allow exhalation. Observe the chest fall to confirm exhalation.
11. During ventilation, check for:
 - Signs of cyanosis
 - Adequacy of ventilation
 - Airway pressure
 - Proper function of all valves
 - Proper function of reservoir and oxygen tubing
12. Should the non-rebreathing valve become contaminated with vomiting, blood, or secretion during ventilation, disconnect the device from the patient and clear the non-rebreathing valve as follows:
 - Rapidly compress the squeeze bag to deliver several sharp breaths through the non-rebreathing valve to expel contaminate. If the contaminate does not clear.
 - Rinse the non-rebreathing valve in water and then rapidly compress the squeeze bag to deliver several sharp breaths through the non-rebreathing valve to expel contaminate. If the contamination still does not clear, discard the resuscitator.
13. When the product working, make sure the patient is closely monitored.

4. Warning/ Caution:

1. Read this instruction manual before use.
2. Do not use the resuscitator set in toxic atmospheres.
3. Do not use if the package is open or damage or contaminated.
4. This product is strictly restricted to single patient use, do not attempt to sterilize.
5. Do not administer supplemental oxygen in the presence of open flames.
6. This device is intended for use by qualified medical and emergency personnel trained in pulmonary ventilation and advanced cardiac life support techniques.
7. Ensure that the mask fits securely and set to prevent air leaks.
8. Proficiency in the assembly, disassembly, and use of this device should be demonstrated before use on a patient.
9. Always monitor airway pressure with a manometer when ventilating a patient.
10. Only qualified personnel trained in the use of Positive End Expiratory Pressure (PEEP) should level the function of the resuscitator before using it on a patient.
11. If overriding the pressure relief valve, great caution must be taken not to allow the pressure in the patient's airways to become too high.
12. Do not attempt to disassemble the pressure relief valve assembly. Disassembly will damage the component.
13. Any serious incident that has occurred in relation to the device should be reported to or its Authorized Representative/competent authorities and the Regulatory Authorities of the country.

5. Disposal:

Dispose of in accordance with local regulations procedures or hospital protocol.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies

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